

Skull Base and Endoscopic Pituitary Surgery






Accessories, High Frequency Cords,
for use with TAKE-APART® bipolar forceps

Bipolar High Frequency Cords

KARL STORZ
Instrument

High Frequency
Surgical Units



	26176 LE	Bipolar High Frequency Cord , for AUTOCON® II 400 SCB system (111, 113, 115, 122, 125), AUTOCON® II 200, AUTOCON® II 80, KARL STORZ Coagulator 26021 B/C/D, 860021 B/C/D, 27810 B/C/D, 28810 B/C/D, AUTOCON® system (50, 200, 350), Erbe-Coagulator, T and ICC series, length 300 cm
	26176 LM	Bipolar High Frequency Cord , for Martin and Berchtold coagulators, length 300 cm
	26176 L	Bipolar High Frequency Cord , with 2 x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Ellman, length 300 cm
	26176 LA	Bipolar High Frequency Cord , with 2 x 4 mm banana plug for KARL STORZ Coagulator 26020 XA/XB and Valleylab, length 300 cm
	26176 LV	Bipolar High Frequency Cord , for KARL STORZ AUTOCON® II 400 SCB system (112, 114, 116, 122, 125) and Valleylab coagulators, length 300 cm

High Frequency Surgical Units AUTOCON® see chapter 12