



Endothelial Punch

This mechanical device gives an excellent donor corneal button punched out from the endothelial side of corneoscleral button.

Features

- Universal - Any size of trephine from 6 mm to 11 mm can be fitted
- Precentered - The trephine and the teflon block are automatically centered for all sized trephines
- Spring Recoil - After completion of the cut, the trephine retracts back with inbuilt spring system
- Sterilization - Endothelial Punch can be sterilized by autoclaving or EO gas
- Each set consists of two plungers




Product Code	Description	Qty. per box
MIPL/N1	Endothelial Punch	1 pc.

Disposable Corneal Trephine

Nanoedge corneal trephines are ready-to-use circular blades designed for accurate cutting of corneal donor buttons to be used in penetrating keratoplasty.

NANOedge
Ophthalmic Blades



- An ultrasharp thin profile cutting edge minimizes damage to the endothelial cells
- Designed for a precise and symmetrical cut of the corneal tissue
- See through design helps in accurate positioning
- Can be used with or without an endothelial punch
- Innovative protective packaging to ensure safety of the cutting edge
- EO Sterile
- CE certified

Product Code	Description	Qty. per box (Sterile)
MIPL/C1/6.00	6.00 mm Diameter	1 pc.
MIPL/C1/6.50	6.50 mm Diameter	1 pc.
MIPL/C1/7.00	7.00 mm Diameter	1 pc.
MIPL/C1/7.25	7.25 mm Diameter	1 pc.
MIPL/C1/7.50	7.50 mm Diameter	1 pc.
MIPL/C1/7.75	7.75 mm Diameter	1 pc.
MIPL/C1/8.00	8.00 mm Diameter	1 pc.
MIPL/C1/8.25	8.25 mm Diameter	1 pc.
MIPL/C1/8.50	8.50 mm Diameter	1 pc.
MIPL/C1/8.75	8.75 mm Diameter	1 pc.
MIPL/C1/9.00	9.00 mm Diameter	1 pc.
MIPL/C1/9.50	9.50 mm Diameter	1 pc.
MIPL/C1/10.00	10.00 mm Diameter	1 pc.
MIPL/C1/10.50	10.50 mm Diameter	1 pc.
MIPL/C1/11.00	11.00 mm Diameter	1 pc.
MIPL/C1/12.00	12.00 mm Diameter	1 pc.
MIPL/C1/14.00	14.00 mm Diameter	1 pc.
MIPL/C1/16.00	16.00 mm Diameter	1 pc.
MIPL/C1/18.00	18.00 mm Diameter	1 pc.

Teflon Block

The corneo sclera tissue is placed over the Teflon block and the donor button then punched out.

- Aspheric curvature to match corneal curvature
- Central hole helps in accurate positioning
- Can be used with or without an endothelial punch
- EO Sterile
- CE certified

Product Code	Description	Qty. per pack
MIPL/D12	Teflon Block (30 mm)	10 pc. (Sterile)

