

minisono

Product Specification

V1.0 Rev03



General Specifications

miniSono L3-12 Physical Specification

- Width: 62.5 mm
- Depth: 150 mm
- Height: 25.5 mm
- Weight: Approx. < 200 g

miniSono C1-6 Physical Specification

- Width: 62.5 mm
- Depth: 157 mm
- Height: 25 mm
- Weight: Approx. < 200 g

Tablet PC Spec(Recommended)

- Windows 10
- Resolution: 1,736 x 1824
- chipset m3
- Integrated SSD
- Thermal printer storage
- Built-in stereo speaker
- System boot-up: 90 sec

Control Panel(Touch UI)

- Alphanumeric keyboard
- Ergonomic key layout
- 8 TGC control levers

System Standard Feature

Architecture

- Digital beamforming
- Multi-beam
- FleXcan™
- MicroFit™
- Optimal Imaging Suite™
- OS: Windows 10

Transducer Type

- Convex
- Linear

Application

- Abdomen
- Spine
- Musculoskeletal
- Anesthesia
- Small Parts
- FAST
- Vascular

Imaging Mode

- 2D (B)
- Color flow (CF)
- Pulse wave Doppler (PW)
- Motion (M)
- Power Doppler (PD)

Complex Mode

- Duplex: B+B, B+PW, B+CF, B+M, B+PD
- Triplex: B+CF+PW

Scanning Parameters

- System frequency range: 1 ~ 20 MHz
- Processing channel: 655,360
- ADC: 10 bits
- Displayed imaging depth: 0 ~ 30cm
- Max. Frame rate (Probe dependent)
 - 2D: 30 FPS
 - Color: 25 FPS
- Dynamic range: Max. 192dB
- Gray scale: 256 levels
- Focus
 - Focal number: Max.2
- Zoom
 - Write zoom : Max. 89.1x (probe dependent)
 - Read zoom : Max. 89.1x (probe dependent)
- CF
 - 300 Hz – 3.3 KHz
- PW
 - 300 Hz – 4.1 KHz

Imaging Presentation

- Single, dual
- Duplex image format
 - 1:1

Display and User Interface

Display Annotation

- Patient ID: Up to 64 characters
- Patient name: First, last and middle
- Gender, age and birth date
- Transducer name
- Institution / hospital name: up to 32 characters
- Date format: 3 types
 - YYYY/MM/DD
 - MM/DD/YYYY
 - DD/MM/YYYY
- Time format: 2 types
 - 12 hours
 - 24 hours

- Acoustic power
 - Mechanical index (MI)
 - Thermal index (TI): TIS, TIB, TIC
- Acoustic power
 - Mechanical index (MI)
- Transducer name
- Preset
- Transducer orientation
- Depth
- Focus
- TGC line
- Frequency
- Frame rate
- Power
- Mode name
- 2D mode
 - Gain
 - Dynamic range
 - Gray scale bar
 - Persistence
- Doppler mode
 - Wall filter
 - PRF
 - Sample volume position, size and angle
 - Velocity or frequency scale
 - Time
- Color flow mode
 - Wall filter
 - ROI box
 - PRF
 - Velocity or frequency scale
 - Color scale bar
- M mode
 - Gain
 - Dynamic range
 - Time
- Arrow
- Indicator
- Text
 - Size: Small, medium, large
- Body pattern (Marker)
- Thumbnail
- Stored image info
- Cine gauge
- Measurement result

Live Mode Display

- 2D
 - Angle steer (Probe dependent): 7 steps (-9 ~ 9°)
 - Dynamic range: 30 ~ 192dB
 - Frequency: 3 Steps
 - SRI: 0 ~ 5
 - Gray map: 0 ~ 14

- Gain: 0 ~ 50
- Harmonic
- Persistence: 4 steps
- Power: 1 ~ 100%
- Reject: 10 steps
- Spatial compound: 3 steps
- M
 - Dynamic range: 30 ~ 150dB
 - Gray map: 0 ~ 13
 - Gain: 0 ~ 50
 - Reject: 10 steps
 - Power: 1 ~ 100%
 - Sweep speed: 6 steps
- PWD
 - Angle: -89° ~ 89°
 - Angle steer (Probe dependent): 7 steps (-15 ~ 15°)
 - Auto angle
 - Base line: 16 steps
 - Frequency: 3 steps
 - Gain: 0 ~ 50
 - Power: 1 ~ 100%
 - Reject: 10 steps
 - Scale (PRF): 0.3 ~ 4.1 KHz
 - SV size: 0.7 ~ 20 mm
 - Sweep speed: 5 steps
 - Wall filter: 9 steps
- CF/PD
 - Angle steer (Probe dependent): 7 steps (-15 ~ 25°)
 - Color map: 9 maps
 - Frequency: 3 steps
 - Gain: 0 ~ 50
 - Power: 1 ~ 100%
 - Scale (PRF): 0.3 ~ 3.3 KHz

Post Processing (on frozen)

- 2D
 - Full SRI
 - TGC
 - Gain
 - Dynamic range
 - Gray map
 - Reject
- PWD
 - Angle correction
 - Auto Calc.
 - Gray map

Technology and Feature

Imaging Technology

- Filtered tissue harmonic imaging (FTHI)
- Full speckle reduction imaging (Full SRI)
- Spatial Compounding Imaging
- Frequency Compounding Imaging (FCI)
- Xpeed™
- Beam steer

Features

- Full SRI
- Xpeed™
- Spatial Compounding
- Frequency Compounding Imaging (FCI)
- Filtered Tissue Harmonic Imaging (FTHI)
- DICOM

Transducer

C1-6

- Convex array
- Frequency: 1.0 - 6.0 MHz
- Field of view: 78°
- Elements: 128
- Radius of curvature: 38.6 mm
- Footprint: 5.9 X 15.8 mm
- Application: Abdomen, Spine, FAST

L3-12

- Linear array
- Frequency: 3.0 - 12.0 MHz
- Field of view: 38.4 mm
- Elements: 128
- Footprint: 44.8 X 7.8 mm
- Application: Musculoskeletal, Anesthesia, Small Parts, FAST, Vascular

Data Management

Image Archive

- E-View™
- Image store: Up to System Storage (Ex: 216,000 images available in 385 GB)
- Image format
 - BMP, JPEG, DICOM, WMV
- Export / Backup / Restore
- Network storage

Cine

- Cine memory: 500 MB, 1,500 frames, Max 5min
- Prospective and retrospective
- Cine review, save
- frame move

Measurement and Calculation

Basic (Caliper)

- 2D
 - Distance
 - Area circumference
 - %Stenosis
 - Volume
 - Angle
 - A/B ratio
 - Disk Volume
 - Histogram
- M
 - Time
 - Slope
 - %Stenosis
 - HR
 - A/B Ratio
 - Distance
- Doppler
 - Velocity
 - Time
 - PI
 - Trace
 - Semi Auto trace
 - Acceleration
 - A/B Ratio
 - VTI
 - RI
 - AC/DC
 - PHT
 - HR- Volume Flow

Abdomen Measurement

- 2D
 - AO
 - Renal L
 - Renal Vol.
 - Renal
 - %Stenosis
 - Volume
 - A/B Ratio
- M
 - HR
- Doppler
 - AO
 - Renal A
 - IVC
 - Hepatic A
 - Hepatic V
 - Mid HV
 - MPV
 - PV
 - Splenic A

- Splenic V
- SMA
- SMV
- IMA
- Iliac A
- Iliac V
- Ves
- Renal V
- HR

MSK Measurement

- 2D
 - Volume
 - Hip (BA)
 - Hip (AB)
- M
 - HR
- Doppler
 - Ves
 - HR

Small Parts

- 2D
 - Thyroid Vol.
 - Testicle Vol.
 - Thyroid
 - Testicle
- M
 - HR
- Doppler
 - Ves
 - HR

Urology

- 2D
 - Renal L
 - Renal Vol.
 - Bladder Vol.
 - Testicle Vol.
 - Prostate Vol.
 - Renal
 - Bladder
 - Testicle
 - Prostate
 - Residual
- M
 - HR
- Doppler
 - Ves
 - HR

Vascular

- 2D
 - CCA
 - ICA
 - Bulb
 - ECA
 - Vert. A
 - Subclav A
 - Ophthalmic A
 - Supratroch A
 - Brachial A
 - Radial A
 - Ulnar A
 - CFA
 - SFA
 - PFA
 - PopA
 - PTA
 - ATA
 - Pero A
 - Arcuate A Sup
 - Arcuate A Med
 - Arcuate A Inf
 - Renal A
 - Kidney
 - Celiac A
 - SMA
 - IMA
 - CIA
 - EIA
 - AbdAo
 - Hepatic A
 - Splenic A
 - Portal Vein
- M
 - HR
- Doppler
 - CCA
 - IC
 - ECA
 - Bulb
 - Vert A
 - Subclav A
 - Ophthalmic A
 - Supratroch A
 - Brachial A
 - Radial A
 - Ulnar A
 - CFA
 - SFA
 - PFA
 - PopA
 - PTA

- ATA
- Pero A
- Arcuate A Sup
- Arcuate A Med
- Arcuate A Inf
- Renal A
- Kidney
- Celiac A
- SMA
- IMA
- AbdAo
- Hepatic A
- Splenic A
- Portal Vein
- CIA
- EIA
- MCA
- ACA
- PCA P1
- PCA P2
- AComA
- PComA
- BA
- Ves

OB

- 2D
 - AC
 - HC
 - FL
 - BPD
 - CLAV
 - CRL
 - GS
 - MSD
 - YS
 - BOD
 - Humerus
 - Radius
 - Ulna
 - Tibia
 - Fibula
 - SL
 - NB
 - EL
 - LL
 - RL
 - Scapular
 - SC
 - FB
 - Foot
 - FK
 - OFD

- OOD
- HA
- MAD
- AD
- AA
- FTA
- TAD
- TTD
- TC
- APTD
- TCD
- APAD
- CTAR
- AxT
- APTD-TTD
- CM
- UmbAD
- UmbVD
- NT
- OV Vol.
- AFI
- EFW
- OV
- Heart
- Teichholz
- CT
- Arteries
- Veins

- M
 - HR
 - Heart
 - Teichholz
- Doppler
 - AO
 - Umbilical A
 - Placenta
 - Des.Aorta
 - Uterine A
 - MCA
 - Ovarian A
 - Ves
 - DV
 - Valves
 - Ventricles
 - Arteries
 - Veins
 - HR

GYN

- 2D
 - Cervix
 - Endo
 - OV Vol.

- UT Vol.
- Bladder Vol.
- OV
- UT
- Follicle
- Bladder
- M
 - HR
- Doppler
 - AO
 - Umbilical A
 - Placenta
 - Des.Aorta
 - Uterine A
 - Mid Cerebral A
 - Ovarian A
 - Ves
 - HR

Connectivity

DICOM

- DICOM 3.0
- Image storage
- Multi-frame image storage
- Secondary image capture storage
- Structured report
- Grayscale print management
- Color print management
- Storage commitment push model
- Modality performed procedure step
- Worklist
- Verification
- Media storage: USB media, DVD, CD-R

Thermal printer

- BW
 - SONYUP-D23MD
 - SonyUP-D897
 - SonyUP-D898MD/X898MD
 - MITSUBISHICP30D/DW(USB)
- Color
 - SonyUP-D25MD
 - MITSUBISHIP95D
 - CanonSELPHYCP800

Safety and Compliance

- EN 60601-1:2006+A11:2011+A1:2013 (IEC 60601-1:2005+A1:2012)
- EN 60601-1-2:2007 (IEC 60601-1-2:2007)
- EN 60601-1-4:1996+A1:1999 (IEC 60601-1-4:1996+A1:1999)

- EN 60601-1-6:2010 (IEC 60601-1-6:2010+A1:2013)
- EN 60601-2-37: 2008+A11:2011 (IEC 60601-2-37:2007)
- EN 55022:2010, Class B (CISPR 22:2008, modified)
- EN 55011:2009+A1:2010, Group 1, Class B
- EN 61000-4-2:2009 (IEC 61000-4-2:2008)
- EN 61000-4-3:2006+A1:2008+A2:2010 (IEC 61000-4-3:2006+A1:2007+A2:2010)
- EN 61000-4-4:2004+A1:2010 (IEC 61000-4-4:2004+A1:2010)
- EN 61000-4-5:2006 (IEC 61000-4-5:2005)
- EN 61000-4-6:2009 (IEC 61000-4-6:2008)
- EN 61000-4-8:2010 (IEC 61000-4-8:2009)
- EN 61000-4-11:2004 (IEC 61000-4-11:2004)
- EN 61000-3-2:2006+A1:2009+A2:2009 (IEC 61000-3-2:2005+A1:2009+A2:2009)
- EN 61000-3-3:2008 (IEC 61000-3-3:2008)
- EN 62304:2006 (IEC 62304:2006)
- EN 62366:2008 (IEC 62366:2014)
- EN ISO 14971:2012 (ISO 14971:2007)
- EN ISO 10993-1:2009 (ISO 10993-1:2009)
- EN ISO 10993-5:2009 (ISO 10993-5:2009)
- EN ISO 10993-10:2013 (ISO 10993-10:2010)
- AIUM/NEMA UD 2:2004 (R2009)
- AIUM/NEMA UD 3:2004 (R2009)
- MDD 93/42/EEC
- CANADA [Regulation] SOR-98-282
- CE
- WEEE
- RoHS

© 2019 ALPINION MEDICAL SYSTEMS – All rights reserved.

2019 ALPINION MEDICAL SYSTEMS reserves the right to make changes in specifications and features shown hereby, or discontinue the product described at any time without notice or obligation. Contact you ALPINION representative or the most current information.

WWW.ALPINION.COM