


AlCath Flux G eXtra

3.5 mm Irrigated GoldTip Ablation Catheter



Ordering Information

Model	Tip	Curve	Spacing	Shaft Diameter	Length	Order Number
AlCath Flux Red G eXtra	3.5 mm irrigated	Red	2-5-2 mm	7 F	110 cm	376379
AlCath Flux Blue G eXtra	3.5 mm irrigated	Blue	2-5-2 mm	7 F	110 cm	376388
AlCath Flux Green G eXtra	3.5 mm irrigated	Green	2-5-2 mm	7 F	110 cm	376384
AlCath Flux Black G eXtra	3.5 mm irrigated	Black	2-5-2 mm	7 F	110 cm	376390
AlCath Flux Cyan G eXtra	3.5 mm irrigated	Cyan	2-5-2 mm	7 F	110 cm	376386

Product Highlights

Irrigated GoldTip ablation

GoldTip ablation results in superior acute efficacy expressed by shorter procedure times in the treatment of Atrial Fibrillation¹

eXtra Design

12 irrigation holes on 3 levels arranged with a 3D irrigation pattern

FullCircle Deflection

Deflectability beyond 180 degrees through FullCircle deflection²

¹ B. Schmidt et al. European Heart Journal Volume 34, Issue suppl 1, 1. August 2013;
DOI: 10.1093/eurheartj/ehz309.P4121.

² VAL11_AlCath_019_Dimensional Test_PG0_Rev-0
VAL11_AlCath_019_Dimensional Test_PG4_Rev-0

AlCath Flux G eXtra

Technical Data

Catheter Specification

Tip electrode length	3.5 mm
Tip electrode type	Irrigated
Tip temperature sensor	Thermocouple
Ring electrode length	1.55 mm
Number of electrodes	4
Electrode spacing	2-5-2 mm
Electrode material	Gold (Au)
Shaft diameter	7 F (2.3 mm)
Curve type	see below*
Curve shape	Steerable – FullCircle
Effective catheter length	110 cm
Connection plug	9-pole connector (REDEL), female, 60° coded Luer-lock (female)

Patient cable PK-109, 111, 112, 142, 143, 147, 150

Patient Cable	Ablation System	Order Number
PK-109	Stockert EP-Shuttle	330079
PK-111	Osypka HAT 300 Smart	330080
PK-112	Medtronic Atakr, Medtronic Atakr II	330081
PK-142	Abbott (SJM) IBI-1500 T7 & T11 Abbott Ampere	362442
PK-143	Boston Scientific Maestro 3000	371597
PK-147	BIOTRONIK Qubic RF, Stockert EP-Shuttle	398853
PK-150	Stockert SmartAblate	402668

Curve Type*	Reach R
Red	48 mm
Blue	57 mm
Green	65 mm
Black	73 mm
Cyan	80 mm

Transport and Storage Conditions

Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight

