

AKKUmed
medical equipment battery

REF 110498-0

LOT 2458

DANGER

- Follow the instructions below to prevent excessive heat generation, ignition, explosion, or leakage.
- Do not throw into fire, heat, disassemble, or modify.
 - Do not drop or apply strong impact.
 - Do not use on unspecified equipment.
 - Do not short-circuit the terminals.
 - Conduct charging on the equipment.

FUKUDA
Denshi



Ni-MH

Nickel-hydride battery pack



Type: 8PH-4/3A3700-H-J18
Voltage: 9.6V
Capacity: 3600mAh

Lot No.: 2012092

Recycle after use
Ni-MH
Beijing Fukuda Denshi
Medical Instruments Co., Ltd.

CE
MADE IN CHINA

(H)
2408

138891498-0

