

# 1. Features of QFT9000

### • Small and Portable

Similar footprint as A4 paper

Verified for ambulance on board use

## • Rapid and Simple

Throughput: ≥150T/h

Test time: 2-15 mins

2 test modes: rapid test + standard test

Swipe RFID card to load calibration curve and set up test


### Good Performance

hs-cTnI detection range: 0.02ng/mL~80.00ng/mL

Intra-batch CV ≤10%

Inter-batch CV ≤15%

Very good correlation against Beckman DxI800





# 2. Test Parameters of Reagent Kit

Product name	Sample type	Test time	Specification
Procalcitonin (PCT) Test Kit (Immunofluorescence Method)	Whole blood, serum, plasma	8 min	20 tests/kit
Whole-process C-Reactive Protein (hs-CRP + conventional	Whole blood, serum, plasma	2 min	20 tests/kit
CRP) Test Kit (Immunofluorescence Method)			
PCT/hs-CRP+conventional CRP Combo Test Kit	Whole blood, serum, plasma	8 min	20 tests/kit
(Immunofluorescence Method)			