

Adhesive Electrodes / EP

Adhesive Electrodes

- Excellent recording performance of the ECG signal.
- Hydrogel ensures excellent and constant adhesion to the patients skin.
- The electrode can be repositioned many times without reducing its adhesion.
- Hydrogel is very gentle so it can be applied on every type of skin, even the most sensitive ones.

Code	Electrode Diameter	Cable Length	Connectors	Box
CDE0201026 (-H)	20 mm	100 cm	 Ø 1,5 mm	40 pack/120 pcs
CDE0201526 (-H)	20 mm	150 cm	 Ø 1,5 mm	40 pack/120 pcs
CDE0241026 (-H)	24 mm	100 cm	 Ø 1,5 mm	40 pack/120 pcs
CDE0241526 (-H)	24 mm	150 cm	 Ø 1,5 mm	40 pack/120 pcs
CDER201026 (-H)	20 mm	100 cm	 Ø 1,5 mm	40 pack/120 pcs
CDER201526 (-H)	20 mm	150 cm	 Ø 1,5 mm	40 pack/120 pcs
CDER241026 (-H)	24 mm	100 cm	 Ø 1,5 mm	40 pack/120 pcs
CDER241526 (-H)	24 mm	150 cm	 Ø 1,5 mm	40 pack/120 pcs
CDES000024 (-H)	24 mm	/	Snap	20 pack/500 pcs

You will receive half the packing unit if you use **-H** in the code.

Tab Adhesive Electrodes

Code	Recording Area	Material	Box
TABSMP000014	22 x 14 mm	Ag/AgCl	25 pack / 2500 pcs
TABSM0000025	22 x 20 mm	Ag/AgCl	25 pack / 2500 pcs
TABSMP000014-H	22 x 14 mm	Ag/AgCl	10 pack / 1000 pcs
TABSM0000025-H	22 x 20 mm	Ag/AgCl	10 pack / 1000 pcs
TABSMP00K100	22 x 14 mm	Ag/AgCl	1 pack / 100 pcs
TABSM000K100	22 x 20 mm	Ag/AgCl	1 pack / 100 pcs