
ラララ DAHLHAUSEN®



Trocar Catheter

The Trocar Catheter is introduced into the pleural cavity after puncture incision. The light metal mandrin is removed and the remaining catheter is connected to the thorax drainage unit.

The Trokar catheter is graduated in steps of 1 cm starting beneath the last lateral eye, for a safe evaluation of the position. The Trocar catheter has a closed atraumatic tip which reduces injuries.

Sizes 10, 12 and 16 are equipped with a funnel-shaped connector.

X-ray contrastable:

X-ray contrast thread

Scale value:

1 cm

Sterile:

Yes

Outer diameter	Length	Number (eyes)
10 Fg	24 cm	. 2
	24 cm	2
4	26 cm	3,
_	40 cm	3
	Outer diameter 10 Fg 12 Fg 16 Fg 20 Fg 24 Fg 28 Fg 32 Fg	10 Fg 24 cm 12 Fg 24 cm 16 Fg 26 cm 20 Fg 40 cm 24 Fg 40 cm 28 Fg 40 cm