

# Disposable Endo Linear Cutter Stapler and Cartridge



- |                    |                     |
|--------------------|---------------------|
| 1. Connector       | 7. Safety           |
| 2. Joint Arrow     | 8. Trigger          |
| 3. Lock Button     | 9. Anvil            |
| 4. Steering Button | 10. Cartridge       |
| 5. Rotate Button   | 11. Closure Cover   |
| 6. Reset Button    | 12. Joint Arrow     |
|                    | 13. Joint Connector |

Stapler No.	Cartridge No.		Cartridge Color	Staple Effective Length (mm)	Staple Height (mm)
	Straight Type	Rotatable Type			
	JQZ-30-2.5	JQD-30-2.5	White	30	2.5
	JQZ-45-2.5	JQD-45-2.5	White	45	2.5
	JQZ-45-3.5	JQD-45-3.5	Blue	45	3.5
JQ-A(Standard)	JQZ-45-4.0	JQD-45-4.0	Golden	45	4.0
JQ-B(Short)	JQZ-45-4.8	JQD-45-4.8	Green	45	4.8
<b>JQ-L(Long)</b>	JQZ-60-2.5	JQD-60-2.5	White	60	2.5
	JQZ-60-3.5	JQD-60-3.5	Blue	60	3.5
	JQZ-60-4.0	JQD-60-4.0	Golden	60	4.0
	JQZ-60-4.8	JQD-60-4.8	Green	60	4.8

- ❖ Standard/short type of endocutter device
- ❖ Straight/rotatable type of cartridge
- ❖ Four staple heights available
  
- ❖ **Scope of Applications:**  
Suitable for tissue cutting under open surgery or endoscopic surgery include lungs, stomach, and intestines
  
- ❖ **Contraindication:**  
Gastrointestinal submucosal edema, over thick muscular layer, poor healing capability, cancerous tissue residue

