



freeflex+ Swabbable Luer Lock Port

- + Luer Lock-compatible injection port prevents needle stick injuries and minimizes risk of bag damage, leakages and contamination of patients or staff
- + Swabbable, reusable injection port ensures safe and convenient disinfection between multiple drug additions
- + Self-sealing septums in both ports help prevent leakages
- + Port accommodates a variety of male Luer Lock syringes



Easy Push & Turn Attachment

- + Only half a turn of a Luer Lock-equipped transfer adaptor or syringe is needed to reach the end position for a safe and tight connection with the bag
- + Minimal force is required to push through the self-sealing septum with a syringe or vial adaptor

Together, we are your ally for infusion delivery solutions.



- **IV Solutions**

A growing portfolio of IV solutions in **free/flex** bags.

- **Pumps**

A portfolio that performs everywhere, from acute care to home care.

- **Sets**

Accessories that support safety and efficiency in every application.

- **Software**

A suite that aligns everyone, from Pharmacy to BioMed.

- **Services**

A team that delivers every day, from initial implementation to ongoing support.

Product Code	NDC Code	Fill Size	Units per Case
232150	65219-218-50	50	60
232200	65219-220-10	100	50
232225	65219-222-25	250	30
232250	65219-224-50	500	20

For more information or to place an order, contact your Sales Representative or call Customer Service at 1.888.386.1300 | www.fresenius-kabi.com/us