

# TECHNICAL DATA SHEET

## AUTOCLAVABLE BIOHAZARD BAG

<b>PRODUCT CODE</b>	:	<b>OT</b>
<b>RELATED STANDARDS</b>	:	Autoclavable Biohazard Bag meets the requirements of the Council Directive 94/62/EC on Packaging.
<b>INTENDED USE</b>	:	Autoclavable Biohazard Bags are used for the sterilization of medical wastes before disposal.
<b>OPERATION CONDITIONS</b>	:	Autoclavable Biohazard Bags are designed for use in steam sterilization methods. The recommended operation temperature is max. 141 °C. Other sterilization conditions should be determined by the end user regarding to material to be sterilized. During sterilization, bag should be left open for an effective sterilization.

### PRODUCT SPECIFICATIONS

PROPERTIES	METHOD	UNIT	VALUE	
THICKNESS (one side)	ASTM D 374	micron	50 ± 5	
TENSILE STRENGTH	ASTM D 882	N/15 mm	MD	22,4
			CD	14,4
ELONGATION AT BREAK	ASTM D 882	%	MD	640
			CD	510

**PACKAGING** : Autoclavable Biohazard Bags are packaged as the listed below.

<u>Product Code</u>	<u>Dimensions</u>	<u>Unit</u>
OT5370	53 cm x 70 cm	300 pcs/box
OT6389	63 cm x 89 cm	200 pcs/box
OT87110	87 cm x 110 cm	120 pcs/box



**PMS Steri-test**

**PMS MEDİKAL Amb. San. Tic Ltd. Şti**  
**pms@pmsmedikal.com**