

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

## Bacterial/Viral Filter, with HME

Bacterial/Viral Filter with HME and CO<sub>2</sub>-port. The purpose of the filter is to protect the patient and to humidify the inhaled air. The filter is placed between the endotracheal tube/tracheal cannula and the Y-piece of the breathing system.

**Shape (Ventilation Filter):** Round

**ISO style (patient end):** Male/Female

**Inner diameter (patient end):** 20/15 mm

**Sterile:** No

**Country availability:** Please note: Not all products are available in

[read more](#)

Article no.	Application area	Inner diameter (machine end)	ISO style (machine end)	Max. tidal volume (Breathing Filter)	minimal quantity
4300906120	Child	13	Male	900	50

<b>Article no.</b>	<b>Application area</b>	<b>Inner diameter (machine end)</b>	<b>ISO style (machine end)</b>	<b>Max. tidal volume (Breathing Filter)</b>	<b>minimal quantity</b>
4300906020	Adult	22/13	Male/Female	1.500	25
4300905708	Adult	22/13	Male/Female	1.500	50