

3D PEMF Therapy Machine

Lumina 3D-e Clinic - With Battery



Advantages and purpose of the PEMF machine

- Option to run two different programs simultaneously
- Time and intensity settings for each program
- Memory for programs + rename option: 200 positions
- Unique programs of the new 3D-e generation

- **Lumina Clinic 3D-e**
Technical parameters

- Technology
- **3D-e patented Biomag® technology**
second generation of 3D pulsed magnetic therapy
- Programs
- **6 innovative 3D-e programs**
with therapeutic effect for a wide range of symptoms
- Display
- **Colour touch**
- Adjustable application time

- **5-90 minutes**
- Outputs
- **4 outputs**
completely independent channels for the latest Biomag® applicators
- Intensity
- **35 mT**
for each output
- Intensity setting
- **1-100%**
- Independent programmes
- Possibility to run up to 2 independent programs
- Renaming programmes
- **200 positions**
for renaming programs
- Frequency
- **1-81 Hz**
for low-frequency magnetic therapy applications
- Repetition of the programme
- **Yes**
auto-repeat can be started
- **Universal**
functional worldwide
- Operation on built-in battery
- **Yes**

3D PEMF therapy

SL60-P - PEMF Therapy Ring

Advantages and parameters of the ring applicator:

- Flat bottom for convenient applications.
- Deep action also on the torso area.
- Diameter 600 mm; depth 340 mm.
- MIMI 1.5 mT.
- Connector 1 x JACK 3.5 mm.
- Weight: 12.65 kg.



3D PEMF therapy

SL-30P - PEMF Therapy Ring Applicator



Advantages and parameters of the ring applicator:

- In-depth applications to the limbs.
- A long time favourite ring.
- Diameter 300 mm; depth 170 mm.
- Connector 1x JACK 3.5 mm.
- MIMI 5,0 mT.
- Weight: 2,55 kg

3D PEMF therapy

A8P - PEMF Therapy Mat Applicator



Advantages and parameters of the mat applicator:

- Deep and surface applications.
- Universal mat applicator for different parts of the body.
- Length 440 mm; width 390 mm; height 40 mm.
- Connector 1x JACK 3.5 mm.
- MIMI 2.5 mT.
- Weight 0,72 kg.

• -----

• -----

3D PEMF therapy

A6P2 - PEMF Therapy Applicator



Advantages and parameters of the applicator:

- Coil switch for localized or extended pulses.
- Intensive local application.
- Length 170 mm; width 130 mm; height 23 mm.
- Connector 1x JACK 3.5 mm.
- MIMI 35,0 mT - SPOT / MIMI 20,0 mT - WIDE.
- Weight: 0,58 kg.
- -----

3D PEMF therapy

A12PM - PEMF Therapy Mat Applicator



of the mat applicator:

- Can be applied to larger parts of the body.
- Universal mat applicator.
- Length 1 400 mm; width 450 mm; height 40 mm.
- Patented 3D technology.
- Possibility of joining into a closed shape or fixation.
- MIMI 3.5 mT.
- Connector 3 x JACK 3.5 mm.
- Weight: 2.60 kg.