

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

## Trocar Catheter

The trocar catheter is inserted into the pleural gap after skin incision and puncture of the pleural cavity. The light metal stylet is removed and the remaining catheter is connected to the connecting tube/stepped connector of a thoracic drainage system.

The sizes CH 10, 12 and 16 have a funnel-shaped connector.

The sizes CH 20, 24, 28 and 32 have a plain connector.

**X-ray contrastable:** X-ray contrast thread

**Sterile:** Yes

[read more](#)

| Article no. | Outer diameter | Length  | Number (eyes) | PZN      | minimal quantity |
|-------------|----------------|---------|---------------|----------|------------------|
| 0200030010  | 10 Fg          | 24 cm   | 2             | 10033868 | 10               |
| 0200030012  | 12 Fg          | 24 cm   | 2             | 10033874 | 10               |
| 0200030016  | 16 Fg          | 26,5 cm | 3             | 10033880 | 10               |
| 0200030020  | 20 Fg          | 40 cm   | 3             | 10033897 | 10               |

| Article no. | Outer diameter | Length | Number (eyes) | PZN      | minimal quantity |
|-------------|----------------|--------|---------------|----------|------------------|
| 0200030024  | 24 Fg          | 40 cm  | 3             | 10033905 | 10               |
| 0200030028  | 28 Fg          | 40 cm  | 3             | 10033911 | 10               |
| 0200030032  | 32 Fg          | 40 cm  | 3             | 10033928 | 10               |