

# AVANTI™ + Sheath Introducer

For atraumatic insertion

Smooth performance  
from start to finish

---




**Cordis**™

# AVANTI™+ Sheath Introducer

## For atraumatic insertion

The proven cannula offers reliable, smooth performance from start to finish



- For smooth, atraumatic insertion
- Kink resistant
- Optimized catheter maneuverability with minimal leakage

### For smooth, atraumatic insertion


- Smooth transitions between cannula, guidewire and dilator for:
  - Effortless vessel entry
  - Increased patient comfort
  - Smallest possible risk for complications

### Kink Resistant

Dual layer material gives AVANTI™+ Sheath Introducer excellent flexibility and the strength needed to prevent kinking - even in tortuous vessels and scar tissue

### Optimized catheter maneuverability with minimal leakage

The patented SLIX haemostatic valve with six-cut spiral design allows free catheter movement without blood reflux or air aspiration



# AVANTI™ + Sheath Introducer

## For atraumatic insertion

### Smooth Performance from Start to Finish

Available with or without mini-guidewire  
(mid-length without wire only)

Available in three lengths  
(usable cannula length)

- 5.5cm **brachial**
- 7.5cm **short length**
- 11cm **standard length**
- 23cm **mid-length**



#### CARDIOLOGY: AVANTI™ + Sheath Introducer

French size	Colour code	5.5cm Cannula		11cm Cannula		23cm Cannula
		With mini-wire	Without mini-wire	With mini-wire	Without mini-wire	Without mini-wire
4	■	402-604P		402-604X	402-604A	
5	■	402-605P		402-605X	402-605A	402-605T
5.5**	■			402-655X		
6	■	402-606P	402-606R	402-606X	402-606A	402-606T
6.5**	■			402-656X		
7	■	402-607P	402-607R	402-607X	402-607A	
8	■			402-608X	402-608A	402-608T
9	■			402-609X	402-609A	402-609T
10	■			402-610X	402-610A	
11	■			402-611X		

\*\*













5F MS = Mid-Size = 5.5F

6F MS = Mid-Size = 6.5F

# AVANTI™+ Sheath Introducer

## For atraumatic insertion

### CARDIOLOGY: AVANTI™+ Sheath Introducer

11 cm Cannula				
French size	Colour code	Mini-wire (inch)	With mini-wire	Without mini-wire
4		.035	504-604X	504-604A
5		.038	504-605X	504-605A
5.5**		.038	504-655X*	
6		.038	504-606X	504-606A
6.5**		.038	504-656X*	504-656A*
7		.038	504-607X	504-607A
7.5**		.038	504-657X*	
8		.038	504-608X	504-608A
8.5**		.038	504-658X*	
9		.038	504-609X	504-609A
10		.038	504-610X	504-610A
11		.038	504-611X	


\* Monitoring sheath: provides slightly larger ID for pressure monitoring or side port aspiration with a catheter in place; they are one-half French size larger in size

\*\* 5F MS = Mid-Size = 5.5F









6F MS = Mid-Size = 6.5F

7F MS = Mid-Size = 7.5F

8F MS = Mid-Size = 8.5F



### CARDIOLOGY: AVANTI™+ Sheath Introducer

23 cm Cannula				
French size	Colour code	Dilator usable (cm)	With obturator	Without obturator
4		30		504-604T
5		30		504-605T
6		30	504-606D	504-606T
6.5**		30		504-656T*
7		30	504-607D	504-607T
8		30	504-608D	504-608T
10		30		504-610T
11		30		504-611T

\*\* 6F MS = Mid-Size = 6.5F




# AVANTI™ + Sheath Introducer

## For atraumatic insertion


### CARDIOLOGY: AVANTI™ + Brachial Sheath Introducer with Mini-Guidewire

5.5 cm Cannula			
French size	Colour code	Mini-wire (inch)	Catalogue number
4	■	.035	504-604P
5	■	.038	504-605P
6	■	.038	504-606P




### CARDIOLOGY: AVANTI™ + Short length Sheath Introducer with mini-guidewire

7.5 cm Cannula			
French size	Colour code	Mini-wire (inch)	Catalogue number
4	■	.021	504-604S
5	■	.021	504-605S
6	■	.021	504-606S



### CARDIOLOGY: AVANTI™ + Transradial Sheath Introducer kit with Mini-Guidewire

		11 cm Cannula		23 cm Cannula	
French size	Colour code	Wire acceptance (inch)	Catalogue number	Wire acceptance (inch)	Catalogue number
4	■	.021	504-614Z	.021 70 cm	504-624Z
5	■	.021	504-615Z	.021 70 cm	504-625Z
6	■	.021	504-616Z	.021 70 cm	504-626Z
7	■	.021	504-617Z		



For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO and AVANTI are trademarks of Cordis and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. © 2024 Cordis. All Rights Reserved. 100545515-5 02/2024

