



COYOTE™ & COYOTE™ ES Balloon Dilatation Catheter Monorail™



COYOTE™ & COYOTE™ ES Balloon Dilatation Catheter Over-the-Wire



Product Matrix

Balloon Diameter	Balloon Length								
	20 mm	30 mm	40 mm	60 mm	80 mm	100 mm	120 mm	150 mm	220 mm
1.5 mm									
2.0 mm									
2.5 mm									
3.0 mm									
3.5 mm									
4.0 mm									

■ COYOTE™ ■ COYOTE™ ES

- Available in balloon lengths up to 220 mm to treat long BTK disease with fewer inflations
- 47 balloon sizes, 156 codes
- OTW and Monorail™ Platforms

20 mm

220 mm

Over-The-Wire Balloon Dilatation Catheter  
COYOTE™ & COYOTE™ ES

Diameter (mm)	Length (mm)	Rated Burst Pressure (ATM)*	Shaft Length		
			90 cm	150 cm	142-145 cm
1.5	20	14	-	-	H74939134152010
1.5	40	14	H74939186150490	H74939186150410	-
2.0	20	14	-	-	H74939134202010
2.0	30	14	-	-	H74939134203010
2.0	40	14	-	-	H74939134204010
2.0	60	14	H74939186200690	H74939186200610	-
2.0	80	14	H74939186200890	H74939186200810	-
2.0	100	14	H74939186201090	H74939186201010	-
2.0	120	14	H74939186201290	H74939186201210	-
2.0	150	14	H74939186201590	H74939186201510	-
2.0	220	14	H74939186202290	H74939186202210	-
2.5	20	14	-	-	H74939134252010
2.5	30	14	-	-	H74939134253010
2.5	40	14	-	-	H74939134254010
2.5	60	14	H74939186250690	H74939186250610	-
2.5	80	14	H74939186250890	H74939186250810	-
2.5	100	14	H74939186251090	H74939186251010	-
2.5	120	14	H74939186251290	H74939186251210	-
2.5	150	14	H74939186251590	H74939186251510	-
2.5	220	14	H74939186252290	H74939186252210	-
3.0	20	14	-	-	H74939134302010
3.0	30	14	-	-	H74939134303010
3.0	40	14	-	-	H74939134304010
3.0	60	14	H74939186300690	H74939186300610	-
3.0	80	14	H74939186300890	H74939186300810	-
3.0	100	14	H74939186301090	H74939186301010	-
3.0	120	14	H74939186301290	H74939186301210	-
3.0	150	14	H74939186301590	H74939186301510	-
3.0	220	14	H74939186302290	H74939186302210	-
3.5	20	12	-	-	H74939134352010
3.5	30	12	-	-	H74939134353010
3.5	40	12	-	-	H74939134354010
3.5	60	14	H74939186350690	H74939186350610	-
3.5	80	14	H74939186350890	H74939186350810	-
3.5	100	14	H74939186351090	H74939186351010	-
3.5	120	14	H74939186351290	H74939186351210	-
3.5	150	14	H74939186351590	H74939186351510	-
3.5	220	14	H74939186352290	H74939186352210	-
4.0	20	12	-	-	H74939134402010
4.0	30	12	-	-	H74939134403010
4.0	40	12	-	-	H74939134404010
4.0	60	14	H74939186400690	H74939186400610	-
4.0	80	14	H74939186400890	H74939186400810	-
4.0	100	14	H74939186401090	H74939186401010	-
4.0	120	14	H74939186401290	H74939186401210	-
4.0	150	14	H74939186401590	H74939186401510	-
4.0	220	14	H74939186402290	H74939186402210	-

\* 14ATM= 1418kPa, 12ATM = 1216kPa

Monorail™ Balloon Dilatation Catheter  
COYOTE™ & COYOTE™ ES

Diameter (mm)	Length (mm)	Rated Burst Pressure (ATM)*	Shaft Length		
			90 cm	150 cm	143-146 cm
1.5	20	14	-	-	H74939135152010
1.5	40	14	H74939185150490	H74939185150410	-
2.0	20	14	-	-	H74939135202010
2.0	30	14	-	-	H74939135203010
2.0	40	14	-	-	H74939135204010
2.0	60	14	H74939185200690	H74939185200610	-
2.0	80	14	H74939185200890	H74939185200810	-
2.0	100	14	H74939185201090	H74939185201010	-
2.0	120	14	H74939185201290	H74939185201210	-
2.0	150	14	H74939185201590	H74939185201510	-
2.0	220	14	H74939185202290	H74939185202210	-
2.5	20	14	-	-	H74939135252010
2.5	30	14	-	-	H74939135253010
2.5	40	14	-	-	H74939135254010
2.5	60	14	H74939185250690	H74939185250610	-
2.5	80	14	H74939185250890	H74939185250810	-
2.5	100	14	H74939185251090	H74939185251010	-
2.5	120	14	H74939185251290	H74939185251210	-
2.5	150	14	H74939185251590	H74939185251510	-
2.5	220	14	H74939185252290	H74939185252210	-
3.0	20	14	-	-	H74939135302010
3.0	30	14	-	-	H74939135303010
3.0	40	14	-	-	H74939135304010
3.0	60	14	H74939185300690	H74939185300610	-
3.0	80	14	H74939185300890	H74939185300810	-
3.0	100	14	H74939185301090	H74939185301010	-
3.0	120	14	H74939185301290	H74939185301210	-
3.0	150	14	H74939185301590	H74939185301510	-
3.0	220	14	H74939185302290	H74939185302210	-
3.5	20	12	-	-	H74939135352010
3.5	30	12	-	-	H74939135353010
3.5	40	12	-	-	H74939135354010
3.5	60	14	H74939185350690	H74939185350610	-
3.5	80	14	H74939185350890	H74939185350810	-
3.5	100	14	H74939185351090	H74939185351010	-
3.5	120	14	H74939185351290	H74939185351210	-
3.5	150	14	H74939185351590	H74939185351510	-
3.5	220	14	H74939185352290	H74939185352210	-
4.0	20	12	-	-	H74939135402010
4.0	30	12	-	-	H74939135403010
4.0	40	12	-	-	H74939135404010
4.0	60	14	H74939185400690	H74939185400610	-
4.0	80	14	H74939185400890	H74939185400810	-
4.0	100	14	H74939185401090	H74939185401010	-
4.0	120	14	H74939185401290	H74939185401210	-
4.0	150	14	H74939185401590	H74939185401510	-
4.0	220	14	H74939185402290	H74939185402210	-

Boston Scientific  
Advancing science for life™

www.bostonscientific-international.com

© 2015 Boston Scientific Corporation or its affiliates. All rights reserved. DINPER4658EA


All cited trademarks are the property of their respective owners. Savvy Long and Sleek are registered trademarks of Cordis Corporation. Amphiron Deep is a registered trademark of Medtronic Inc. Ultraverse is a registered trademark of CR Bard inc. Armada is a registered trademark of Abbott Laboratories. Nanocross is a registered trademark of ev3. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Illustrations are property of Boston Scientific Corporation. Illustrations are not necessarily to scale.

Information not intended for use or distribution in France.

PI-19303-AB MAR2015 Printed in Germany by medicalvision.

COYOTE™ & COYOTE™ ES PTA Balloon Dilatation Catheter  
GET TO. GET THROUGH. GET GOING.



You'll want to get your hands on  
**COYOTE™** Balloon Dilatation Catheter.

**COYOTE™** Balloon Dilatation Catheter puts the numbers in your favour.

Introducing **COYOTE**, the go-to below-the-knee (BTK) peripheral balloon. Confidently take on BTK procedures with the ultimate in deliverability, crossability, and ultra-low profile.




**Get to**  
 With four catheter configurations optimized for maximum push and excellent track



**Get through**  
 With an ultra-low lesion entry profile (0.017"/0.43 mm) and best-in-class crossing profile (0.030"/0.76 mm)



**Get going**  
 With deflation times of less than 10 seconds\* and balloon lengths up to 220 mm



Balloon dilatation of the dorsalis pedis artery carried out using 2,5 x 220 mm COYOTE™ Over-the-Wire (OTW) balloon Dilatation Catheter at 14 ATM (1419kpa)




Balloon dilatation of the posterior tibial artery carried out using 3,0 x 220 mm COYOTE™ Over-the-Wire (OTW) balloon Dilatation Catheter at 14 ATM (1419kpa)



Balloon dilatation of the anterior tibial artery carried out using 3,0 x 220 mm COYOTE™ Over-the-Wire (OTW) balloon Dilatation Catheter at 14 ATM (1419 kpa)

**<10**  
 seconds\*  
 Best-in-class deflation time


**COYOTE Catheter offers rapid deflation time**



**4x**

Available in four configurations: Over-the-Wire and Monorail™ platforms on 90 cm and 150 cm shaft lengths

**COYOTE Catheter is optimized for maximum push and excellent track\***



**220 mm**

Available in lengths up to 220 mm

**0.017"**  
 (0.43 mm)\*

Ultra-low lesion entry profile

**COYOTE Catheter offers ultra-low 0.017" (0.43 mm) lesion entry profile with balloon lengths up to 220 mm\***

