


RT[®]ED
SUTURE



山东海迪科医用制品有限公司

SHANDONG HAIDIKE MEDICAL PRODUCTS CO.,LTD.

Add: Tianfu Road, Dongcheng District, Shan County, 274300

Heze City, Shandong Province, China.

Tel: +86-530-4660062

Mobile: +86-13105409308

E-mail: info@suturescn.com

Web: www.suturescn.com



SHANDONG HAIDIKE MEDICAL PRODUCTS CO.,LTD.

山东海迪科医用制品有限公司

Shandong Haidike Medical Products Co., Ltd. is a leading and world famous medical devices manufacturer with ISO13485, ISO14001, ANVISA GMP, CE, FDA, TGA, FSC and NMPA certificates. "RTMED" and "Haidike" are registered trademark and protected under the Madrid Agreement. Our medical devices have been registered in more than 20 countries and exported to more than 150 countries and owning good reputation all over the world!

- ◆ Absorbable and nonabsorbable Surgical sutures
- ◆ Lifting threads for plastic surgery
- ◆ Hernia repair meshes
- ◆ Orthopedic casting tapes & splints
- ◆ Electric Wheelchair & Manual Wheelchair
- ◆ Veterinary surgical supplies
- ◆ Medical protection products etc.



ENTERPRISE BUSINESS LICENSE



MEDICAL DEVICE MANUFACTURING PERMIT



GOLD PLUS SUPPLIER SGS ASSESSMENT CERTIFICATE



FSC BY EUROPEAN UNION



CE CERTIFICATE



ISO13485 CERTIFICATE



ISO14001 CERTIFICATE



SURGICAL SUTURES 510K



BGMP





FREE SALE CERTIFICATE



RTMED MADRID TRADEMARK



Needle Shapes



Advantages

- 1 Strong hardness
- 2 Outstanding sharpness
- 3 Excellent tenacity and ductility

- Round bodied (heavy)
- ▲ Curved cutting (heavy)
- ▼ Reverse cutting (heavy)
- ⊗ Taper cutting
- ▼ Microfine spatula
- ⊙ Blunt point
- Straight round bodied
- ▲ Straight curved cutting



PGA POLYGLYCOLIC ACID

Properties:

Structure: Multifilament, Braided

Colour: Violet, Undyed

Chemical composition: Polyglycolic acid

Coating: Polycaprolactone + Calcium Stearate

Sizes: USP 8/0 - USP 4


Type of absorption: Hydrolysis

Sterilization: Ethylene Oxide

Absorption: 60 - 90 days

Advantages:

- High initial tensile strength
- Excellent knotting ability
- Reliable absorbability
- Easy to handle
- Outstanding holding power through the critical wound periode
- Smooth surface, slides smoothly through the tissue



PGLA POLYGLACTIN 910

Properties:

Structure: Multifilament, Braided

Colour: Violet, Undyed

Chemical composition: Poly(glycolide-co-l-lactide 90/10)

Coating: Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate

Sizes: USP 8/0 - USP 2


Type of absorption: Hydrolysis

Sterilization: Ethylene Oxide

Absorption: 56 - 70 days

Advantages:

- High tensile strength
- Smooth tissue passage
- Easy knotting and tie down
- Excellent knot security
- Reliable absorbability



PDO POLYDIOXANONE

Properties:

Structure: Monofilament

Colour: Violet, Undyed


Chemical composition: Poly (p-dioxanone)

Sizes: USP 7/0 - USP 2

Type of absorption: Hydrolysis

Sterilization: Ethylene Oxide

Absorption: 180 - 210 days



Advantages:

- High initial tensile strength
- Prolonged wound support
- Minimal tissue reaction
- Lower incidence of infection and trauma due to smooth surface

PGCL POLIGLECAPRONE 25

Properties:

Structure: Monofilament

Colour: Violet, Undyed


Chemical composition: Poly (glycolide-co-caprolactone)

Sizes: USP 6/0 - USP 2

Type of absorption: Hydrolysis

Sterilization: Ethylene Oxide

Absorption: 90 - 120 days



Advantages:

- High initial tensile strength
- Easy handling
- Smooth passage through tissue
- Excellent knotting ability


CHROMIC CATGUT

Properties:

Colour: Brown
Composition: Purified connective tissue
Sizes: USP 7/0 - USP 4
Type of absorption: Proteolysis
Sterilization: Gamma radiation
Absorption: 90 - 100 days

Advantages:

- High tensile strength
- Excellent handling properties
- Minimal tissue reaction
- Uniform absorption




PLAIN CATGUT

Properties:

Colour: Beige
Composition: Purified connective tissue
Sizes: USP 7/0 - USP 4
Type of absorption: Proteolysis
Sterilization: Gamma radiation
Absorption: 60 - 70 days

Advantages:

- High tensile strength
- Smooth tissue passage
- Easy knotting and tie down
- Excellent knot security
- Reliable absorbability




SILK

Properties:

Structure: Braided
Colour: Black, White
Composition: Natural silk filaments
Coating: Silicon / wax
Sizes: USP 8/0 - USP 4
Sterilization: Ethylene Oxide

Advantages:

- Excellent handling properties
- High tensile strength
- Highly smooth and flexible
- Good knot security




POLYPROPYLENE

Properties:

Structure: Monofilament
Colour: Blue
Chemical composition: Polypropylene
Sizes: USP 7/0 - USP 3
Sterilization: Ethylene Oxide

Advantages:

- High tensile strength
- Smooth tissue passage
- Easy and secure knotting
- Minimal acute inflammatory reaction in tissues
- Excellent and permanent tissue support



NYLON

Properties:

Structure: Monofilament

Colour: Blue, Black


Chemical composition: polyamide 6 - 6.6

Sizes: USP 12/0 - USP 4

Sterilization: Ethylene Oxide

Advantages:

- High tensile strength
- Smooth tissue passage
- Easy and secure knotting
- Excellent and permanent tissue support



POLYESTER

Properties:

Structure: Braided

Colour: Green, White

Chemical composition: Polyethylene terephthalate


Coating: Silicon

Sizes: USP 6/0 - USP 4

Sterilization: Ethylene Oxide

Advantages:

- High tensile strength
- Excellent handling properties
- Minimal tissue reaction
- Good knot security



Ultra High Molecular Weight Polyethylene (UHMWPE)

Properties:

Structure: Braided

Colour: White, white+blue, white+black, white+green


Chemical composition: Polyethylene

Sizes: USP 7/0 - USP 5

Sterilization: Ethylene Oxide

Advantages:

- Extremely strong, flexible and smooth
- Outstanding tensile strength
- Minimal inflammatory reaction



PTFE

Properties:

Structure: Monofilament

Colour: White

Chemical composition: polytetrafluoroethylene

Sizes: USP 2/0 - 6/0

Sterilization: Ethylene Oxide

Advantages:

- High tensile strength; does not weaken when moistened.
- Decreased inflammation and faster wound repair.
- Soft, comfortable for patients.
- Great handling properties and knot security.



THREADS		NEEDLES	
THREADS TYPE		THREADS DIAMETER	THREADS LENGTH
ABSORBABLE	POLYGLYCOLIC ACID (size: 5-8/0)	PGA	
	POLYGLYCOLIC ACID RAPID (size: 5-8/0)	PGAR	
	POLYGLYCOLIC ACID PLUS (size: 5-8/0)	PGA+	
	POLYGLACTIN 910 (size: 2-8/0)	PGLA	
	POLYGLYGLACTIN 910 RAPID (size: 2-8/0)	PGLAR	
	POLYGLYGLACTIN 910 PLUS (size: 2-8/0)	PGLA+	
	POLYDIOXANONE (size: 2-7/0)	PDO	
	POLYDIOXANONE RAPID (size: 2-7/0)	PDOR	
	POLYDIOXANONE PLUS (size: 2-7/0)	PDO+	
	CHROMIC CATGUT (size: 4-7/0)	CC	
	PLAIN CATGUT (size: 4-7/0)	PC	
	POLIGLECAPRONE 25 (size: 2-6/0)	PGCL	
	POLIGLECAPRONE 25 RAPID (size: 2-6/0)	PGCLR	
	POLIGLECAPRONE 25 PLUS (size: 2-6/0)	PGCL+	
NON-ABSORBABLE	SILK (size: 4-6/0)	SK	
	NYLON (size: 4-12/0)	NL	
	POLYPROPYLENE (size: 3-7/0)	PP	
	POLYESTER (size: 4-7/0)	PB	
	Ultra High Molecular Weight Polyethylene (UHMWPE) (size: 5-7/0)	UHMWPE	
	PTFE (size: 2/0-6/0)	PTFE	

PGA

20

F

01

-2

PGA20F0126-2

Polyglycolic Acid 2/0 75mm 1/2 circle round bodied double needles 26mm

Disposable Skin Stapler and Remover

Application:
Used for wound closure and tissue suture

- Features and Advantages:**
- Single use, avoid cross-infection
 - Biologically compatible
 - Significantly decrease wound closure time and reduce procedure & anaesthesia time
 - Suture easily and efficiency as well as easily removed
 - With better cosmetic results.



Disposable Trocars

Application:
It is used for laparoscopy surgery,the puncture of the abdominal wall of the human body is used to establish a working channel for abdominal surgery.

- Features and Advantages:**
- Unique air-tight design to prevent gas leakage.
 - Anti-slip grain design, increasing the friction between cannular and abdominal wall, to avoid cannular fall off.
 - Universal positioning of puncture cone without directional insertion.



Disposable Endo Cutter Stapler and The Cartridge

Application:
Suitable for endoscopic or open surgical tissue removal in abdominal surgery, gynecology, pediatrics, thoracic surgery.

- Features and Advantages:**
- One-handed control ensures complete tissue cutting and reduces tissue damage.
 - Better tissue thickness adaptability.
 - Larger rotation angle for easier operation.




SPECIALTY SUTURES PORTFOLIO

Dental sutures




Optimal performance in DENTAL SURGERY




Ophthalmic sutures



NON-ABSORBABLE SUTURE
(11-0, 10-0, 9-0, 8-0, 7-0)



Hernia Repair Mesh



Large mesh aperture

Promotes liquid circulation and smoothes tissue growth.

Less micropores


Greatly reduces the possibility of bacteria hiding and is not easy to infect.

Light weight


Reduce postoperative pain and foreign body sensation.

Products Description	Size
Polypropylene Mesh	6x11cm
Polypropylene Mesh	7.6x15cm
Polypropylene Mesh	10x15cm
Polypropylene Mesh	15x15cm
Polypropylene Mesh	30x30cm

Hernia Example



After Mesh Repair





Cassette Sutures

Length	25m 50m 75m 100m
Absorbable	Chromic catgut, Plain catgut PGA, PGLA, PDO
Non-absorbable	Silk, Nylon , Polyester, Polypropylene



Suture Materials

Threads	Sizes
Silk	USP4-6/0
Nylon	USP4-8/0
Polyester	USP5-6/0
Polypropylene	USP2-6/0



Self-adhesive Elastic Bandage

Standard	Features
5.0cm×450cm	Self-adhesive
7.5cm×450cm	High elasticity
10cm×450cm	Waterproof
15cm×450cm	Easy to tear

I.V. Cannula with Wings



specification

22G 24G 26G



Dosing mouth shape





Pen type



Butterfly type

Veterinary catheter



Lifting Threads Needle

SHARP

Basic type, generally used for linear filling and simple lifting, Beginners can also operate effortlessly.

W TYPE

Circular head, reducing the damage to the skin tissue. After the injection, the needle can be taken out then correct the position and depth of the needle, which improves the lack of a stitch of the other types.

L TYPE

L type blunt needle, reduce the damage of blood vessels, bruising and hematoma. The pain was significantly reduced.

U-COG

Double needle is used for eyebrow lifting and other parts requiring strong lifting. U-COG is a special structure line with strong lifting ability. The local lifting effect is almost as good as the little facelift.

PDO/PCL Lifting Thread

Mono

- Needle** Sharp Needle
- Parts** Face, Body
- Effects** Tightening, Filling, Wrinkles Removing, Whitening

Tornado Screw

- Needle** Sharp Needle
- Parts** Eye, Neck, Crow's Feet, Cheek
- Effects** Tightening, Filling, Whitening, Wrinkles Removing

Nasal COG

- Needle** L type Blunt Needle
- Parts** Nasal columella, Nasal bridge
- Effects** Rhinoplasty, Nose Shaping

COG Lifting

- Needle** L/W type Blunt Needle
- Parts** Lifting Eyebrow, Apple Zone, Cheek, Chest
- Effects** Face Lifting, V Face, Chest Lifting

Multi Threads

- Needle** L type Blunt Needle
- Parts** Face, Neck, Eye Bag, Crow's Feet
- Effects** Tightening, Filling, Whitening, Wrinkles Removing

Molding COG

- Needle** W type Blunt Needle
- Parts** Face, Neck, Crow's Feet, Eye Bag
- Effects** Tightening, Filling, Whitening, Wrinkles Removing



Forehead

Multi Threads / Tornado Screw / Mono Screw

Crow's feet / Eye bags / Eye Corner's Lifting

Blunt Needle For Eye 30G-25mm

Apple Zone

Multi Threads / Tornado Screw / Mono Screw

Nasolabial folds / Mouth Corner's Lifting

Multi Threads / Tornado Screw / Mono Screw



Face Lifting

18G-100mm / 19G-100mm

Nose Dorsum / Nose Bridge


19G-60mm

Nose Columella

19G-38mm

Cervical Stripe / Neck Lifting

23G-38mm / 18G-100mm



Orthopedic Casting Tape



USE INSTRUCTION

- Apply 1~2 layers stockinet over the affected body part;
- Wear surgical gloves;
- Choose proper size orthopedic casting tape, then open the package;
- Immerse the casting tape in normal temperature(20℃~30℃) water for 5~6 seconds and squeeze 2~3 times. Then remove from the water and squeeze out the excess water. The curing time of orthopedic casting tape is directly related to the temperature of water. The curing time is shorter in high temperature water than in low temperature water. If the water temperature is over 27℃, the curing time is too short to increase the difficulty of the operation.
- Wrap the casting tape spirally around the affected body part, overlapping the previous layer by one-half or two-thirds the width of the tape. The tension is moderate. If it is too loose, it will weaken the adhesion of the orthopedic casting tape. If it is too tight, it will affect the blood circulation.
- The intensity of the casting tape is related to the number of overlapped layers. Overlap 3~4 layers for the non-weight bearing parts, while overlap 5~6 layers for the weight bearing parts. If more layers are needed, please ensure the adhesivity between the layers.
- To increase the adhesion between layers, the casting tape will be shaped after wrapped. The operating time should be in 3~5 minutes. The adhesion and shaping effect cannot reach the best if the operating time is longer or the immersed time is too long. The patients do not allow to move freely in 10 minutes before the orthopedic casting tape are completely hardened.
- The electric plaster saw is needed if the affected parts need to open a window or remove the orthopedic casting tape.

Orthopedic Casting Tape				
Product Description	Sizes	Specification(W×L)	Packaging	Operating Part
Fiberglass Orthopedic Casting Tape	RT002	5.0cm × 360cm	10bags/box 12boxes/ctn	Wrists/Childrens
Fiberglass Orthopedic Casting Tape	RT003	7.5cm × 360cm	10bags/box 12boxes/ctn	Arms
Fiberglass Orthopedic Casting Tape	RT004	10.0cm × 360cm	10bags/box 12boxes/ctn	Arms / Legs
Fiberglass Orthopedic Casting Tape	RT005	12.5cm × 360cm	10bags/box 9boxes/ctn	Legs
Fiberglass Orthopedic Casting Tape	RT006	15.0cm × 360cm	10bags/box 9boxes/ctn	Legs

Polymer Splint



Wrist Fixation
Recommended sizes:
312 (7.5cmX30cm)
415 (10cmX40cm)




Thumb Fixation
Recommended sizes:
312 (7.5cmX30cm)
415 (10cmX40cm)



The Fourth and Fifth Finger Fixation
Recommended sizes:
312 (7.5cmX30cm)
415 (10cmX40cm)



Elbow Fixation
Recommended sizes:
312 (7.5cmX30cm)
415 (10cmX40cm)



Forearm Fixation
Recommended sizes:
335 (7.5cmX90cm)
430 (10cmX75cm)



Horseshoe Fixation
Recommended sizes:
335 (7.5cmX90cm)
430 (10cmX75cm)



Ankle Fixation
Recommended sizes:
430 (10cmX75cm)
530 (12.5cmX75cm)



Legs Fixation
Recommended sizes:
335 (7.5cmX90cm)
545 (12.5cmX115cm)

Polymer splint				
Product Description	Sizes	Specification(WxL)	Packaging	Operating part
Medical Polymer Splint	312	7.50cmX30cm	20bags/box 6boxes/ctn	wrists/Childrens
Medical Polymer Splint	335	7.50cmX90cm	10bags/box 6boxes/ctn	arms/legs
Medical Polymer Splint	415	10.0cmX40cm	15bags/box 6boxes/ctn	arms/legs
Medical Polymer Splint	430	10.0cmX75cm	10bags/box 6boxes/ctn	legs
Medical Polymer Splint	530	12.5cmX75cm	10bags/box 6boxes/ctn	legs
Medical Polymer Splint	545	12.5cmX115cm	5bags/box 6boxes/ctn	legs
Medical Polymer Splint	630	15.0cmX75cm	10bags/box 6boxes/ctn	legs
Medical Polymer Splint	645	15.0cmX115cm	5bags/box 6boxes/ctn	legs

Powered Wheelchair



Front wheel size: 10 inches	Net weight: 47kg (actually prevail)
Rear wheel size: 16 inches	Maximum speed: 6km/h
Dimensions: 1180×630×1360mm	Theoretical travel: 15km
Safe load: 100kg	Motor (motor): DC24V250W*2pcs
Frame material: Carbon steel	Gradeability (gradient): 6 degrees
Battery: Standard	



Front wheel size: 8 inches	Net weight: 21.5kg (actually prevail)
Rear wheel size: 12 inches	Maximum speed: 6km/h
Dimensions: 990×590×950mm	Theoretical travel: 15km
Safe load: 100kg	Motor (motor): DC24V250W*2pcs
Frame material: Aluminum alloy	Gradeability (gradient): 6 degrees
Battery: 12AH/20AH lithium battery	Folding size: 765×590×350mm

Manual Wheelchair



Front wheel size: 7 inches	Net weight: 18kg (actually prevail)
Rear wheel size: 24 inches	Armrest: blow molding/leather
Dimensions: 950×680×900mm	Seat width: 460mm
Safe load: 100kg	Seat height: 520mm
Frame material: Aluminum alloy	
Tire: solid tire	



Front wheel size: 6 inches	Net weight: 7.2kg (actually prevail)
Rear wheel size: 8 inches	Armrest: blow molding/leather
Dimensions: 1000×500×910mm	Seat width: 340mm
Safe load: 100kg	Seat height: 490mm
Frame material: Aluminum alloy	Backrest height: 410mm
Folding size: 220×500×700mm	

Seniors Shopping Cart

A good helper for parents in life



Tie-rod brake design



Comfortable seat



Wear-resistant reinforced brake rear wheel



New Generation Multi-claw Crutches

Material upgrade



Smooth all the way

- Anti-stumbling design
- Steady walk without hindrance
- The inside and outside width design is balanced and stable.



Feel comfortable

- High wear-resistant handrail
- Surface cracked anti-skid treatment, anti-skid and wear-resistant, Ergonomic design, Fits to the hand without getting tired from holding for a long time.



Firmer and stronger

- 1.2mm wall thickness Thicken reinforcement
- Lightweight and load-bearing



Firm and convenient

- Annular clavicle clamp
- Fixed main pole Three simple steps



Medical Grade Crutches

Stable and non-slip height adjustable



PE environmentally friendly foam material arm rest

Thick PE foam material, high resilience, non-slip, soft and comfortable, strong abrasion resistance, not tired after long-term use.



Leg height adjustment in 9 positions

The height of 9 levels is adjustable to meet the needs of people of different heights. Strengthen the aluminum alloy to provide super good support.




Grip three-position adjustment

The grip has three gears adjustment, the adjustment part is locked firmly and stably, and it is convenient to use and feel comfortable and soft



Non-slip durable rubber pad

The rubber pad is made of rubber, which is thicker, stable and non-slip, and the built-in metal gasket is wear-resistant and durable



Nursing Bed

Bent leg angle: -45°~25°

Back angle: 0°~80°

Rollover angle: 0°~55°

Toilet hole size: 280×215mm


Product Size: 2100×900×550mm

Bed size: 1970×900mm


Bed material: Thickened steel pipe

For people: Elderly, paralyzed patients, etc

Product accessories: Wash basin, infusion stand, mattress, dining table, guardrail, castor, bedpan



Medical Disposable Protective Clothing



Head



Neck



Waist

Executive Standard GB 19082-2009



Specifications: One-piece Sizes:

□ 160 □ 165 □ 170 □ 175 □ 180 □ 185

Intended use: It is used to provide barrier and protection for the potentially infectious patients' blood, body fluids, secretions, etc. that clinical medical personnel come into contact with at work. It cannot be used to isolate places with strict microbiological index control requirements such as intensive care unit (room).

Expiration date:

2 years after production

Surgical Gowns



Collar tape



Waist



Knit elastic cuffs

CE BS EN 13795-1:2019



Model: □S □M □L □XL □XXL □XXXL □XXXXL

Specifications:

□High performance type □Standard type

Intended use: For the surgical staff to wear, reduce the source of infection to the patient's surgical wound to prevent postoperative wound infection.

Storage:

Store in an environment with relative humidity not exceeding 80%, no corrosive gas, cool, dry, well ventilated and clean.

Expiration date:

3 years after production

Isolation Gown



- Good air permeability
- Strong barrier performance
- Non-sterile
- Do not reuse

Isolation Gown TYPE A/TYPE B



Isolation gowns - type A/type B

Specifications:

120cm×120cm/One-size-fits-all

Material:100% Non-woven fabric


Product description:

Non-woven fabrics are usually used as the main raw materials, which are made by cutting and sewing. Non-sterile, single-use

Expiration date:

3 years after production

Medical Scrubs Uniforms



Model: type A: split type; type B: gown type.



Specifications:

□S □M □L □XL □XXL

Product Description: Usually made of cotton fiber or non-woven fabric. The medical scrubs uniforms is short-sleeve or long-sleeve style that have the function of blocking dandruff, and do not have liquid barrier properties. It is non-sterile supply and can be reused. It should be sterilized before using.

Expiration date:

3 years after production

Medical Goggle



Model specifications

□White □Black □Transparent
□Large □Medium □Small

Product description

It is composed of a protective cover made of polymer material, a foam strip and a fixing device. Non-sterile, single use.

Expiration date

3 years after production

Storage

This product should be stored in a dry, ventilated, non-corrosive gas warehouse.

Medical Isolation Shoe Cover



size:

□Long □Short

style:

□Non-woven fabric □PVC □PP+PE

Medical Cap



Model specifications

Round cap. φ25.5cm FZ/T64005-2011

Product description

Non-woven fabric is usually used as the main raw material, which is made by cutting and sewing. Non-sterile, single-use.

Material

□Non-woven fabric □PP+PE

Colour

Blue

Storage

This product should be stored in a dry, ventilated, non-corrosive gas warehouse.

Medical Isolation Mask



Model: Ear hook type / Headband type

EN166:2001
EN167:2002



Product description

It is composed of a protective cover made of polymer material, a foam strip and a fixing device. Non-sterile, single use.

Expiration date

3 years after production

Storage

This product should be stored in a dry, ventilated, non-corrosive gas warehouse.

Medical Bed Sheet



Materials: 70% Polypropylene, 30% Polyethylene non-woven fabric

Product description

It is made by compounding or sewing non-woven fabric and plastic film. Non-sterile, single use.

Expiration date

3 years after production

Storage

This product should be stored in a dry, ventilated, non-corrosive gas warehouse.

Disposable Virus Specimen Collection Tube



Disposable sampler

Model: Inactivated and non-inactivated

Product description

Disposable virus specimen collection tube is composed of swabs and/or collection cups and tubes containing storage solution. Non-sterile.

Expiration date

12 months

Storage

The product should be stored in a non-corrosive gas and well-ventilated warehouse at room temperature.

Medical examination gloves



Model/Specification: According to the different materials used, it is divided into: polyvinyl chloride, butyronitrile; according to the size, it is divided into: XS, S, M, L, XL.

Product Description: Usually made of materials such as polyvinyl chloride. Have sufficient strength and barrier properties. Non-sterile supply, single use.

Intendad use: It will be wore by doctor to do checking for the patients or lab testing.

Medical rubber gloves



Model/Specification: XS, S, M, L, XL.

Product Description: Usually made of materials such as polyvinyl chloride and rubber. Have sufficient strength and barrier properties. Non-sterile supply, single use.

Intendad use: It will be wore by doctor to do checking for the patients or lab testing.

Disposable Surgical Face Mask



3 Ply Flat type	Ear Loop	EN 14683:2019	Type IIR	NON STERILE
Size	17.5 cm x 9.5 cm			
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, metal wire with plastic covering nose clip, polyester/nylon spandex blend ear loops			
Product performance:	The bacterial filtration efficiency of the mask should be not less than 98%.			
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.			
Expiration date	2 years after production			
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.			



3 Ply Flat type	Ear Loop	EN 14683:2019	Type IIR	NON STERILE
Size	14.5 cm x 9.5 cm			
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, metal wire with plastic covering nose clip, polyester/nylon spandex blend ear loops			
Product performance:	The bacterial filtration efficiency of the mask should be not less than 98%.			
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.			
Expiration date	2 years after production			
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.			

Disposable Medical Face Mask



3 Ply Flat type	Ear Loop	EN 14683:2019	Type I	NON STERILE
Size	17.5 cm x 9.5 cm			
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, metal wire with plastic covering nose clip, polyester/nylon spandex blend ear loops			
Product performance:	The bacterial filtration efficiency of the mask should be not less than 95%.			
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.			
Expiration date	2 years after production			
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.			



3 Ply Flat type	Ear Loop	EN 14683:2019	Type IIR	NON STERILE
Size	17.5 cm x 9.5 cm			
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, metal wire with plastic covering nose clip, polyester/nylon spandex blend ear loops			
Product performance:	The bacterial filtration efficiency of the mask should be not less than 98%.			
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.			
Expiration date	2 years after production			
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.			



3 Ply Flat type	Ear Loop	EN 14683:2019	Type I	NON STERILE
Size	14.5 cm x 9.5 cm			
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, metal wire with plastic covering nose clip, polyester/nylon spandex blend ear loops			
Product performance:	The bacterial filtration efficiency of the mask should be not less than 95%.			
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.			
Expiration date	2 years after production			
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.			




5 Ply Folding type	EN 14683:2019	Type IIR	NON STERILE
Size	11.5 cm x 15 cm		
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, nose clip, polyester/nylon spandex blend ear loops		
Product performance:	The bacterial filtration efficiency of the mask should be not less than 98%.		
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.		
Expiration date	2 years after production		
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.		



3 Ply Cup type	Head Loop	EN14683:2019	Type IIR	NON STERILE
Size	13.5 cm x 13.5 cm			
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, nose clip, polyester/nylon spandex blend head loops			
Product performance:	The bacterial filtration efficiency of the mask should be not less than 98%.			
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.			
Expiration date	2 years after production			
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.			



5 Ply Folding type	Head Loop	N95-V1	STERILE EO
Size	11.5 cm x 15 cm		
Material	Polypropylene spunbond nonwoven fabric, polypropylene melt-blown nonwoven, nose clip, polyester/nylon spandex blend head loops		
Product performance:	The bacterial filtration efficiency of the mask should be not less than 95%.		
Application	It is worn by clinical staff during non-invasive procedures, covering the user's mouth, nose, and jaw, providing a certain physical barrier to prevent direct penetration of pathogenic microorganisms, particulate matter, in patients, etc.		
Expiration date	11 months after production		
Storage	Store at room temperature, the relative humidity should not exceed 80%, non-corrosive gas and well-ventilated room. Transportation means should be clean and isolated from fire sources.		

**Tissue Scissors****Stitch Scissors****Umbilical Scissors****Hemostatic Forceps****Tissue Forceps****Needle Forceps****Edison Forceps****Dissecting Forceps****Dressing Forceps****Surgical Blade****Scalpel 3#****Scalpel 4#****Scalpel 7#****Bone Hook****Skin Retractor****Surgery Suction Pipe**



Infusion set



Blood transfusion set



Disposable blood bag



Syringes

