

 MEDICI



CYSTOSCOPY



ENDOMED





Cystoscopy



ENDOMED



|                                                                                |             |
|--------------------------------------------------------------------------------|-------------|
| <b>Urology Scopes</b>                                                          | <b>page</b> |
| Cystoscopes                                                                    | 06          |
| Fiber Light Cables                                                             | 07          |
| <b>Cistoscopy-Urethrotomy</b>                                                  | <b>page</b> |
| Sheaths                                                                        | 08-09       |
| Working Sheaths - Bridges - Albarran - Sealing Cup                             | 10          |
| Laser Sheaths                                                                  | 11          |
| Monopolar Flexible Electrodes - Monopolar cables / adapter                     | 12          |
| Bipolar Flexible Electrodes, Bipolar cables / adapters                         | 13          |
| Flexible/Semirigid Forceps-Scissors, WL 400mm                                  | 14-19       |
| Optical Forceps                                                                | 20          |
| Flexible Stone Extractors                                                      | 21          |
| <b>Lithotripsy</b>                                                             | <b>page</b> |
| Optical Lithotripsy Forceps/Sheaths                                            | 22-23       |
| MAUERMAYER/Optical Lithotripsy Forceps                                         | 24-25       |
| <b>Monopolar Resectoscopy</b>                                                  | <b>page</b> |
| Working Elements                                                               | 27          |
| Sheaths                                                                        | 28-29       |
| Obturators - Cables / adapters                                                 | 30          |
| HF Electrodes / Cold Preparation Knives                                        | 31          |
| Slim Working Elements - 19Ch Slim HF Electrodes                                | 32          |
| Slim Sheaths                                                                   | 33          |
| <b>Laser Resectoscopy</b>                                                      | <b>page</b> |
| Laser Working Elements - Optical forceps for tissue removing after enucleation | 34          |
| Optical Forceps for tissue/sheath                                              | 35          |
| <b>Bipolar Resectoscopy</b>                                                    | <b>page</b> |
| Bipolar Working Elements                                                       | 37          |
| Sheaths                                                                        | 38-39       |
| Obturators - Cables                                                            | 40          |
| 24Ch Bipolar HF Electrodes / Cold Preparation Knives                           | 41          |
| Mini Working Elements                                                          | 42          |
| Mini Sheaths - Obturators / 17Ch Mini HF Electrodes - Cables                   | 43          |

# Cystoscopy

## Content

### **Urethrotomy** *page*

|                                  |    |
|----------------------------------|----|
| <i>Working Elements</i>          | 44 |
| <i>Sheaths - Working bridges</i> | 45 |
| <i>Obturators</i>                | 46 |
| <i>Cold Knives</i>               | 46 |
| <i>OTIS Urethrotome</i>          | 47 |
| <i>Evacuators - Accessories</i>  | 48 |

### **Nephroscopy** *page*

|                                     |    |
|-------------------------------------|----|
| <i>Nephroscopes</i>                 | 49 |
| <i>Nephroscopes Forceps</i>         | 50 |
| <i>Nephroscope Sheaths</i>          | 51 |
| <i>Mini Nefroscope complete SET</i> | 52 |
| <i>Mini Nefroscope - Forceps</i>    | 53 |

### **Uretherorenoscopy** *page*

|                                                    |         |
|----------------------------------------------------|---------|
| <i>Urethero-renoscopes - Lithotripter Adaptors</i> | 54 - 55 |
| <i>Forceps and Scissors</i>                        | 56 - 61 |

### **Accessories** *page*


|                             |    |
|-----------------------------|----|
| <i>Catheter introducers</i> | 62 |
| <i>Bougies Dilatings</i>    | 63 |
| <i>Metal Catheters</i>      | 63 |

### **Pediatric Urology** *page*

|                                  |    |
|----------------------------------|----|
| <i>Endoscopes</i>                | 64 |
| <i>Cystoscope complete SET</i>   | 65 |
| <i>Urethrotome Complete SET</i>  | 66 |
| <i>Resectoscope complete SET</i> | 67 |

### Cystoscopes

High-Definition, wide angle  
Autoclavable, max.134°C/2,3bar/5min



Medici Endoscope-line consists of high quality telescopes which have been developed according to the most modern technology and production standards. A new generation optic system assures a unique color, contrast and highly defined details reproduction, both in the center of the image and at the edges. A particularly true color reproduction and a markable depth of focus guarantee a natural picture. Distal and proximal sapphire windows provide an extremely high level of scratches resistance.

#### Diameter 2,9mm

| Ordering n°  | Ø     | Viewing Angle | Working length |
|--------------|-------|---------------|----------------|
| 40 302 LA.HD | 2,9mm | 0°            | 302mm          |
| 40 302 LB.HD | 2,9mm | 12°           | 302mm          |
| 40 302 LE.HD | 2,9mm | 30°           | 302mm          |

#### Diameter 4,0mm

| Ordering n°  | Ø     | Viewing Angle | Working length |
|--------------|-------|---------------|----------------|
| 40 302 NA.HD | 4,0mm | 0°            | 302mm          |
| 40 302 NB.HD | 4,0mm | 12°           | 302mm          |
| 40 302 ND.HD | 4,0mm | 25°           | 302mm          |
| 40 302 NE.HD | 4,0mm | 30°           | 302mm          |
| 40 302 NG.HD | 4,0mm | 70°           | 302mm          |



Olympus  
Cystoscope Connection



Wolf Connection



Olympus  
Arthroscop Connection



Stryker Connection

# Cystoscopy

## Fiber Light Cables

### Fiber Light Cables, Xenon Safe, Led White Light

Medici Optical Fiber Cable is made to resist and to shield against very high thermic energies especially in the use of **Xenon** light sources due to the fusion of light fiber in the cable itself.

This guarantees **20%** light more in the operating field.



**New High Performances White Line-LED Light Cables**  
New «WhiteLine» LED series Optical Fibers Light Cables have been designed and made up to be used with cutting edge LED Light Sources. By enhancing the trasmission effectiveness of an already well-established and first-choice optical fibers light cable, light output and illumination, when operating with a LED light source togheter with a modern HD telescope, can be visibly and considerably increased.

### 3,5mm diameter , xenon safe, for Ø 2,9 - 4,0mm scopes

| Xenon    | LED       | Working length | Xenon    | LED       | Working length |
|----------|-----------|----------------|----------|-----------|----------------|
| 49 180 X | 49 180 XL | 1800mm         | 49 300 X | 49 300 XL | 3000mm         |
| 49 230 X | 49 230 XL | 2300mm         | 49 400 X | 49 400 XL | 4000mm         |
| 49 250 X | 49 250 XL | 2500mm         | 49 500 X | 49 500 XL | 5000mm         |

### 4,9mm diameter , xenon safe, for Ø 10,00mm scopes

| Xenon    | LED       | Working length | Xenon    | LED       | Working length |
|----------|-----------|----------------|----------|-----------|----------------|
| 49 180 P | 49 180 PL | 1800mm         | 49 300 P | 49 300 PL | 3000mm         |
| 49 230 P | 49 230 PL | 2300mm         | 49 400 P | 49 400 PL | 4000mm         |
| 49 250 P | 49 250 PL | 2500mm         | 49 500 P | 49 500 PL | 5000mm         |

### Adaptors for optic fiber light cables

| Ordering n° | Projector end                                                                       | Description  | Ordering n° | Instrument end                                                                        |
|-------------|-------------------------------------------------------------------------------------|--------------|-------------|---------------------------------------------------------------------------------------|
| 49 110 A    |  | Medici/Storz | 49 010 A    |  |
| 49 110 B    |  | Acmi/Stryker | 49 010 B    |  |
| 49 110 C    |  | Wolf/Dyonics | 49 010 C    |  |
| 49 110 D    |  | Olympus/Acme |             |                                                                                       |
| 49 110 F    |  | Olympus new  | 49 010 F    |  |

# ENDOMED

# Cystoscopy

## Cysto-Urethroscope Sheaths



Cysto-Urethroscope sheaths for adults, with 2 fixed stopcocks,  
for Ø 2,9/4,0 - w/302mm Scopes, continuous flow.

| System  | Colour | Ø     | Sheath only | Obturator only | Sheath QL | Obturator QL | Endoscope |
|---------|--------|-------|-------------|----------------|-----------|--------------|-----------|
| 4104109 | Black  | 15 Ch | 4104009     | 4114009        |           |              | Ø 2,9 mm  |
| 4110108 | Black  | 16 Ch |             |                | 4101008   | 4111008      | Ø 2,9 mm  |
| 4103110 | Yellow | 17 Ch | 4103010     | 4113010        |           |              | Ø 4,0 mm  |
| 4110110 | Yellow | 17 Ch |             |                | 4101010   | 4111010      | Ø 4,0 mm  |
| 4103111 | Green  | 19 Ch | 4103011     | 4113011        |           |              | Ø 4,0 mm  |
| 4110111 | Green  | 19 Ch |             |                | 4101011   | 4111011      | Ø 4,0 mm  |
| 4103112 | Red    | 21 Ch | 4103012     | 4113012        |           |              | Ø 4,0 mm  |
| 4110112 | Red    | 21 Ch |             |                | 4101012   | 4111012      | Ø 4,0 mm  |
| 4103113 | Blue   | 23 Ch | 4103013     | 4113013        |           |              | Ø 4,0 mm  |
| 4110113 | Blue   | 23 Ch |             |                | 4101013   | 4111013      | Ø 4,0 mm  |
| 4103114 | White  | 25 Ch | 4103014     | 4113014        |           |              | Ø 4,0 mm  |
| 4110114 | White  | 25 Ch |             |                | 4101014   | 4111014      | Ø 4,0 mm  |



Visual Obturator

| Ordering n° | Ordering QL n° | Colour | Ø     |
|-------------|----------------|--------|-------|
| 4114039     | 4111039        | Black  | 15 Ch |
| 4113040     | 4111040        | Yellow | 17 Ch |
| 4113041     | 4111041        | Green  | 19 Ch |
| 4113042     | 4111042        | Red    | 21 Ch |
| 4113043     | 4111043        | Blue   | 23 Ch |
| 4113044     | 4111044        | White  | 25 Ch |



| Ordering n° | Description                                |
|-------------|--------------------------------------------|
| 4133060     | Cystoscope adapter for female urethroscope |

# ENDOMED



# Cystoscopy

## Cysto-Urethroscope Sheaths



Cysto-Urethroscope sheaths for adults, with central valve, for Ø4,0mm WL302mm Scopes, continuous flow.

| System  | Colour | Ø     | Sheath only | Obturator only | Sheath QL | Obturator QL |
|---------|--------|-------|-------------|----------------|-----------|--------------|
| 4105115 | Yellow | 17 Ch | 4105015     | 4113010        |           |              |
| 4110115 | Yellow | 17 Ch |             |                | 4102015   | 4111010      |
| 4105116 | Green  | 19 Ch | 4105016     | 4113011        |           |              |
| 4110116 | Green  | 19 Ch |             |                | 4102016   | 4111011      |
| 4105117 | Red    | 21 Ch | 4105017     | 4113012        |           |              |
| 4110117 | Red    | 21 Ch |             |                | 4102017   | 4111012      |
| 4105118 | Blue   | 23 Ch | 4105018     | 4113013        |           |              |
| 4110118 | Blue   | 23 Ch |             |                | 4102018   | 4111013      |
| 4105119 | White  | 25 Ch | 4105019     | 4113014        |           |              |
| 4110119 | White  | 25 Ch |             |                | 4102019   | 4111014      |



Visual Obturator

| Ordering n° | Ordering QL n° | Colour | Ø     |
|-------------|----------------|--------|-------|
| 4113040     | 4111040        | Yellow | 17 Ch |
| 4113041     | 4111041        | Green  | 19 Ch |
| 4113042     | 4111042        | Red    | 21 Ch |
| 4113043     | 4111043        | Blue   | 23 Ch |
| 4113044     | 4111044        | White  | 25 Ch |

## Cysto-Urethroscope Working-Sheath Bridges



| Ordering n° | Ordering QL n° | Description                                         |
|-------------|----------------|-----------------------------------------------------|
| 4115020     | 4111020        | Diagnostic Endoscope bridge                         |
| 4115021     | 4111021        | Endoscope bridge, with 1 Instrument channel, 12 Ch  |
| 4115022     | 4111022        | Endoscope bridge, with 2 Instrument channels, 12 Ch |



| Ordering n° | Ordering QL n° | Description                                                                                        |
|-------------|----------------|----------------------------------------------------------------------------------------------------|
| 4115030     | 4111030        | Albarran Working Element with deflector and 1 Ø 12 Ch Instrument channel, with notches             |
| 4115035     |                | Albarran Working Element with deflector and 1 Ø 12 Ch Instrument channel, with continuous movement |
| 4115032     |                | Protection Rod for Albarran                                                                        |




| Ordering n° | Ordering QL n° | Description                                                                                         |
|-------------|----------------|-----------------------------------------------------------------------------------------------------|
| 4115031     | 4111031        | Albarran Working Element with deflector and 2 Ø 12 Ch Instrument channels, with notches             |
| 4115036     |                | Albarran Working Element with deflector and 2 Ø 12 Ch Instrument channels, with continuous movement |
| 4115032     |                | Protection Rod for Albarran                                                                         |

| Ordering n° | Description |
|-------------|-------------|
| 090 -100    | Sealing Cup |



# Cystoscopy

## Laser Cysto-Urethroscope Sheaths



| System  | Colour | Ø     | Sheath only | Obturator only |
|---------|--------|-------|-------------|----------------|
| 4106100 | Blue   | 23 Ch | 4106000     | 4116003        |



Visual obturator, for laser cystoscope-urethroscope sheath

| Ordering n° | Colour | Ø     |
|-------------|--------|-------|
| 4116001     | Blue   | 23 Ch |



Working Insert for continuous flow Laser Cystoscope Urethroscope-Sheath with 8Ch working channel, colour code blue

| Ordering n° | Colour | Ø     |
|-------------|--------|-------|
| 4116002     | Blue   | 23 Ch |

### Albarran deflector for laser cystoscope

| Ordering n° | Colour                                                             | Ø |
|-------------|--------------------------------------------------------------------|---|
| 4116004     | Albarran Working Element with deflector and 1 Instrument channels. |   |
| 4116005     | Albarran Working Element with deflector and 2 Instrument channels. |   |

## Monopolar Flexible Electrodes



Flexible loop monopolar electrode



Flexible needle monopolar electrode



Flexible button monopolar electrode



| WL 470mm | Ø    | WL 670mm | Ø    |
|----------|------|----------|------|
| 4323090  | 5 Ch | 4323090L | 5 Ch |
| 4323091  | 7 Ch | 4323091L | 7 Ch |
| 4323092  | 9 Ch | 4323092L | 9 Ch |

| WL 470mm | Ø    | WL 670mm | Ø    |
|----------|------|----------|------|
| 4323082  | 3 Ch | 4323082L | 3 Ch |
| 4323088  | 5 Ch | 4323088L | 5 Ch |
| 4323089  | 7 Ch | 4323089L | 7 Ch |

| WL 470mm | Ø     | WL 670mm | Ø     |
|----------|-------|----------|-------|
| 4323083  | 3 Ch  | 4323083L | 3 Ch  |
| 4323084  | 4 Ch  | 4323084L | 4 Ch  |
| 4323085  | 5 Ch  | 4323085L | 5 Ch  |
| 4323086  | 6 Ch  | 4323086L | 6 Ch  |
| 4323087  | 7 Ch  | 4323087L | 7 Ch  |
| 4323094  | 8 Ch  | 4323094L | 8 Ch  |
| 4323093  | 9 Ch  | 4323093L | 9 Ch  |
| 4323095  | 10 Ch | 4323095L | 10 Ch |

### Monopolar Cables with 2mm HF connection



Ø 6mm Medici



Ø 4mm Martin/Berchthold



Ø 5mm Erbe/STORZ



Ø 8mm Valleylab/Bovie

| Length | Ordering n° | Ordering n° | Ordering n° | Ordering n° |
|--------|-------------|-------------|-------------|-------------|
| 3mt    | 2373068     | 2373066     | 2373069     | 2373067     |
| 5mt    | 2375068     | 2375066     | 2375069     | 2375067     |

| Adapter | Ordering n° |
|---------|-------------|
|---------|-------------|

Monopolar adapter  
Medici / Ø 2-8mm



2370016



# Cystoscopy

## Bipolar Flexible Electrodes



| Ordering n° | Ø    | WL    |
|-------------|------|-------|
| 4323237     | 5 Ch | 360mm |



Flexible straight needle bipolar electrode

| Ordering n° | Ø    | WL    |
|-------------|------|-------|
| 4323236     | 5 Ch | 360mm |



Flexible needle bipolar electrode

| Ordering n° | Ø    | WL    |
|-------------|------|-------|
| 4323235     | 5 Ch | 360mm |



Flexible button bipolar electrode

### Bipolar Cables



Martin

| Ordering n° | Side plug, generator | Length |
|-------------|----------------------|--------|
| 2373003     | Martin               | 3mt    |
| 2375003     | Martin               | 5mt    |
| 2373001     | 2 Banana             | 3mt    |
| 2375001     | 2 Banana             | 5mt    |



2 Banana Ø 4mm



Valleylab

| Ordering n° | Side plug, generator | Length |
|-------------|----------------------|--------|
| 2373004     | Valleylab            | 3mt    |
| 2375004     | Valleylab            | 5mt    |
| 2373002     | Erbe                 | 3mt    |
| 2375002     | Erbe                 | 5mt    |



Erbe



Medici

| Ordering n° | Side plug, generator | Length |
|-------------|----------------------|--------|
| 2373093     | Medici               | 3mt    |
| 2375093     | Medici               | 5mt    |

### Bipolar adapters



| Ordering n° | Side plug, generator |
|-------------|----------------------|
| 2370017     | Medici/Martin        |
| 2370018     | Medici/Erbe          |
| 2370019     | 2 Banana 4mm         |




# ENDOMED

# Cystoscopy

## Cysto-Urethroscope Forceps, w/400mm



Fix handle system  
With spring




Biopsy Oval Cup Sharp




Grasping forceps, interlocking teeth, concave, short, for stone street



Alligator



3-Prongs Retriever



| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42004008 | 42004008S  | 2,5  | 0,8  |
| 42004010 | 42004010S  | 3,0  | 1,0  |
| 42004012 | 42004012S  | 3,5  | 1,2  |
| 42004014 | 42004014S  | 4,0  | 1,4  |
| 42004016 | 42004016S  | 5,0  | 1,6  |
| 42004021 | 42004021S  | 6,0  | 2,1  |
| 42004023 | 42004023S  | 7,0  | 2,3  |
| 42004030 | 42004030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42074012 | 42074012S  | 3,5  | 1,2  |
| 42074016 | 42074016S  | 5,0  | 1,6  |
| 42074023 | 42074023S  | 7,0  | 2,3  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42044008 | 42044008S  | 2,5  | 0,8  |
| 42044010 | 42044010S  | 3,0  | 1,0  |
| 42044012 | 42044012S  | 3,5  | 1,2  |
| 42044014 | 42044014S  | 4,0  | 1,4  |
| 42044016 | 42044016S  | 5,0  | 1,6  |
| 42044023 | 42044023S  | 7,0  | 2,3  |
| 42044030 | 42044030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42094010 | 42094010S  | 3,0  | 1,0  |
| 42094016 | 42094016S  | 5,0  | 1,6  |
| 42094023 | 42094023S  | 7,0  | 2,3  |

# ENDOMED

# Cystoscopy

## Cysto-Urethroscope Forceps, w/400mm




Fix handle system  
With spring

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42084012 | 42084012S  | 3,5  | 1,2  |
| 42084016 | 42084016S  | 5,0  | 1,6  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42254012 | 42254012S  | 3,5  | 1,2  |
| 42254016 | 42254016S  | 5,0  | 1,6  |
| 42254023 | 42254023S  | 7,0  | 2,3  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42014010 | 42014010S  | 3,0  | 1,0  |
| 42014012 | 42014012S  | 4,0  | 1,3  |
| 42014016 | 42014016S  | 5,0  | 1,6  |
| 42014021 | 42014021S  | 6,0  | 2,1  |
| 42014023 | 42014023S  | 7,0  | 2,3  |
| 42014030 | 42014030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 42024023 | 42024023S  | 7,0  | 2,3  |
| 42024030 | 42024030S  | 9,0  | 3,0  |




Grasping forceps, long, for stone street




Stone fragments  
Grasping Forceps




Single action Scissors



Single action  
Hook Scissors



## Optical Instruments



Optical Instruments with Coated Handle  
for Ø4,0mm-wl302mm Scopes

|                                                                                     | Ordering n° | QUICK LOCK n° | Description                                                         | for Scope                   |
|-------------------------------------------------------------------------------------|-------------|---------------|---------------------------------------------------------------------|-----------------------------|
|    | 4133050     | 4131050       | Optical biopsy spoon forceps, 4,3mm, double action jaws             | Ø 4,0 mm, 0/12°, wl 302 mm  |
|    | 4133051     | 4131051       | Optical grasping forceps, serrated, double action jaws              | Ø 4,0 mm, 0/12°, wl 302 mm  |
|    | 4133052     | 4131052       | Optical biopsy punch forceps, cutting, double action jaws           | Ø 4,0 mm, 0/12°, wl 302 mm  |
|   | 4133053     | 4131053       | Optical biopsy forceps with short jaws 2,3mm<br>double action jaws  | Ø 4,0 mm, 70°, wl 302 mm    |
|  | 4133054     | 4131054       | Optical biopsy forceps with short jaws 2,3mm,<br>double action jaws | Ø 4,0 mm, 12/30°, wl 302 mm |
|  | 4133055     | 4131055       | Optical scissor, hook shaped, double action jaws                    | Ø 4,0 mm, 0/12°, wl 302 mm  |
|  | 4133065     | 4131065       | Grasping forceps for Stents, double action jaws                     | Ø 4,0 mm, 0/12°, wl 302 mm  |
|  | 4133058     | 4131058       | Optical Forceps Adapter,<br>for outer 26Ch resectoscope-sheath      |                             |
|  | 4133059     | 4131059       | Optical Forceps Adapter,<br>for resectoscope-sheath                 |                             |

for sheaths check cystoscope pages

# Cystoscopy

## Flexible Loops Stone Extractor/Injection cannula



| Ordering n° | Description                                         | Ordering n° | Description                    |
|-------------|-----------------------------------------------------|-------------|--------------------------------|
| 4133552     | Stone extractor with 4 wires w/o handle, 5Ch, 600mm | 4360559     | Teflon Tubing                  |
| 4133553     | Stone extractor with 4 wires w/o handle, 7Ch, 600mm | 4360561     | Teflon 3Ch tubing              |
| 4133554     | Stone extractor with 5 wires w/o handle, 5Ch, 600mm | 4360550.3   | Handle for 3Ch stone extractor |
| 4133554A    | Stone extractor with 5 wires w/o handle, 7Ch, 600mm | 4360550.5   | Handle for 5Ch stone extractor |
| 4133555     | Stone extractor with 6 wires w/o handle, 5Ch, 600mm | 4360550.7   | Handle for 7Ch stone extractor |
| 4133555A    | Stone extractor with 6 wires w/o handle, 7Ch, 600mm |             |                                |
| 4133560     | Stone extractor with 4 wires w/o handle, 3Ch, 600mm |             |                                |
| 4133560A    | Stone extractor with 4 wires w/o handle, 3Ch 900mm  |             |                                |
| 4133560B    | Stone extractor with 4 wires w/o handle, 3Ch 1200mm |             |                                |



| Ordering n° | Description                                                                                                    |
|-------------|----------------------------------------------------------------------------------------------------------------|
| 4133556     | Flexible Injection Cannula with luer lock connection for Ø 6Ch (2,0mm) working channels, working lenght 500 mm |

ENDOMED



## Optical Lithotripsy Forceps



Optical stone grasping forceps for lithotripsy of small stones or capturing of stone fragments after therapy with lithotripter, quick lock connection, for Ø4,0mm wL302mm 12/30° scopes.

| Bayonet n° | Quick Lock n° | Description                                                     |
|------------|---------------|-----------------------------------------------------------------|
| 4133056    | 4131056       | Optical Stone Punch, for removal of small stones fragments only |



for WL302mm Scopes, continuous flow

| System  | Colour | Ø     | Sheath only | Obturator only |
|---------|--------|-------|-------------|----------------|
| 4103114 | White  | 25 Ch | 4103014     | 4113014        |




Cysto-Urethroscope sheaths for adults, with central valve, for Ø4,0mm WL302mm Scopes, continuous flow

| System  | Colour | Ø     | Sheath only | Obturator only |
|---------|--------|-------|-------------|----------------|
| 4105119 | White  | 25 Ch | 4105019     | 4113014        |


# Cystoscopy

## Sheath for Lithotripsy Forceps



Continuous Flow outer sheath, double sheath system, quick-lock, with 2 stopcocks, with obturator, 24/26Ch

| Ordering n° | Description                                                  |
|-------------|--------------------------------------------------------------|
| 4417122A    | <b>Fixed</b> outer sheath, all around perforations, 26Ch     |
| 4417127A    | <b>Rotatable</b> outer sheath, all around perforations, 26Ch |



Standard sheath with central valve, quick-lock, **oblique sheath end**

| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4417153 | 4417152 | 27Ch | 4417123   |



Standard sheath, with 1 stopcock, **oblique sheath end**, quick-lock

| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4417118 | 4417116 | 27Ch | 4417123   |



Adapter for Optical Punch, with resectoscope sheath


| Ordering n° | Description                                                   |
|-------------|---------------------------------------------------------------|
| 4133059     | Adapter for Optical Stone Punch, with resectoscope sheath     |
| 4131059     | Adapter for Optical Stone Punch, QL, with resectoscope sheath |

Cystoscopy


Lithotripsy Stone Punch



| Ordering n° | Description                                                            |
|-------------|------------------------------------------------------------------------|
| 4201077     | Working Element for Punch Lithotripsy, (MAUERMAYER), quick lock device |



| System  | Colour | Ø     | Sheath only | Obturator only |
|---------|--------|-------|-------------|----------------|
| 4211181 | Blue   | 23 Ch | 4211081     | 4211080        |




| System  | Colour | Ø     | Sheath only | Obturator only |
|---------|--------|-------|-------------|----------------|
| 4211178 | Blue   | 23 Ch | 4211078     | 4211080        |



| Ordering n° | Description                                                                                     |
|-------------|-------------------------------------------------------------------------------------------------|
| 4211079     | Telescope adapter with 7Ch Instrument channel, for flexible instruments only, quick lock device |


# Cystoscopy

## Lithotripsy Stone Forceps



Cystoscopy Lithotripter/Stone crusher  
with handle and screw mechanism

| Ordering n° | Description           |
|-------------|-----------------------|
| 4133076     | Optical Stone Forceps |



DETACHABLE


All sheaths and working elements can be ordered with bayonet connection or Quick Lock connection



Bayonet connection sheaths



Quick-Lock connection sheaths



Bayonet connection working element



Quick-Lock connection working element

# Cystoscopy

## Monopolar Resectoscope Working Elements for Ø4,0mm Scopes



| Ordering n° | Description                                                                               | Scope                      |
|-------------|-------------------------------------------------------------------------------------------|----------------------------|
| 4403100     | Working Element, passive action, closed handle, 24-26Ch / 27-28,5Ch, Bayonet connection   | Ø 4mm, WI302mm, 0°-12°-30° |
| 4401100     | Working Element, passive action, closed handle, 24-26Ch /27-28,5Ch, Quick-Lock connection | Ø 4mm, WI302mm, 0°-12°-30° |



| Ordering n° | Description                                                                    | Scope                      |
|-------------|--------------------------------------------------------------------------------|----------------------------|
| 4403110     | Working Element, active action, closed handle, 24-26Ch / 27/28,5Ch, Bayonet    | Ø 4mm, WI302mm, 0°-12°-30° |
| 4401110     | Working Element, active action, closed handle, 24-26Ch / 27/28,5Ch, Quick-Lock | Ø 4mm, WI302mm, 0°-12°-30° |




# Cystoscopy

## Resectoscope Sheaths



Standard sheath, with 1 stopcock, **oblique sheath end**, with obturator

| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4413117 | 4413115 | 24Ch | 4413123   |
| 4413118 | 4413116 | 27Ch | 4413123   |



Standard sheath with central valve, **oblique sheath end**, with obturator

| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4413151 | 4413150 | 24Ch | 4413123   |
| 4413153 | 4413152 | 27Ch | 4413123   |


# Cystoscopy

## Resectoscope Sheaths



Continuous Flow sheath, double sheath system, with 2 stopcocks, **ROTATABLE**, oblique shaft 30° end, with obturator, 24/26Ch

| System    | Inner sheath | Inner Tip | Outer sheath | Outer sheath Model                        |                                                                                     |
|-----------|--------------|-----------|--------------|-------------------------------------------|-------------------------------------------------------------------------------------|
| 4413125   | 4413126      | Obblique  | 4413127A     | all around perforations                   |  |
| 4413125A  | 4413126A     | Straight  | 4413127B     | long holes all around, perforations       |  |
|           |              |           | 4413127C     | long holes upper perforations             |                                                                                     |
|           |              |           | 4413127D     | upper and lower perforations              |                                                                                     |
|           |              |           | 4413127E     | upper perforations                        |                                                                                     |
|           |              |           | 4413127F     | extra small holes all around perforations |  |
| Obturator |              | 4413123   |              |                                           |                                                                                     |



Continuous Flow sheath, double sheath system, with 2 stopcocks, **FIXED**, oblique shaft 30° end, with obturator, 24/26Ch

| System    | Inner sheath | Inner Tip | Outer sheath | Outer sheath Model                  |                                                                                       |
|-----------|--------------|-----------|--------------|-------------------------------------|---------------------------------------------------------------------------------------|
| 4413120   | 4413121      | Obblique  | 4413122A     | all around perforations             |  |
| 4413120A  | 4413121A     | Straight  | 4413122B     | all around, long holes perforations |                                                                                       |
|           |              |           | 4413122C     | upper/lower perforations            |                                                                                       |
| Obturator |              | 4413123   |              |                                     |                                                                                       |

For Quick-Lock sheaths refer to bipolar system

ENDOMED

## Resectoscope Obturators



| Ordering n° | Description                           |
|-------------|---------------------------------------|
| 4413130     | Deflecting tip 24Ch obturator, yellow |
| 4413131     | Deflecting tip 27Ch obturator, yellow |




| Ordering n° | Description                                               |
|-------------|-----------------------------------------------------------|
| 4413137     | 24Ch Visual Obturator with 6Ch Instrument channel, yellow |
| 4413138     | 27Ch Visual Obturator with 9Ch Instrument channel, black  |



| Ordering n° | Description                   | Ordering n° | Description                  |
|-------------|-------------------------------|-------------|------------------------------|
| 4413139     | 24Ch Visual Obturator, yellow | 4413134     | 27Ch Visual Obturator, black |

### Monopolar Cables with 2 mm HF connection



| Length | Ordering n° | Ordering n° | Ordering n° | Ordering n° |
|--------|-------------|-------------|-------------|-------------|
| 3mt    | -           | 2373066     | 2373069     | 2373067     |
| 5mt    | 2375068     | 2375066     | 2375069     | 2375067     |

| Adapter                               | Ordering n° |
|---------------------------------------|-------------|
| Monopolar adapter<br>Medici / Ø 2-8mm | 2370016     |



# Cystoscopy

## Monopolar Resectoscope HF Electrodes



Monopolar Resectoscope HF Electrodes

| Description                                                     | 24Ch Yellow | 27Ch Brown |                                                                                       |
|-----------------------------------------------------------------|-------------|------------|---------------------------------------------------------------------------------------|
| HF 30° angled Loop Electrode, for 30°scope, cutting wire 0,35mm | 4423200     | 4423201    |    |
| HF Loop Electrode, straight cutting wire 0,35mm                 | 4423202     | 4423212    |    |
| HF Loop Electrode, 90°angled, for 0° scope, cutting wire 0,35mm | 4423213     | 4423214    |    |
| HF Knife Electrode                                              | 4423203     | 4423203A   |    |
| HF Ball Electrode, 3mm                                          | 4423204     | 4423204A   |    |
| HF Ball Electrode, 5mm                                          | 4423205     | 4423205A   |    |
| HF coagulating cylinder Electrode, 3mm                          | 4423206     | 4423206A   |   |
| HF coagulating cylinder Electrode, 5mm                          | 4423207     | 4423207A   |  |
| HF Cone shaped Electrode                                        | 4423208     |            |  |
| HF Spiked Vapor electrode, 3mm                                  | 4423215     | 4423215A   |  |
| HF Spiked Vapor electrode, 5mm                                  | 4423216     | 4423216A   |  |
| HF Vaporisation Electrode, 3mm                                  | 4423217     | 4423217A   |  |
| HF Vaporisation Electrode, 5mm                                  | 4423218     | 4423218A   |  |
| HF Needle Electrode                                             | 4423225     |            |  |
| HF smooth band-Electrode for 0°scope                            | 4423226     |            |  |
| HF smooth band Electrode for 30°scope                           | 4423230     |            |  |
| HF Band-Vaporisation Electrode                                  | 4423228     |            |  |
| Curette, sharp                                                  | 4423220     |            |                                                                                       |
| Curette, blunt                                                  | 4423227     |            |  |
| Preparation Cold Knife, Straight                                | 4321341     |            |  |
| Preparation Cold Knife, 90°                                     | 4321342     |            |  |
| Preparation Cold Knife H, 90°                                   | 4321345     |            |  |



Sterilization stainless steel tube for Electrodes and cold knives 4423990

ENDOMED

## *Slim Monopolar Resectoscope Working Elements for Ø2,9mm Scopes/Slim Resectoscope HF Electrodes*



| Ordering n° | Description                                                   | Scope                 |
|-------------|---------------------------------------------------------------|-----------------------|
| 4402100     | Working Element, passive action, closed handle, SLIM, 19-22Ch | Ø 2,8mm, WI302mm, 12° |
| 4402110     | Working Element, active action, closed handle, SLIM, 19-22Ch  | Ø 2,8mm, WI302mm, 12° |




### Monopolar Resectoscope HF Electrodes

|                                                                                    | Description                                                     | 19Ch    |
|------------------------------------------------------------------------------------|-----------------------------------------------------------------|---------|
|  | HF 30° angled Loop Electrode, for 30°scope, cutting wire 0,35mm | 4422200 |
|  | HF Loop Electrode, straight, cutting wire 0,35mm                | 4422202 |
|  | HF Knife Electrode                                              | 4422203 |
|  | HF 3mm Ball Electrode                                           | 4422204 |
|  | HF Cylinder Electrode                                           | 4422206 |
|  | HF 90° angled Loop Electrode, for 0°scope, cutting wire 0,35mm  | 4422213 |
|  | HF Vaporisation Electrode                                       | 4422217 |


# Cystoscopy

## Slim Resectoscope Sheaths



Continuous Flow sheath, double sheath system, with 2 stopcocks, **ROTATABLE**, oblique shaft end, without obturator, 19/22Ch.

| System  | Inner sheath | Outer sheath | Model                   | Obturator |
|---------|--------------|--------------|-------------------------|-----------|
| 4412125 | 4412126      | 4412127      | all around perforations | 4412123   |



Continuous Flow sheath, double sheath system, with 2 stopcocks, **FIXED**, oblique shaft end, without obturator, 19/22Ch.

| System  | Inner sheath | Outer sheath | Model                   | Obturator |
|---------|--------------|--------------|-------------------------|-----------|
| 4412120 | 4412121      | 4412122      | all around perforations | 4412123   |



| Sheath  | Description                                                                                | Obturator |
|---------|--------------------------------------------------------------------------------------------|-----------|
| 4412115 | Standard Resectoscope sheath, with 1 stopcock, oblique shaft end, without obturator, 19Ch. | 4412123   |




| Sheath  | Description                                    |
|---------|------------------------------------------------|
| 4412135 | Visual Obturator for 19Ch resectoscope sheath. |




# Cystoscopy

## Laser Resectoscope



| Ordering n° | Description                                                                                    |
|-------------|------------------------------------------------------------------------------------------------|
| 4404100/0,8 | Laser Working Element, passive action, for Ø4,0mm scopes for flexible laser probe until 0,8 mm |
| 4404100/1,2 | Laser Working Element, passive action, for Ø4,0mm scopes for flexible laser probe until 1,2 mm |
| 4404110/0,8 | Laser Working Element, passive action, for Ø2,9mm scopes for flexible laser probe until 0,8 mm |
| 4404110/1,2 | Laser Working Element, passive action, for Ø2,9mm scopes for flexible laser probe until 1,2 mm |
| 4404111     | Luer Lock Connector                                                                            |
| 4404111L    | Long Luer Lock Connector                                                                       |



| Ordering n° | Description                                                              |
|-------------|--------------------------------------------------------------------------|
| 4404100     | Laser Working Element, passive action, without button, for Ø4,0mm scopes |
| 4404050     | Laser Working sheath, for Ø4,0mm scopes, for Laser probes (up to 0,7mm)  |
| 4404100A    | Laser probe guiding channel for Laser probes (up to 0,8mm)               |
| 4404100B    | Laser probe guiding channel for Laser probes (up to 1,0mm)               |
| 4404100C    | Laser probe guiding channel for Laser probes (up to 1,2mm)               |
| 4404100D    | Laser probe guiding channel for Laser probes (up to 1,5mm)               |
| 4404118     | Fiber stop, T-B adapter with 3 silicon seals Long Luer Lock Connector    |

# ENDOMED


Cystoscopy

Optical Forceps for tissue/Sheath




Optical forceps for tissue removing after enucleation, to be used with following sheaths, for Ø4,0mm w/302mm 12/30° Scopes

| Ordering n° | Description                                           |
|-------------|-------------------------------------------------------|
| 4413000     | Optical forceps for tissue removing after enucleation |
| 4413120     | Courette for tissue removing Ø24Ch                    |



Continuous Flow outer sheath, double sheath system, quick-lock, with 2 stopcocks, **FIXED/ROTATABLE**, with obturator, 26Ch

|                                                              | Description |
|--------------------------------------------------------------|-------------|
| <b>Fixed</b> outer sheath, all around perforations, 26Ch     | 4417122A    |
| <b>Rotatable</b> outer sheath, all around perforations, 26Ch | 4417127A    |



Standard sheath with central valve, quick-lock, **oblique sheath end**

|      |         |         |
|------|---------|---------|
| 27Ch | 4417152 | 4417123 |
|------|---------|---------|




Standard sheath, with 1 stopcock, **oblique sheath end**, quick-lock

|      |         |         |
|------|---------|---------|
| 27Ch | 4417116 | 4417123 |
|------|---------|---------|

All sheaths and working elements can be ordered with bayonet or quick lock connection




Bayonet connection sheaths



Bayonet connection  
working element




Quick-Lock connection  
sheaths



Quick-Lock connection  
working element

# Cystoscopy

## Bipolar Resectoscope Working Elements for Ø4,0mm Scopes



| Ordering n° | Description                                                                                                   | Scope                      |
|-------------|---------------------------------------------------------------------------------------------------------------|----------------------------|
| 4407100     | Bipolar Working Element, passive action, closed handle, quick lock connection, for adult, 24-26Ch / 27-28,5Ch | Ø 4mm, Wl302mm, 0°-12°-30° |
| 4407100STO  | Bipolar Working Element, passive action, closed handle, Bayonet connection, for adult, 24-26Ch / 27-28,5Ch    | Ø 4mm, Wl302mm, 0°-12°-30° |




| Ordering n° | Description                                                                          | Scope                      |
|-------------|--------------------------------------------------------------------------------------|----------------------------|
| 4407110     | Bipolar Working Element, active action, closed handle, quick lock, 24/26Ch 27/28,5Ch | Ø 4mm, Wl302mm, 0°-12°-30° |
| 4407110STO  | Bipolar Working Element, active action, closed handle, bayonet, 24/26Ch 27/28,5Ch    | Ø 4mm, Wl302mm, 0°-12°-30° |

## Resectoscope Quick-Lock Sheaths



Continuous Flow sheath, double sheath system, quick-lock, with 2 stopcocks, **FIXED**, oblique shaft end, with obturator, 24/26Ch

| System    | Inner sheath | Inner Tip | Outer sheath | Outer sheath Model                  |                                                                                     |
|-----------|--------------|-----------|--------------|-------------------------------------|-------------------------------------------------------------------------------------|
| 4417120   | 4417121      | Obblique  | 4417122A     | all around perforations             |  |
| 4417120A  | 4417121A     | Straight  | 4417122D     | upper/lower perforations            |  |
|           |              |           | 4417122B     | all around, long holes perforations |  |
| Obturator |              |           | 4417123      |                                     |                                                                                     |



Continuous Flow sheath, double sheath system, quick-lock, with 2 stopcocks, **ROTATABLE**, oblique shaft end, with obturator, 24/26Ch.


| System    | Inner sheath | Inner Tip | Outer sheath | Outer sheath Model                        |                                                                                       |
|-----------|--------------|-----------|--------------|-------------------------------------------|---------------------------------------------------------------------------------------|
| 4417125   | 4417126      | Obblique  | 4417127A     | all around perforations                   |  |
| 4417125A  | 4417126A     | Straight  | 4417127B     | all around, long holes perforations       |                                                                                       |
|           |              |           | 4417127C     | long holes upper perforations             |  |
|           |              |           | 4417127D     | upper and lower perforations              |  |
|           |              |           | 4417127E     | upper perforations                        |                                                                                       |
|           |              |           | 4417127F     | all around/extra small holes perforations |  |
| Obturator |              |           | 4417123      |                                           |                                                                                       |

For Bayonet connector sheaths refer to Monopolar system

# ENDOMED

# Cystoscopy

## Resectoscope Sheaths, Quick-Lock



Standard sheath with central valve, quick-lock, **oblique sheath end**

| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4417151 | 4417150 | 24Ch | 4417123   |
| 4417153 | 4417152 | 27Ch | 4417123   |



Standard sheath, with 1 stopcock, **oblique sheath end**, quick-lock

| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4417117 | 4417115 | 24Ch | 4417123   |
| 4417118 | 4417116 | 27Ch | 4417123   |



Cystoscopy

Resectoscope Obturators



| Description                           | Quick-Lock | Bayonet |
|---------------------------------------|------------|---------|
| Deflecting tip 24Ch Obturator, yellow | 4417130    | 4413130 |
| Deflecting tip 27Ch Obturator, black  | 4417131    | 4413131 |



| Description                                               | Quick-Lock | Bayonet |
|-----------------------------------------------------------|------------|---------|
| 24Ch Visual Obturator with 6Ch Instrument channel, yellow | 4417137    | 4413137 |
| 27Ch Visual Obturator with 9Ch Instrument channel, black  | 4417138    | 4413138 |




| Description                   | Quick-Lock | Bayonet |
|-------------------------------|------------|---------|
| 24Ch Visual Obturator, yellow | 4417139    | 4413139 |
| 27Ch Visual Obturator, black  | 4417134    | 4413134 |



Resectoscope Bipolar Cables

| Ordering n° | Description               | Ordering n° | Description                  |
|-------------|---------------------------|-------------|------------------------------|
| 2375079     | Medici Bipolar cable, 5mt | 2373080     | Valleylab Bipolar cable, 3mt |



# Cystoscopy

## Bipolar Resectoscope HF Electrodes



Bipolar Resectoscope HF Electrodes

| 24Ch TC    | 24Ch PT | 24Ch STO   | 22Ch WF | 24Ch OLY | Description                                                                 |
|------------|---------|------------|---------|----------|-----------------------------------------------------------------------------|
| 4427200    | 4427220 | 4427200STO | 4427240 | 4427230M | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,35mm, MEDIUM size |
| 4427200/B  |         |            |         |          | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,30mm, MEDIUM size |
| 4427200/04 |         |            |         |          | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,40mm, MEDIUM size |
| 4427202    | 4427222 | 4427202STO |         | 4427232  | HF Loop Electrode, straight, cutting wire 0,35mm                            |
| 4427213    | 4427226 | 4427213STO | 4427241 | 4427235  | HF Loop Electrode, 90°angled, for 0° scope, cutting wire 0,35mm             |
| 4427203    |         | 4427203STO | 4427243 | 4427233  | HF Knife Electrode                                                          |
| 4427204    |         | 4427204STO |         | 4427234  | HF Ball Electrode, 3mm                                                      |
| 4427206    |         | 4427206STO | 4427246 | 4427236  | HF Vaporitation Electrode                                                   |



| 24Ch OLY | Description                                                                  |
|----------|------------------------------------------------------------------------------|
| 4427230M | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,35mm, MEDIUM size  |
| 4427230S | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,35mm SMALL size    |
| 4427230L | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,35mm, LARGE size   |
| 4427231M | HF 12° angled, Loop Electrode for 12° scope, cutting wire 0,35mm MEDIUM size |
| 4427231S | HF 12° angled, Loop Electrode for 12° scope, cutting wire 0,35mm SMALL size  |
| 4427231L | HF 12° angled, Loop Electrode for 12° scope, cutting wire 0,35mm, LARGE size |



| 24Ch OLY | Description                                                                  |
|----------|------------------------------------------------------------------------------|
| 4427237M | HF 45° angled, Loop Electrode for 30° scope, cutting wire 0,35mm MEDIUM size |
| 4427237S | HF 45° angled, Loop Electrode for 30° scope, cutting wire 0,35mm SMALL size  |
| 4427237L | HF 45° angled, Loop Electrode for 30° scope, cutting wire 0,35mm, LARGE size |


| 24Ch GYR | Description                                                                                                |
|----------|------------------------------------------------------------------------------------------------------------|
| 4427250  | HF 30°angled, Loop Electrode for 30°scope, cutting wire 0,35mm, Use with Gyrus super Pulse Generator only. |




| Description                                                       | Ordering n° |
|-------------------------------------------------------------------|-------------|
| Sterilization stainless steel tube for Electrodes and cold knives | 4423990     |

# ENDOMED

Mini Bipolar Resectoscope Working Elements for Ø2,9mm Scopes



| Ordering n° | Description                                                                                                                                |
|-------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| 4408100     | Bipolar Working Element, passive action, closed handle, quick-lock, for 17,0Ch - 17,5/18,5Ch. Sheaths, Ø 2,9mm, WI 302mm, 0°-12°-30° Scope |



Continuous Flow sheath, double sheath system, with 2 stopcocks, oblique shaft end

| ROTATABLE System | Inner sheath | Inner Tip | Outer sheath | Obturator |
|------------------|--------------|-----------|--------------|-----------|
| 4418125          | 4418126      | Oblique   | 4418127      | 4418123   |

| FIX System | Inner sheath | Inner Tip | Outer sheath | Obturator |
|------------|--------------|-----------|--------------|-----------|
| 4418120    | 4418121      | Oblique   | 4418122      | 4418123   |

# Cystoscopy

## Resectoscope Sheaths



| System  | Sheath  | Ø    | Obturator |
|---------|---------|------|-----------|
| 4418117 | 4418115 | 17Ch | 4418123   |



### Bipolar Resectoscope HF Electrodes


| Description                                                     | 24Ch Yellow |
|-----------------------------------------------------------------|-------------|
| HF 30° angled Loop Electrode, for 30°scope, cutting wire 0,35mm | 4428200     |
| HF 30° angled Loop Electrode, for 30°scope, cutting wire 0,30mm | 4428200/0,3 |
| HF Loop Electrode, straight cutting wire 0,35mm                 | 4428202     |
| HF 90° angled Loop Electrode, for 0°scope, cutting wire 0,35mm  | 4428213     |
| HF Knife Electrode                                              | 4428203     |
| HF Ball Electrode, 3mm                                          | 4428204     |



### Resectoscope Bipolar Cables



| Ordering n° | Description               | Ordering n° | Description                  |
|-------------|---------------------------|-------------|------------------------------|
| 2375079     | Medici Bipolar cable, 5mt | 2373080     | Valleylab Bipolar cable, 3mt |



### Bipolar adapters



| Ordering n° | Side plug, generator |
|-------------|----------------------|
| 2370017     | Medici/Martin        |
| 2370018     | Medici/Erbe          |
| 2370019     | 2 Banana 4mm         |



## Urethrotome Working Element for Ø4,0mm Scopes



Sachse Urethrotome

Do not use with HF energy!!

| Ordering n° | Description                                                                                | Scope                |
|-------------|--------------------------------------------------------------------------------------------|----------------------|
| 4403100     | Working Element, passive action, closed handle, with guide tube for cold knives, for adult | Ø 4mm, WI302mm, 0°   |
| 4403101     | Working Element, passive action, open handle, with guide tube for could knives, for adult  | Ø 4mm, WI302mm, 0°   |
| 4303100     | Working element, passive, for 15,5Ch Urethrotome, without HF connector                     | Ø 2,9mm, WI302mm, 0° |



| Ordering n° | Description                                                                                           | Scope                 |
|-------------|-------------------------------------------------------------------------------------------------------|-----------------------|
| 4311400     | 21Ch Urethrotome sheath with lateral catheter insertion, 5Ch, w/o Obturator, for adult                | Ø 4,0mm, WL 302mm, 0° |
| 4313450     | 15,5Ch Urethrotome sheath with lateral catheter insertion, 3Ch, w/o Obturator to be used with 4303100 | Ø 2,9mm, WL 302mm, 0° |

# Cystoscopy

## Urethrotome Sheaths



| Ordering n° | Description                                                                                          | Scope                 |
|-------------|------------------------------------------------------------------------------------------------------|-----------------------|
| 4311403     | Outer Sheath, half round for placing balloon catheter, slides to be used with 21 Ch sheath (4311400) | Ø 4,0mm, WI 302mm, 0° |



| Ordering n° | Description                                                                          | Scope                 |
|-------------|--------------------------------------------------------------------------------------|-----------------------|
| 4311404     | Continuous flow 21Ch outer sheath, optional for Continuous flow, slides over 4311400 | Ø 4,0mm, WI 302mm, 0° |



| Ordering n° | Description                                                                                          | Scope                 |
|-------------|------------------------------------------------------------------------------------------------------|-----------------------|
| 4311402     | Telescope adapter with 1 Instrument channel for 5Ch flexible Instruments to be used with 21Ch sheath | Ø 4,0mm, WI 302mm, 0° |
| 4312402     | Telescope adapter with 1 Instrument channel for 5Ch flexible Instruments to be used with 17Ch sheath | Ø 2,9mm, WI 302mm, 0° |
| 4313454     | Telescope adapter 1 Instrument channel for 5Ch flexible Instruments to be used with 15,5Ch sheath    | Ø 2,9mm, WI 302mm, 0° |
| 090-100     | Sealing Cup 10pcs                                                                                    |                       |

## Urethrotome Obturators-Knives



| Ordering n° | Description                                   |
|-------------|-----------------------------------------------|
| 4311405     | Standard Obturator                            |
| 4312401     | Standard obturator for 7Ch Urethrotome sheath |
| 4313451     | Standard obturator for 5Ch Urethrotome sheath |



| Ordering n° | Description      | Scope              |
|-------------|------------------|--------------------|
| 4311406     | Visual Obturator | Ø 4mm, WI302mm, 0° |



| Ordering n° | Description                                                              | Scope              |
|-------------|--------------------------------------------------------------------------|--------------------|
| 4311401     | Hollow obturator with central channel, for insertion of guiding catheter | Ø 4mm, WI302mm, 0° |




| Description                   | 21Ch    | 17/15Ch |
|-------------------------------|---------|---------|
| Cold Knife, Straight          | 4321410 | 4323410 |
| Cold Knife, Round             | 4321411 | 4323411 |
| Cold Knife, Straight Serrated | 4321412 | 4323412 |
| Cold Knife, ring-like         | 4321413 | 4323413 |
| Cold Knife, Hook-like         | 4321414 | 4323414 |
| Cold Knife, Sickle-like       | 4321415 |         |



# Cystoscopy

## Urethrotome and Accessories



| Ordering n° | Description                                                               |
|-------------|---------------------------------------------------------------------------|
| 166-1702-00 | OTIS Urethrotome Dilator SET, 2 Knives, 3 different Tips in a wooden case |



| Ordering n° | Description      |
|-------------|------------------|
| 166-1702-01 | OTIS Urethrotome |



| Ordering n° | Description                      |
|-------------|----------------------------------|
| 166-1702-02 | Ball Tip for OTIS Urethrotome    |
| 166-1702-03 | Flat Tip for OTIS Urethrotome    |
| 166-1702-04 | Conical Tip for OTIS Urethrotome |



| Ordering n° | Description                |
|-------------|----------------------------|
| 166-1702-99 | Knife for OTIS Urethrotome |

# Cystoscopy

## Accessories



| Ordering n° | Description                                               |
|-------------|-----------------------------------------------------------|
| 4133500     | Evacuator (Glass bulb + Rubber bulb) w/o sheath connector |
| 4133504     | Glass Bulb                                                |
| 4133505     | Rubber Bulb                                               |



| Ordering n° | Description                                     |
|-------------|-------------------------------------------------|
| 4133502     | Sheath connector, simple with tubing            |
| 4133503     | Sheath connector with locking device and tubing |
| 4133503A    | Sheath connector with tubing, QUICK-LOCK device |
| 4133058     | Sheath connector, QUICK-LOCK device to BAYONET  |



| Ordering n° | Description                                  |
|-------------|----------------------------------------------|
| 4133507     | Bladder Syringe 50ml with connector          |
| 4133507A    | Bladder Syringe 75ml with connector          |
| 4133507B    | Bladder Syringe 100ml with connector         |
| 4133507C    | Bladder Syringe 150ml with connector         |
| 4133507D    | Bladder Syringe 200ml with connector         |
| 4133507Q    | Bladder Syringe 50ml with quick-lock device  |
| 4133507AQ   | Bladder Syringe 75ml with quick-lock device  |
| 4133507BQ   | Bladder Syringe 100ml with quick-lock device |
| 4133507CQ   | Bladder Syringe 150ml with quick-lock device |
| 4133507DQ   | Bladder Syringe 200ml with quick-lock device |


# Cystoscopy

## Nephroscopes



*High Definition, Autoclavable,  
max.134°C/2,2bar/5min*

| Ordering n° | Description                                                                                                                        |
|-------------|------------------------------------------------------------------------------------------------------------------------------------|
| 45200YO     | Medici Nephroscope, 45° angled ocular, 6° vision angle, WL 200mm, Ø5 mm working channel, 80° field of view (+/-5° ), autoclavable  |
| 45200YC     | Medici Nephroscope, 45° angled ocular, 15° vision angle, WL 200mm, Ø5 mm working channel, 80° field of view (+/-5° ), autoclavable |
| 45250YO     | Medici Nephroscope, 45° angled ocular, 6° vision angle, WL 250mm, Ø5 mm working channel, 70° field of view (+/-5° ), autoclavable  |
| 45250YC     | Medici Nephroscope, 45° angled ocular, 15° vision angle, WL 250mm, Ø5 mm working channel, 70° field of view (+/-5° ), autoclavable |



| Ordering n° | Description                                                                                                                        |
|-------------|------------------------------------------------------------------------------------------------------------------------------------|
| 45201YO     | Medici Nephroscope, 90° angled ocular, 6° vision angle, WL 200mm, Ø5 mm working channel, 70° field of view (+/-5° ), autoclavable  |
| 45201YC     | Medici Nephroscope, 90° angled ocular, 15° vision angle, WL 200mm, Ø5 mm working channel, 70° field of view (+/-5° ), autoclavable |
| 45251YO     | Medici Nephroscope, 90° angled ocular, 6° vision angle, WL 250mm, Ø5 mm working channel, 70° field of view (+/-5° ), autoclavable  |
| 45251YC     | Medici Nephroscope, 90° angled ocular, 15° vision angle, WL 250mm, Ø5 mm working channel, 70° field of view (+/-5° ), autoclavable |

| Ordering n° | Description                                                                      |
|-------------|----------------------------------------------------------------------------------|
| 45251YY     | Irrigation adaptor M18x1 Ø5,0 with valve for simultaneous working and irrigation |



| Ordering n° | Description              |
|-------------|--------------------------|
| 4133131     | Stopcock for Nephroscope |

# ENDOMED

## Nephroscope Forceps, Semirigid, Ø3,5mm

Fix handle system  
with luerlock port for easy  
and optimal cleaning



360°rotatable  
With or without spring



| WL 320mm | WL 450mm | WL 320mm -<br>rotatable | WL 450mm -<br>rotatable |
|----------|----------|-------------------------|-------------------------|
| 45013235 | 45014535 | 45013235R               | 45014535R               |

Stone grasping forceps for nephroscope, Alligator teeth



| WL 320mm | WL 450mm | WL 320mm -<br>rotatable | WL 450mm -<br>rotatable |
|----------|----------|-------------------------|-------------------------|
| 45033235 | 45034535 | 45033235R               | 45034535R               |

Stone grasping forceps, for nephroscope, with fenestrated jaws




| WL 320mm | WL 450mm | WL 320mm -<br>rotatable | WL 450mm -<br>rotatable |
|----------|----------|-------------------------|-------------------------|
| 45023235 | 45024535 | 45023235R               | 45024535R               |

Biopsy forceps for nephroscope, cup, sharp, single action jaws



| WL 320mm | WL 450mm | WL 320mm -<br>rotatable | WL 450mm -<br>rotatable |
|----------|----------|-------------------------|-------------------------|
| 45043235 | 45044535 | 45043235R               | 45044535R               |

Micro-scissors, straight, single action jaws



| WL 320mm | WL 450mm | WL 320mm -<br>rotatable | WL 450mm -<br>rotatable |
|----------|----------|-------------------------|-------------------------|
| 45053235 | 45054535 | 45053235R               | 45054535R               |


3-prongs Retriever forceps



| Ordering n° | Description                                                                                                                           |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------|
| 45400530    | Three-prong Stone grasper, for 10Ch channel (3mm), Working length 400mm, with luer lock for irrigation, U-spring handle, self closing |

# Cystoscopy

## Nephroscope Sheaths



| 26CH    | 24CH    | Description                                      |
|---------|---------|--------------------------------------------------|
| 4511124 | 4511100 | Nephroscope-sheath, 1 fixed stopcock, WL200mm    |
| 4511126 | 4511102 | Nephroscope-sheath, 1 rotating stopcock, WL200mm |
| 4511121 | 4511101 | Obturator for Nephroscope sheath, WL200mm        |
| 4501126 | 4501102 | Nephroscope-sheath, 1 rotating stopcock, WL250mm |
| 4501121 | 4501101 | Obturator for Nephroscope sheath, WL250mm        |
| 4501124 | 4501100 | Nephroscope-sheath, 1 fixed stopcock, WL250mm    |




| 26CH    | Description                                                               |
|---------|---------------------------------------------------------------------------|
| 4511122 | Nephroscope sheath, 26Ch, 1 fixed stopcock, WL200mm, distal irrigation    |
| 4511123 | Nephroscope sheath, 26Ch, 1 rotating stopcock, WL200mm, distal irrigation |
| 4511121 | Obturator for 26Ch Nephroscope sheath                                     |
| 4501123 | Nephroscope sheath, 26Ch, 1 rotating stopcock, WL250mm, distal irrigation |
| 4501121 | Obturator for 26Ch Nephroscope sheath, WL250mm                            |
| 4501122 | Nephroscope sheath, 26Ch, 1 fixed stopcock, WL250mm, distal irrigation    |

# Cystoscopy


## Mini-Nephroscope

High Definition

Autoclavable, max.134°C/2,2bar/5min



| Ording n° | Description                                                                                                          |
|-----------|----------------------------------------------------------------------------------------------------------------------|
| 45220NA   | Nephroscope with 30° angled ocular, Ø 12Ch, viewing angle 12°, WL 220mm, Ø 7,0Ch working channel, 1x5,0mm or 2x3,0mm |




| Ording n° | Description                                                          |
|-----------|----------------------------------------------------------------------|
| 4511250   | Operative sheath, 15/18Ch, for Continuous flow /Mini-Nephroscopy     |
| 4511251   | Operative sheath, 16,5/19,5Ch, for Continuous flow /Mini-Nephroscopy |
| 4511252   | Operative sheath, 21/24Ch, for Continuous flow / Mini-Nephroscopy    |
| 4511255   | Applicator for seal, to be used with sheath + plunger                |
| 4511256   | Dilator with central channel for wire guiding, for 15/18Ch sheath    |
| 4511257   | Dilator, double channel for wire guiding, for 16,5/19,5Ch sheath     |
| 4511258   | Dilator, double channel for wire guiding, for 21/24Ch sheath         |

# ENDOMED

# Cystoscopy

## Mini-Nephroscope Forceps, Semirigid, Ø1,8mm




| WL 320mm                                                       | WL 320mm - rotatable |                                                                                       |
|----------------------------------------------------------------|----------------------|---------------------------------------------------------------------------------------|
| 45013218                                                       | 45013218R            |                                                                                       |
| Stone grasping forceps for nephroscope, Alligator teeth        |                      |    |
| WL 320mm                                                       | WL 320mm - rotatable |                                                                                       |
| 45023218                                                       | 45023218R            |                                                                                       |
| Biopsy forceps for nephroscope, cup, sharp, single action jaws |                      |   |
| WL 320mm                                                       | WL 320mm - rotatable |                                                                                       |
| 45043218                                                       | 45043218R            |                                                                                       |
| Micro-scissors, straight, single action jaws                   |                      |  |



# Cystoscopy


## Ureterorenoscopes



*High Definition 30.000 PIXEL*  
*Autoclavable, max.134°C/2,2bar/5min*

| Ordering n° | MARBERGER Ø 4,5/6,5 Fr                                                                                                                                                              |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 46315JA1    | Semirigid Ureterorenoscope, 10.000 PIXEL, Marberger ocular, Ø steps by 4,5/6,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 4,0 Ch for 1x3ch.            |
| 46430JA1    | Semirigid Ureterorenoscope, 10.000 PIXEL, Marberger ocular, Ø steps by 4,5/6,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 4,0 Ch for 1x3ch.            |
| Ordering n° | MARBERGER Ø 6,0/7,5 Fr                                                                                                                                                              |
| 46315RA1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.   |
| 46430RA1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.   |
| 46315RB1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.  |
| 46430RB1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.  |
| Ordering n° | MARBERGER Ø 8,0/9,8 Fr                                                                                                                                                              |
| 46315TB1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.  |
| 46430TB1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.  |
| Ordering n° | MARBERGER Ø 8,5/11,5 Fr                                                                                                                                                             |
| 46315VB1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 8,5/11,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x6 or 2X4,0ch. |
| 46430VB1    | Semirigid Ureterorenoscope, 30.000 PIXEL, Marberger ocular, Ø steps by 8,5/11,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x6 or 2X4,0ch. |


| Ordering n | Description                      |
|------------|----------------------------------|
| 46ACC2     | Instrument Bridge                |
| 46ACC4     | Adapter for EMS Lithotripter LC1 |
| 46ACC6     | Adapter for EMS Lithoreipter LC2 |
| 46ACC8     | Stopcock 6mm, Luer Lock          |



# ENDOMED

# Cystoscopy

## Ureterorenoscopes



*High Definition 30.000 PIXEL  
Autoclavable, max.134°C/2,2bar/5min*


| Ordering n° |                                                                                                                                                                                          | BICHLER Ø 4,5/6,5 Fr  |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|
| 46315JA2    | Semirigid Ureterorenoscope, 10.000 PIXEL, with 30° angled ocular, Ø steps by 4,5/6,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 4,0 Ch for 1x3ch.           |                       |
| 46430JA2    | Semirigid Ureterorenoscope, 10.000 PIXEL, with 30° angled ocular, Ø steps by 4,5/6,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 4,0 Ch for 1x3ch.           |                       |
| Ordering n° |                                                                                                                                                                                          | BICHLER Ø 6,0/7,5 Fr  |
| 46315RA2    | Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.  |                       |
| 46430RA2    | Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.  |                       |
| 46315RB2    | Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch. |                       |
| 46430RB2    | Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch. |                       |
| Ordering n° |                                                                                                                                                                                          | BICHLER Ø 8,0/9,8 Fr  |
| 46315TB2    | Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch. |                       |
| 46430TB2    | Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch. |                       |
| Ordering n° |                                                                                                                                                                                          | BICHLER Ø 8,5/11,5 Fr |
| 46315VB2    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,5/11,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x6 or 2X4,0ch.       |                       |
| 46430VB2    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,5/11,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x6 or 2X4,0ch.       |                       |
| Ordering n° |                                                                                                                                                                                          | STRAIGHT Ø 6,0/7,5 Fr |
| 46315RA3    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.         |                       |
| 46430RA3    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.         |                       |
| 46315RB3    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.        |                       |
| 46430RB3    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.        |                       |
| Ordering n° |                                                                                                                                                                                          | STRAIGHT Ø 8,0/9,8 Fr |
| 46315TB3    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.        |                       |
| 46430TB3    | Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.        |                       |

# ENDOMED

## Ureterorenoscope Forceps, wl600mm




Biopsy Oval Cup Sharp



Grasping forceps, cup, serrated



Alligator



Grasping forceps,  
for stone street



| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46006008 | 46006008S  | 2,5  | 0,8  |
| 46006010 | 46006010S  | 3,0  | 1,0  |
| 46006012 | 46006012S  | 4,0  | 1,4  |
| 46006016 | 46006016S  | 5,0  | 1,6  |
| 46006021 | 46006021S  | 6,0  | 2,0  |
| 46006023 | 46006023S  | 7,0  | 2,3  |
| 46006030 | 46006030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46076012 | 46076012S  | 4,0  | 1,4  |
| 46076016 | 46076016S  | 5,0  | 1,6  |
| 46076023 | 46076023S  | 7,0  | 2,3  |
| 46076030 | 46076030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46046008 | 46046008S  | 2,5  | 0,8  |
| 46046010 | 46046010S  | 3,0  | 1,0  |
| 46046012 | 46046012S  | 4,0  | 1,4  |
| 46046016 | 46046016S  | 5,0  | 1,6  |
| 46046023 | 46046023S  | 7,0  | 2,3  |
| 46046030 | 46046030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46086012 | 46086012S  | 4,0  | 1,4  |
| 46086016 | 46086016S  | 5,0  | 1,6  |

# Cystoscopy

## Ureterorenoscope Scissors, Forceps, wl600mm



Fix handle system  
with spring

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46016010 | 46016010S  | 3,0  | 1,0  |
| 46016012 | 46016012S  | 4,0  | 1,2  |
| 46016016 | 46016016S  | 5,0  | 1,6  |
| 46016023 | 46016023S  | 7,0  | 2,3  |
| 46016030 | 46016030S  | 9,0  | 3,0  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46056023 | 46056023S  | 7,0  | 2,3  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46626012 | 46626012S  | 4,0  | 1,4  |
| 46626016 | 46626016S  | 5,0  | 1,6  |

| Flexible | Semi Rigid | Ø Ch | Ø mm |
|----------|------------|------|------|
| 46256012 | 46256012S  | 4,0  | 1,4  |
| 46256016 | 46256016S  | 5,0  | 1,6  |
| 46256023 | 46256023S  | 7,0  | 2,3  |




Single action Scissors




Double action Scissors



Uretero-Slitting  
retrograde Scissors



Stone fragments  
Grasping Forceps




## Catheter introducers-Clamps-Forceps




| Ordering n° | Description                          | WL     |
|-------------|--------------------------------------|--------|
| 166-1804-39 | Catheter introducer, GUYON, straight | 390 mm |




| Ordering n° | Description                        | WL     |
|-------------|------------------------------------|--------|
| 166-1805-35 | Catheter introducer, GUYON, curved | 350 mm |



| Ordering n° | Description           | WL     |
|-------------|-----------------------|--------|
| 166-1707-11 | Meatus clamp, strauss | 110 mm |



| Ordering n° | Description             | WL    |
|-------------|-------------------------|-------|
| 166-1705-08 | Meatus clamp, stockmann | 80 mm |



| Ordering n° | Description                   | WL     |
|-------------|-------------------------------|--------|
| 166-1608-20 | Foreign body forceps, MATHIEU | 200 mm |
| 166-1608-28 | Foreign body forceps, MATHIEU | 280 mm |

# Cystoscopy

## Dilating Bougies-Metal Catheters



|       | Dilating Bougie<br>van BUREN | Dilating Bougie<br>GUYON | Metal Catheter<br>for man | Metal Catheter<br>for women | Coxeter<br>Prostate Catheter |
|-------|------------------------------|--------------------------|---------------------------|-----------------------------|------------------------------|
| 6 Ch  | 166-2001-06                  | 166-2005-06              |                           |                             |                              |
| 8 Ch  | 166-2001-08                  | 166-2005-08              | 166-2101-08               | 166-2103-08                 | 166-2109-08                  |
| 10 Ch | 166-2001-10                  | 166-2005-10              | 166-2101-10               | 166-2103-10                 | 166-2109-10                  |
| 12 Ch | 166-2001-12                  | 166-2005-12              | 166-2101-12               | 166-2103-12                 | 166-2109-12                  |
| 14 Ch | 166-2001-14                  | 166-2005-14              | 166-2101-14               | 166-2103-14                 | 166-2109-14                  |
| 16 Ch | 166-2001-16                  | 166-2005-16              | 166-2101-16               | 166-2103-16                 | 166-2109-16                  |
| 18 Ch | 166-2001-18                  | 166-2005-18              | 166-2101-18               | 166-2103-18                 | 166-2109-18                  |
| 20 Ch | 166-2001-20                  | 166-2005-20              | 166-2101-20               | 166-2103-20                 | 166-2109-20                  |
| 22 Ch | 166-2001-22                  | 166-2005-22              | 166-2101-22               | 166-2103-22                 | 166-2109-22                  |
| 24 Ch | 166-2001-24                  | 166-2005-24              | 166-2101-24               | 166-2103-24                 | 166-2109-24                  |
| 26 Ch | 166-2001-26                  | 166-2005-26              |                           |                             |                              |
| 28 Ch | 166-2001-28                  | 166-2005-28              |                           |                             |                              |
| 30 Ch | 166-2001-30                  | 166-2005-30              |                           |                             |                              |
| 32 Ch | 166-2001-32                  | 166-2005-32              |                           |                             |                              |
| 34 Ch | 166-2001-34                  | 166-2005-34              |                           |                             |                              |
| 36 Ch | 166-2001-36                  | 166-2005-36              |                           |                             |                              |
| 38 Ch | 166-2001-38                  | 166-2005-38              |                           |                             |                              |
| 40 Ch | 166-2001-40                  | 166-2005-40              |                           |                             |                              |

# Cystoscopy

## Pediatric Endoscopes




High Definition, wide angle,  
autoclavable, max. 134° C / 2,2 bar / 5 min

### Diameter 2,0

| Ordering n°  | Ø     | Viewing Angle | Working length |
|--------------|-------|---------------|----------------|
| 40 188 CA.HD | 2,0mm | 0°            | 188mm          |
| 40 188 CE.HD | 2,0mm | 30°           | 188mm          |

### Diameter 2,7

| Ordering n°  | Ø     | Viewing Angle | Working length |
|--------------|-------|---------------|----------------|
| 40 188 GA.HD | 2,7mm | 0°            | 188mm          |
| 40 188 GE.HD | 2,7mm | 30°           | 188mm          |
| 40 188GG.HD  | 2,7mm | 70°           | 188mm          |




| Ordering n° | Ø     | Description                          |
|-------------|-------|--------------------------------------|
| 4117020     | 13 Ch | Examination insert, telescope bridge |




# Cystoscopy

## Pediatric Cystoscope Sheaths



Cystoscope-urethroscope sheath,  
with two fixed stopcocks, 3 Ch instrument channel

| System   | Ø     | Sheath   | Obturator | Scope           |
|----------|-------|----------|-----------|-----------------|
| 4107113  | 13 Ch | 4107013  | 4117013   | Ø 2,7mm WL187mm |
| 4108113  | 11 Ch | 4108013  | 4118013   | Ø 2,0mm WL187mm |
| 4109113  | 09 Ch | 4109013  | 4119013   | Ø2,0mm WL187mm  |
| 4109113A | 09 Ch | 4109013A | 4119013A  | Ø1,9mm WL120mm  |
| 4109113B | 7,5Ch | 4109013B | 4119013B  | Ø1,9mm WL117mm  |



Working insert with ALBARRAN lever

| Ordering n° | Ø     | Scope           |
|-------------|-------|-----------------|
| 4117015     | 13 Ch | Ø 2,7mm WL187mm |


## Pediatric Urethrotome



Pediatric Working Element, Passive

| Ordering n° | Ø     |
|-------------|-------|
| 4405100     | 13 Ch |
| 4406100     | 11 Ch |

Urethrotome Sheath, with two rotatable stopcocks



| Ordering n° | Ø     | Obturator |
|-------------|-------|-----------|
| 4314400     | 13 Ch | 4314401   |
| 4315400     | 11 Ch | 4315401   |



| Ordering n° | Description                                          |
|-------------|------------------------------------------------------|
| 4314402     | Examination insert, with 1 instrument channel, 13 Ch |
| 4315402     | Examination insert, with 1 instrument channel, 11 Ch |


|  | 13 Ch   | 11 Ch   | Description              |
|--|---------|---------|--------------------------|
|  | 4324410 | 4325410 | Cold Knife, straight     |
|  | 4324411 | 4325411 | Cold Knife, round        |
|  | 4324412 | 4325412 | Cold Knife, hook shaped  |
|  | 4324413 | 4325412 | Cold Knife, sickle knife |

Cystoscopy

Pediatric Resectoscope

Pediatric Working Element, Passive

| Ordering n° | Ø     |
|-------------|-------|
| 4405100     | 13 Ch |
| 4406100     | 11 Ch |



Resectoscope Sheath, with rotatable stopcock

| System  | Description                                           | Obturator |
|---------|-------------------------------------------------------|-----------|
| 4415115 | Resectoscope Sheath, with rotatable stopcock, Ø 13 Ch | 4415116   |
| 4416115 | Resectoscope Sheath, with rotatable stopcock, Ø 11 Ch | 4416116   |

| Ordering n° | Description                                            |
|-------------|--------------------------------------------------------|
| 4314402     | Examination insert, with 1 instrument channel, Ø 13 Ch |
| 4315402     | Examination insert, with 1 instrument channel, Ø 11 Ch |



| 13 Ch   | 11 Ch   | Description                   |
|---------|---------|-------------------------------|
| 4425200 | 4426200 | HF Loop Electrode, Gr.1 small |
| 4425201 | 4426201 | HF Loop Electrode, Gr.2 large |
| 4425203 | 4426203 | HF Knife Electrode            |
| 4425204 | 4426204 | HF Ball Electrode             |
| 4425204 | 4425205 | HF Hook Electrode             |



Cystoscopy

*Memo*

ENDOMED







# MEDICI

**Company:**



**ENDOMED**

web: [www.endomedsrl.com](http://www.endomedsrl.com)

tel. 0771.726355

tel. 0823.1656213

fax. 0823.932337

mail: [info@endomedsrl.com](mailto:info@endomedsrl.com)

**Distributed By:**