



Detro Enzym

Enzymatic Cleaner for Medical Instruments and Endoscopes

Detro Enzym is a high performance concentrated cleaning solution for medical instruments and equipment specifically designed for your professional cleaning needs. It has been specially developed to ensure effective cleaning of surgical instruments, medical devices and endoscopes. The special aldehyde

and phenol-free formula is a versatile product that can be used in a variety of areas.

The product provides a more effective cleaning by removing organic wastes (blood, oil, protein, starch, etc.) from surfaces and dissolving residues on instruments thanks to the enzymes it contains. It can be used compatible with different materials such as glass, rubber, wood, ceramic, silicone, plastic and stainless steel. Removal of organic contamination contributes to the reduction of microbial load. It can be used in harmony with various materials and minimizes the risk of corrosion. This helps you to safely clean your instruments and equipment. It offers effective results within 5 minutes in both manual and ultrasonic cleaning processes. The concentrated nature of the product allows you to clean effectively with a small amount of product.

USAGE AREAS:

Pre-cleaning of all kinds of surgical instruments, medical devices and endoscopes Organic wastes such as blood, fat, protein, starch on the endoscope surface Glass, rubber, wood, ceramic, silicone, plastic and stainless steel materials.

Detro Matic NS

Neutralization Solution for Medical Instrument Washing Machines

Detro Matic NS is a neutralization solution specially designed for automatic washing machines. This product is used for neutralization after alkaline cleaning of instruments made of various materials such as metal, glass, ceramic, porcelain, and plastic. This solution containing citric acid is used for

neutralization of alkaline residues remaining at the end of the

main washing process, elimination of lime deposits that may form in the machine, and provides effective neutralization of alkaline residues. Detro Matic NS, which is compatible with all washer-disinfector devices, effectively removes alkaline residues thanks to its strong neutralization capacity while ensuring that machines and instruments are not damaged. Containing corrosion inhibitors, this product is compatible with all washing machine brands and offers economic advantages with its low usage concentration.

USAGE AREAS:

Rinsing after alkaline cleaning in automatic washing machines

