



27005 FA



HOPKINS® Telescope 12°,
enlarged view, diameter 4 mm, length 30 cm, autoclavable, fiber optic light transmission incorporated,
color code: black



495 GO




Adaptor,
to connect KARL STORZ HOPKINS® telescopes to Olympus/Winter & Ibe light cable, with attached lens



27076 A



Lithotrite,
incorporating both handle and turning screw action, 24 Fr.



27293 BD




Operating Sheath,
26 Fr., for continuous irrigation and suction, with LUER-Lock stopcock, rotating,
color code: black-red



27293 CD



Operating Sheath,
24 Fr., for continuous irrigation and suction, with LUER-Lock stopcock, rotating,
color code: white-red



Percutaneous Nephroscope with Continuous Irrigation and Suction

Length 25 cm

Sheaths for continuous irrigation and suction



27293 BD/CD

27293 BD **Operating Sheath, 26 Fr.,**
for continuous irrigation and suction,
with LUER-Lock stopcock, rotating,
color code: black-red

27293 CD **Operating Sheath, 24 Fr.,**
for continuous irrigation and suction,
with LUER-Lock stopcock, rotating,
color code: white-red

27293 BL **Hollow Obturator and Fascial Dilator,**
color code: black-red

27293 CL **Hollow Obturator and Fascial Dilator,**
color code: white-red



27093 GM **Adaptor,** for use with PCN sheaths,
for introducing flexible cysto-urethroscopes
through a PCN sheath



27093 GN **Adaptor,** for use with resectoscope outer sheaths
and PCN sheaths in conjunction with evacuators
and bladder syringes



27040 LB **Adaptor,** for use with Telescopes 27293 AA
and 27292 AMA with resectoscope outer sheaths
of 27040 SL/SD and 27050 SL





27040 SC **Adaptor,** for use with Telescopes 27293 AA
and 27292 AMA with resectoscope outer sheaths
of 27050 SC/SD



27293 AA



HOPKINS® Wide-Angle Straight Forward Telescope 6°,
with angled eyepiece, autoclavable, fiber optic light transmission incorporated with working channel, with LUER-Lock connection for inflow,
color code: green-red




27093 GM



Adaptor,
for use with PCN sheaths, for introducing flexible cysto-urethrosopes through a PCN sheath




6127390



Sealing Cap,
(50/2.6)



6127290



Sealing Cap,
(50/2.2)



Bladder Syringes



27220



27211 LO

27211 LO

TOOMEY Syringe, 50 cc

27212 LO

TOOMEY Syringe, 100 cc



27215 LO

27215 LO

REINER-ALEXANDER Syringe, 75 cc

27215 EL

Spare Barrel, for use with Syringe 27215 LO

27216 LO

REINER-ALEXANDER Syringe, 100 cc

27216 EL

Spare Barrel, for use with Syringe 27216 LO

27218 LO

REINER-ALEXANDER Syringe, 150 cc

27218 EL


Spare Barrel, for use with Syringe 27218 LO

27220

Catheter Adaptor, for use with Syringes 27211 LO – 27218 LO

27213

Adaptor, to connect Syringe 27211 LO to miniature cystoscope-urethroscope sheaths





Locking device "LO"

Syringes 27211 – 27218 and Evacuator 27224 are supplied with locking device "LO" to secure watertight closure between the syringe and the sheath.

3-874



27220



Catheter Adaptor,
for use with Syringes 27211 LO - 27218 LO

