

Certificate of Analysis
RSID™-SALIVA Detection Kits

Product Name: RSID™-Saliva Lab Kits with Dual Buffers or Universal Buffer

Weight: 1.0 pound

Catalog Number: 0130, 0130-LARGE, 0100

Shipping Conditions: Ambient Temperature

Storage Conditions: RSID™-Saliva cassettes should be stored at room temperature. Buffer should be removed from kit and stored at 2-8°C.

Shelf Life: Twenty-one months from production of kit components. Specific expiration dates are printed on cassette pouches and buffer bottle. Buffers have a 2 year shelf life from date of manufacture.

Description: RSID™-Saliva is a lateral flow strip test that detects human Saliva in forensic stains on evidence from extracted swabs or fabric cuttings.

Quality Control: RSID™-Saliva Lab Kit components (cassettes and buffer) are tested for sensitivity, specificity, and reproducibility and each batch of kit components must meet pre-determined quality standards prior to release.

Components:

1. 25 RSID™-Saliva Cassettes – Lot#022124A1 Exp. DEC 2025
2. RSID™ Universal Buffer (30 ml) Lot#UB2403 Exp: MAR 2026
OR RSID™ Dual Buffers – Extraction and Running Buffers
3. One RSID™-Saliva Technical Information and Protocol Sheet

SYNTHETIC REAGENT FOR FORENSIC STAIN ANALYSIS

DOES NOT CONTAIN HUMAN OR ANIMAL PRODUCTS

DOES NOT CONTAIN INFECTIOUS AGENTS OR ARTICLES


DOES NOT CONTAIN BIOLOGICALLY ACTIVE MATERIALS

ALL PRODUCTS ARE FOR "RESEARCH USE ONLY"

NOT FOR IN VITRO DIAGNOSTIC USE

NOT USED FOR HUMAN HEALTH, MEDICAL DIAGNOSTICS, VETERINARY SCIENCE

Approved by:



Karl Reich, Ph.D.

Chief Scientific Officer

Independent Forensics

500 Waters Edge, Suite 210

Lombard IL 60148 USA

p (1) 708.234.1200 f 708.978.5115 skype: karl.reich

www.ifi-test.com

May 10, 2024

MANUFACTURERS DECLARATION

The kits contain test strips and buffers for forensic stain identification.

Product Name: RSID™-Saliva Lab Kit (*for the detection of human saliva*)

Weight: 1.0 pound

Catalog Number: 0100, 0130

Shipping Conditions: Ambient Temperature

Storage Conditions: RSID™-Saliva cassettes should be stored at room temperature. Buffers should be removed from kit and stored at 2-8°C.

Cassette Shelf Life: Twenty-one months from production of kit components. Specific expiration dates are printed on cassette pouches and buffer bottles. Buffer expiration is 24 months from manufacture date.

Quality Control: RSID™-Saliva Lab Kit components (cassettes and buffers) are tested for sensitivity, specificity, and reproducibility and each batch of kit components must meet pre-determined quality standards prior to release.

Components:

- | | | |
|-----------------------------------|--------------|---------------|
| 1. 25 RSID™-Saliva Cassettes | Lot#022124A1 | Exp: DEC 2025 |
| 2. RSID™ Universal Buffer (30 ml) | LOT#UB2403 | Exp: MAR 2026 |

OR RSID™-Saliva Extraction Buffer (25 ml) and RSID™-Saliva Running Buffer (5 ml)

3. One RSID™-Saliva Technical Information and Protocol Sheet

SYNTHETIC REAGENT FOR FORENSIC STAIN ANALYSIS

DOES NOT CONTAIN HUMAN OR ANIMAL PRODUCTS

DOES NOT CONTAIN INFECTIOUS AGENTS OR ARTICLES

DOES NOT CONTAIN BIOLOGICALLY ACTIVE MATERIALS

ALL PRODUCTS ARE FOR "RESEARCH USE ONLY"

NOT FOR IN VITRO DIAGNOSTIC USE

NOT USED FOR HUMAN HEALTH, MEDICAL DIAGNOSTICS, VETERINARY SCIENCE

Approved by:



Karl Reich, Ph.D.

Chief Scientific Officer

Independent Forensics

500 Waters Edge, Suite 210

Lombard IL 60148 USA

p (1) 708.234.1200 f 708.978.5115 skype: karl.reich

www.ifi-test.com
