





# Reliability. Safety. Quality.

Since early 2000, Antmed high pressure syringes and pressure connecting tubes have been used by medical facilities large and small. With seamless compatibility, reliable performance and uncompromised operational safety, cardiologist, radiologists, technicians and nurses all come to trust and prefer Antmed quality products. When our customers count on our products for accurate diagnosis and effective treatments, we make sure our high pressure syringes and connecting tubes are of the highest quality in order to meet their needs and ensure patient safety.

#### **Compatible For Medrad**



- Medrad MCT Plus, Vistron CT, EnVision CT, For CT
- Volume: 200ml



- Medrad Vistron CT, EnVision CT,
  Medrad Mark V, For CathLab For CT
- Volume: 125ml



- Volume: 150ml



- Medrad Mark V, For CathLab
- Volume: 200ml



- Medrad Mark V, For CathLab
- Volume: 60ml



- Medrad Mark IV, For CathLab
- Volume: 130ml



- Medrad Spectris, Spectris Solaris EP, For MR
- Volume: 65ml



- · Medrad Spectris, Spectris Solaris EP, For MR
- Volume: 115ml



- · Imaxeon Biotel, For CT
- Volume: 140ml



- · Imaxeon Visimax, For CT
- Volume: 150ml



#### **Compatible For LF**



 LF CT9000, CT9000ADV, OptiVantage DH, For CT

Volume: 200ml



• LF Angiomat 6000, For CathLab

Volume: 150ml



• LF Angiomat 3000, For CathLab

Volume: 100ml



• LF Angiomat Illumena, For CathLab

Volume: 150ml

#### **Compatible For LF**



• LF Optistar MR, For MR

Volume: 60ml

### Compatible For Medtron



 Medtron Accutron CT, Accutron CT-D, For CT, Accutron HP-D, For CathLab

Volume: 200ml

#### Compatible For EZEM



 EZEM EmpowerCT, EmpowerCTA, For CT

Volume: 200ml



EZEM EmpowerMR, For MR

Volume: 100ml

#### **Compatible For Nemoto**



 Nemoto A-25, A-60, Dual Shot, For CT

Volume: 100ml



Nemoto A-60, Dual Shot, For CT

• Volume: 200ml



Nemoto 120S, For CathLab

Volume: 125ml



 Nemoto Dual Shot, For CT, Sonic Shot 50, For MR

Volume: 60ml



### **High Pressure Syringes**

- Ensure Safety Tight seal, high pressure resistence, stringent bio-compatibility testing
- · Enhance Efficiency High dimensional accuracy, rapid injector identification, easy connect and disconnect of syringes and tubes
- · Meet Customer Needs Full range of syringes, compatible with major injector brands, meet medical facility requirements

| Compatible Brand | Model                  | Compatible Injector                                      |
|------------------|------------------------|--|
| Medrad           | 125ml/200ml            | MCT plus, Vistron CT,<br>EnVision CT                     |
|                  | 60ml/130ml/150ml/200ml | Mark IV, Mark V, Mark V Provis                           |
|                  | 65ml/115ml             | Spectris, Spectris Solaris EP                            |
|                  | 140ml/150ml            | Imaxeon Biotel, Imaxeon Visimax                          |
| LF               | 150ml/200ml            | CT9000 & CT9000 ADV,<br>OptiVantage DH,<br>Angiomat 6000 |
|                  | 100ml/150ml            | Angiomat 3000,<br>Angiomat 6000,<br>Angiomat Illumena    |
|                  | 60ml                   | OptiStar   |
| Nemoto           | 60ml/100ml/200ml       | A-25,<br>A-60,<br>Dual Shot                              |
|                  | 125ml 120S             |  |
|                  | 60ml                   | Sonic Shot 50  |
| EZEM             | 200ml                  | EmpowerCT, EmpowerCTA                                    |
|                  | 100ml                  | EmpowerMR  |
| Medtron          | 200ml                  | Accutron CT, Accutron CT-D                               |
|                  | 200ml                  | Accutron HP-D  |
|                  | 130ml/200ml            | Injektron 82 MRT   |



# Shenzhen Antmed Co., Ltd.

18 Jinhui Ave., Pingshan New District, Shenzhen 518122, China Tel: +86 755 8606 0992 Fax: +86 755 8606 0991 E-mail: info@antmed.com http://www.antmed.com







P/N: EN-HPS-201806







### **Pressure Connecting Tubes**

- Ensure Patient Safety DEHP Free, comprehensive clinical validation and stringent bio-compatibility testing
- Enhance Operational Efficiency Excellent in temperature and pressure resistant, the deformation of tube is smaller, yields great effect for angiography
- A variety of styles of tubing to best fit your needs Lengthened, coiled, straight, T-tube, one-way checkvalve and backflow prevention device. Customization requests can be accommodated



• Tubes for CT (straight & coiled)



Tubes for MR (Y & T type)



• Tubes for CathLab (600 psi)



• Tubes for CathLab (1200 psi)

| Connecting<br>Tube Type    | Length      | Pressure Range<br>MPa(psi) | Spectifications  |
|----------------------------|-------------|----------------------------|--|
| CT Connecting Tube         | 5cm ~ 250cm | 2.4(350)                   | Y & T straight, coiled tube with checkvalve(s), stopcock & spike |
| MR Connecting Tube         | 5cm ~ 250cm | 2.4(350)                   | Y & T straight, coiled tube with checkvalve(s), stopcock & spike |
| CathLab<br>Connecting Tube | 5cm ~ 250cm | 4.2(600)                   | Straight line with rotating male Luer-Lock                       |
|                            | 5cm ~ 250cm | 8.3(1200)                  | Braided line with rotating male Luer-Lock                        |

# antmed

# Shenzhen Antmed Co., Ltd.

18 Jinhui Ave., Pingshan New District, Shenzhen 518122, China Tel: +86 755 8606 0992 Fax: +86 755 8606 0991 E-mail: info@antmed.com http://www.antmed.com







P/N: EN-CT-201806