

Chapter 14 Replace the Accessory

**Warning:**

- This anesthesia machine parts are of latex free.
- The parts which contact with human body are all comply with the biological compatibility of ISO 10993-1, and will not create adverse reaction.
- Only use the accessories specified in this manual otherwise may lead to danger or incorrect measurements.
- Read the packing instruction before use the accessories. The disposable accessories can only be used once, otherwise may lead to cross infection or incorrect performance. The reusable accessories must be cleaned and disinfected before usage to avoid the cross infection.
- To dispose the accessories according to the local laws, regulations or the hospital’s waste disposal regulations to avoid the environment pollution.

Replaced accessories list:

Accessory description	Notes
corrugated pipe	
Silicon corrugated pipe (adult, 100 cm)	
Silicon corrugated pipe (pediatric, 100 cm)	
Adult corrugated pipe sets(including corrugated pipe、Y-connector,L-connector,filter, manual breathing bag)	
Pediatric corrugated pipe set(including corrugated pipe,Y-connector, L-connector, filter, manual breathing bag)	
Manual breathing bag	
Silicon bag (1L)	
silicon bag (2L)	
silicon bag (3L)	
rubber bag (latex free, 1L)	
rubber bag (latex free, 2L)	
rubber bag (latex free, 3L)	
mask	



soft siliconmask (pediatric, 2#)	
siliconmask (pediatric big size, 3#)	
siliconmask (adult, 4#)	
siliconmask (adult big size, 5#)	
inflatable anaestheticmask (pediatric, 3#, disposable)	
inflatable anaestheticmask (pediatric big size, 4#, disposable)	
inflatable anaestheticmask (adult, 5#, disposable)	
soft siliconmask (pediatric, 2#)	
siliconmask (pediatric big size, 3#)	
siliconmask (adult, 4#)	
siliconmask (adult big size, 5#)	
inflatable anaestheticmask (pediatric, 3#, disposable)	
inflatable anaestheticmask (pediatric big size, 4#, disposable)	
inflatable anaestheticmask (adult, 5#, disposable)	
inflatable anaestheticmask (adult big size, 6#, disposable)	
connector	
High temperature resistance L-connector	
High temperature resistance Y-connector	
Flow sensor	
Breathing flow sensor set inspiratory	
Breathing flow sensor set expiratory	
Oxygen sensor	
Oxygen sensor	
Oxygen sensor cable	
soda lime canister	
soda lime canister	
soda lime canister sponge	
Side stream CO ₂ and AG module accessories	
Phase in Nomolinesampling pipe (disposable)	
Main stream CO ₂ and AG module accessories	
Airway adapter (adult, disposable)	
Airway adapter (pediatric, disposable)	
BIS accessories	
BIS patient cable (adult/pediatric)	
BIS sensor (adult)	



BIS sensor (pediatric)	
Module	
Main stream CO ₂ module (IRMA CO ₂)	
Side stream CO ₂ module (ISA CO ₂)	
Main stream AG module (IRMA AX+)	
Side stream AG module (ISA AX+)	
Power cable	
Power cable (0A, 125V/5m)	
battery	
Li-ion battery, 4800mAh 11.1V)	



Chapter 25. Accessories



Warning

- Use the accessories specified in the Operator's Manual; using other accessories may damage the monitor, or cannot reach the performance claimed in this Manual.
- The operating and storage environment of the monitor should meet the requirements of the accessories. Please refer to the manual of the accessories for these requirements.
- Disposable accessories can only be used once, because repeated use can cause performance degradation.
- If the packaging or accessories have any sign of damage, do not use such accessories.

25.1 Standard Accessory List

Standard accessories are as follows:

No.	Name	Specifications & characteristics	Standard quantity
1	5-lead ECG cable	Adult 5-lead all-in-one clip	1
2	Electrode slice	50mm×55mm, adult, disposable	1
3	Oxygen probe	Adult finger clip	1
4	Blood pressure cuff	Adult repetitive single-tube cuff	1
	Blood pressure gas nozzle	Metal	1
	Blood pressure extension tube	3M gray	1
5	2.25K-type temperature probe	Adult body surface	1
6	2600mAh lithium battery	Detachable three-cell lithium battery (14.8V 2600mAh)	1
7	Ground wire	V0020A	1



No.	Name	Specifications & characteristics	Standard quantity
8	Power cord	3×0.75mm ² 10A 250V L=2M	1

25.2 ECG

12-pin integrated 5-lead left/right clamp ECG cable

Supported lead wires	Symbol standard	Type	Applicable objects	Length (main cable)	Length (5-lead ECG cable)
3/5-lead wire	IEC	Defibrillation type	Adult	2.7m	0.9m

25.3 SpO2


Adult finger clip oxygen probe

Applicable objects	Applicable parts	Type	Length
Adult	Fingers, toes	Reusable	3m

25.4 NIBP

1. Blood pressure cuff

Applicable objects	Applicable parts	Type	Remark
Adult	Arm or thigh	Reusable	Single-tube



2. Blood pressure gas nozzle

Operating temperature	Applicable parts	Mechanical life (connecting and disconnecting)	Remark
-30°C~+60°C	Blood pressure cuff	> 1000 times	Pull-off force of gas connector and extension tube is greater than 5kg

3. Blood pressure extension tube

Applicable objects	Type	Length	Remark
Adult, child	Reusable	3m	Used together with cuff

25.5 Temp

2.25K-type temperature probe

Applicable objects	Applicable parts	Type	Length
Adult	Body surface	Reusable	3m

