


# A Sampling of Cannula from our Catalog



Anesthesia		Description	Item #
<b>Retrobulbar (Sharp Point)</b> 10 per box.		.50 x 38 mm (25G x 1 1/2 in)	581275
		.40 x 32 mm (27G x 1 1/4 in)	585109
<b>Retrobulbar/Peribulbar [Atkinson]</b> 10 per box.		.50 x 32 mm (25G x 1 1/4 in)	585108
		.60 x 32 mm (23G x 1 1/4 in)	585001
<b>Retrobulbar [Atkinson]</b> 10 per box.		.50 x 38 mm (25G x 1 1/2 in)	581637
		.60 x 38 mm (23G x 1 1/2 in)	581638
<b>Sub-Tenon's Anesthesia Cannula</b> 10 per box.		1.10 x 25 mm (19G x 1 in)	585176

Irrigation		Description	Item #
<b>Anterior Chamber Cannula [Rycroft]</b> Angled 45 degrees, 4 mm from end. 10 per box.		.30 x 22 mm (30G x 7/8 in)	581273
		.40 x 22 mm (27G x 7/8 in)	581280
		.50 x 22 mm (25G x 7/8 in)	581639
		.60 x 22 mm (23G x 7/8 in)	585025
		.80 x 22 mm (21G x 7/8 in)	581281
<b>Anterior Chamber Cannula [Rycroft]</b> Angled 45 degrees, 6 mm from end. 10 per box.		.30 x 22 mm (30G x 7/8 in)	581285
		.40 x 22 mm (27G x 7/8 in)	581284
		.50 x 22 mm (25G x 7/8 in)	581642
<b>Viscoflow® Cannula</b> Angled 45 degrees, 9 mm from end. 10 per box.		.30 x 22 mm (30G x 7/8 in)	585046
		.40 x 22 mm (27G x 7/8 in)	585026
		.50 x 22 mm (25G x 7/8 in)	585045
		.60 x 22 mm (23G x 7/8 in)	585009

Capsule Polishers		Description	Item #
<b>Flexible Tip Cannula</b> Angled 45 degrees, 7 mm from end, 10 mm sleeve. 10 per box.		.40 x 22 mm (27G x 7/8 in)	585005
		.50 x 22 mm (25G x 7/8 in)	585017
<b>Olive Tip Polisher</b> .33 mm end port, 30 degrees, 10 mm from end. 10 per box.		.60 x 25 mm (23G x 1 in)	585143
<b>Capsule Polisher [Kratz]</b> .30 mm top port, 35 degrees, 9 mm from end. 10 per box.		.60 x 22 mm (23G x 7/8 in)	581605

For a complete list of our cannulas, please refer to our full product catalog.















Hydrodissection		Description	Item #	
<b>Nucleus Hydrodissector [Pearce]</b> Angled 35 degrees, 8 mm from end. 10 per box.			.40 x 22 mm (27G x 7/8 in)	585036
			.50 x 22 mm (25G x 7/8 in)	585037
<b>Nucleus Hydrodissector [Sauter]</b> Tip is flattened on the vertical plane. 10 per box.			.40 x 22 mm (27G x 7/8 in)	585099
<b>Cortical Cleaving Hydrodissector (Flat End)</b> 45 degrees, 11 mm from end. 10 per box.			.40 x 22 mm (27G x 7/8 in)	585158
<b>Nucleus Hydrodissector</b> Tip is flattened on the horizontal plane. 10 per box.			.40 x 22 mm (27G x 7/8 in)	585155
<b>Micro J Shape Cannula</b> .25 mm open end port. 10 per box.			.40 x 22 mm (27G x 7/8 in)	585018
			.50 x 22 mm (25G x 7/8 in)	585003
<b>Nucleus Hydrodissector [Chang]</b> 90 degree bend, 1 mm flattened tip. 10 per box.			.40 x 18 mm (27G x 7/8 in) bevel angle 33 degrees.	585294
			.40 x 18 mm (27G x 7/8 in) bevel angle 47 degrees.	585295
<b>Nucleus Hydrodissector [Akahoshi Style]</b> Tapered tip, thin wall. 10 per box.			.40 x 22 mm (27G x 7/8 in)	585198

Cystotomes		Description	Item #	
<b>Irrigating Cystotome</b> Formed cystotomes are bent to conform with the convexity of the crystalline lens. 10 per box.			.40 x 16 mm (27G x 5/8 in)	581618
			.50 x 16 mm (25G x 5/8 in) [Pearce]	581610
<b>Irrigating Cystotome</b> Formed 8.5 mm from bend to tip. 10 per box.			.40 x 16 mm (27G x 5/8 in)	585008
			.50 x 16 mm (25G x 5/8 in)	585004
<b>Irrigating Cystotome [Blumenthal]</b> Formed 8.5 mm from bend to tip. 10 per box.			.40 x 16 mm (27G x 5/8 in)	585042
<b>Irrigating Cystotome [Simcoe]</b> 45 degrees, 10.5 mm from tip. 10 per box.			.40 x 12.5 mm (27G x 1/2 in)	585093
<b>Irrigating Cystotome [Joyce]</b> Tip .50 mm long. 10 per box.			.40 x 16 mm (27G x 5/8 in)	585039
<b>Irrigating Cystotome</b> 10 per box.			.50 x 16 mm (25G x 5/8 in)	581617









For a complete list of our cannulas, please refer to our full product catalog.









Irrigation & Aspiration		Description	Item #
<b>I and A Cannula [Simcoe]</b> Thin wall cannula. 5 per box			.60 mm/.60 mm (23G/23G) 585047
<b>Bimanual Aspiration Cannula</b> Curved 45 degrees, .40 mm top aspiration port. 10 per box.			.60 x 17 mm (23G x 5/8 in) 585222
<b>Bimanual Irrigation Cannula</b> Thin wall cannula, curved 45 degrees, .51 mm open end irrigation port. 10 per box.			.80 x 17 mm (21G x 5/8 in) 585223

Lacrimal		Description	Item #
<b>Lacrimal Cannula</b> Curved, 28 mm length. 10 per box.			.45 x 28 mm (26G x 1 1/8 in) 581276
<b>Lacrimal Cannula</b> Curved, 35 mm length. 10 per box.			.80 x 35 mm (21G x 1 3/8 in) 581615
<b>Lacrimal Cannula</b> Straight, 13 mm length. 10 per box.			.50 x 13 mm (25G x 1/2 in) 585068

NEW

Refractive		Description	Item #
<b>LASIK Irrigating Cannula</b> 25 degrees, formed 7 mm from end. 10 per box.			.40 x 22 mm (27G x 7/8 in) 585278
<b>LASIK Irrigating Cannula [Solomon]</b> Thin wall cannula. Angled tip features three .25 mm ports on 8 mm extension. 10 per box.			.50 x 22 mm (25G x 7/8 in) 585280
<b>LASIK Irrigating Cannula [Lindstrom]</b> Features three .25 mm ports located on 7 mm extension. 10 per box.			.50 x 22 mm (25G x 7/8 in) 585250
<b>LASIK Flap Spatula/Irrigator [Hersh]</b> 45 degrees, 9 mm from spatulated tip. 10 per box.			.60 x 22 mm (23G x 7/8 in) 585252

Retinal		Description	Item #
<b>Subretinal Fluid Cannula [Grizzard]</b> Sleeve extends 6 mm beyond end. Overall length with tubing 35 mm (1 3/8 in). 5 per box.			.90 x 35 mm (20G x 1 3/8 in) 585021
<b>Soft Tip Aspirator</b> Flexible tip extends 3 mm beyond end. Overall length with tubing 36 mm. 10 per box			.60 x 36 mm (23G x 1 4/10 in) 585226
<b>Soft Tip Aspirator</b> Flexible tip extends 3 mm beyond end. Overall length with tubing 36 mm. 10 per box.			.50 x 36 mm (25G x 1 4/10 in) 585229

NEW

NEW




## *Instrucție de utilizare*

Rev 01 2010-11

• Atenție! Utilizate instrumente ascuțite contaminate și pot purta de agenți patogeni ai bolilor. După utilizare, aruncați-le imediat într-un recipient rezistent la înțepături în conformitate cu standardele acceptate.

• Nu utilizați dacă ambalajul este deteriorat sau deschisă



 Beaver-Visitec International Ltd.  
Waterloo Industrial Estate  
Bidford on Avon Warwickshire,  
B504JHUK  
Beaver-Visitec International, Inc.  
Waltham, MA 02452 USA.

US tel: [866.906.8080](tel:866.906.8080) sau [781.906.8080](tel:781.906.8080)  
[www.Beaaver-Visitec.com](http://www.Beaaver-Visitec.com)

Reprezentant Autorizat  
Lenses Grup SRL  
str. Fierarilor 2,  
Chisinau MD 2001  
Republica Moldova

