





Company JIANGSU HUIDA Profile

Jiangsu Huida Medical Instruments Co.,Ltd is located in the National Environmental Protection Industrial Park-Jiangsu Yancheng Environmental Protection Industrial Park, where it is committed to the research and production of laboratory and medical consumables. It is a large-scale professional manufacturer, producing Microscope slides, Cover glass, plastic labwares, laboratory glassware, sterile surgical atraumatic suture needle, sterile surgical blades, blood lancets, etc.

Established in August 2003, the company has 12,000 square meters of standard workshops and 1,000 square meters of one hundred thousand grade of purifying rooms, designed as the National Industry Rules of Medical Instruments required. We have a team with strong professional skills and the most advanced manufacturing and inspection equipment. We also passed ISO9001, ISO13485, got CE certified and FDA registered. Of our products, 98% is exported to the USA, Canada, the European Union, Latin America, the Middle East, etc. Guided by the market, surviving on the basis of quality and developed on science and technology, Huida aspires to supply the best products and service, becoming the most reliable partners to all customers.

Qualification Certificate



















Product Catalog

| - | T 4" | | C11 1 | - 0 | ~1 |
|-------------------|-------|---------|--------|---------|--------------|
| JSHD | IVI1C | roscone | Slides | & Cover | Glass(01-11) |
| The second second | TATTO | 1000000 | DITUOD | | Oldoblol II |

| Plain Slides | 01 |
|------------------------|-------|
| Frosted Slides | 02-03 |
| Color Frosted Slides | 04 |
| Adhesive Slides | 05-08 |
| Microscope Cover Glass | 09-11 |

Laboratory Plasticware(12-34)

| Embedding Cassettes and Accessories | 13-16 |
|---|-------|
| Slide Storage | 17-18 |
| Petri Dish & Culture Plate | 19-21 |
| Inoculation Loop | 22 |
| Cell Spreader & Transport Swab | 23 |
| Centrifuge Tube | 24-25 |
| Cryovial Tubes \ Urine Collection Tube \ Culture Tube | 26 |
| Test Tube & Test Tube Stopper | 27-28 |
| Racks for Centrifuge Tube, Test Tube and Cryotubes | 28-29 |
| Transfer Pipette | 30 |
| Pipette Tips & Tip Box | 31-32 |
| Slot for Reagent | 32 |
| Specimen Containers | 33-34 |

*Capillary Tubes

Laboratory Glassware(35-60)





Huida microscope slides are manufactured from superior glass, including soda lime glass and super white glass. All slides are pre-cleaned and ready to use. Huida can also produce slides based on your special requirements and needs. A variety of colors, corner styles and edges are available.

Huida 's cover glass is made of high quality transparent glass with no blister or flaws. Advanced technology and equipment make sure that the edges, thickness and dimension of our cover glass is uniform and accurate. All cover glass is specially treated to eliminate sticking. Recommended for automated and manual coverslipping.

Huida cover glass is available in eight standard sizes and four thicknesses. Other unique sizes and circular shapes are also available.

Standard Microscope Slides

Standard Microscope Slides







| Code no. | Edges | Corner | Material | Dimension | Thickness | Packaging |
|------------|---------------|----------------------------------|--|--|------------------|----------------------------|
| 7101 | Ground edges | ☑ 90° ☑ 45° | | | | |
| 7102 | Cut edges | ☑ 90° □ 45° □ 90° ☑ 45° | ☑Soda lime glass ☐Super white glass | ☑ 76 x 26mm ☑ 1" x3" mm ☑ 25 x75mm | ☑1.0mm ☑1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS016-1 | Beveled edges | | | | | |
| 7101A | Ground edges | ☑ 90° ☑ 45° | | | | |
| 7102A | Cut edges | ☑ 90° □ 45° | ☐ Soda lime glass ☑ Super white glass | ☑ 76 x 26mm ☑ 1" x 3" mm ☑ 25 x 75mm | ☑1.0mm ☑1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS016-1A | Beveled edges | □ 90° ☑ 45° | | | | 1 |

Single frosted


| d | Code no. | Edges | Corner | Material | Dimension | Thickness | Packaging |
|---|------------|---------------|----------------------------------|--|--|--------------------|----------------------------|
| ì | 7105 | Ground edges | ☑ 90° ☑ 45° | | | | |
| ۱ | 7105-1 | Cut edges | ☑ 90° □ 45° | ☑Soda lime glass ☐Super white glass | ☑ 76 x 26mm ☑ 1" x3" mm ☑ 25 x75mm | ☑ 1.0mm | ☑ 50pcs/box ☑ 72pcs/box |
| | HDAS016-5 | Beveled edges | □ 90° ☑ 45° | | | | |
| | 7105A | Ground edges | ☑ 90° ☑ 45° | | | | |
| 4 | 7105-1A | Cut edges | ☑ 90° □ 45° □ 90° ☑ 45° | ☐ Soda lime glass ☑ Super white glass | | ☑ 1.0mm ☑ 1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| H | HDAS016-5A | Beveled edges | | | | | |



www.chinahuida.cn

Standard Microscope Slides

Color Frosted Microscope Slides



Double frosted

| Code no. | Edges | Corner | Material | Dimension | Thickness | Packaging |
|------------|---------------|----------------|--|--|--------------------|----------------------------|
| 7107 | Ground edges | ☑ 90° ☑ 45° | | | | |
| 7107-1 | Cut edges | ☑ 90° □ 45° | ☑Soda lime glass ☐Super white glass | ☑ 76 x 26mm ☑ 1" x3" mm ☑ 25 x75mm | ☑ 1.0mm ☑ 1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS016-7 | Beveled edges | □ 90° ☑ 45° | | | | |
| 7107A | Ground edges | ☑ 90° ☑ 45° | | | | |
| 7107-1A | Cut edges | ☑ 90° □ 45° | ☐ Soda lime glass ☑ Super white glass | ☑ 76 x 26mm ☑ 1" x3" mm ☑ 25 x75mm | ☑1.0mm ☑1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS016-7A | Beveled edges | □ 90° ☑ 45° | | | | |

Concaves slides

| Code no. | Concaves | Corner | Material | Edges | Dimension | Thickness | Packaging |
|----------|-----------------|----------------|-------------------|-------------|---------------------------|-----------|-------------|
| 7103 | Single concave | ☑ 90° ☑ 45° | Soda lime glass | Ground edge | ☑76 x 26mm | ☑1.0mm | ☑ 50pcs/box |
| 7104 | Double concaves | ☑ 90° ☑ 45° | Super white glass | Ground edge | ☑ 1" x3" mm ☑ 25 x75mm | ☑1.1mm | ☑ 72pcs/box |
| 7103A | Single concave | ☑ 90° ☑ 45° | Soda lime glass | Ground edge | ☑ 76 x 26mm | ☑ 1.0mm | ☑ 50pcs/box |
| 7104A | Double concaves | ☑ 90° ☑ 45° | Super white glass | Ground edge | ☑ 1" x3" mm ☑ 25 x75mm | ☑1.1mm | ☑ 72pcs/box |




Huida Color Frosted Microscope Slides come with a light opaque coating in six standard colors; white, blue, pink, yellow, green and orange, resistant to commonly used laboratory solvents. Its painted background never change color in routine H&E staining.



- ◆ Frosted one end, one side
- ◆ Suitable for marking with inkjet and thermal transfer printers and permanent markers

| Code no. | Colors | Edges | Material | Corner | Dimension | Thickness | Packaging |
|----------|--------|-------------|---------------------|-----------------------------|--|------------------|----------------------------|
| 7109-W | White | | | | | | |
| 7109-0 | Orange | Ground edge | | | | | ☑ 50pcs/box ☑ 72pcs/box |
| 7109-G | Green | | ☑Soda lime glass | ime glass ☑ 90° ☑ 76 x 26mm | | ☑ 1.0mm | |
| 7109-P | Pink | | ☐ Super white glass | ☑ 45° | ☑ 45 ☑ 1" x3" mm ☑ 25 x75mm | ☑1.1mm | |
| 7109-B | Blue | | | ≥ 25 x/5mm | | | |
| 7109-Y | Yellow | | | | | | |
| 7109A-W | White | | | | ☑ 76 x 26mm ☑ 1" x 3" mm ☑ 25 x 75mm | ☑1.0mm ☑1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| 7109A-O | Orange | | | | | | |
| 7109A-G | Green | Ground edge | ☐ Soda lime glass | ☑ 90° | | | |
| 7109A-P | Pink | | ✓Super white glass | ☑ 30 ☑ 45° | | | |
| 7109A-B | Blue | | | E 40 | | | |
| 7109A-Y | Yellow | | | | | | |






Adhesive Microscope Slides



These Color-Frosted (+) glass slides are made by a new process which offers the following features: They place a permanent positive charge in the microscope slide.


- 1). They electrostatically attracts frozen tissue sections and cytology preparations, binding them to the slide.
- 2). They form a bridge so that covalent bonds develop between formalin fixed sections and the glass.
- 3). Tissue sections and cytological preparations adhere better to the Plus glass slides without the need for special adhesives or protein coatings.



- ◆ Recommended for routine H&E stains, IHC, ISH, frozen sections, cytology smear.
- Suitable for marking with inkjet and thermal transfer printers and permanent markers.

| Code no. | Colors | Edges | Material | Corner | Dimension | Thickness | Packaging |
|------------|--------|-------------|---------------------|-----------------------|---------------------------|---|----------------------------|
| HDAS001-W | White | | | acs ⋈ 90° ⋈ 76 x 26mm | | | |
| HDAS001-0 | Orange | | | | | ☑1.0mm ☑1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS001-G | Green | 4 | ☑Soda lime glass | | ☑ 76 x 26mm | | |
| HDAS001-P | Pink | Ground edge | ☐ Super white glass | M 45° | ☑ 1" x3" mm ☑ 25 x75mm | | |
| HDAS001-B | Blue | | | | | | |
| HDAS001-Y | Yellow | | | | | | |
| HDAS001A-W | White | | | | | | |
| HDAS001A-O | Orange | | | | | | |
| HDAS001A-G | Green | Ground edge | ☐ Soda lime glass | ☑ 90° | ☑76 x 26mm | ₩ 1.0mm | ☑ 50pcs/box |
| HDAS001A-P | Pink | | Super white glass | ☑ 30 ☑ 45° | ☑1" x3" mm | ₫1" x3" mm ☑ 1.0mm ₫25 x75mm ☑ 1.1mm | ☑ 72pcs/box |
| HDAS001A-B | Blue | | | E 40 | ≥ 25 x/5mm | | |
| HDAS001A-Y | Yellow | | | | | | |

Our premium glass slides with a frosted end are now available pre-coated with Polysine, These slides improve the adhesion of tissues to the slide.



- ◆ Recommended for routine H&E stains, IHC, ISH, frozen sections,
- Suitable for marking with inkjet and thermal transfer printers and permanent markers.

| Code no. | Colors | Edges | Material | Corner | Dimension | Thickness | Packaging |
|------------|--------|-------------|---------------------|--|-------------|------------------|----------------------------|
| HDAS002-W | White | | | | | | |
| HDAS002-0 | Orange | Ground edge | ☑Soda lime glass | ☑ 90° ☑ 76 × 26mm ☑ 1° x3" mm ☑ 25 x75mm | | ☑ 1.0mm | |
| HDAS002-G | Green | | | | | | ☑ 50pcs/box |
| HDAS002-P | Pink | | ☐ Super white glass | | ☑ 1.1mm | ☑ 72pcs/box | |
| HDAS002-B | Blue | | | | ≥ 20 X/5mm | | |
| HDAS002-Y | Yellow | | | | | | |
| HDAS002A-W | White | | | | | ☑1.0mm ☑1.1mm | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS002A-O | Orange | | | | | | |
| HDAS002A-G | Green | Ground edge | ☐ Soda lime glass | ☑ 90° | ☑76 x 26mm | | |
| HDAS002A-P | Pink | | ✓ Super white glass | ☑ 30 ☑ 45° | ☑ 1" x3" mm | | |
| HDAS002A-B | Blue | | | E 40 | ☑ 25 x75mm | | |
| HDAS002A-Y | Yellow | | | | | | |

www.chinahuida.cn








This kind of slides is prepared by Silane, to enhance adhesion of histological and plastic sections to Microscope Slides.




- ◆Recommended for routine H&E stains, IHC,ISH, frozen sections.
- ◆Suitable for marking with inkjet and thermal transfer printers and permanent markers.

| Code no. | Colors | Edges | Material | Corner | Dimension | Thickness | Packaging |
|------------|--------|-------------|---------------------|-------------------------|-------------------|------------------|----------------------------|
| HDAS003-W | White | | | | | | |
| HDAS003-0 | Orange | Ground edge | | nlass ☑ 90° ☑ 76 x 26mm | | | ☑ 50pcs/box ☑ 72pcs/box |
| HDAS003-G | Green | | ☑Soda lime glass | | | ☑1.0mm ☑1.1mm | |
| HDAS003-P | Pink | | ☐ Super white glass | ☑ 45° | 7/45° ☑ 1″ x3″ mm | | |
| HDAS003-B | Blue | | | | ☑ 25 x75mm | | |
| HDAS003-Y | Yellow | | | | | | |
| HDAS003A-W | White | | | | | | |
| HDAS003A-O | Orange | | | | | | |
| HDAS003A-G | Green | Ground edge | ☐ Soda lime glass | ☑ 90° | ☑76 x 26mm | ☑ 1.0mm | ☑ 50pcs/box |
| HDAS003A-P | Pink | | ✓ Super white glass | | ☑1" x3" mm | ☑ 1.1mm | ☑ 72pcs/box |
| HDAS003A-B | Blue | | | E 40 | 25 x75mm | | |
| HDAS003A-Y | Yellow | | | | | | |



end, Hydrophilic surface.
Suitable for marking with inkjet and thermal transfer printers and permanent markers.

| Perm | ancn | marke | J. | | | | |
|------------|--------|-------------|---------------------|--------|-------------|-----------|-------------|
| Code no. | Colors | Edges | Material | Corner | Dimension | Thickness | Packaging |
| HDAS005-W | White | | | | | | |
| HDAS005-0 | Orange | | | | | | |
| HDAS005-G | Green | | ☑Soda lime glass | ☑ 90° | ☑ 76 x 26mm | ☑ 1.0mm | ☑ 50pcs/box |
| HDAS005-P | Pink | Ground edge | ☐ Super white glass | ☑ 45° | ☑1" x3" mm | ☑ 1.1mm | ☑ 72pcs/box |
| HDAS005-B | Blue | | | | ☑ 25 x75mm | | |
| HDAS005-Y | Yellow | | | | | | |
| HDAS005A-W | White | | | | | | |
| HDAS005A-O | Orange | | | | | | |
| HDAS005A-G | Green | | ☐ Soda lime glass | ☑ 90° | ☑76 x 26mm | ☑ 1.0mm | ☑ 50pcs/box |
| HDAS005A-P | Pink | Ground edge | ✓Super white glass | ☑ 45° | ☑1" x3" mm | ☑ 1.1mm | ☑ 72pcs/box |
| HDAS005A-B | Blue | | | _ 10 | ☑ 25 x75mm | | |
| HDAS005A-Y | Yellow | | | | | | |



Microscope Cover Glass



Standard sizes of Square and Rectangular Huida cover glass

18 x 18mm 20 x 20mm 22 x 22mm 24 x 24mm

24 x 32mm 24 x 40mm 24 x 50mm 24 x 60mm Thickness:

No. 0: 0.09mm – 0.13mm No. 1: 0.13mm – 0.16mm No. 1 1/2: 0.16mm – 0.19mm No. 2: 0.19mm – 0.23mm

Different packaging is available for square and rectangular cover glass.

1.Two parts box containing 100pcs, 10boxes/vacuum pack

| Dimension | 18x18 | 20x20 | 22x22 | 24x24 |
|-----------|------------|------------|------------|------------|
| Packaging | 100pcs/box | 100pcs/box | 100pcs/box | 100pcs/box |
| Dimension | 24x32 | 24x40 | 24x50 | 24x60 |
| Packaging | 100pcs/box | 100pcs/box | 100pcs/box | 100pcs/box |



2. Hinged box containing 100pcs or 200pc, 5boxes/vacuum pack

| Dimension | 18x18 | 20x20 | 22x22 | 24x24 |
|-----------|--------------------------|--------------------------|--------------------------|--------------------------|
| Packaging | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box |
| Dimension | 24x32 | 24x40 | 24x50 | 24x60 |
| Packaging | 100pcs/box | 100pcs/box | 100pcs/box | 100pcs/box |



3 Hinged box containing 100pcs or 200pcs, 1box/vacuum pack, 10boxes/dispenser box

| Dimension | 18x18 | 20x20 | 22x22 | 24x24 |
|-----------|--------------------------|--------------------------|--------------------------|--------------------------|
| Packaging | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box |
| Dimension | 24x32 | 24x40 | 24x50 | 24x60 |
| Packaging | 100pcs/box | 100pcs/box | 100pcs/box | 100pcs/box |



4.Packing by weight: Hinged box containing 1OZ or 2OZ, 10boxes/dispenser box

| - | | | | |
|-----------|--------------------|--------------------|--------------------|--------------------|
| Dimension | 18x18 | 20x20 | 22x22 | 24x24 |
| Packaging | 1OZ/box 2OZ/box | 1OZ/box 2OZ/box | 1OZ/box 2OZ/box | 1OZ/box 2OZ/box |
| Dimension | 24x32 | 24x40 | 24x50 | 24x60 |
| Packaging | 1OZ/box 2OZ/box | 1OZ/box 2OZ/box | 1OZ/box 2OZ/box | 1OZ/box 2OZ/box |



www.chinahuida.cn

Microscope Cover Glass

5. Outside buckled hinged-lid box containing 100pcs or 200pcs

| Dimension | 18x18 | 20x20 | 22x22 | 24x24 |
|-----------|--------------------------|--------------------------|--------------------------|--------------------------|
| Packaging | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box | 100pcs/box 200pcs/box |
| Dimension | 24x32 | 24x40 | 24x50 | 24x60 |
| Packaging | 100pcs/box | 100pcs/box | 100pcs/box | 100pcs/box |




Cover Glass Circles

Packaging for circular cover glass: Hinged box containing 100pcs





| Dimension | dia.8 | dia.9 | dia.10 | dia.11 | dia.12 | dia.14 | dia.16 | dia.18 | dia.20 | dia.22 | |
|-----------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|--|
| Packaging | 100pcs/box | |



Huida offers different kinds of plastic labwares, which are used in research application for the medical and life sciences. You can find them in laboratories as well as hospitals. Huida disposable plasticwares are manufactured from high-quality raw material such as PP, PS, PE, ABS, AS etc. Sterile and non-sterile products are both available.

Embedding Cassette


Huida produce many kinds of embedding cassettes. Most models can be used with automated labeling machines.






Code no. HP20281

Special hinge hold the base and cover together. The cassettes can be opened as often as possible. 45 °angled writing surface.



Snap-latch and one-piece with attached lid allows for one-hand operation. The lids can be opened and closed as often as possible. 35 °angled writing surface.








Code no.HP2030

Detachable lid. It can be opened and closed as often as possible. 45 °angled writing surface.



Code no.HP20283

Front-hinged design. It can be opened and closed as often as possible. 45 °angled writing surface.

Code no.HP20284

Front-hinged design with four square compartments. It can be opened and closed as often as possible. 45 °angled writing surface.




Code no.HP2028

Snap-latch and one-piece with attached lid allows for one-hand operation. The lids can be opened and closed as often as possible. 45 °angled writing surface.





www.chinahuida.cn



Code no.HP2032

Embedding Rings
Designed for holding and identifying tissue sample blocks and fit well in microtome chuck adapters.




Code no.HP2031

Round holes, without lid. Use with metal lid. 45 °angled writing surface.




Code no.HP20291

Special hinge hold the base and cover together. The cassettes can be opened as often as possible. 45 °angled writing surface.



Code no.HP20292

Snap-latch and one-piece with attached lid allows for one-hand operation. The lids can be opened and closed as often as possible. 35 °angled writing surface.



Snap-latch and one-piece with attached lid allows for one-hand operation. The lids can be opened and closed as often as possible. With fine square pore. 45 °angled writing surface.




Code no.HP2027

Snap-latch and one-piece with attached lid allows for one-hand operation. The lids can be opened and closed as often as possible. With square pore 45 °angled writing surface.



In addition to bulk packaging, taped and sleeved package is also available.








Code NO. HP 2033 HP 2034 HP 2035 HP 2036 HP 2037 Well Dimension Well Dimension Well Dimension Well Dimension Well Dimension 15×15×9mm 24×24×9mm 30×24×9mm 34×24×9mm



Reusable lid for cassette HP2031. Stainless steel Code no.: HPCSC

Biopsy Foam pads



Code no. HPBF-01

Hold biopsies in place and prevent them from being lost during processing.

Cardboard slide mailer

Designed for holding Microscope slides.


| Code NO. | Description | Specification |
|----------|---------------|--|
| HPCM-01 | Slide mailers | Holds 1 piece of slides |
| HPCM-02 | Slide mailers | Holds 2 pieces of slides |
| HPCM-03 | Slide mailers | Holds 3 pieces of slides |
| HPCM-04 | Slide mailers | Holds 4 pieces of slides |
| HPCM-05 | Slide mailers | Holds 5 pieces of slides |
| HPCM-06 | Slide mailers | Holds 6 pieces of slides |
| HPCM-09 | Slide mailers | Holds 9 pieces of slides |
| HPCM-10 | Slide mailers | Holds 10 pieces of slides |
| HPCM-20 | Slide mailers | Holds 20 pieces of slides |
| HPCM-20N | Slide mailers | Holds 20 pieces of slides without cover |
| HPCM-24 | Slide mailers | Holds 24 pieces of slides |
| HPCM-32 | Slide mailers | Holds 32 pieces of slides |



Plastic slide tray

| ı | Code No. | Description | Specification | Material | Packing | |
|---|----------|-----------------------|------------------------|----------|--------------------|--|
| | HP4001 | Staining Jar and Rack | For 24pieces of slide | PP | 100pcs/case | |
| 1 | HP4002 | Slide tray | For 20 pieces of slide | ABS | 1pc/bag,50pcs/case | |





Plastic slide mailer

| Code No. | Description | Specification | Material | Packing |
|----------|---------------|------------------------------|----------|-----------------------------------|
| HP4003 | Slide mailers | For 1 piece of slide, clear | PP | 100pcs/polybag, 1000pcs/case |
| HP4004 | Slide mailers | For 1 piece of slide, blue | PP | 100pcs/polybag, 1000pcs/case |
| HP4005 | Slide mailers | For 2 pieces of slide, clear | PP | 100pcs/polybag, 1000pcs/case |
| HP4006 | Slide mailers | For 2 pieces of slide, blue | PP | 100pcs/polybag, 1000pcs/case |
| HP4006-1 | Slide mailers | For 3 pieces of slide, blue | PP | 50pcs/polybag, 1000pcs/case |
| HP4007 | Slide mailers | For 5 pieces of slide, clear | PP | 100pcs/dispenser box, 800pcs/case |











Plastic slide storage box


| ı | Code No. | Description | Specification | Material | Packing |
|---|----------|-------------------|-------------------------|----------|----------------------------|
| | HP4008 | Slide storage box | For 12 pieces of slide | ABS/PS | 4pcs/polybag, 240 pcs/case |
| | HP4009 | Slide storage box | For 25 pieces of slide | ABS/PS | 2pcs/polybag, 200 pcs/case |
| | HP4010 | Slide storage box | For 50 pieces of slide | ABS/PS | 2pcs/polybag,100 pcs/case |
| ı | HP4011 | Slide storage box | For 100 pieces of slide | ABS/PS | 1pcs/polybag,50 pcs/case |
| | | | | | |

Petri Dish



| Code No. | Specification | Sterile | Material | Packing |
|----------|---------------|---------|----------|-----------------------------|
| HP0001 | 35 × 15mm | EO | PS | 10pcs/polybag, 2000pcs/case |
| HP0002 | 60 × 15mm | EO | PS | 10pcs/polybag, 1000pcs/case |







| Code No. | Specification | Sterile | Material | Packing | |
|----------|---------------------------|---------|----------|----------------------------|--|
| HP0005 | 90 × 15mm for machine use | EO | PS | 10pcs/polybag, 500pcs/case | |




| Code No. | Specification | Sterile | Material | Packing |
|----------|------------------------|---------|----------|----------------------------|
| HP0006 | 90 × 15mm, one room | EO | PS | 10pcs/polybag, 500pcs/case |
| HP0007 | 90 × 15mm, two rooms | EO | PS | 10pcs/polybag, 500pcs/case |
| HP0008 | 90 × 15mm, three rooms | EO | PS | 10pcs/polybag, 500pcs/case |
| HP0009 | 90 × 15mm, four rooms | EO | PS | 10pcs/polybag, 500pcs/case |
| | | | | |




| Code No. | Specification | Sterile | Material | Packing |
|----------|---------------|---------|----------|----------------------------|
| HP0010 | 90×20mm | EO | PS | 10pcs/polybag, 500pcs/case |




| Code No. | Specification | Sterile | Material | Packing |
|----------|---------------|---------|----------|----------------------------|
| HP0011 | 150 × 15mm | EO | PS | 10pcs/polybag, 200pcs/case |





Petri Dish

Inoculation Loop



Packing

10pcs/polybag, 300pcs/case

| Code No. | Description | Specification | Sterile | Material | Packing |
|----------|---------------|----------------------------|---------|----------|--------------------------|
| HP4037 | Culture plate | 96wells, U-shape bottom | Gamma | PS | 1pc/polybag, 100pcs/case |
| HP4038 | Culture plate | 96wells, V-shape bottom | Gamma | PS | 1pc/polybag, 100pcs/case |
| HP4039 | Culture plate | 96wells, Flat-shape bottom | Gamma | PS | 1pc/polybag, 100pcs/case |
| HP4040 | Elisa plate | 96 wells detachable | Gamma | PS | 1pc/polybag, 100pcs/case |







Specification

130*130mm

Code No.

HP0012





| Code No. | Specification | Sterile | Material | Packing |
|----------|---------------------------|----------|----------|---------------------------|
| HP40321 | 10ul with needle,flexible | EO/Gamma | ABS | 1pc/pack, 5000pcs/case |
| HP40322 | 10ul with needle,flexible | EO/Gamma | ABS | 5pcs/pack, 5000pcs /case |
| HP40323 | 10ul with needle,flexible | EO/Gamma | ABS | 10pcs/pack,10000pcs/case |
| HP40324 | 10ul with needle,flexible | EO/Gamma | ABS | 20pcs/pack,10000pcs/case |
| HP40331 | 1ul with needle,flexible | EO/Gamma | ABS | 1pc/pack, 5000pcs/case |
| HP40332 | 1ul with needle,flexible | EO/Gamma | ABS | 5pcs/pack, 5000pcs /case |
| HP40333 | 1ul with needle,flexible | EO/Gamma | ABS | 10pcs/pack, 10000pcs/case |
| HP40334 | 1ul with needle,flexible | EO/Gamma | ABS | 20pcs/pack, 10000pcs/case |
| | | | | |

| Code No. | Specification | Sterile | Material | Packing |
|----------|-------------------------|----------|----------|--------------------------|
| HP40341 | 1ul+10ul, rigid | EO/Gamma | AS | 1pc/pack, 5000pcs/case |
| HP40342 | 1ul+10ul, rigid | EO/Gamma | AS | 5pcs/pack, 5000pcs /case |
| HP40343 | 1ul+10ul, rigid | EO/Gamma | AS | 10pcs/pack,10000pcs/cas |
| HP40344 | 1ul+10ul, rigid | EO/Gamma | AS | 20pcs/pack,10000pcs/cas |
| HP40351 | 1ul with needle, rigid | EO/Gamma | AS | 1pc/pack, 5000pcs/case |
| HP40352 | 1ul with needle, rigid | EO/Gamma | AS | 5pcs/pack, 5000pcs /case |
| HP40353 | 1ul with needle, rigid | EO/Gamma | AS | 10pcs/pack,10000pcs/cas |
| HP40354 | 1ul with needle, rigid | EO/Gamma | AS | 20pcs/pack,10000pcs/cas |
| HP40361 | 10ul with needle, rigid | EO/Gamma | AS | 1pc/pack, 5000pcs/case |
| HP40362 | 10ul with needle, rigid | EO/Gamma | AS | 5pcs/pack, 5000pcs /case |
| HP40363 | 10ul with needle, rigid | EO/Gamma | AS | 10pcs/pack,10000pcs/cas |
| HP40364 | 10ul with needle, rigid | EO/Gamma | AS | 20pcs/pack,10000pcs/cas |
| | | | | |

| Code No. | Specification | Sterile | Material | Packing |
|----------|---------------------------|----------|----------|--------------------------|
| HP50011 | 10ul with needle,flexible | EO/Gamma | ABS | 1pc/pack, 5000pcs/case |
| HP50012 | 10ul with needle,flexible | EO/Gamma | ABS | 5pcs/pack, 5000pcs /case |
| HP50013 | 10ul with needle,flexible | EO/Gamma | ABS | 10pcs/pack,10000pcs/case |
| HP50014 | 10ul with needle,flexible | EO/Gamma | ABS | 20pcs/pack,10000pcs/case |
| HP50021 | 1ul with needle,flexible | EO/Gamma | ABS | 1pc/pack, 5000pcs/case |
| HP50022 | 1ul with needle,flexible | EO/Gamma | ABS | 5pcs/pack, 5000pcs /case |
| HP50023 | 1ul with needle,flexible | EO/Gamma | ABS | 10pcs/pack,10000pcs/case |
| HP50024 | 1ul with needle,flexible | EO/Gamma | ABS | 20pcs/pack,10000pcs/case |









Cell Spreader / Transport Swab

Centrifuge Tube


Cell Spreader

| Code No. | Specification | Sterile | Material | Packing |
|----------|------------------|-----------|----------|-------------------------|
| HP50031 | L-shape, sterile | E.O/Gamma | PP | 1pc/pack, 2000pcs/case |
| HP50032 | L-shape, sterile | E.O/Gamma | PP | 5pcs/pack,3000pcs /case |
| HP50033 | L-shape, sterile | E.O/Gamma | PP | 10pcs/pack,3000pcs/case |
| HP50041 | T-shape, sterile | E.O/Gamma | PP | 1pc/pack, 2000pcs/case |
| HP50042 | T-shape, sterile | E.O/Gamma | PP | 5pcs/pack,3000pcs /case |
| HP50043 | T-shape, sterile | E.O/Gamma | PP | 10pcs/pack,3000pcs/case |



Transport swab

| Code No. | Stick Material | Head Material | Tube material | Sterile | Packing |
|----------|----------------|---------------|---------------|-----------|-----------------------|
| HP5005 | wooden | Cotton | pp | E.O/Gamma | In bulk, 2000pcs/case |
| HP5006 | wooden | Viscose | pp | E.O/Gamma | In bulk, 2000pcs/case |
| HP5007 | snapple PS | Cotton | рр | E.O/Gamma | In bulk, 2000pcs/case |
| HP5008 | snapple PS | Viscose | pp | E.O/Gamma | In bulk, 2000pcs/case |



| Code No. | Stick Material | Head Material | Tube material | Sterile | Packing |
|----------|----------------|---------------|---------------|-----------|---------------------------|
| HP5009 | wooden | Cotton | pp | E.O/Gamma | 1pc/polybag, 1000pcs/case |
| HP5010 | wooden | Viscose | рр | E.O/Gamma | 1pc/polybag, 1000pcs/case |
| HP5011 | snapple PS | Cotton | pp | E.O/Gamma | 1pc/polybag, 1000pcs/case |
| HP5012 | snapple PS | Viscose | pp | E.O/Gamma | 1pc/polybag, 1000pcs/case |
| | | | | | |


| | onappio i o | 0011011 | PP | | .babaijaagi isoobaaiaaa |
|----------|----------------|---------------|---------------|-----------|-----------------------------|
| HP5012 | snapple PS | Viscose | pp | E.O/Gamma | 1pc/polybag, 1000pcs/case |
| Code No. | Stick Material | Head Material | Tube material | Sterile | Packing |
| HP5013 | wooden | Cotton | pp | E.O/Gamma | 1pc/peel pack, 1000pcs/case |
| HP5014 | wooden | Viscose | pp | E.O/Gamma | 1pc/peel pack, 1000pcs/case |
| HP5015 | snapple PS | Cotton | pp | E.O/Gamma | 1pc/peel pack, 1000pcs/case |
| HP5016 | snapple PS | Viscose | DD | E.O/Gamma | 1pc/peel pack, 1000pcs/case |



| Code No. | Stick Material | Head Material | Tube material | Medium | Sterile | Packing |
|----------|----------------|---------------|---------------|-----------------|---------|-----------------------------|
| HP5017 | snapple PS | Viscose | рр | Amies | Gamma | 1pc/peel pack, 1000pcs/case |
| HP5018 | snapple PS | Viscose | pp | Amies+charcoal | Gamma | 1pc/peel pack, 1000pcs/case |
| HP5019 | snapple PS | Viscose | pp | Cary Blair | Gamma | 1pc/peel pack, 1000pcs/case |
| HP5020 | snapple PS | Viscose | рр | Stuart | Gamma | 1pc/peel pack, 1000pcs/case |
| HP5021 | snapple PS | Viscose | pp | Stuart+charcoal | Gamma | 1pc/peel pack, 1000pcs/case |
| | | | | | | |




| A Code No. | Specification | Sterile | Material | Packing |
|------------|--|---------|----------|----------------------------|
| HP1001 | 50ml self-standing bottom with screw cap | Non | PP | 25pcs/polybag, 500pcs/case |
| HP10011 | 50ml self-standing bottom with screw cap | EO | PP | 25pcs/polybag, 500pcs/case |
| HP10012 | 50ml self-standing bottom with screw cap | EO | PP | 1pc/polybag, 500pcs/case |
| HP10014 | 50ml self-standing bottom with screw cap | EO | PP | 25pcs/rack, 500pcs/case |



| (| Code No. | Specification | Sterile | Material | Packing |
|---|----------|------------------------------------|---------|----------|-----------------------------|
| | HP1004 | 15ml conical bottom with screw cap | Non | PP | 50pcs/polybag, 1000pcs/case |
| | HP10041 | 15ml conical bottom with screw cap | EO | PP | 50pcs/polybag, 1000pcs/case |
| | HP10042 | 15ml conical bottom with screw cap | EO | PP | 1pc/polybag, 1000pcs/case |
| | HP10044 | 15ml conical bottom with screw cap | EO | PP | 50pcs/rack, 1000pcs/case |
| | | | | | |

| Code No. | Specification | Sterile | Material | Packing |
|----------|------------------------------------|---------|----------|-----------------------------|
| HP1005 | 10ml conical bottom with screw cap | Non | PP | 50pcs/polybag, 1000pcs/case |
| HP10051 | 10ml conical bottom with screw cap | EO | PP | 50pcs/polybag, 1000pcs/case |
| HP10052 | 10ml conical bottom with screw cap | EO | PP | 1pc/polybag, 1000pcs/case |


| Code No. | Specification | Material | Packing |
|----------|---------------|----------|------------------------------|
| HP1006 | 10ml | PP | 200pcs/polybag, 3200pcs/case |
| HP1007 | 7ml | PP | 200pcs/polybag, 4000pcs/case |
| HP1008 | 5ml | PP | 300pcs/polybag, 6000pcs/case |

















| Code No. | Specification | Material | Packing |
|----------|-----------------------|----------|--------------------------------|
| HP1009 | 0.2ml | PP | 1000pcs/polybag, 70000pcs/case |
| HP1010 | 0.5ml | PP | 1000pcs/polybag, 25000pcs/case |
| HP1011 | 1.5ml | PP | 500pcs/polybag, 10000pcs/case |
| HP1012 | 1.5ml with shield lid | PP | 500pcs/polybag, 10000pcs/case |
| HP1013 | 2ml | PP | 500pcs/polybag, 10000pcs/case |

| CodeNo. | Specification | Material | Packing |
|---------|----------------------|----------|------------------------------|
| HP1014 | 1.5ml with screw cap | PP | 500pcs/polybag,10000pcs/case |
| HP1015 | 1.5ml with screw cap | PP | 500pcs/polybaq.12000pcs/case |





| CodeNo. | Description | Specificaion |
|---------|-------------|------------------|
| HP4058 | PCR Tube | 8 strips, 0.2ml |
| HP4059 | PCR Tube | 12 strips, 0.2ml |




Cryovial Tube







| Code No. | Specification | Material | Packing |
|----------|---------------|----------|------------------------------|
| HP1016 | 5ml | PP | 200pcs/polybag,3200pcs/case |
| HP1017 | 1.8ml | PP | 500pcs/polybag,10000pcs/case |
| HP10171 | 10ml | PP | 100pcs/polybag,10000pcs/case |





Urine collection tube


| Code No. | Specification | Material | Packing |
|----------|---------------|----------|------------------------------------|
| HP1019 | 12ml | PS | 250pcs/dispenser box, 2000pcs/case |



Culture tube

| Code No. | Specification | Material | Packing |
|----------|---------------|----------|------------------------------|
| HP1020 | 15×65mm | PP | 500pcs/polybag, 5000pcs/case |






www.chinahuida.cn



| Code No. | Specification | Material | Packing |
|-------------------|---------------|----------|------------------------------|
| A HP1021 | ø12×60mm | PP | 500pcs/polybag, 5000pcs/case |
| B HP1022 | ø12×60mm | PS | 500pcs/polybag, 5000pcs/case |
| ● HP1023 | ø12×75mm | PP | 500pcs/polybag, 5000pcs/case |
| D HP1024 | ø12 × 75mm | PS | 500pcs/polybag, 5000pcs/case |
| B HP1025 | ø13 × 75mm | PP | 500pcs/polybag, 5000pcs/case |
| P HP1026 | ø13 × 75mm | PS | 500pcs/polybag, 5000pcs/case |
| G HP1027 | ø13 × 100mm | PP | 150pcs/polybag, 3000pcs/case |
| HP1028 | ø13 × 100mm | PS | 150pcs/polybag, 3000pcs/case |
| I HP1029 | ø15 × 100mm | PP | 150pcs/polybag, 2500pcs/case |
| U HP1030 | ø15 × 100mm | PS | 150pcs/polybag, 2500pcs/case |
| () HP1031 | ø16 × 100mm | PP | 500pcs/polybag, 2500pcs/case |
| HP1032 | ø16 × 100mm | PS | 500pcs/polybag, 2500pcs/case |
| | | | |

| Code No. | Description | Specification | Material | Packing |
|----------|-------------------|---------------|----------|------------------------------|
| HP1033 | Conical Test Tube | ø16 × 102mm | PP | 150pcs/polybag, 2400pcs/case |
| HP1034 | Conical Test Tube | ø16 × 102mm | PS | 150pcs/polybag, 2400pcs/case |



Test Tube Stopper

| Code No. | Specification | Material | Packing |
|----------|---------------|----------|--------------------------------|
| HP1035 | ø12mm | PE | 1000pcs/polybag, 25000pcs/case |
| HP1036 | ø13mm | PE | 1000pcs/polybag, 25000pcs/case |
| HP10361 | ø13mm | PE | 1000pcs/polybag, 25000pcs/case |
| HP1037 | ø16mm | PE | 500pcs/polybag, 16000pcs/case |
| HP1038 | ø12mm | PE | 1000pcs/polybag, 25000pcs/case |
| HP1039 | ø13mm | PE | 1000pcs/polybag, 25000pcs/case |
| HP1040 | ø16mm | PE | 1000pcs/polybag, 25000pcs/case |






Centrifuge tube box







| Code No. | Specification | Material | Packing | |
|----------|------------------------------------|----------|-------------|--|
| HP4012 | 72wells for 0.5ml Centrifuge tube | PP | 84pcs/case | |
| HP4013 | 72wells for 1.5ml Centrifuge tube | PP | 84pcs/case | |
| HP4014 | 96wells for 0.2ml Centrifuge tube | PP | 250pcs/case | |
| HP40131 | 100wells for 1.5ml Centrifuge tube | PP | 84pcs/case | |
| | | | | |




Centrifuge tube rack


| Code No. | Specification | Material | Packing |
|----------|---|----------|-------------|
| HP4015 | 60wells for 0.5ml/1.5ml centrifuge tube, double use | PP | 200pcs/case |
| HP4016 | 96wells for 0.5ml/1.5ml centrifuge tube, double use | PP | 100pcs/case |

| Code No. | Specification | Material | Packing |
|----------|---------------------------|----------|-------------|
| HP40151 | 15ml/18wells+50ml/10wells | ABS | 200pcs/case |
| HP40161 | 15ml/30wells+50ml/20wells | ABS | 100pcs/case |









Freezing box





| HP4020 | |
|--------|--|

| Code No. | Specification | Material | Packing |
|----------|-----------------------------|----------|-------------|
| HP4018 | 50wells for 1.8ml Cryotubes | PP | 150pcs/case |
| HP4019 | 81well for 1.8ml Cryotubes | PP | 100pcs/case |
| HP4020 | 100well for 1.8ml Cryotubes | PP | 90pcs/case |
| HP4021 | 100well for 1.8ml Cryotubes | PP | 90pcs/case |



Multi-purpose tube racks

| Code No. | Specification | Material | Packing |
|----------|------------------|----------|------------|
| HP4022 | ø10-18mm/50wells | ABS | 50pcs/case |



Detachable test tube racks

| Code No. | Material | Specification | Packing |
|----------|----------|---------------|------------|
| HP4023 | PS | ø13/50wells | 50pcs/case |
| HP4024 | PS | ø15/50wells | 50pcs/case |
| HP4025 | PS | ø18/40wells | 50pcs/case |




Test tube racks


| Code No. | Material | Specification | Packing |
|----------|----------|------------------|---------------------------|
| HP40231 | PP | ø12-13mm/40wells | 5pcs/poly bag, 50pcs/case |
| HP40241 | PP | ø15-17mm/60wells | 5pcs/poly bag, 50pcs/case |
| HP40251 | PP | ø20-21mm/90wells | 5pcs/poly bag, 50pcs/case |



Transfer Pipette

| | Code No. | Specification | Material | Serile | Packing |
|---|-------------------|--------------------|----------|------------|---|
| A | HP2002 | 65mm | PE | Non | In bulk,10000pcs/case |
| | HP2006 | 1ml 145mm | PE | Non | In bulk, 5000pcs/case |
| | HP20061 | 1ml 145mm, sterile | PE | E.O | 1pc/peel pack, 3000pcs/case |
| 0 | HP20062 | 1ml 145mm, sterile | PE | E.O | 1pc/polybag, 3000pcs/case |
| 0 | HP20063 | 1ml 145mm, sterile | PE | E.O | 5pcs/peel pack, 3000pcs/cas |
| | HP20064 | 1ml 145mm, sterile | PE | E.O | 10pcs/peel pack, 3000pcs/cas |
| | HP20065 | 1ml 145mm, sterile | PE | E.O | 20pcs/peel pack, 3000pcs/cas |
| | 1100007 | 4-1400 | DE | 100 | L- L- III - FAAA/ |
| | HP2007 | 1ml 160mm | PE | Non | In bulk, 5000pcs/case |
| | HP20071 | 1ml 160mm, sterile | PE | E.O | 1pc/peel pack, 3000pcs/case |
| 0 | HP20072 | 1ml 160mm, sterile | PE | E.O | 1pc/polybag, 3000pcs/case |
| | HP20073 | 1ml 160mm, sterile | PE | E.0 | 5pcs/peel pack, 3000pcs/cas |
| | HP20074 | 1ml 160mm, sterile | PE | E.O | 10pcs/peel pack, 3000pcs/cas |
| | HP20075 | 1ml 160mm, sterile | PE | E.O | 20pcs/peel pack, 3000pcs/cas |
| | HP2008 | 2ml 150mm | PE | Non | In bulk, 5000pcs/case |
| | HP20081 | 2ml 150mm, sterile | PE | E.O | 1pc/peel pack, 3000pcs/case |
| | HP20082 | 2ml 150mm, sterile | PE | E.O | 1pc/polybag, 3000pcs/case |
| 0 | HP20083 | 2ml 150mm, sterile | PE | E.O | 5pcs/peel pack, 3000pcs/cas |
| | HP20084 | 2ml 150mm, sterile | PE | E.O | 10pcs/peel pack, 3000pcs/cas |
| | HP20085 | 2ml 150mm, sterile | PE | E.O | 20pcs/peel pack, 3000pcs/cas |
| | LIDOOOO | 3ml 160mm | PE | Non | In halls F000/ |
| | HP2009 HP20091 | 3ml 160mm, sterile | PE | Non E.O | In bulk, 5000pcs/case |
| | HP20091 | 3ml 160mm, sterile | PE | E.O | 1pc/peel pack, 3000pcs/case |
| 0 | HP20092 | 3ml 160mm, sterile | PE | E.O | 1pc/polybag, 3000pcs/case |
| | HP20093 | 3ml 160mm, sterile | PE | E.O | 5pcs/peel pack, 3000pcs/cas 10pcs/peel pack, 3000pcs/cas |
| | HP20095 | 3ml 160mm, sterile | PE | E.0 | 20pcs/peel pack, 3000pcs/cas |
| | | 200,000,000 | | | |
| | HP2010 | 3ml 160mm | PE | Non | In bulk, 5000pcs/case |
| | HP20101 | 3ml 160mm, sterile | PE | E.O | 1pc/peel pack, 3000pcs/case |
| a | HP20102 | 3ml 160mm, sterile | PE | E.O | 1pc/polybag, 3000pcs/case |
| • | HP20103 | 3ml 160mm, sterile | PE | E.O | 5pcs/peel pack, 3000pcs/cas |
| | HP20104 | 3ml 160mm, sterile | PE | E.O | 10pcs/peel pack, 3000pcs/cas |
| | HP20105 | 3ml 160mm, sterile | PE | E.O | 20pcs/peel pack, 3000pcs/cas |
| | HP5030 | 3ml 178mm | PE | Non | In bulk, 5000pcs/case |
| | HP50301 | 3ml 178mm, sterile | PE | E.O | 1pc/peel pack, 3000pcs/case |
| 0 | HP50302 | 3ml 178mm, sterile | PE | E.O | 1pc/polybag, 3000pcs/case |
| 0 | HP50303 | 3ml 178mm, sterile | PE | E.O | 5pcs/peel pack, 3000pcs/cas |
| | HP50304 | 3ml 178mm, sterile | PE | E.O | 10pcs/peel pack, 3000pcs/cas |
| | | | | | |














www.chinahuida.cn



| | Code No. | Description | Specification | Material | Packing |
|---|----------|---|---------------|----------|--------------------------------|
| 2 | HP2011 | Tips for Eppendorf Pipette | 10ul White | PP | 1000pcs/polybag, 50000pcs/case |
| | HP2012 | Tips for Eppendorf Pipette | 200ul Yellow | PP | 1000pcs/polybag, 50000pcs/case |
| | HP2013 | Tips for Eppendorf Pipette | 1000ul Blue | PP | 500pcs/polybag, 15000pcs/case |
| | Code No. | Description | Specification | Material | Packing |
| | HP2014 | Tips for Gilson Pipette,universal | 10ul White | PP | 1000pcs/polybag,140000pcs/case |
| | HP2015 | Tips for Gilson Pipette,universal | 200ul Yellow | PP | 1000pcs/polybag, 50000pcs/case |
| | HP2016 | Tips for Gilson Pipette,universal | 1000ul Blue | PP | 500pcs/polybag,15000pcs/case |
| | HP20141 | Tips for Gilson Pipette,universal,Graduated | 10ul White | PP | 1000pcs/polybag,140000pcs/case |
| | HP20151 | Tips for Gilson Pipette,universal,Graduated | 200ul Yellow | PP | 1000pcs/polybag, 50000pcs/case |
| | HP20161 | Tips for Gilson Pipette,universal,Graduated | 1000ul Blue | PP | 500pcs/polybag,15000pcs/case |
| | HP2017 | Tips for Finland Pipette,universal | 10ul White | PP | 1000pcs/polybag,140000pcs/case |
| | HP2018 | Tips for Finland Pipette,universal | 200ul White | PP | 1000pcs/polybag, 50000pcs/case |
| | HP2019 | Tips for Finland Pipette,universal | 1000ul White | PP | 500pcs/polybag,15000pcs/case |
| | HP2020 | Tips for Gilson Pipette | 5000ul White | PP | 300pcs/polybag,7500pcs/case |
| | HP2021 | Tips for Finland Pipette | 5000ul White | PP | 300pcs/polybag,4800pcs/case |
| | HP2022 | Tips for Dalong Pipette | 5000ul White | PP | 300pcs/polybag,4000pcs/case |
| | HP2023 | Filter tips for Gilson | 10ul White | PP | 1000pcs/polybag,140000pcs/case |
| | HP2024 | Filter tips for Gilson | 200ul Yellow | PP | 1000pcs/polybag, 50000pcs/case |
| | HP2025 | Filter tips for Finland | 200ul White | PP | 1000pcs/polybag, 50000pcs/case |
| | HP2026 | Filter tips for Gilson | 1000ul Blue | PP | 500pcs/polybag,15000pcs/case |
| | | | | | |
















Tip box





| Code No. | Description | Specification | Material | Packing |
|----------|-------------|--------------------------|----------|-------------|
| HP4027 | Tip box | 96wells for 10ul tips | PP | 170pcs/case |
| HP4028 | Tip box | 96wells for 200ul tips | PP | 120pcs/case |
| HP4029 | Tip box | 96wells for 200ul tips | PP | 100pcs/case |
| HP4030 | Tip box | 60wells for 1000ul tips | PP | 100pcs/case |
| HP4031 | Tip box | 100wells for 1000ul tips | PP | 70pcs/case |

Slot for reagent

| Code No. | Specification | Material | Sterile | Packing |
|----------|---------------|----------|---------|----------------------------|
| HP4026 | 50ml | PP | Non | 20pcs/polybag, 800pcs/case |
| HP40261 | 50ml | PP | EO | 1pc/polybag, 200pcs/case |






Specimen Containers



Urine container


| Code No. | Specification | Sterile | Material | Packing |
|------------------|---------------|---------|----------|---------------------------|
| HP3001 | 120ml | Non | PP | 200pcs/case |
| A HP30011 | 120ml | Non | PP | 2*250lids+500cups/case |
| HP30012 | 120ml | E.O | PP | 1pc/polybag, 200pcs/case |
| HP3002 | 100ml | Non | PP | 200pcs/case |
| B HP30021 | 100ml | Non | PP | 2*250lids+500cups/case |
| HP30022 | 100ml | E.O | PP | 1pc/polybag, 200pcs/case |
| HP3003 | 80ml | Non | PP | 500pcs/case |
| HP30031 | 80ml | Non | PP | 4*250lids+1000cups/case |
| HP30032 | 80ml | E.O | PP | 1pc/polybag, 500pcs/case |
| HP3004 | 60ml | Non | PP | 500pcs/case |
| HP30041 | 60ml | E.O | PP | 1pc/polybag, 500pcs/case |
| 3 HP3005 | 50ml | Non | PP | 500pcs/case |
| HP30051 | 50ml | E.O | PP | 1pc/polybag, 500pcs/case |
| HP3006 | 40ml | Non | PP | 800pcs/case |
| P HP30061 | 40ml | Non | PP | 6*250lids+1500cups/case |
| HP30062 | 40ml | E.O | PP | 1pc/polybag, 400pcs/case |
| HP3007 | 30ml | Non | PS | 50pcs/polybag, 500pcs/cas |
| HP30071 | 30ml | E.O | PS | 1pc/polybag, 500pcs/case |
| HP30072 | 30ml | Non | PP | 50pcs/polybag, 500pcs/cas |
| HP30073 | 30ml | E.O | PP | 1pc/polybag, 500pcs/case |





0


0





Sputum container

| Code No. | Specification | Sterile | Material | Packing |
|-----------|---------------|---------|----------|--------------------------|
| HP3021 | 20ml | Non | PP | 1000pcs/case |
| A HP30211 | 20ml | Non | PP | 4*250lids+1000cups/case |
| HP30212 | 20ml | E.O | PP | 1pc/polybag, 500pcs/case |
| HP3020 | 30ml | Non | PP | 1000pcs/case |
| B HP30201 | 30ml | Non | PP | 4*250lids+1000cups/case |
| HP30202 | 30ml | E.O | PP | 1pc/polybag, 500pcs/case |
| HP3011 | 40ml | Non | PP | 800pcs/case |
| ● HP30111 | 40ml | Non | PP | 6*250lids+1500cups/case |
| HP30112 | 40ml | E.O | PP | 1pc/polybag,400pcs/case |
| | | | | |





| Code No. | Specification | Material | Packing |
|----------|---------------|----------|------------------------------|
| HP3012 | 30ml | PP | 100pcs/polybag, 5000pcs/case |
| HP3022 | 30ml | PP | 100pcs/polybag, 5000pcs/case |
| HP3021 | 90ml | PP | 1440pcs/case |

| Code No. | Description | Specification | Material | Packing |
|----------|-------------|---------------|----------|----------------------------|
| HP3013 | Hitachi Cup | 16×38mm | PS | 500pcs/bag, 5000pcs/case |
| HP3014 | Beckman Cup | 13×24mm | PS | 1000pcs/bag, 10000pcs/case |
| HP3015 | 700Cup | 14×25mm | PS | 1000pcs/bag, 10000pcs/case |

| | 115010 | 700Cup | 14^2011111 F3 | 1000 | pcs/bag, 10000pcs/case |
|---|----------|-----------------|----------------------------------|----------|--------------------------|
| 1 | Code No. | Description | Specification | Material | Packing |
| | HP3016 | Urine container | With needle, W/lid, labeled | PP | 200pcs/case |
| | HP3017 | Urine tube | Vacuum Urine tube conical bottom | PET | 100pcs/rack, 800pcs/case |







*Capillary Tubes

- 1). Hematocrite capillary tubes without heparine
- 2). Hematocrite capillary tubes with sodium heparine 80IU/mL

Vacuum Urine tube round bottom

Technical specifications and dimensions:

- 1,1 1,2mm internal diameter
- 1,5 1,6mm external diameter



Product Presentation:

Our micro-haematocrit capillary tubes are made of glass and are available in different lengths and diameters. The standard length is 75 mm. The standard inside diameter is 1.1-1.2 mm and the standard outside diameter is 1.5-1.6 mm. 100 capillaries are packed in glass vials, 10 vials are in a carton. The capillaries are available without heparine and sodium-heparinized. HDA micro-haematocrit tubes are quality products that conform to international standard ISO 12772. The tubes are supplied fire polished at the color-coded end, Supplied with color-coding on the packaging and on the tubes in plain (blue color), or with uniform sodium heparinisation (red colour)

100pcs/rack, 800pcs/case



www.chinahuida.cn

Stool container



| IIUI | eaker Low form with spout and prin | nted graduations |
|--------------|------------------------------------|------------------|
| Capacity(ml) | O.D.(mm) | Heigh(mm) |
| 5 | 22 | 30 |
| 10 | 26 | 35 |
| 25 | 34 | 50 |
| 50 | 42 | 60 |
| 100 | 51 | 70 |
| 150 | 60 | 80 |
| 200 | 65 | 88 |
| 250 | 70 | 95 |
| 300 | 80 | 110 |
| 400 | 80 | 110 |
| 500 | 87 | 118 |
| 600 | 90 | 125 |
| 800 | 100 | 135 |
| 1000 | 106 | 145 |
| 2000 | 130 | 185 |
| 3000 | 150 | 210 |
| 5000 | 170 | 270 |
| 10000 | 217 | 350 |

Laboratory Glassware

Huida produces a variety of laboratory glasswares with different specifications, which are transparent, heat-resistant, and easy to customize. Our laboratory glassware mainly work in chemistry and biology

laboratories. They can be used for volumetric measuring, holding, mixing or storing chemicals, heating,

In 20°C

600 ml


| 1102 | Beaker Tall form with spout and prin | ted graduations |
|--------------|--------------------------------------|-----------------|
| Capacity(ml) | O.D.(mm) | Height(mm) |
| 25 | 30 | 55 |
| 50 | 38 | 70 |
| 100 | 48 | 80 |
| 150 | 54 | 95 |
| 250 | 60 | 120 |
| 400 | 70 | 130 |
| 500 | 75 | 140 |
| 600 | 80 | 150 |
| 800 | 90 | 175 |
| 1000 | 95 | 185 |
| 2000 | 120 | 240 |
| 3000 | 135 | 280 |
| | | |






www.chinahuida.cn

cooling, distillation, extraction etc.














| 1115-1 | Boiling Flask Round bottom | n, Wide neck with beaded r | rim |
|--------------|----------------------------|----------------------------|------------|
| Capacity(ml) | O.D. of neck(mm) | O.D. of bulb(mm) | Height(mm) |
| 50 | 34 | 51 | 105 |
| 100 | 34 | 64 | 110 |
| 250 | 50 | 85 | 145 |
| 500 | 50 | 105 | 168 |
| 1000 | 50 | 131 | 210 |
| 2000 | 76 | 166 | 260 |

| Capacity(ml) | O.D. of neck(mm) | O.D. of bulb(mm) | Height(mm) |
|--------------|------------------|------------------|------------|
| 25 | 32 | 42 | 70 |
| | | | |
| 50 | 34 | 51 | 85 |
| 100 | 34 | 64 | 105 |
| 200 | 50 | 79 | 131 |
| 250 | 50 | 85 | 140 |
| 300 | 50 | 87 | 156 |
| 500 | 50 | 105 | 175 |
| 1000 | 50 | 131 | 220 |
| 2000 | 72 | 153 | 275 |
| 3000 | 106 | / 190 | 285 |


| 1121 | Conical flask, (Erlenmeyer F | Flask), Narrow neck, with g | raduations |
|--------------|------------------------------|-----------------------------|------------|
| Capacity(ml) | O.D. of neck(mm) | O.D. of bulb(mm) | Height(mm) |
| 5 | 14 | 25 | 50 |
| 10 | 16 | 30 | 55 |
| 25 | 22 | 42 | 75 |
| 50 | 22 | 51 | 90 |
| 100 | 22 | 64 | 105 |
| 150 | 28 | 71 | 120 |
| 200 | 34 | 79 | 135 |
| 250 | 34 | 85 | 145 |
| 300 | 34 | 87 | 160 |
| 500 | 34 | 105 | 180 |
| 1000 | 42 | 131 | 220 |
| 2000 | 50 | 166 | 280 |
| 3000 | 50 | 187 | 310 |
| 5000 | 50 | 220 | 365 |














| 1122 | Conical flask With ground-in glass stopper | |
|--------------|--|------------|
| Capacity(ml) | O.D.of bulb(mm) | Height(mm) |
| 50 | 50 | 90 |
| 100 | 60 | 108 |
| 250 | 83 | 140 |
| 500 | 98 | 185 |
| 1000 | 130 | 230 |



| 1131 | Distilling Flask Wi | th side tube | | |
|--------------|---------------------|------------------|----------------------|------------|
| Capacity(ml) | O.D. of neck(mm) | O.D. of bulb(mm) | O.D.of Side tube(mm) | Height(mm) |
| 25 | 22 | 41 | 110 | 46 |
| 50 | 22 | 51 | 120 | 46 |
| 100 | 22 | 64 | 150 | 58 |
| 125 | 23 | 70 | 150 | 58 |
| 250 | 34 | 85 | 200 | 76 |
| 500 | 34 | 105 | 250 | 96 |
| 1000 | 34 | 131 | 300 | 112 |
| 2000 | 46 | 175 | 410 | 120 |



| 1161 | Retort With ground-in glass stopper | |
|--------------|-------------------------------------|-------------------------|
| Capacity(ml) | O.D. of bulb(mm) | Length of Side neck(mm) |
| 250 | 86 | 225 |
| 500 | 105 | 300 |
| 1000 | 130 | 340 |



Laboratory Glassware

| 1171 Evapora | ting Dish Flat bottom with spout |
|----------------------|----------------------------------|
| O.D.of top(mm) | Height(mm) |
| 60 | 30 |
| 90 | 45 |
| 120 | 60 |
| 150 | 75 |
| 1172 Evapora | ting Dish Flat bottom with spout |
| O.D.of top(mm) | Height(mm) |
| 60 | 30 |
| 90 | 45 |
| | |
| 120 | 60 |
| O.D.of top(mm) 60 | Height(mm) 30 45 |











| 0.88 | |
|-----------------|----------------------|
| 1173 Crystalliz | ring Dish with spout |
| O.D.of top(mm) | Height(mm) |
| 50 | 38 |
| 60 | 40 |
| 70 | 42 |
| 80 | 43 |
| 90 | 45 |
| 100 | 50 |
| 125 | 63 |
| 150 | 75 |
| 180 | 90 |
| 200 | 110 |
| 230 | 115 |
| | |

| 1175 | Watch Glas |
|----------|------------|
| Diar | meter(mm) |
| A | 45 |
| 17 | 50 |
| | 60 |
| 1 | 70 |
| 1 | 80 |
| 1 | 90 |
| 4 | 100 |
| 16 | 120 |
| | 150 |
| | 180 |
| | |

| 1177 Pe | tri Culture Dish | |
|----------------|------------------|------------|
| Cover dia(mm) | Bottom dia.(mm) | Height(mm) |
| 65 | 60 | 15 |
| 82 | 75 | 15 |
| 98 | 90 | 18 |
| 107 | 100 | 20 |
| 130 | 120 | 25 |
| 160 | 150 | 30 |

www.chinahuida.cn



| 1213 | Condenser Liebig with fused i | nner tube | |
|----------------------|-------------------------------|-------------------|------------|
| Length of jacket(mm) | O.D.of.Upper tube(mm) | Cooling area(dm2) | Length(mm) |
| 200 | 18 | 0.7 | 350 |
| 250 | 18 | 1 | 420 |
| 300 | 20 | 1.1 | 480 |
| 400 | 22 | 1.6 | 590 |
| 500 | 24 | 2.2 | 710 |
| 600 | 24 | 2.8 | 810 |
| 800 | 28 | 4.3 | 1050 |
| 1000 | 28 | 5.7 | 1250 |



| 1214 co | ondenser Allihn with bulbed | inner tube | |
|----------------------|-----------------------------|-------------------|------------|
| Length of jacket(mm) | O.D.of.Upper tube(mm) | Cooling area(dm2) | Length(mm) |
| 200 | 18 | 1 | 350 |
| 250 | 18 | 1.5 | 420 |
| 300 | 20 | 1.8 | 480 |
| 400 | 22 | 2.7 | 590 |
| 500 | 24 | 3.7 | 710 |
| 600 | 24 | 4.9 | 800 |
| 800 | 28 | 7.6 | 1050 |
| 1000 | 28 | 10 | 1250 |




| 1215 co | ondenser Graham with coile | d inner tube | |
|----------------------|----------------------------|-------------------|------------|
| Length of jacket(mm) | O.D.of.Upper tube(mm) | Cooling area(dm2) | Length(mm) |
| 200 | 18 | 2.3 | 350 |
| 250 | 20 | 2.7 | 420 |
| 300 | 20 | 3.3 | 480 |
| 400 | 22 | 4.7 | 590 |
| 500 | 24 | 7 | 710 |
| 600 | 24 | 8.2 | 810 |






| 1230 Test cultu | re tube with screw cap | 1231 Test tub | e without rim plain |
|------------------------|------------------------|----------------------|---------------------|
| O.D.of. tube(mm) | Length (mm) | O.D.of. tube(mm) | Length (mm) |
| 12 | 75 | 12 | 75 |
| 12 | 100 | 12 | 100 |
| 13 | 100 | 13 | 100 |
| 15 | 100 | 15 | 100 |
| 15 | 150 | 15 | 150 |
| 16 | 100 | 16 | 100 |
| 16 | 150 | 16 | 150 |
| 16 | 160 | 16 | 160 |
| 18 | 150 | 18 | 150 |
| 18 | 180 | 18 | 180 |
| 20 | 125 | 20 | 125 |
| 20 | 150 | 20 | 150 |
| 25 | 200 | 25 | 125 |
| | | | |













| 1236 | Test Tube With graduations | | |
|--------------|----------------------------|------------------|------------|
| Capacity(ml) | Graduation(ml) | O.D.of. tube(mm) | Length(mm) |
| 5 | 0.1 | 11 | 110 |
| 10 | 0.2 | 14 | 130 |
| 15 | 0.2 | 15 | 150 |
| 20 | 0.5 | 17 | 165 |
| 25 | 0.5 | 18 | 180 |
| 30 | 0.5 | 19 | 180 |
| 50 | 1 | 23 | 200 |



www.chinahuida.cn

Laboratory Glassware

| 1241 | Centrifuge Tube Conical | |
|--------------|-------------------------|------------|
| Capacity(ml) | O.D.of. tube(mm) | Length(mm) |
| 5 | 14 | 90 |
| 10 | 17 | 110 |
| 15 | 19 | 120 |
| 25 | 23 | 140 |
| 50 | 33 | 150 |



| 1242 | Centrifuge Tube Conical, with graduation | | |
|--------------|--|------------------|------------|
| Capacity(ml) | Graduation(ml) | O.D.of. tube(mm) | Length(mm) |
| 5 | 0.1 | 13 | 90 |
| 10 | 0.1 | 17 | 110 |
| 15 | 0.1 | 19 | 120 |
| 25 | 0.2 | 23 | 140 |
| 50 | 0.2 | 33 | 150 |











| 1301 Weighting Bottle tall form O.D.(mm) Height(mm) 25 25 30 20 25 25 30 25 25 25 40 35 25 25 25 60 30 40 20 20 30 30 40 25 30 30 40 25 30 30 30 40 30 30 30 40 30 3 |
|--|
| O.D.(mm) Height(mm) 25 |
| 25 |
| 25 30 30 25 25 40 35 25 25 50 35 70 25 60 40 20 30 30 40 25 30 40 30 30 50 40 70 30 60 45 25 30 70 45 30 40 50 50 30 45 70 60 30 45 70 60 30 45 70 60 35 45 70 60 35 50 50 330 70 40 70 40 70 60 35 60 35 60 35 60 40 70 60 35 60 35 60 40 70 60 35 60 35 60 40 70 60 40 60 35 60 40 70 60 40 60 35 60 40 70 60 40 60 40 60 40 70 60 40 60 35 60 40 70 60 40 |
| 25 |
| 25 50 35 70 25 60 40 20 30 30 40 25 30 40 40 30 30 50 40 70 30 60 45 25 30 70 45 30 40 50 50 30 45 70 60 30 45 70 60 35 45 70 60 35 45 70 60 35 45 70 30 60 40 45 35 50 30 46 35 50 30 47 70 70 70 70 70 48 70 70 70 70 70 49 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 50 70 70 70 70 70 70 50 70 70 70 70 70 70 70 50 70 70 70 70 70 70 70 70 70 70 70 70 70 |
| 25 60 40 20 30 30 40 25 30 40 30 30 30 30 40 40 30 30 30 30 40 40 30 30 30 30 40 40 70 30 45 25 30 40 50 30 45 30 30 45 70 60 30 60 35 60 40 40 30 30 30 30 30 30 30 30 30 30 30 30 30 |
| 30 |
| 30 |
| 30 50 40 70 30 60 45 25 30 70 45 30 40 50 50 30 45 70 60 30 60 35 60 35 60 40 60 35 60 40 60 35 60 40 60 40 60 35 60 40 70 30 70 30 70 35 500 60 40 |
| 30 60 45 25 30 40 40 50 40 50 30 45 70 60 30 60 35 60 40 40 70 30 70 35 500 60 400 70 40 |
| 30 70 45 30 40 50 30 40 50 30 45 70 60 30 60 35 60 40 70 30 70 35 500 60 40 70 40 |
| 40 50 50 30 30 45 70 60 30 60 35 60 35 60 35 60 40 70 30 70 35 500 60 40 70 40 |
| 45 70 60 30 60 35 60 35 60 40 60 35 60 35 60 40 70 30 70 30 70 35 70 35 70 40 |
| 60 35 60 40 60 35 60 40 70 30 70 35 500 60 40 70 35 70 40 |
| 1311 Gas Drying Tower apacity(ml) O.D.of body(mm) Height(mm) 250 50 330 500 60 400 60 40 70 30 70 35 70 40 |
| apacity(ml) O.D.of body(mm) Height(mm) 250 50 330 500 60 400 70 30 70 35 70 35 70 40 |
| 250 50 330 70 35 500 60 400 70 40 |
| 500 60 400 70 40 |
| 30×60 50×30 |
| 50x30 |
| |
| 1311 |
| 1351 Desiccator with porcelain plate, clear |
| I.D.(mm) Height of plate(mm) Total height(mm) |
| 120 180 110 |
| 150 220 135 |
| 180 260 160 |
| 210 275 180 |
| 240 325 210 |
| 300 325 235 |
| 350 400 255 |
| 400 460 300 |





| 1352 | Desiccator with porcelain plate, amber | | |
|----------|--|------------------|--|
| I.D.(mm) | Height of plate(mm) | Total height(mm) | |
| 120 | 180 | 110 | |
| 150 | 220 | 135 | |
| 180 | 260 | 160 | |
| 210 | 275 | 180 | |
| 240 | 325 | 210 | |
| 300 | 355 | 235 | |
| 350 | 400 | 235 | |
| 400 | 460 | 300 | |

| 1354 Vacuu | m Desiccator with ground-in stoped | ock and porcelain plate, clear |
|-------------------|------------------------------------|--------------------------------|
| I.D.(mm) | Height of plate(mm) | Total height(mm) |
| 150 | 220 | 135 |
| 180 | 260 | 160 |
| 210 | 275 | 180 |
| 240 | 325 | 210 |
| 300 | 355 | 235 |
| 350 | 400 | 235 |
| 400 | 460 | 300 |

| 1355 Vacuum D | esiccator with ground-in stopco | ck and porcelain plate, amber |
|----------------------|---------------------------------|-------------------------------|
| I.D.(mm) | Height of plate(mm) | Total height(mm) |
| 150 | 220 | 135 |
| 180 | 260 | 160 |
| 210 | 275 | 180 |
| 240 | 325 | 210 |
| 300 | 355 | 235 |
| 350 | 400 | 255 |
| 400 | 460 | 300 |



1381 Alcohol Lamp with plastic cap Capacity(ml) Height(mm) 150 118 250 130





| 1401 | Reagent Bottle clear, narrow mo | outh, With groud-in glass stop | per or plastic stopper |
|--------------|---------------------------------|--------------------------------|------------------------|
| Capacity(ml) | Dia.of bottle(mm) | Dia.of mouth(mm) | Height(mm) |
| 30 | 40 | 18 | 76 |
| 60 | 46 | 22 | 85 |
| 125 | 57 | 24 | 110 |
| 250 | 70 | 27 | 135 |
| 500 | 85 | 33 | 172 |
| 1000 | 106 | 38 | 202 |
| 2500 | 145 | 48 | 270 |
| 5000 | 185 | 58 | 338 |
| 10000 | 225 | 68 | 430 |
| 20000 | 265 | 82 | 520 |



| 1402 | Reagent Bottle amber, narrow me | outh, With groud-in glass sto | pper or plastic stopper |
|--------------|---------------------------------|-------------------------------|-------------------------|
| Capacity(ml) | Dia.of bottle(mm) | Dia.of mouth(mm) | Height(mm) |
| 30 | 40 | 18 | 76 |
| 60 | 46 | 22 | 85 |
| 125 | 57 | 24 | 110 |
| 250 | 70 | 27 | 135 |
| 500 | 85 | 33 | 172 |
| 1000 | 106 | 38 | 202 |
| 2500 | 145 | 48 | 270 |
| 5000 | 185 | 58 | 338 |
| 10000 | 225 | 68 | 430 |
| | | | |



| 1403 | Reagent Bottle clear, wide mouth, With groud-in glass stopper or plastic | | |
|--------------|--|------------------|------------|
| Capacity(ml) | Dia.of bottle(mm) | Dia.of mouth(mm) | Height(mm) |
| 30 | 40 | 25 | 72 |
| 60 | 46 | 30 | 80 |
| 125 | 57 | 38 | 108 |
| 250 | 70 | 50 | 130 |
| 500 | 85 | 58 | 165 |
| 1000 | 106 | 65 | 188 |
| 2500 | 145 | 90 | 260 |
| 5000 | 185 | 110 | 330 |
| 10000 | 225 | 125 | 415 |
| 20000 | 288 | 122 | 492 |



| Capacity(ml) | Dia.of bottle(mm) | Dia.of mouth(mm) | Height(mm) |
|--------------|-------------------|------------------|------------|
| 30 | 40 | 25 | 72 |
| 60 | 46 | 30 | 80 |
| 125 | 57 | 38 | 108 |
| 250 | 70 | 50 | 130 |
| 500 | 85 | 58 | 165 |
| 1000 | 106 | 65 | 188 |
| 2500 | 145 | 90 | 260 |
| 5000 | 185 | 110 | 330 |
| 10000 | 225 | 125 | 415 |
| | | | |



| 1408 | Reagent bottle(Media bott | tle) with plastic blue screw | cap, clear |
|--------------|---------------------------|------------------------------|------------|
| Capacity(ml) | O.D.of body(mm) | Height(mm) | GL |
| 100 | 56 | 100 | 45 |
| 250 | 70 | 138 | 45 |
| 500 | 86 | 176 | 45 |
| 1000 | 101 | 225 | 45 |
| 2000 | 136 | 260 | 45 |
| 5000 | 186 | 330 | 45 |
| 10000 | 234 | 410 | 45 |
| 20000 | 299 | 505 | 45 |



| 1408-1 | Reagent bottle(Media bo | ottle) with plastic blue scre | w cap, amber |
|--------------|-------------------------|-------------------------------|--------------|
| Capacity(ml) | O.D.of body(mm) | Height(mm) | GL |
| 100 | 56 | 100 | 45 |
| 250 | 70 | 138 | 45 |
| 500 | 86 | 176 | 45 |
| 1000 | 101 | 225 | 45 |
| 2000 | 136 | 260 | 45 |
| 5000 | 186 | 330 | 45 |
| 10000 | 234 | 410 | 45 |
| 20000 | 299 | 505 | 45 |



| 1421 | Aspirator Bottle | | |
|--------------|------------------|-------------------------|------------|
| Capacity(ml) | O.D.of neck(mm) | O.D. of lower mouth(mm) | Height(mm) |
| 2500 | 48 | 30 | 270 |
| 5000 | 58 | 32 | 345 |
| 10000 | 68 | 35 | 420 |
| 20000 | 82 | 38 | 500 |
| | | | |



| Capacity(ml) | O.D.of side tube(mm) | Height(mm) |
|--------------|----------------------|------------|
| 125 | 9 | 120 |
| 250 | 9 | 160 |
| 500 | 9 | 210 |
| 1000 | 9 | 245 |
| 2500 | 11 | 290 |
| 5000 | 11 | 370 |
| 10000 | 13 | 450 |
| 15000 | 15 | 500 |



| 1451 | Dropping Bottle with ground-in pipette and latex rubber n | | |
|--------------|---|------------|--|
| Capacity(ml) | O.D. of tube(mm) | Height(mm) | |
| 30 | 40 | 76 | |
| 60 | 46 | 85 | |
| 125 | 57 | 110 | |




| 1452 | Dropping Bottle Amber glass with ground-in pip | pette and latex rubber ni |
|--------------|--|---------------------------|
| Capacity(ml) | O.D. of tube(mm) | Height(mm) |
| 30 | 40 | 76 |
| 60 | 46 | 85 |
| 125 | 57 | 110 |

1478 Dyeing Jar square form, For 5 pcs











| Diameter(mm) | O.D. of stem(mm) | Length(mm) |
|--------------|------------------|------------|
| 40 | 7 | 50 |
| 50 | 7 | 50 |
| 60 | 8 | 60 |
| 75 | 9 | 75 |
| 90 | 11 | 90 |
| 100 | 13 | 100 |
| 120 | 15 | 120 |
| 150 | 17 | 150 |
| 180 | 18 | 150 |
| 210 | 21 | 160 |
| 240 | 24 | 180 |
| 300 | 26 | 200 |



| 1543/8543 | Separatory Funnel Squib pear shape with ground - in glass sto | opper and GLASS/PTFE stopcock |
|--------------|---|-------------------------------|
| Capacity(ml) | Dia. of stem × Length (mm) | Length(mm) |
| 10 | 8×60 | 126 |
| 20 | 8×60 | 135 |
| 60 | 10×70 | 220 |
| 125 | 10×70 | 260 |
| 250 | 11×70 | 290 |
| 500 | 12×70 | 330 |
| 1000 | 12×70 | 400 |



| 1601H | Measuring Cylinder with glass | s hexagonal base, with graduation | |
|--------------|-------------------------------|-----------------------------------|------------|
| Capacity(ml) | graduation(ml) | Capacity tolerance(±ml) | Height(mm) |
| 5 | 0.1 | 0.1 | 115 |
| 10 | 0.2 | 0.2 | 135 |
| 25 | 0.5 | 0.5 | 160 |
| 50 | 1 | 0.5 | 195 |
| 100 | 1 | 1 | 240 |
| 250 | 2 | 2 | 320 |
| 500 | 5 | 5 | 380 |
| 1000 | 10 | 10 | 465 |
| 2000 | 20 | 20 | 565 |




| 1601 | Measurin | g Cylinder with spout and g | graduations, with glass round base | |
|-------|-----------|-----------------------------|------------------------------------|-------------|
| Capac | city (ml) | graduation (ml) | Capacity tolerance(± ml) | Height (mm) |
| 5 | 5 | 0.1 | 0.1 | 110 |
| 1 | 0 | 0.2 | 0.2 | 135 |
| 2 | 5 | 0.5 | 0.5 | 160 |
| 5 | 0 | 1 | 0.5 | 195 |
| 10 | 00 | 1 | 1 | 250 |
| 25 | 50 | 2or5 | 2 | 300 |
| 50 | 00 | 5 | 5 | 350 |
| 10 | 00 | 10 | 10 | 430 |
| 20 | 00 | 20 | 20 | 500 |
| | | | | |



| 1601P | Measurin | g Cylinder with plastic hexa | gonal base, with graduation | |
|-------|----------|------------------------------|-----------------------------|------------|
| Capa | city(ml) | graduation(ml) | Capacity tolerance(±ml) | Height(mm) |
| | 5 | 0.1 | 0.1 | 110 |
| 1 | 0 | 0.2 | 0.2 | 135 |
| 2 | 5 | 0.5 | 0.5 | 160 |
| 5 | 60 | 1 | 0.5 | 195 |
| 10 | 00 | 1 | 1 | 250 |
| 25 | 50 | 2or5 | 2 | 300 |
| 50 | 00 | 5 | 5 | 350 |
| 10 | 000 | 10 | 10 | 430 |
| 20 | 000 | 20 | 20 | 500 |
| | | | | |









| 1611 | Graduate with spout, short lines | | |
|--------------|----------------------------------|--------------------------|------------|
| Capacity(ml) | graduation(ml) | Capacity tolerance(± ml) | Height(mm) |
| 5 | 0.5 | 0.3 | 80 |
| 10 | 1 | 0.5 | 95 |
| 20 | 2 | 0.6 | 115 |
| 25 | 2 | 0.6 | 118 |
| 50 | 5 | 1 | 135 |
| 100 | 10 | 1.5 | 165 |
| 250 | 25 | 3 | 205 |
| 500 | 50 | 6 | 250 |
| 1000 | 50 | 10 | 310 |
| | | | |





Laboratory Glassware



| 1621A | Volumetric Flask grade A with ground | in glass stopper or plastic st | opper, clear |
|--------------|--------------------------------------|--------------------------------|--------------|
| Capacity(ml) | Capacity tolerance(± ml) | Ground mouth | Height(mm) |
| 5 | 0.02 | 7/11 | 74 |
| 10 | 0.02 | 7/11 | 90 |
| 25 | 0.03 | 10/13 | 110 |
| 50 | 0.05 | 10/13 | 140 |
| 100 | 0.1 | 12/14 | 170 |
| 200 | 0.15 | 14/15 | 210 |
| 250 | 0.15 | 14/15 | 220 |
| 500 | 0.25 | 16/16 | 260 |
| 1000 | 0.4 | 19/17 | 310 |
| 2000 | 0.6 | 24/20 | 370 |



| 1622A | Volumetric Flask amber, grade A, with | ground-in glass stopper or p | lastic stopper |
|--------------|---------------------------------------|------------------------------|----------------|
| Capacity(ml) | Capacity tolerance(± ml) | Ground mouth | Height(mm) |
| 10 | 0.02 | 7/11 | 90 |
| 25 | 0.03 | 10/13 | 110 |
| 50 | 0.05 | 10/13 | 140 |
| 100 | 0.1 | 12/14 | 170 |
| 200 | 0.15 | 14/15 | 210 |
| 250 | 0.15 | 14/15 | 220 |
| 500 | 0.25 | 16/16 | 260 |
| 1000 | 0.4 | 19/17 | 310 |



| 1630A | Measuring Pipette | | | |
|--------------|-------------------|--------------------------|--------------|------------|
| Capacity(ml) | graduation(ml) | Capacity tolerance(± ml) | Flow rate(s) | Length(mm) |
| 0.1 | 0.01 | 0.006 | 2-5 | 320 |
| 0.2 | 0.01 | 0.006 | 2-5 | 320 |
| 0.5 | 0.02 | 0.006 | 2-5 | 320 |
| 1 | 0.01 | 0.008 | 4-10 | 350 |
| 2 | 0.02 | 0.012 | 4-12 | 350 |
| 5 | 0.05 | 0.025 | 6-14 | 350 |
| 10 | 0.1 | 0.05 | 7-17 | 350 |
| 20 | 0.1 | 0.1 | 11-21 | 470 |
| 25 | 0.1 | 0.01 | 11-21 | 470 |
| 50 | 0.2 | 0.01 | 15-25 | 480 |



| 1633 | 633 Volumetric Pipette with one mark | | |
|--------------|--------------------------------------|--------------|------------|
| Capacity(ml) | Capacity tolerance(±ml) | Flow rate(s) | Length(mm) |
| 1 | 0.007 | 39641 | 330 |
| 2 | 0.01 | 39641 | 340 |
| 3 | 0.015 | 15-25 | 360 |
| 5 | 0.015 | 15-25 | 400 |
| 10 | 0.02 | 20-30 | 450 |
| 15 | 0.025 | 20-30 | 520 |
| 20 | 0.03 | 25-35 | 560 |
| 25 | 0.03 | 25-35 | 570 |
| 50 | 0.05 | 30-40 | 600 |
| 100 | 0.08 | 35-45 | 640 |








| 7/87 | Water Determination Apparatus with ground-in glass joints | | |
|---------|---|-----------------------------------|--|
| Art No. | Spare parts | | |
| 1781-1 | 500ml round | boiling flask | |
| 1781-2 | 10ml receiving | ng tube | |
| 1781-3 | 400mm strai | ght condenser | |
| | | | |
| 1902 | Filter Funnel (Bu | uchner funnel), with fritted disc | |
| | Capacity(ml) | Porosity | |
| | 30 | 3 | |
| | 30 | 4 | |
| | 70 | 3 | |
| | 70 | 4 | |
| | 125 | 3 | |
| | 125 | 4 | |
| | 250 | 3 | |
| | 250 | 4 | |
| | 500 | 3 | |
| | 500 | 4 | |



| 424 Burette with stra | aight stopcock | | |
|-----------------------|----------------|-------------------------|------------|
| Capacity(ml) | graduation(ml) | Capacity tolerance(±ml) | Length(mm) |
| 10 | 0.05 | 0.025 | 600 |
| 25 | 0.1 | 0.05 | 630 |
| 50 | 0.1 | 0.05 | 860 |
| 100 | 0.2 | 0.1 | 860 |


www.chinahuida.cn

| 1908 | 1908 Filtration apparatus (Filtration Assembly, Vacuum filtration system) aluminum clamp or stainless steel clamp, with fritted disc. | |
|------|---|--------------------|
| | Funnel(ml) | Conical flask (ml) |
| | 10 | 125 |
| | 15 | 125 |
| | 25 | 125 |
| | 300 | 250 |
| | 300 | 500 |
| | 300 | 1000 |
| | 300 | 2000 |
| | 300 | 3000 |
| | 300 | 5000 |
| | | |













| 5004 | Flask Oblique shape, with | three necks, Standard | Ground Mouth |
|--------------|---------------------------|-----------------------|--------------|
| Capacity(ml) | Middle ground mouth | | Height(mm) |
| 250 | 24/29 | 19/26 | 139 |
| 250 | 24/29 | 24/29 | 139 |
| 250 | 29/32 | 14/23 | 145 |
| 150 | 29/32 | 19/26 | 145 |
| 500 | 19/26 | 14/23 | 158 |
| 500 | 19/26 | 19/26 | 158 |
| 500 | 24/29 | 14/23 | 161 |
| 500 | 24/29 | 19/26 | 161 |
| 500 | 24/29 | 24/29 | 161 |
| 500 | 29/32 | 14/23 | 167 |
| 500 | 29/32 | 19/26 | 167 |
| 1000 | 19/26 | 14/23 | 188 |
| 1000 | 19/26 | 19/26 | 180 |
| 1000 | 24/29 | 14/23 | 196 |
| 1000 | 24/29 | 19/26 | 196 |
| 1000 | 24/29 | 24/29 | 196 |
| 1000 | 29/32 | 14/23 | 200 |
| 1000 | 29/32 | 19/26 | 200 |
| 1000 | 29/32 | 24/29 | 200 |
| 1000 | 29/32 | 29/32 | 200 |
| 2000 | 24/29 | 24/29 | 231 |
| 2000 | 29/32 | 19/26 | 234 |
| 2000 | 29/32 | 24/29 | 234 |
| 2000 | 29/32 | 29/32 | 234 |
| 3000 | 29/32 | 24/29 | 281 |
| 3000 | 34/35 | 24/29 | 284 |
| 5000 | 40/38 | 24/29 | 333 |
| 5000 | 50/42 | 24/29 | 337 |



| <i>5007</i> | Boiling Flask flat bo | ttom, short neck Standa | rd ground mouth |
|--------------|-----------------------|-------------------------|-----------------|
| Capacity(ml) | Ground mouth | O.D.of bulb(mm) | Height(mm) |
| 25 | 14/23 | 41 | 70 |
| 50 | 14/23 | 53 | 88 |
| 50 | 19/26 | 53 | 94 |
| 100 | 14/23 | 65 | 100 |
| 100 | 19/26 | 65 | 106 |
| 100 | 24/29 | 65 | 109 |
| 150 | 19/26 | 76 | 117 |
| 150 | 24/29 | 76 | 120 |
| 250 | 19/26 | 88 | 134 |
| 250 | 24/29 | 88 | 139 |
| 500 | 19/26 | 110 | 158 |
| 500 | 24/29 | 110 | 161 |
| 500 | 29/32 | 110 | 167 |
| 1000 | 24/29 | 140 | 196 |
| 1000 | 29/32 | 140 | 199 |
| 2000 | 29/32 | 175 | 234 |
| 2000 | 34/35 | 175 | 237 |
| 5000 | 34/35 | 235 | 320 |



| 5009 | Conical Flask St | andard ground mouth | |
|--------------|------------------|---------------------|------------|
| Capacity(ml) | Ground mouth | O.D.of bottle(mm) | Height(mm) |
| 5 | 10/19 | 25 | 50 |
| 5 | 14/23 | 25 | 50 |
| 10 | 10/19 | 32 | 60 |
| 10 | 14/23 | 32 | 60 |
| 10 | 19/26 | 32 | 60 |
| 25 | 10/19 | 42 | 65 |
| 25 | 14/23 | 42 | 70 |
| 25 | 19/26 | 42 | 75 |
| 50 | 14/23 | 51 | 85 |
| 50 | 19/26 | 51 | 85 |
| 50 | 24/29 | 51 | 90 |
| 100 | 14/23 | 64 | 105 |
| 100 | 19/26 | 64 | 105 |
| 100 | 24/29 | 64 | 105 |
| 100 | 29/32 | 64 | 105 |
| 150 | 19/26 | 71 | 115 |
| 150 | 24/29 | 71 | 115 |
| 250 | 19/26 | 85 | 135 |
| 250 | 24/29 | 85 | 135 |
| 250 | 29/32 | 85 | 128 |
| 250 | 34/35 | 85 | 128 |
| 500 | 19/26 | 105 | 170 |
| 500 | 24/29 | 105 | 170 |
| 500 | 29/32 | 105 | 165 |
| 500 | 34/35 | 105 | 160 |
| 1000 | 24/29 | 131 | 210 |
| 1000 | 29/32 | 131 | 210 |
| 1000 | 34/35 | 131 | 210 |
| 2000 | 29/35 | 166 | 275 |
| 2000 | 34/35 | 166 | 275 |







| | 198 | | | |
|--------------|---------------------------------|-------------------------|-----------------|------------|
| 5010 | Filtering Flask Standard groun | d mouth | | |
| Capacity(ml) | Ground mouth | O.D.of bottle(mm) | Height(mm) | |
| 10 | 10/19 | 32 | 60 | |
| 10 | 14/23 | 32 | 60 | |
| 25 | 10/19 | 42 | 70 | |
| 25 | 14/23 | 42 | 70 | |
| 50 | 14/23 | 51 | 85 | |
| 50 | 19/26 | 51 | 85 | |
| 100 | 19/26 | 64 | 105 | |
| 100 | 24/29 | 64 | 105 | |
| 250 | 24/29 | 85 | 135 | |
| 250 | 29/32 | 85 | 135 | |
| 500 | 24/29 | 105 | 170 | |
| 500 | 29/32 | 105 | 170 | |
| 1000 | 24/29 | 131 | 210 | |
| 1000 | 29/32 | 131 | 210 | |
| 5013 | Flask pear shape Standard gro | und mouth | | |
| Capacity(ml) | Ground mouth | O.D.of bulb(mm) | Height(mm) | |
| 10 | 19/10 | 30 | 76 | |
| 10 | 14/23 | 30 | 79 | |
| 25 | 10/19 | 40 | 90 | |
| 25 | 14/23 | 40 | 93 | |
| 25 | 19/26 | 40 | 95 | |
| 50 | 14/23 | 50 | 108 | |
| 50 | 19/26 | 50 | 110 | |
| 100 | 14/23 | 60 | 118 | |
| 100 | 19/26 | 60 | 121 | |
| 5014 | Flask pear shape, with two neck | s Standard ground mouth | | |
| Capacity(ml) | Middle ground mouth | Side ground mouth | O.D.of bulb(mm) | Height(mm) |
| 25 | 10/19 | 10/19 | 40 | 102 |
| 25 | 14/23 | 14/23 | 40 | 108 |
| 25 | 19/26 | 14/23 | 40 | 110 |
| 50 | 10/19 | 10/19 | 50 | 112 |
| 50 | 14/23 | 14/23 | 50 | 118 |
| 50 | 19/26 | 14/23 | 50 | 120 |
| 100 | 14/23 | 14/23 | 60 | 128 |
| 100 | 19/26 | 14/23 | 60 | 132 |
| 100 | 19/26 | 19/26 | 60 | 132 |




| 5021 Co | ndenser straight sha | ape Standard ground mou | ıth |
|---------------------|----------------------|-------------------------|------------|
| ength of jacket(mm) | Cround mouth | Cooling area(dm2) | length(mm) |
| 80 | 10/19 | 0.16 | 165 |
| 80 | 14/23 | 0.16 | 175 |
| 150 | 14/23 | 0.6 | 254 |
| 200 | 14/23 | 1 | 305 |
| 200 | 19/26 | 1 | 332 |
| 200 | 24/29 | 1 | 350 |
| 200 | 29/32 | 1 | 360 |
| 300 | 19/26 | 1.6 | 432 |
| 300 | 24/29 | 1.6 | 450 |
| 300 | 29/32 | 1.6 | 460 |
| 400 | 19/26 | 2.1 | 535 |
| 400 | 24/29 | 2.1 | 555 |
| 400 | 29/32 | 2.1 | 560 |
| 500 | 19/26 | 2.7 | 632 |
| 500 | 24/29 | 2.7 | 650 |
| 500 | 29/32 | 2.7 | 660 |
| 600 | 19/26 | 3.2 | 735 |
| 600 | 24/29 | 3.2 | 735 |
| 600 | 29/32 | 3.2 | 760 |




| 5022 Co | ondenser with bulbe | ed inner tube Standard | ground mouth |
|----------------------|---------------------|--------------------------------|------------------|
| Length of jacket(mm) | Cround mouth | Cooling area(dm ²) | Total length(mm) |
| 150 | 14/23 | 0.6 | 235 |
| 200 | 14/23 | 1.1 | 295 |
| 200 | 19/26 | 1.1 | 330 |
| 200 | 24/29 | 1.1 | 350 |
| 200 | 29/32 | 1.1 | 360 |
| 300 | 14/23 | 1.6 | 405 |
| 300 | 19/26 | 1.6 | 430 |
| 300 | 24/29 | 1.6 | 450 |
| 300 | 29/32 | 1.6 | 460 |
| 400 | 19/26 | 2.2 | 530 |
| 400 | 24/29 | 2.2 | 550 |
| 400 | 29/32 | 2.2 | 560 |
| 500 | 19/26 | 3 | 630 |
| 500 | 24/29 | 3 | 650 |
| 500 | 29/32 | 3 | 660 |
| 600 | 19/26 | 3.7 | 730 |
| 600 | 26/29 | 3.7 | 750 |
| 600 | 29/32 | 3.7 | 760 |













58

www.chinahuida.cn

Laboratory Glassware



| 5054B | Receiving Tube 105° | bent Standard ground mouth | |
|--------------|---------------------|----------------------------|--|
| Ground mouth | | Ground stopper | |
| 10/19 | | 10/19 | |
| 14/2 | 23 | 14/23 | |
| 14/23 | | 19/26 | |
| 19/26 | | 19/26 | |
| 19/26 | | 24/29 | |
| 24/29 | | 24/29 | |
| 24/ | 29 | 29/32 | |
| 29/3 | 32 | 29/32 | |

| 5054B | A |
|-------|---|

| <i>5055</i> | Adapter 105° | bent Standard ground mouth | |
|--------------|--------------|----------------------------|--|
| Ground mouth | | Ground stopper | |
| 10/19 | | 110 | |
| 14/23 | | 120 | |
| 19/26 | | 130 | |
| 24/29 | | 140 | |
| 29/32 | | 150 | |
| | | | |



| <i>5056</i> | Suction Adapter | Standard ground mouth |
|--------------|-----------------|-----------------------|
| Ground mouth | | Ground stopper |
| 10/19 | | 10/19 |
| 14/23 | | 14/23 |
| 14/23 | | 19/26 |
| 19/26 | | 19/26 |
| 19/26 | | 24/29 |
| 24/29 | | 24/29 |
| 29/32 | | 29/32 |
| 29/32 | | 29/32 |



| 5057 Receiving Test Tube Standard ground mouth | | | | |
|--|----------------------|--|--|--|
| Ground mouth | O.D.of side tube(mm) | | | |
| 14/23 | 10 | | | |
| 19/26 | 12 | | | |
| 24/29 | 12 | | | |



