



SCI - SFI: 300

Cortical Screw 3.5 mm

Profile

Thread Diameter 3.5 mm
Drill Bit for Threaded Hole 2.5 mm
Drill Bit for Gliding Hole 3.5 mm
Core Diameter 2.4 mm
Hexagonal Socket 2.5 mm
Head Diameter 6.0 mm

Cat No.	Length	Cat No.	Length
SCI-300-10	10mm	SCI-300-32	32mm
SCI-300-12	12mm	SCI-300-34	34mm
SCI-300-14	14mm	SCI-300-36	36mm
SCI-300-16	16mm	SCI-300-38	38mm
SCI-300-18	18mm	SCI-300-40	40mm
SCI-300-20	20mm	SCI-300-42	42mm
SCI-300-22	22mm	SCI-300-44	44mm
SCI-300-24	24mm	SCI-300-46	46mm
SCI-300-26	26mm	SCI-300-48	48mm
SCI-300-28	28mm	SCI-300-50	50mm
SCI-300-30	30mm		




SCI - SFI: 301

Cancellous Screw 3.5 mm

Profile

Thread Diameter 3.5 mm
Drill Bit for Threaded Hole 2.5 mm
Drill Bit for Gliding Hole 2.5 mm
Core Diameter 1.9 mm
Hexagonal Socket 2.5 mm
Head Diameter 6.0 mm

SCI-301-10	10mm	SCI-301-28	28mm
SCI-301-12	12mm	SCI-301-30	30mm
SCI-301-14	14mm	SCI-301-35	35mm
SCI-301-16	16mm	SCI-301-40	40mm
SCI-301-18	18mm	SCI-301-45	45mm
SCI-301-20	20mm	SCI-301-50	50mm
SCI-301-22	22mm	SCI-301-55	55mm
SCI-301-24	24mm	SCI-301-60	60mm
SCI-301-26	26mm		



Small Fragment Implants



SCI - SFI : 9


One Third Tubular Plate

Used as buttress plate for fibula, metatarsals and metacarpals. Fixation with 3.5 & 4.0 mm screws. Compression effect can be achieved by eccentric insertion of remote screws.

Profile

Thickness 1.0 mm
Width 9.0 mm
Hole Spacing 12.0 and 16.0 mm
(One-third tube 12 12.0 and 16.0 mm)

Cat No	Holes	Cat No	Holes
SCI-109-02	2	SCI-109-07	7
SCI-109-03	3	SCI-109-08	8
SCI-109-04	4	SCI-109-09	9
SCI-109-05	5	SCI-109-10	10
SCI-109-06	6	SCI-109-12	12



SCI - SFI : 10

One Third Tubular Plate with Collar

Used as buttress plate for fibula, metatarsals and metacarpals. Fixation with 3.5 & 4.0 mm screws. Compression effect can be achieved by eccentric insertion of remote screws.

Profile

Thickness 1.0 mm
Width 9.0 mm
Hole Spacing 12.0 and 16.0 mm
(One-third tube 12 mm dia)

Cat No	Holes	Cat No	Holes
SCI-110-02	2	SCI-110-07	7
SCI-110-03	3	SCI-110-08	8
SCI-110-04	4	SCI-110-09	9
SCI-110-05	5	SCI-110-10	10
SCI-110-06	6	SCI-110-12	12



SCI - SFI : 11

Small DCP 3.5 mm

For Ulna and Radius

Profile

Thickness 3.5 mm
Width 10.0 mm
Hole Spacing 12.0 and 16.0 mm

Cat No	Holes	Cat No	Holes
SCI-111-02	2	SCI-111-07	7
SCI-111-03	3	SCI-111-08	8
SCI-111-04	4	SCI-111-09	9
SCI-111-05	5	SCI-111-10	10
SCI-111-06	6	SCI-111-12	12



SCI - SFI : 12

Straight Reconstruction Plate 3.5 mm

Note: Do not Bend More than 15

Profile

Thickness 2.8 mm
Width 10.0 mm
Hole Spacing 12.0 mm

Cat No	Holes	Cat No	Holes
SCI-112-05	5	SCI-112-12	12
SCI-112-06	6	SCI-112-14	14
SCI-112-07	7	SCI-112-16	16
SCI-112-08	8	SCI-112-18	18
SCI-112-09	9	SCI-112-20	20
SCI-112-10	10	SCI-112-22	22

