

DETACHABLE ELECTRODE IN SOFT RUBBER

The innovative structure of this electrode allows a brief assembly on the cap with elastic tubes in elastomer and an easy position.

The surface in contact with the patient is in soft rubber. The electrode works as fluid-electrode so the internal space must be filled up with conducting gel, using a syringe before the registration.

The connecting cable is fixed permanently to the electrode so the structure results more compact.

The electrodes can be required also with a identifier mark on the two extremities.

The standard length is about 1300 mm with connection and safe plug.

CODE for the request:: **CUPTMM1**



Detachable electrode model **CUPTMM1**