




**GO DISTAL...  
MORE TRACKABILITY  
FOR MORE SUPPORT<sup>a 1,2</sup>**



**CERENOVUS**

PART OF THE **Johnson & Johnson** FAMILY OF COMPANIES

**CEREBASE™ DA**  
guide sheath

# ENGINEERED TO SECURE DISTAL ACCESS FOR GEOMETRIC ANCHORING<sup>b1,2,4</sup>

8F Short Sheath  
Compatible<sup>2</sup>

Available in 80, 90,  
and 95cm lengths

## MORE SUPPORT<sup>b1,2</sup>

ENGINEERED TO SECURE DISTAL ACCESS FOR GEOMETRIC  
ANCHORING, RESULTING IN MORE PROCEDURAL SUPPORT<sup>b1,2,4</sup>

### FEATURES:

- More proximal shaft stiffness designed for arch support<sup>2</sup>
- More procedural support<sup>b1,2,4</sup>



Dexterous (DEX) Tip  
facilitates tight turns<sup>1,2,4</sup>

Large Inner  
Diameter ID  
(0.90")<sup>1,2</sup>

## MORE TRACKABILITY<sup>a1,2</sup>


DEX TIP FACILITATES TIGHT TURNS (SMALL TURNING RADIUS),  
PROVIDING TRACKABILITY TO DISTAL ANATOMY<sup>1,2,3</sup>

### FEATURES:

- DEX Tip - short, flexible and soft<sup>a1,2,3</sup>
- Soft and compliant, rounded distal edges<sup>2</sup>
- Flexibility of the CEREBASE™ DEX Tip is similar to the flexibility of an aspiration catheter<sup>1,2,3</sup>
- Smooth transition zones for pushability<sup>2</sup>
- Excellent kink resistance<sup>a2</sup>

# CEREBASE™ DA

## guide sheath



### CEREBASE™ DA Guide Sheath

| Product Code | ID     | OD (compatible short sheath) | Usable Length |
|--------------|--------|------------------------------|---------------|
| GS9080SD     | 0.090" | 8F                           | 80cm          |
| GS9090SD     | 0.090" | 8F                           | 90cm          |
| GS9095SD     | 0.090" | 8F                           | 95cm          |

<sup>a</sup> CEREBASE™ Sheath has demonstrated more trackability and more support when compared to the four leading selling competitor long guide sheaths. CEREBASE™ has demonstrated more support when compared to one leading selling competitor long guide sheath.

<sup>b</sup> CEREBASE™ Sheath has demonstrated more support when compared to five leading selling long guide sheath competitors when placed in a distal cervical position.

<sup>c</sup> CEREBASE™ Sheath has demonstrated more shaft stiffness when compared to five leading selling competitors long guide sheaths.

<sup>d</sup> CEREBASE™ Sheath has demonstrated that its distal tip is shorter, more flexible and softer than five leading selling competitor long guide sheaths.

<sup>e</sup> CEREBASE™ Sheath has demonstrated more kink resistance than five leading selling competitor long guide sheaths.

<sup>1</sup> CERENOVUS, Competitive Test Report #103618595 Rev2, April 2020, Data on file

<sup>2</sup> CERENOVUS, Competitive Test Report #10364599 Rev1, April 2020, Data on file


<sup>3</sup> CERENOVUS, Competitive Test Report #103628641 October 2019 (Aspiration Catheter), Data on file

<sup>4</sup> Turk A, Manzoor MU, Nyberg EM, et al. Initial experience with distal guide catheter placement in the treatment of cerebrovascular disease: clinical safety and efficacy. J Neurointerv Surg. 2013 May;5(3):247-52. doi: 10.1136/neurintsurg-2011-010256. Epub 2012 Mar 13

**Johnson & Johnson Medical Limited**  
Baird House, 4 Lower Gilmore Bank,  
Edinburgh,  
EH3 9QP  
United Kingdom

**Legal Manufacturer:**  
Medos International SARL  
Chemin-Blanc 38  
CH-2400 Le Locle  
Switzerland

**Important information:** Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, warnings and precautions.  
This publication is not intended for distribution outside of EMEA  
© CERENOVUS 2021. 161748-201204 EMEA



**CERENOVUS**  
PART OF THE **Johnson & Johnson** FAMILY OF COMPANIES