

DOWNLOAD APP



CORONARY

PERIPHERAL

MICROCATHETERS

NEUROVASCULAR

ANGIOGRAPHIC

ALL PRODUCTS

## ASAHI CHIKAI® 008

FRONTLINE/FINESSE

Precision Cases: Tiny Vessels

Soft Tip and Body for Distal Anatomy




TIP RADIOGRAPHY

9 cm

COATING

Hydrophilic, 200cm  
Coated with SLIP-COAT® coating



[ACT ONE >](#)

## Ordering Information

To order, please call us toll-free at 855-286-9473, or email us at [customersupport@asahi-intecc.com](mailto:customersupport@asahi-intecc.com). To see a complete listing of all available ASAHI products, go to our [Full Product Listing](#) page.

PRODUCT	CATALOG NO.	TIP SHAPE	USABLE LENGTH (CM)	DIAMETER (INCH(MM))	RADIOPAQUE SEGMENT (CM)	SPRING COIL LENGTH (CM)	COIL DIAMETER (CM)
ASAHI CHIKAI 008	WAIN-CKI-008-200	Straight	200	0.008 / 0.010 (0.20 / 0.26)	9	9	1.0

\*Coated with SLIP-COAT® coating



## Indications for Use

The ASAHI CHIKAI Neurovascular Guide Wires are intended to be used in the neurovasculature to facilitate the placement and exchange of therapeutic devices such as cerebral catheters during intravascular therapy. These guide wires are intended for use only in the neurovasculature.

[510k Letter >](#)

[< Back to Neurovascular](#)

[BACK TO TOP ^](#)


## ASAHI FULL PRODUCT LISTING

Select from our comprehensive portfolio of purpose-driven tools to master all your interventional cases.

FOUR CORE TECHNOLOGIES  
FULL PRODUCT LISTING  
CUSTOMER SERVICE  
QUALITY POLICY

DOWNLOAD APP

PRODUCT LISTING [>](#)



CONTACT ASAHI INTECC

Toll-free: 855-286-9473  
[customersupport@asahi-intecc.com](mailto:customersupport@asahi-intecc.com)

**CONTACT FORM [>](#)**

Copyright © 2022 ASAHI INTECC USA, INC. All Rights Reserved | [PRIVACY POLICY](#) | [COOKIE POLICY](#)

Please read the terms of use for further details. Unless otherwise specified, all product names appearing in this website are trademarks or registered trademarks owned by or licensed to

ASAHI. All other registered trademarks or trademarks are the property of their respective owners. No use of any ASAHI trademark, trade name, or trade dress in this website may be made without the prior written authorization of ASAHI, except to identify the product or services of the company.