



Successfully designed, our Ares Balloon Inflation Device is intended to inflate the balloon with the desired fluid, monitor the pressure during the procedure and deflate the balloon when procedure is completed.

Remarkably ergonomic design of Ares Balloon Inflation Device assures easy manipulation from beginning to end. The rotating handle is particular for comfortable grip and handling.

Available with a **secure system**, Ares Balloon Inflation Device offers **one step locking and unlocking mechanism** and facilitates to the overall success of the procedure.

Display: Analog

Dial Range: 0-30 ATM (0-400 PSI)

Accuracy: FS ± 1.6%

Stopcock: 3-Way 500 PSI

Extension Tubing: 25 cm Braided with Male

Rotating Luer

Barrel Volume: 30 mL

Y-Connector: Push-Pull Type

Inner Lumen: 9F

Insertion Tool: Blunt Tip - 18G 10 cm

Torquer: For guidewires 0.010" (0.26 mm)

to 0.025" (0.63 mm)

The crystal clear barrel ensure perfect visibility of the fluid allowing the user to observe the volume and content. The 30 mL barrel is suitable for fast and easy inflation and deflation even for balloons with bigger volume.

Thanks to the ergonomic rotating handle, high pressures (up to 30 atm) are achievable with only 3 rotations of the handle.

With its easy-to-handle push and pull head, the y-connector ensures perfect hemostasis. The large lumen allows effective insertion of devices while the accessories, torquer and insertion tool serves as perfect helpers to manipulate.

REF CODE: BID6004



