

ERIC[®] Retrieval Device

Ischemic Stroke Devices

ERIC[®] Retrieval Device Embolus Retriever with Interlinked Cages

1 per box and includes introducer sheathPu

Reference No.	Product Name	# of Spheres	Diameter (mm)	Working Length (mm)	Distal Tip Length (mm)	Min. Catheter ID (inch)	Pusher Length (cm)
ER173015	ERIC 3	3	3.0	15	5.0	0.017	203
ER173020	ERIC 3	4	3.0	20	5.0	0.017	203
ER174024	ERIC 4	4	4.0	24	5.0	0.017	203
ER174030	ERIC 4	5	4.0	30	5.0	0.017	203
ER176035	ERIC 6	4	6.0	35	5.0	0.017	203
ER176044	ERIC 6	5	6.0	44	5.0	0.017	203

INDICATIONS FOR USE:

The ERIC[®] Retrieval Device is intended for use in the revascularization of acute ischemic stroke caused by the intracranial occlusive vessels of patients who are not eligible for intravenous tissue plasminogen activator, IV tPA, or who fail IV tPA therapy.