

direction for use

Sterile insulin syringes for single use

2.1.1 The raw material used in the product, which is producing at the present time, is PP, PE, PVC, SUS304, rubber. The types of product contain only with needle and without needle. The former is made up of needle, needle cap, piston, plunger, barrel etc.

2.1.2 The product has been manufactured since 2005s in the company. The product intended for use the aspiration of insulin or for the injection of fluids immediately after filling.

2.2 ensuring product free from bacteria

The product is sterile after being sterilized by EO gas. It is safer and more convenient for clinical using which can protect the users and the patients against cross-infections.

2.3 Cautions

2.3.1 Only for single use. Use as soon as pouch opened. Discard after single use.

2.3.2 Guarantee contents non-contaminated and package easy to open before expiry date.

2.3.3 Do not use if previously opened or if pouch damaged.

2.3.4 Protect product against being humidified, heated and contaminated during storage and transport.

direcția de utilizare

Seringi sterile pentru insulină de unică folosință

2.1.1 materia primă utilizată în produs, care produce în prezent, este PP, PE, PVC, SUS304, cauciuc. Tipurile de produse conțin numai cu ac și fără ac. Primul este alcătuit din AC, capac de ac, piston, Piston, butoi etc.

2.1.2 produsul a fost fabricat din anii 2005 în cadrul companiei. Produsul destinat utilizării aspirației de insulină sau injectării de lichide imediat după umplere.

2.2 asigurarea produsului fără bacterii

Produsul este steril după ce a fost sterilizat cu gaz EO. Este mai sigur și mai convenabil pentru utilizarea clinică, care poate proteja utilizatorii și pacienții împotriva infecțiilor încrucișate.

2.3 Precauții

2.3.1 numai pentru o singură utilizare. Utilizați imediat ce s-a deschis punga. A se arunca după o singură utilizare.

2.3.2 garanție conținut necontaminat și ambalaj ușor de deschis înainte de data expirării.

2.3.3 nu utilizați dacă s-a deschis anterior sau dacă punga este deteriorată.

2.3.4 protejați produsul împotriva umidificării, încălzirii și contaminării în timpul depozitării și transportului.

direction for use

Sterile syringes for single use

2.1 The product is composed of barrel, plunger, piston and so on.

2.2 The sterile syringe for single use that our factory produce is used to injection of muscles and dosage.

2.3 Cautions

2.3.1 Only for single use. Use as soon as pouch opened. Discard after single use.

2.3.2 Do Not use if previously opened or if pouch damaged.

2.3.3 Guarantee contents non-contaminated and package easy to open before expiry date.

2.3.4 Protect product against being humidified, heated and contaminated during storage and transport.