

# Conventional

Conventional syringes are available with and without an attached needle. Syringe-needle combinations feature a clear barrel with a bold scale.


# Auto disable

Auto disable syringes are distinguished by the fact that they cannot be reused: following injection, the syringe is automatically blocked, thus becoming unusable.




3 ml	Ordinary Syringe	23G x 1 1/4
5 ml	Ordinary Syringe	21G x 1 1/2
10 ml	Ordinary Syringe	18G x 1 1/2
20 ml	Ordinary Syringe	18G x 1 1/2
50 ml	Ordinary Syringe	18G x 1 1/2


Other needle sizes are available



**Luer slip tip**  
Our luer slip tip provides a frictionfit connection that requires you to push and twist the syringe tip into the needle hub.



**Luer Lock tip**  
Our Luer-Lok tip is generally used for injections requiring a secure connection of the syringe to another device.



**Eccentric luer slip tip**  
The eccentric luer slip tip is generally used for venipunctures and aspiration of fluids, allowing close proximity to the skin.

