

# Absorbable Sutures & Surgical Needles

## **Violet Braided & Monofilament Absorbable Sutures**





### **SECURE**

The biggest advantage of using absorbable suture is that they eliminate the need for a second invasive procedure. The superior coating we use ensures first-throw holding, easy knot-repositioning and ultimately knot security you can rely on.

#### **STRONG**

The ideal absorbable suture is one that would be strong enough to provide sufficient tissue support while degrading. Our advanced braiding technology guarantees excellent tensile strength.

#### SOFT

Poly caprolactone calcium stearate coating ensures a smooth & drag free passage through tissues.

NEEDLE	SUTURE	SUTURE	ADDITIONS		SUTURE SIZES			
NEEDLE	MATERIAL	LENGTH	APPLICATIONS	5-0	6-0	7-0	8-0	10-0
Advanced Micro point Spatula AS - 090 - 1308 Diameter - 0.33mm Length - 8.0mm 1/4 circle 90°	Poly Glycolic Acid Violet Braided	45cm (18")	Lacrimal Surgery, Lid & Socket Reconstruction, Retinal Detachment, Vitreous Surgery, Squint Correction, Pterygium, Orbit & Oculoplasty and Membranectomy	A5108				
Advanced Micro point Spatula AS - 090 - 1308 Diameter - 0.33mm Length - 8.0mm	Poly Glycolic Acid Violet Braided	45cm (18")	Lacrimal Surgery, Lid & Socket Reconstruction, Retinal Detachment, Vitreous Surgery, Squint Correction, Pterygium, Orbit & Oculoplasty and Membranectomy		A6108			
Conventional cutting CC - 135 - 2020 Diameter - 0.51 mm Length - 20.0 mm 3/8 circle 135°	Poly Glycolic Acid Violet Braided	45cm (18")	Orbit & Oculoplasty Surgery		A6409			
Advanced Micro point Spatula AS - 140 - 0907 Diameter - 0.230mm Length - 6.5mm 3/8 circle 140°	Poly Glycolic Acid Violet Braided	30cm (12")	Scleratomies, Vitrectomy, Lensectomy and Retinal Detachment			A7102		
Advanced Micro point Spatula AS - 140 - 0806 Diameter - 0.20mm Length - 6.0mm 3/8 circle 140°	Poly Glycolic Acid Violet Braided	30cm (12")	Trabeculectomy, and Conjunctival suturing				A8104	
Advanced Micro point Spatula AS - 140 - 0606 Diameter - 0.15mm Length - 6.0mm 3/8 circle 140°	Poly Glycolic Acid Violet monofilament	10cm (4")	Conjunctival suturing & Cataract Surgeries					A0101

Retains perfect tensile strength during the crucial wound Violet Braided & Monofilament Absorbable Sutures go healing period: Approximately 65% of original tensile hand-in-glove with Aurolab's Ophthalmic Suture strength at 14 days & 35% of original tensile strength at 21 Needles that deliver flexible product features to meet days. Predictable absorption (Loss of mass) in 60-90 days. every surgeon's needs. Call us to know more.

"PGA suture product of sizes from 2-0 to 8-0 are USFDA Approved".



Information published in this catalogue is subject to change without notification

No 1, Sivagangai Main Road, Veerapanjan, Madurai 625 020, India.

Phone: 91 452 3096 100 Fax: 91 452 2446200

E-mail: info@aurolab.com Web: www.aurolab.com