AVANTI[™]+ Sheath Introducer For atraumatic insertion

The proven cannula offers reliable, smooth performance from start to finish



- For smooth, atraumatic insertion
- Kink resistant
- Optimized catheter maneuverability with minimal leakage

For smooth, atraumatic insertion

• Smooth transitions between cannula, guidewire and dilator for:

- Effortless vessel entry
- Increased patient comfort
- Smallest possible risk for complications

Kink Resistant

Dual layer material gives AVANTI[™]+ Sheath Introducer excellent flexibility and the strength needed to prevent kinking - even in tortuous vessels and scar tissue

Optimized catheter maneuverability with minimal leakage

The patented SLIX haemostatic valve with six-cut spiral design allows free catheter movement without blood reflux or air aspiration



Smooth Performance from Start to Finish

Available in 4 lengths with or without mini-guidewire (mid-length without wire only).

- ▶ 5.5cm = brachial
- 7.5cm = short length
- 11cm = standard length
- 23cm = mid-length



ENDOVASCULAR

		5.5cm Cannula		11cm Cannula		23cm Cannula
French size	Colour code	With mini-wire	Without mini-wire	With mini-wire	Without mini-wire	Without mini-wire
4		402-604P		402-604X	402-604A	
5		402-605P		402-605X	402-605A	402-605T
5.5**				402-655X		
6		402-606P	402-606R	402-606X	402-606A	402-606T
6.5**				402-656X		
7		402-607P	402-607R	402-607X	402-607A	
8				402-608X	402-608A	402-608T
9	•			402-609X	402-609A	402-609T
10				402-610X	402-610A	
11				402-611X		



AVANTI[™]+ Standard Sheath Introducer CARDIOLOGY: 11 cm Cannula French Colour **Mini-wire** With Without code (inch) mini-wire mini-wire size .035 504-604A 4 504-604X 5 .038 504-605X 504-605A .038 5.5** 504-655X* .038 504-606A 6 504-606X 6.5** 504-656A* .038 504-656X* 7 .038 504-607X 504-607A 7.5** .038 504-657X* .038 8 504-608X 504-608A 8.5** .038 504-658X* 9 .038 504-609X 504-609A 10 .038 504-610X 504-610A 11 .038 504-611X * Monitoring sheath: provides slightly larger ID for pressure monitoring or side port aspiration with a catheter in place; they are one-half French size larger in size ** 5F MS = Mid-Size = 5.5F 6F MS = Mid-Size = 6.5F 7F MS = Mid-Size = 7.5F 8F MS = Mid-Size = 8.5F

with mini-quidewire

without mini-guidewire



CARDIOLOGY:

AVANTI™+ Mid-length Sheath Introducer

			23 cm Cannula	
French size	Colour code	Dilator usable (cm)	With obturator	Without obturator
4		30		504-604T
5		30		504-605T
6		30	504-606D	504-606T
6.5**		30		504-656T*
7		30	504-607D	504-607T
8		30	504-608D	504-608T
10		30		504-610T
11		30		504-611T
6F MS= Mid-Size=6.5F				


Cordis AVANTI[™] + Sheath Introducer

CARDIOLOGY: AVANTI [™] + Brachial Sheath Introducer with mini-guidewire				
		5.5 cm Cannula		
French size	Colour code	Mini-wire (inch)	Catalogue number	
4		.035	504-604P	
5		.038	504-605P	
6		.038	504-606P	

CARDIOLOGY:

AVANTI[™]+ Short length Sheath Introducer with mini-guidewire

		7.5 cm Cannula		
French size	Colour code	Mini-wire (inch)	Catalogue number	
4		.021	504-604S	
5		.021	504-605S	
6		.021	504-606S	



For Healthcare Professionals Only. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy. Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and AVANTI are trademarks of Cordis and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. © 2023 Cordis. All Rights Reserved. 100545515-4 08/2023

