

Instruments

Pituitary Instruments

Pituitary



PAPAVERO-CASPAR

FF589R 80 x 11 mm

FF590R 90 x 13 mm

FF591R 100 x 15 mm

Slim profile and lightweighted specula
for transsphenoidal surgery with key TE749R



FF590B 90 x 13 mm


Noir®

Transsphenoidal specula



TE749R

Key



COTTLE

OK089R-OK090R

Self-retaining nasal speculum,
suitable as introducer for FF589R, FF590R, FF591R
135 mm, 5 1/4"

OK089R

76 x 7 mm

OK090R


90 x 7 mm



PAPAVERO-HARTMANN

FF595R

Spreader for specula
FF589R, FF590R, FF591R
140 mm, 5 1/2"




backwards
cutting


FA076R

Antrum punch for removal of posterior nasal septum
Kerrison bone punches and rongeurs see page 4-28.

Micro Suture Tying Forceps with Round Handles




with plateau at the tips




FD280R-FD281R

120 mm, 4 3/4"



FD284R-FD285R


150 mm, 6"



FD286R-FD287R

180 mm, 7"

DUROGRIP® Micro Needle Holders



smooth



0.2



0.4



0.5

- Aesculap needle holders have jaws with tungsten carbide inserts and according to the use required the resp. structure can be chosen. A smooth jaw surface as well as 3 different serrations ensure a sure „grip“ of the resp. needle.

- We recommend tungsten carbide inserts

with smooth jaw surface for suture material	9/0 - 11/0
with 0.2 mm pitch of serration for suture material	6/0 - 10/0
0.4 mm	4/0 - 6/0
0.5 mm	up to 3/0

with
smooth
jaw**BM001R**

straight

1/1
0.2**BM002R**

straight

CASTROVIEJO

BM001R-BM002Rwith ratchet
145 mm, 5 3/4"

Instruments

The Long Micro Forceps

Micro

0.3 mm



1/1

MÜLLER
FM033R

160 mm, 6 1/4"

0.5 mm



1 x 2 teeth



1/1

MÜLLER
FM034R

160 mm, 6 1/4"



Diam.
of ring:
1.2 mm



1/1

MÜLLER
FM035R

160 mm, 6 1/4"

0.4 mm



with plateau



1/1

MÜLLER
FM036R

Micro suturing forceps
160 mm, 6 1/4"

Instruments

Pituitary Instruments

Pituitary

140 mm, 5 1/2"

FF616R Diam. 6.5 mm
cutting to the right

FF617R Diam. 6.5 mm

FF618R Diam. 6.5 mm





NICOLA
FF616R-FF618R
Curettes, malleable
260 mm, 10 1/4"

1/1

215 mm, 8 1/2"

NICOLA
FF619R
Raspatory
215 mm, 8 1/2"


1/1



Instruments

Pituitary Instruments

Pituitary



FA072R

straight



FA073R

left curved



FA074R

right curved




FA075R

angular


FA072R-FA075R

Scissors



FA069R

straight



FA070R

right curved





FA071R

left curved

FA069R-FA071R


Forceps



Instruments

Pituitary Instruments

Pituitary



FF640R Diam. 2.2 mm
Raspatory, rigid



FF641R Diam. 3.2 mm
Raspatory, rigid



FF642R Diam. 2.5 mm
Curette, rigid





FF643R Diam. 1.7 mm
Blunt hook, rigid



FF644R Diam. 2.0 mm
Dissector, rigid

LANDOLT-REULEN
FF640R-FF644R
260 mm, 10 1/4"






PLESTER

Art. No.	Diam.
GF939R	1.5 mm
GF940R	2.0 mm
GF941R	2.5 mm
GF942R	3.0 mm

with finger cut-off and stylet, thumb piece
parallel to curve, with LUER hub
195 mm, 7 3/4"



FERGUSON

Working length	110 mm, 4 3/8"	130 mm, 5 1/8"	160 mm, 6 1/4"	200 mm, 8"
Total length	195 mm, 7 3/4"	215 mm, 8 1/2"	245 mm, 9 3/4"	280 mm, 11"
Diam.				


1.5 mm	GF350R	GF360R	GF370R	GF380R
2.0 mm	GF351R	GF361R	GF371R	GF381R
2.5 mm	GF352R	GF362R	GF372R	GF382R
3.0 mm	GF353R	GF363R	GF373R	GF383R
4.0 mm	GF354R	GF364R	GF374R	GF384R
5.0 mm	GF355R	GF365R	GF375R	GF385R



MÜLLER
FM021R
with ratchet
145 mm, 5 3/4"



MÜLLER
FM061R
with ratchet
180 mm, 7"




FD120R
without ratchet, bayonet-shaped
195 mm, 7 3/4"

Instruments

Micro Scissors, Straight

Micro



MÜLLER

FM010R

120 mm, 4 3/4"



MÜLLER

FM011R

120 mm, 4 3/4"



MÜLLER


FM012R

120 mm, 4 3/4"

Instruments

KERRISON Bone Punches

Bone



Jaw position 130°, upwards opening

Shaft length	Width (mm)	Footplate	Non detachable, without ejector	Detachable	Ejector	Noir®, detachable	Ejector	Jaw opening (mm)
200 mm, 8"	1.5	standard	FF649R	FK966R	–	FK966B	–	9
	2.0	standard	FF737R	FK913R	✓	FK913B	✓	9
	2.5	standard	FF657R	FK967R	✓	FK967B	✓	10
	3.0	standard	FF738R	FK914R	✓	FK914B	✓	10
	4.0	standard	FF727R	FK915R	✓	FK915B	✓	12
	5.0	standard	FF728R	FK916R	✓	FK916B	✓	12
	6.0	standard	FF729R	FK917R	✓	FK917B	✓	12
	1.0	thin	FF661R	FK961R	–	FK961B	–	8
	1.5	thin	FF658R	FK968R	–	FK968B	–	9
	2.0	thin	FF662R	FK962R	✓	FK962B	✓	9
	2.5	thin	FF659R	FK969R	✓	FK969B	✓	10
	3.0	thin	FF663R	FK963R	✓	FK963B	✓	10
	4.0	thin	FF664R	FK964R	✓	FK964B	✓	12
	5.0	thin	FF665R	FK965R	✓	FK965B	✓	12
	2.0	standard	FF660R	FK970R	✓	FK970B	✓	15
	3.0	standard	FF739R	FK982R	✓	FK982B	✓	15
	4.0	standard	FF745R	FK983R	✓	FK983B	✓	15
	5.0	standard	FF751R	FK984R	✓	FK984B	✓	15
	6.0	standard	FF753R	FK985R	✓	FK985B	✓	15
230 mm, 9"	2.0	standard	FF682R	FK919R	✓	–	–	9
	3.0	standard	FF683R	FK920R	✓	–	–	10
	4.0	standard	FF684R	FK921R	✓	–	–	12
	5.0	standard	FF685R	FK922R	✓	–	–	12
	2.0	thin	FF746R	FK925R	✓	FK925B	✓	9
	3.0	thin	FF747R	FK926R	✓	FK926B	✓	10
	4.0	thin	FF748R	FK927R	✓	FK927B	✓	12
	5.0	thin	FF749R	FK928R	✓	FK928B	✓	12
280 mm, 11"	2.0	thin	FF666R	FK976R	✓	FK976B	✓	9
	3.0	thin	FF667R	FK977R	✓	FK977B	✓	10
	4.0	thin	FF668R	FK978R	✓	FK978B	✓	12
	5.0	thin	FF669R	FK979R	✓	FK979B	✓	12

Width (mm)

1.0

1.5

2.0

2.5

3.0

4.0

5.0

6.0



MÜLLER
FM001R-FM003R
110 mm, 4 $\frac{3}{8}$ "



MÜLLER
FM004R
110 mm, 4 $\frac{3}{8}$ "



MÜLLER
FM005R
110 mm, 4 $\frac{3}{8}$ "



MÜLLER
FM006R
110 mm, 4 $\frac{3}{8}$ "

Storages for Instruments Racks for Micro Instruments



MÜLLER

Art. No.	Dimensions L/W/H
----------	---------------------

FM024R 270 x 124 x 23 mm

Rack with silicone pads and compartments for accessories provides organized and protective storage for the MÜLLER micro instruments FM001R–FM021R, applying forceps FD568R or FD569R and various micro vessel clips and approximators. (Instruments not included)




MÜLLER

Art. No.	Dimensions L/W/H
----------	---------------------

FM026R 270 x 124 x 23 mm

Rack with silicone pads and compartments for accessories provides organized and protective storage for the MÜLLER micro instruments FM033R–FM061R, applying forceps FD568R or FD569R and various micro vessel clips and approximators. (Instruments not included)



Art. No.	Dimensions L/W/H
----------	---------------------

JG381R 273 x 176 x 30 mm

Perforated basket with instrument racks and lid for storage of AL MEFTY skull base instruments

PRODUCT DATA SHEET

Micro Dressing Forceps, straight, bayonet-shaped, delicate, smooth, 185 mm (7 1/4") (FD110R)



Attribute Name	Attribute Value
Product search type	Micro Dressing Forceps
Product family	Instruments
Tungsten Carbide	no
Geometry	straight
Form	delicate, bayonet-shaped
Length	185.00 mm
Length (Inch)	7 1/4"
Working length (mm)	90.0 mm
Working length (Inch)	3 1/2"
Jaw profile	smooth
Jaw/Blade width	0.60 mm
Handle	flat handle
Discipline	Neurosurgery, Spine Surgery
UDI Direct Marking	Yes
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	reusable
Quality standard	Classic
Shape working end	straight
Design	bayonet

Attribute Name	Attribute Value
Product category	micro

PRODUCT DATA SHEET

Micro Dressing Forceps, straight, bayonet-shaped, broad, smooth, 185 mm (7 1/4") (FD111R)



Attribute Name	Attribute Value
Product search type	Micro Dressing Forceps
Product family	Instruments
Tungsten Carbide	no
Geometry	straight
Form	broad, bayonet-shaped
Length	185.00 mm
Length (Inch)	7 1/4"
Working length (mm)	90.0 mm
Working length (Inch)	3 1/2"
Jaw profile	smooth
Jaw/Blade width	0.70 mm
Handle	flat handle
Discipline	Neurosurgery, Spine Surgery
UDI Direct Marking	Yes
Product contains Latex	No
Sterility	non-sterile
Disposable / Reusable	reusable
Quality standard	Classic
Shape working end	straight
Design	bayonet

Attribute Name	Attribute Value
Product category	micro

PRODUCT DATA SHEET

KILLIAN SEPTUM SPECUL.ASEPTIC 75X7MM #3 (OK083R)



Attribute Name	Attribute Value
Product search type	Nasal Specula
Author	KILLIAN
Aesculap Product Groups	General Surgical Instruments, ENT
Length	135.00 mm
Length (Inch)	5 1/4"
Dimenions	75 x 7 mm
Figure / Size	Fig. 3
Indication	Rhinoplasty, Rhinology
Surgical Step	Surgical Access
UDI Direct Marking	Yes
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No

Aesculap[®] Detachable Rongeurs

Simply clean



Aesculap Neurosurgery

Aesculap® Detachable Rongeurs

Simply clean

Ever increasing demands in the reprocessing of surgical instruments require professional and efficient processes, together with intelligent product solutions, which meet modern hygiene standards. Our new detachable rongeurs fulfill these demands.

Simply clean; even without pre-cleaning.

Einfach sauber

Die stetig wachsenden Anforderungen an die Aufbereitung chirurgischer Instrumente erfordern, neben professionellen und effizienten Aufbereitungsprozessen, intelligente Produktlösungen, die den Ansprüchen der Hygiene gerecht werden.

Unsere neuen zerlegbaren Rongeurs tun dies.

Einfach sauber; auch ohne Vorreinigung.

NOIR® coating


No irritating light reflections

NOIR®-Beschichtung


Keine störenden Lichtreflexionen



straight / gerade



curved / aufgebogen



standard / standard



semi-serrated / halb-gezahnt

The choice is yours

86 models with working lengths from 150 mm to 280 mm, jaw widths from 1.5 mm to 5 mm, standard and semi-serrated jaws in straight and 150° upward curved. All models come with the proven NOIR® coating for the prevention of distracting light reflections.

Sie haben die Wahl

86 Modelle mit Arbeitslängen von 150 mm bis 280 mm, Maulbreiten zwischen 1,5 mm bis 5 mm, gezahnte und ungezahnte Arbeitsenden in gerade und 150° aufgebogen. Alle Modelle haben die bewährte NOIR®-Beschichtung, um störende Lichtreflexionen zu verhindern.

Intuitive handling

The golden button indicates the detachability and enables an intuitive, easy and quick assembly and disassembly. In the open cleaning position all parts remain connected. This prevents parts getting lost during reprocessing.

Intuitives Handling

Der goldene Knopf weist auf die Zerlegbarkeit hin und ermöglicht eine intuitive, einfache und schnelle Montage und Demontage. In der offenen Reinigungsposition bleiben alle Teile miteinander verbunden. Dies verhindert, dass Einzelteile während der Aufbereitung verloren gehen.



Easy identification

- Working length
- Jaw width

Einfache Identifizierung

- Arbeitslänge
- Maulbreite

Simple disassembly / Einfaches Zerlegen


- 1 Push the golden button
Drücken Sie den goldenen Knopf
- 2 Pull down the front handle
Ziehen Sie den vorderen Griff nach unten
- 3 Lift the slider
Heben Sie den Schieber an

Aesculap® Detachable Rongeurs



Rongeurs standard / Rongeur standard


Working length Arbeitslänge	Jaw Maul	1.5 mm	2 mm	3 mm	4 mm	5 mm
150 mm		FH801B	FH802B	FH803B	FH804B	FH805B
150 mm		FH901B	FH902B	FH903B	FH904B	FH905B
180 mm		FH811B	FH812B	FH813B	FH814B	FH815B
180 mm		FH911B	FH912B	FH913B	FH914B	FH915B
200 mm		FH821B	FH822B	FH823B	FH824B	FH825B
200 mm		FH921B	FH922B	FH923B	FH924B	FH925B
230 mm		FH891B	FH892B	FH893B	FH894B	FH895B
230 mm		FH931B	FH932B	FH933B	FH934B	FH935B
280 mm		-	-	FH837B	FH838B	FH839B



Rongeurs semi-serrated / Rongeurs halb-gezahnt

Working length Arbeitslänge	Jaw Maul	1.5 mm	2 mm	3 mm	4 mm	5 mm
150 mm		FH851B	FH852B	FH853B	FH854B	FH855B
150 mm		FH951B	FH952B	FH953B	FH954B	FH955B
180 mm		FH861B	FH862B	FH863B	FH864B	FH865B
180 mm		FH961B	FH962B	FH963B	FH964B	FH965B
200 mm		FH871B	FH872B	FH873B	FH874B	FH875B
200 mm		FH971B	FH972B	FH973B	FH974B	FH975B
230 mm		FH881B	FH882B	FH883B	FH884B	FH885B
230 mm		FH981B	FH982B	FH983B	FH984B	FH985B
280 mm		-	-	FH847B	FH848B	FH849B

Aesculap® Instrument Organisation System



Aesculap Instrument Organisation System

Surgical instruments are valuable assets for every hospital, which is why they should be processed as efficiently and gently as possible.

Therefore, protecting the instruments during sterilisation, as well as in transport to sterile supply stores and to the operating room, is of crucial importance.

The new Aesculap organisation system for instrument trays facilitates the storage especially of delicate instruments and small items.

Aesculap Instrument Organisation System

Chirurgische Instrumente stellen ein wertvolles Kapital für jedes Krankenhaus dar und sollten daher so effizient und schonend wie möglich aufbereitet werden.


Der Schutz der Instrumente während der Sterilisation und beim Transport in das Sterilgutlager und in den OP sind dabei von zentraler Bedeutung.

Das neue Aesculap Organisationssystem für Instrumentensiebe erleichtert die Lagerung insbesondere von empfindlichen Instrumenten und Kleinteilen.



For further information on the Aesculap Instrument Organisation System, please refer to the brochure C81902.

Weitere Informationen zum Aesculap-Instrumenten-Organisationssystem finden Sie in unserem Prospekt C81901.



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

Aesculap – a B. Braun company


The main product trademark "Aesculap" and the product trademarks "NOIR", are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. C96011

0115/2.0/1

Pos No.	Item Code	Description / Comment	Weight (g)	Qty.
	10 FF589R	PAPAVERO-CASPAR SPEC SLIM-PROF 80X11MM  Customer Reference: 129.1	216	1
	20 FF590R	PAPAVERO-CASPAR SPEC SLIM-PROF 90X13MM  Customer Reference: 129.2	224	1
	30 FF591R	PAPAVERO-CASPAR SPEC SLIM-PROF 100X15MM  Customer Reference: 129.3	227	1
	40 FF595R	PAPAVERO-HARTMAN SPREADER F/FF589R-91R Customer Reference: 129.4	89	1

Pos No.	Item Code	Description / Comment	Weight (g)	Qty.
50	OC485R	CASTROVIEJO Micro Scissors (Eye Scissors), angled to side, blunt/blunt, 95 mm (3 3/4")  Customer Reference: 129.10	37	1
60	OL153R	MASING Elevator, 200 mm (7 7/8"), double ended, graduated, blunt/blunt, non-sterile, reusable  Customer Reference: 129.14	27	1
70	JK444	Bottom, without base perforation, Standard 1/1, outer length: 592 mm, outer width: 274 mm, outer height: 187 mm, inner length: 544 mm, inner width: 258 mm, inner height: 172 mm  Customer Reference: 129.37	3,009	1
80	JK489	Inner Lid, basis, Standard 1/1, silver, outer length: 582 mm, outer width: 291 mm, outer height: 36 mm  Customer Reference: 129.37	1,269	1

Pos No.	Item Code	Description / Comment	Weight (g)	Qty.
	90 JF932	Positioning mat, Standard DIN, 470 x 230 x 32 mm, trimmed and perforated, slip-free storage of instruments, good drying characterist	832	2
				
		Customer Reference: 129.37		
	100 JF223R	Perforated basket, Standard 1/1, 540 x 253 x 76 mm, with feet	1,099	2
				
		Customer Reference: 129.37		
	110 JG785B	Accessories Container, identification label, standard, aluminium, red, outer length: 56 mm, outer width: 18 mm	2	1
				
		Customer Reference: 129.37		
	120 JK090	Reusable Filter, package of 10 pieces, diam. 190 mm	32.2	1
				
		Customer Reference: 129.37		