

Product Data Sheet

ActiFast Two Way Stretch Tubular Retention Bandage

1. Composition of the product

ActiFast Two Way Stretch Tubular Retention Bandage consists of:

- 52% Viscose
- 48% Elastane (Spandex)

This product data sheet covers the following items:

REF 285-6490 ActiFast Tubular Retention Bandage	Red	3.5cm x 1m
REF 285-6508 ActiFast Tubular Retention Bandage	Red	3.5cm x 10m
REF 285-6516 ActiFast Tubular Retention Bandage	Green	5cm x 1m
REF 285-6524 ActiFast Tubular Retention Bandage	Green	5cm x 3m
REF 285-6532 ActiFast Tubular Retention Bandage	Green	5cm x 5m
REF 285-6540 ActiFast Tubular Retention Bandage	Green	5cm x 10m
REF 285-6557 ActiFast Tubular Retention Bandage	Blue	7.5cm x 1m
REF 285-6565 ActiFast Tubular Retention Bandage	Blue	7.5cm x 3m
REF 285-6573 ActiFast Tubular Retention Bandage	Blue	7.5cm x 5m
REF 285-6581 ActiFast Tubular Retention Bandage	Blue	7.5cm x 10m
REF 285-6599 ActiFast Tubular Retention Bandage	Yellow	10.75cm x 1m
REF 285-6615 ActiFast Tubular Retention Bandage	Yellow	10.75cm x 3m
REF 285-6623 ActiFast Tubular Retention Bandage	Yellow	10.75cm x 5m
REF 285-6631 ActiFast Tubular Retention Bandage	Yellow	10.75cm x 10m
REF 292-4298 ActiFast Tubular Retention Bandage	Beige	17.5cm x 1m
REF 290-5750 ActiFast Tubular Retention Bandage	Beige	17.5cm x 10m
REF 369-0302 ActiFast Tubular Retention Bandage	Purple	20cm x 1m
REF 369-0310 ActiFast Tubular Retention Bandage	Purple	20cm x 5m
REF 369-0328 ActiFast Tubular Retention Bandage	Purple	20cm x 10m

Not on FP10/GP10

2. Availability

FP10 & GP10

Products not available on FP10/GP10 are available in Hospitals.

3. Indications for use

Dressing retention, wet wrapping, stocking and bandage liner, skin protection and also helps penetration of emollients.

4. Warnings and precautions

Not suitable for use as a compression bandage.

5. Aftercare

ActiFast cannot be washed, for single use only. Cool Iron. Do not bleach.

6. Limb/Trunk Measurement Guide

Red – Small limbs

Green – Small/medium limbs

Blue – Large limbs

Yellow – Extra-large limbs/heads and children's trunks

Beige – Small/medium adult's trunks

Purple – Large/extra-large adult's trunks

7. Extra information

Latex free unique self-supporting, conforming tubular bandage. ActiFast tape measure is available. ActiFast is very easy to apply and is colour coded to denote size.

ActiFast can also be cut to any length.

8. Storage

Keep ActiFast dry.

9. Shelf Life

5 years.

10. Wear Time

As advised by Healthcare professional.