


DAHLHAUSEN®

10 Jul 38

**CPAP Masks**

CPAP mask with cuff. The mask cuff can be adapted to the shape of the head by the amount of air administered into the mask.
The mask has a double port. The smaller port with a 22 mm diameter is used to attach it to a tube for connection to the CPAP generator. The larger port with 30 mm is used to attach it the PEEP-valve required.

Art.no.	Size
5160005521	Small/medium
<u>5160005520</u>	<u>Large</u>

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Str. 53
50996 Köln

Tel.: (02236) 39 13 - 0
Fax: (02236) 39 13 - 109