
ACCESSORIES

FILTER ACCESSORIES




JP055

Mounting handle for filter insert, PRIMELine® and PRIMELine® PRO lids.

CONSUMABLES

REUSABLE FILTERS






AESCULAP® reusable filter, round for use in up to 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 10






AESCULAP® reusable filter, square for legacy standard sterile containers, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 10






AESCULAP® reusable filter for mini, dental and scope sterile containers with filter retention plate JK098, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 2





AESCULAP® reusable filter for dental sterile containers without filter retention plate JK098, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 2

MARKING PENS





Marking pen for AESCULAP® reusable filters

Order unit (PAK): Pack of 3







Hi-Line XS handpiece, angled 1 ring code; spray nozzle GB761SU (not included)



GB757R



Hi-Line XS handpiece, angled 2 ring code; spray nozzle GB762SU (not included)



CR7COD



Hi-Line XS handpiece, angled 3 ring code; spray nozzle GB763SU (not included)



GB740R

GB614R

Hi-Line XS craniotome and multi purpose handpiece

to be used with dura guards, drill depth guard and holding sleeve GB741R - GB748R, spray nozzle GB767SU (not included)



GB741R GB615R

Hi-Line XS dura guard

fixed;

1 ring code;

to be used with craniotome burr **GE420R, GE420SU, GE429SU**



Hi-Line XS dura guard

fixed;

2 ring code;

to be used with craniotome burr GE520R, GE520SU, GE529SU



GB743R

Hi-Line XS dura guard

fixed;

3 ring code;

to be used with craniotome burr GE620R, GE620SU, GE629SU



Hi-Line XS dura guard

steerable;

1 ring code;

to be used with craniotome burr GE420R, GE420SU, GE429SU



Hi-Line XS dura guard

steerable;

2 ring code;

to be used with craniotome burr GE520R, GE520SU, GE529SU