

Sclerotherapy Injection Needle

Used for endoscopic injection into esophageal varices above esophageal junction




Art. No.	Needle Gauge	Needle Length (mm)	Working Length (cm)	Working Channel Ø (mm)
SGN-2319	19G with lock	7	170	2.8
SGN-2319-C	19G with lock	7	230	2.8
SGN-2320	20G with lock	7	170	2.8
SGN-2320-C	20G with lock	7	230	2.8
SGN-2321	21G with lock	7	170	2.8
SGN-2321-C	21G with lock	7	230	2.8
SGN-2322	22G with lock	7	170	2.8
SGN-2322-C	22G with lock	7	230	2.8
SGN-2323	23G with lock	7	170	2.8
SGN-2323-C	23G with lock	7	230	2.8
SGN-2324	24G with lock	4	170	2.8
SGN-2324-C	24G with lock	4	230	2.8
SGN-2325	25G with lock	4	170	2.8
SGN-2325-C	25G with lock	4	230	2.8

C = Colono

Class IIa device – sterile – for single use only

Clip Applicator – single

Used to mechanically bind blood vessels



Art. No.	Working Length (cm)	Sheath Ø (mm)	Jaw angle	Opening width (mm)
S-CD-16	165	2.5	135°	6
S-CD-21	210	2.5	135°	6

Class IIa device – sterile – for single use only