











< 1/5 >

<











>

## BLUESTAR® Tablets - 8 Applications

\$147.00 USD

Quantity

 - +

Add to cart

### Description

Very compact and multi-use product providing ability to make 1 to 8 applications of BLUESTAR® chemiluminescent reagent. Will not harm DNA.

### BLUESTAR® TABLETS - 8 APPLICATIONS:

- 2 tubes of tablets
- Red Tube holds 8 beige tablets (reagent)
- White Tube hold 8 white tablets (catalyst)
- User guide

### PREPARE THE WORKING SOLUTION:

Place 1 white and 1 beige tablet in 4 oz. spray bottle of distilled or purified water. Dissolve tablets and spray on target area. \*Bottle not included.

### SUGGESTED STORAGE:

Store BLUESTAR® products in a dry environment with temperatures ranging 39° - 86° F.


### DOCUMENTATION:

- User's manual
- Safety Data Sheets :
  - + Beige tablet
  - + White tablet
  - + Working solution

### PRODUCT ORDER CODE:

**BL-FOR-TAB8**

## You may also like



BLUESTAR® MINI KIT <b>\$69.00 USD</b>	BLUESTAR® Tablets - 4 applications <b>\$94.00 USD</b>	BLUESTAR® TRAINING TABLETS - 8 Applications <b>\$61.00 USD</b>	BLUESTAR® TRAINING TABLETS - 4 Applications <b>\$33.00 USD</b>
--	---	---	---



### BLUESTAR USA

P.O Box 1446 - Greer, SC - 29651 - [alicia.conners@bluestar-forensic.com](mailto:alicia.conners@bluestar-forensic.com) - US Phone number: 1-803-394-3014

[FAQs](#)

[Search](#)

[Return Policy](#)

[Privacy Policy](#)

[Terms of Service](#)