



Test report No. 055023id

EVALUATION OF FUNGICIDAL OR YEASTICIDAL ACTIVITY FOR INSTRUMENTS
USED IN MEDICAL AREA (EN 14562)

Name of the product: STERISEPT WIPES

Batch number: 14151022

Date of test report: 04.03.2023

Client, representative:
Chemi-Pharm Ltd.
Tänassilma tee 11
Tänassilma küla
Saku vald 76406
ESTONIA

Test report No. 055023id

EVALUATION OF FUNGICIDAL OR YEASTICIDAL ACTIVITY FOR INSTRUMENTS
USED IN MEDICAL AREA (EN 14562)

Name of the product: STERISEPT WIPES

Batch number: 14151022

Order number: 20135

Manufacturer: Chemi-Pharm Ltd.

Client, representative: Chemi-Pharm Ltd., Tännassilma tee 11, Tännassilma küla, Saku vald, 76406, ESTONIA, Siimu Rom

Date of delivery: 16.03.2023

Test material conditions: No specific features, sample in the manufacturers tare

Storage conditions: At room temperature, dark

Active substance – conc.: N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine 0.45%;
Didecyl-Dimethyl-Ammonium Chloride (DDAC)0.45%

Appearance of the product: Transparent liquid

Test concentration: Ready to use

Contact time: 3 min

Interfering substance: 3g/l bovine albumin + 3ml/l sheep blood erythrocytes (dirty conditions)

Rinsing liquid: -

Neutralizer: Polysorbate 80 30g/l; saponin 30 g/l, lecithin 3 g/l

Test organisms: *Candida albicans* ATCC 10231

Testing method: EVS-EN 14562:2006
Quantitative carrier test for the evaluation of fungicidal or yeasticidal activity for instruments used in medical area (phase 2, step 2)

Testing period: 27.03.2023-30.03.2023

Results: look appendix 1-2



Allar Laaneleht

Microbiology chief specialist

Date of test report: 04.03.2023

Appendix 1

TEST RESULTS (yeastocidal carrier test)

EVS-EN 14562; Phase 2, step 2

Dilution – neutralization method

Rinsing liquid: -

Test organism: *Candida albicans* ATCC 102;

Test temperature: +20° C; Incubation temperature: +30° C

Interfering substance: 3g/l bovine albumin + 3ml/l sheep blood erythrocytes

Nordic Tersus Laboratory LLC.

Date of test: 27.03.2023

Responsible person: Allar Laaneleht

Validation and controls

Validation suspension (N_{vo})		Experimental Conditions control (A)		Neutralizer control (B)		Method validation (C)	
V_{C1}	V_{C2}	V_{C1}	V_{C2}	V_{C1}	V_{C2}	V_{C1}	V_{C2}
82	79	96	75	82	87	80	99
$\bar{x} = 60.5$		$\bar{x} = 85.5$		$\bar{x} = 84.5$		$\bar{x} = 89.5$	
$30 \leq \bar{x} \text{ of } N_{vo} \leq 160?$ yes <input checked="" type="checkbox"/> ; no <input type="checkbox"/>		$\bar{x} A \text{ is } \geq 0,5 \bar{x} \text{ of } N_{vo}?$ yes <input checked="" type="checkbox"/> ; no <input type="checkbox"/>		$\bar{x} B \text{ is } \geq 0,5 \bar{x} \text{ of } N_{vo}?$ yes <input checked="" type="checkbox"/> ; no <input type="checkbox"/>		$\bar{x} C \text{ is } \geq 0,5 \bar{x} \text{ of } N_{vo}?$ yes <input checked="" type="checkbox"/> ; no <input type="checkbox"/>	

Test suspension

Test suspension (N):	N	V_{C1}	V_{C2}	$\bar{x} \text{ wm} = 3.31 \times 10^8$; $\lg N = 8.68$ $8.17 \leq \lg N \leq 8.7$ yes <input checked="" type="checkbox"/> ; no <input type="checkbox"/>
	10^{-6}	>330	>330	
	10^{-7}	23	34	

Water control

Water control (N_w):	N_w	V_{C1}	V_{C2}	$\bar{x} \times 10 = 1.6 \times 10^6$ $6.15 \leq \lg N_w = 6.38 \leq (\lg N - 1.3)?$ yes <input checked="" type="checkbox"/> ; no <input type="checkbox"/>
	10^{-4}	21	17	
	10^{-5}	2	2	

Test

Conc. of the product %	Dilution step	V_{C1}	V_{C2}	$Lg Na = Lg (\bar{x} \text{ or } \bar{x}_{wm}) + 1$	$Lg R (lg N_w = 6.18)$	Contact time
100%	10^0	14	16	2.20	4.28	3 min
	10^{-1}	<14	<14			
	10^{-2}	<14	<14			
	10^{-3}	<14	<14			
50%	10^0	>330	>330	4.96	1.42	3 min
	10^{-1}	>330	>330			
	10^{-2}	85	99			
	10^{-3}	>14	>14			
25%	10^0	>330	>330	5.76	0.62	3 min
	10^{-1}	>330	>330			
	10^{-2}	>330	>330			
	10^{-3}	66	50			

Explanations

V_c = count per ml (one plate or more)

\bar{x} = average of V_{C1} and V_{C2}

R = reduction ($lg R = lg N_w - lg N_a$)

If $N_a < 140$, $lg R = > [lg N_w - 2,15]$

Appendix 2

Interpretation

The product for instrument disinfection **STERISEPT WIPES** (batch no. 14151022) was tested according to the test method EVS-EN 14562:2006. The test was performed at 20 ± 1 °C, under dirty conditions with the contact time of 3 min. The dilution – neutralization method was used for testing the products' effectiveness against the reference strain *Candida albicans* ATCC 10231. Under dirty conditions the product was effective against the reference strain within 3 min.

Conclusion

The surviving count of the reference strain showed at least 4 lg reduction meaning that according to the standard EVS-EN 14562:2006, **under dirty conditions the ready to use product STERISEPT WIPES has a yeasticidal effect within 3 min.**



Allar Laaneleht

Microbiology chief specialist

04.03.2023