

## Technical Data Sheet


**Model No.:N7601**  
**BUDI-DI: 697442484L33QL**

- **Product Name:** Universal Pack
- **Classification:** Class Is, according to Annex IX, MDD 93/42/EEC
- **Description:** Universal pack is mainly made of AS-Only SMMS by cutting, punching and processing with certain protective characteristics, including 6 accessories: table cover, Mayo stand cover, side drape, hand towel, etc.
- **Intended use:** Universal pack is used in various departments of medical institutions can be used alone or paired with other surgical packages..
- **Material:** 45g AS-Only SMMS
- **Approvals:** CE, ISO 13485, EN13795-1
- **Manufacturer:** Leboo Healthcare Products Limited
- **Country of origin:** China
- **Specification:** Detail see specification as below.

Component	Description	Color	Quantity
Table Cover	150x190cm (75x190cm)	Blue	1
Mayo Stand Cover	80x145cm (48x60cm)	Blue	1
Top Drape	270x150 cm (60x30cm) with adhesive	Blue	1
Bottom Drape	190x195 cm (60x30cm) with adhesive	Blue	1
Tape Strip (op-tape)	10x50cm	Blue	1
Hand towel	30x40cm	White	2

### ● Packaging:

- (1) Carton Requirements: Five layers corrugated
- (2) Packing Quantity: 1pc/header pouch, 12pcs/ctn  
5 Layers Carton(Paper)



### ●Storage:

- (1) Store in dry, clean conditions in original packaging.
- (2) Store away from direct sunlight, source of high temperature and solvent vapors.
- (3) Store with the temperature range -5°C to +45°C and with relative humidity below 80%.

### ●Shelf Life:

Shelf life is 36 months from the date of manufacture when stored as stated above.

Design:

Check:

Approve:

Date: 2022-9-13

**●Cautions and disposal:**

- (1) Do not contact with heavy oils, sparks or flames;
- (2) Do not use in the environment with conditions of excessive heat;
- (3) Replace if damaged, heavily contaminated or in accordance with local practice;
- (4) Hand and dispose of contaminated products with care and in accordance with national regulations.

CE 0197 CE Symbol

MD Medical Device



Keep away from sunlight



Do not reuse



Keep dry



Don't use when packing damaged  
and consult instructions for use



Temperature limit

STERILE EO EO Sterile