

## 3D PEMF Therapy Machine

# Lumina 3D-e Clinic - With Battery



## Advantages and purpose of the PEMF machine

- Option to run two different programs simultaneously
- Time and intensity settings for each program
- Memory for programs + rename option: 200 positions
- Unique programs of the new 3D-e generation

## • Lumina Clinic 3D-e Technical parameters

- Technology
- **3D-e patented Biomag® technology**  
second generation of 3D pulsed magnetic therapy
- Programs
- **6 innovative 3D-e programs**  
with therapeutic effect for a wide range of symptoms
- Display
- **Colour touch**
- Adjustable application time

- **5-90 minutes**
- Outputs
- **4 outputs**  
completely independent channels for the latest Biomag® applicators
- Intensity
- **35 mT**  
for each output
- Intensity setting
- **1-100%**
- Independent programmes
- Possibility to run up to 2 independent programs
- Renaming programmes
- **200 positions**  
for renaming programs
- Frequency
- **1-81 Hz**  
for low-frequency magnetic therapy applications
- Repetition of the programme
- **Yes**  
auto-repeat can be started
- **Universal**  
functional worldwide
- Operation on built-in battery
- **Yes**

• -----

- Lumina clinic 3D-e

**3D PEMF therapy**

# **A11P - PEMF Therapy Mat Applicator**



### Advantages and parameters of the mat applicator:

- Deep and surface applications on multiple parts of the body.
- The most demanded and most equipped Biomag mat applicator.
- Length 1 170 mm; width 420 mm; height 40 mm.
- Patented 3D technology.
- Possibility of joining into a closed shape or fixation.
- MIMI 3.0 mT.
- Connector 3 x JACK 3.5 mm.
- Weight: 2 kg.
- -----

### 3D PEMF therapy

## A9P - PEMF Therapy Applicator



### Advantages and parameters of the applicator

- Intensive topical application.
- Shaped loical head applicator.
- Length 800 mm; width 160 mm; height 550 mm.
- Specific use in therapeutic practice.
- Connector 1 x JACK 3.5 mm.
- Weight 0,58 kg.
- -----