


Manufacturer site: No.600, Sec.1, Dong Shan Rd., Dong Shan Township, Yilan County, Taiwan TEL: +886 3 959 4777 FAX: +886 3 959 4707

Taipei office: Rm. 1420, 14F., No. 206, Sec. 1, Keelung Rd., Xinyi Dist., Taipei City, Taiwan








#### **Unique water trap design**

- Circuit auto-seal mechanism, care-free water pouring
  - Less ventilator alarma
  - Less opportunistic infection<sup>b</sup>
- a. Compared with non auto-seal water trap model
- b. Respiratory Care May 2019, 64 (5) 545-549



#### **Patented auto-refill chamber**

- Automatically keep water in a steady level<sup>a</sup>
- 100% tested during production, hassle free, load and forget<sup>a</sup>
- PP materil, zero BPA, zero interference to human body, safe for all ages
- a. DC Medical, production and inspection standard operation process



No water

Adding Water: Pausing water safely



#### **Optimal wire length**

- Efficient and stable heating
- Superior heating function, 21% faster in reaching 37 degree celcius<sup>a</sup>
- Less condensation inside circuit, 20% reduction in the volume of water<sup>b</sup>
- a. Compared with economic alternative
- b. Respiratory Care May 2019, 64 (5) 545-549

# **Unique combination** of material

- 21%-44% less expansion against pressure
- Air-tight without adhesive/sealant, zero leakage under 60mmH2Ob
- a. Compared with common PP circuit
- b. DC Medical product master file

#### **ADULT**

1 Heated expiratory limb: 150cm

2 Non-heated expiratory limb: 90 cm + 60cm

3 Inspiratory limb: 150 cm, 22ID, MR850 system chamber adaptor

4 Chamber limb: 60 cm


5 55ml Water trap 22 mmOD

6 Auto-fill chamber connection ports: 22 mmOD

Connector kit: straight connector, 22nnOD; elbow connector, 22 mmOD, 22 mmID

**13** Temperature port: 7.6 mm/D


Patient wye connector :22/15 mm





### **NEONATE**

- 1 Heated expiratoryl imb: 140 cm
- 2 Non-heated expiratory limb: 90+60 cm
- 3 Inspiratory limb: 110+30 cm, 22 mmID, MR 850 style humidifier adaptor
  - 4 Chamber limb: 60 cm
  - **5**55ml Water trap 10 mmOD
  - 6 Auto-fill chamber connection ports: 22 mmOD
    - Connector kit
    - 8 Temperature port: 7.6 mmlD
    - 9 Patient wye connector: 15 mmID
      - 10 Pressure port



## **ADULT**

REF	Heated Wire	Description	Package	Accessory
3410	Insp.	Single Heated Wire Circuit	1 set/bag, 10 sets/carton	Chamber included
3420	Insp. & Exp.	Dual Heated Wire Circuit	1 set/bag, 10 sets/carton	Chamber included
3411	Insp.	Single Heated Wire Circuit	1 set/bag, 10 sets/carton	
3421	Insp. & Exp.	Dual Heated Wire Circuit	1 set/bag, 10 sets/carton	
3413	Insp.	Inspiratory Heated Wire Circuit	1 set/bag, 15 sets/carton	Chamber included
3490		Auto-filling Humidification Chamber	1 set/bag, 40 sets/carton	
3493		Manual-filling Humidification Chamber	1 set/bag, 40 sets/carton	



## **NEONATE**

REF	Heated Wire	Description	Package	Accessory
3440	Insp.	Single Heated Wire Circuit	1 set/bag, 20 sets/carton	Chamber, Connector Kit
3450	Insp. & Exp.	Dual Heated Wire Circuit	1 set/bag, 20 sets/carton	Chamber, Connector Kit
3441	Insp.	Single Heated Wire Circuit	1 set/bag, 20 sets/carton	Connector Kit
3451	Insp. & Exp.	Dual Heated Wire Circuit	1 set/bag, 20 sets/carton	Connector Kit
3445	Insp.	Single-heated Wire Circuit with Chamber (Auto-filling) for SLE5000	1 set/bag, 20 sets/carton	Chamber, Connector Kit
3449	Insp.	Single-heated Wire Circuit with Chamber (Manual-filling) for SLE5000	1 set/bag, 20 sets/carton	Chamber, Connector Kit
3446	Insp.	Single Heated Wire Circuit for SLE5000	1 set/bag, 20 sets/carton	Connector Kit
3443	Insp.	Inspiratory Heated Wire Circuit with Chamber (Auto-filling)	1 set/bag, 20 sets/carton	Chamber, Connector Kit
3484		SLE5000 Connector (Restrictor)	1000 pcs/bag, 1 bag/carton	Bulk package
3495		Auto-filling Humidification Chamber	1 set/bag, 40 sets/carton	
3498		Manual-filling Humidification Chamber	1 set/bag, 40 sets/carton	