

## Product Summary Sheet


Catalog/Product/Reference#: **C855-5308**

Description: **Disposable Nemoto 200mL/100mL CT Syringe Dual Pack with 2 Spikes**

### Components / Materials

Part		Component		Material
Name	Qty/Unit	Description	Qty/Ea	
200mL Syringe	1	Barrel	1	Natural PET
		Plunger	1	White Polypropylene (PP)
		Plunger Jacket	1	Black Synthetic Rubber
		Dust Cap	1	White Polyethylene (PE)
100mL Syringe	1	Barrel	1	Natural PET
		Plunger	1	White PP
		Plunger Jacket	1	Black Synthetic Rubber
		Dust Cap	1	White PE
Y-Set w/Coiled Extension	1	1.77mm (0.070") ID x3.25mm (0.128") OD Tubing	1	1500mm (60") Clear PVC
			1	76.2mm (3") [or 101.6mm (4")] Clear PVC
			1	177.8mm (7") [or 127mm (5")] Clear PVC w/Blue Stripe
		Male Luer	1	Clear PVC
		Female Luer	2	Clear Polycarbonate (PC)
		Female Dust Cap	1	White PE
		Male Dust Cap	2	(1) Blue & (1) Natural PP
		Check Valve	1	Translucent Purple PC and Blue Silicone
Y-Connector	1	Clear PVC or ABS		
		Fill Spike	2	Body
Cap	2	Green PE		

### Photo



### Label Representation

Left Panel Right Panel

Distributed by **Nemoto**

**COEUR**



**200mL / 100mL Disposable CT Dual Syringe Pack**  
 Includes (1) Y-Set / (1) Check Valve / (2) Spikes

Distributed by **Nemoto**  
 Nemoto Kiyomoto Co., Ltd.  
 2-27-30 Ichige, Yamaguchi  
 Tokyo 113-0033 Japan

REF. C855-5308  
 200mL / 100mL CT Dual Syringe Pack with Spikes

Contents are sterile and non-pyrogenic unless package has been damaged.  
 内容物は無菌・無芽胞・無熱原です。パッケージに損傷がある場合は使用しないでください。  
 Il contenuto è sterile e apirogenico, a meno che il contenitore non sia danneggiato.  
 El contenido es estéril y apirogénico, a menos que el envase no esté dañado o estropeado.  
 Use immediately and aseptically, before the packaging is deconstructed.  
 直ちに無菌状態で使用し、包装を破棄する前に個人防護具を装着してください。  
 内容物は無菌・無芽胞・無熱原です。パッケージに損傷がある場合は使用しないでください。  
 除菌・無菌状態で使用し、包装を破棄する前に個人防護具を装着してください。


USE BY 有效期 

DATE OF MANUFACTURE 製造日 

LOT NUMBER 批号 

 SINGLE USE ONLY 一次使用  DO NOT RESTERILISE 不要重新消毒

 ATTENTION SEE INSTRUCTIONS 注意 请参考使用说明书

 CE  MDC  ISO

Manufactured by **COEUR**  
 COEUR Inc., Lebanon, TN 37090, USA  
 www.coeurinc.com

**MDC** MDC Medical Supplies  
 Torrance, CA 90503, USA

745-918-1326 Fax: 82

### Packaging / Sterilization

Attribute	Material/Detail
Lid	White Tyvek
Rigid Tray	Natural Polystyrene
Quantity per Case	20 Dual Packs
Approx. Weight per Case	6.6 kg (14.5 lbs)
Shipper Dimensions	56.18cm x 34.92cm x 31cm (22 1/8" x 13 3/4" x 12 13/16")
Sterilant	ETO

### Technical Characteristics

Attribute	Detail
Pressure Rating	2.1MPa (300psi 21kgf/cm2)
Intended Use	To deliver, in low-pressure injection systems, contrast media to the vascular system for diagnostic purposes
Injector with Which Product is to be Used	Nemoto Dual Shot Alpha Injector

This summary sheet has been developed based on an approved manufacturing specification.

Title: VP, Quality & Regulatory Affairs

Signature: 

Date: 02/18/11