


SELF ADHESIVE TYPE 4 STEAM INDICATOR

- PRODUCT CODE** : 4AC4I
- RELATED STANDARDS** : 4A Medical Steam Indicator, indicator is clasified as Type 4 according to ISO 11140-1 standard.
- INTENDED USE** : 4A Medical Steam Indicators Type 4 are allow the end user to identify incorrect packaging, loading problems and sterilization failures. It is easy for documentation due to its self adhesive properties.
- OPERATION CONDITIONS** : 4A Medical Steam Indicators Type 4 are designed to operate at steam cycle.
- CHEMICAL COMPOSITION** : The sterilization indicator ink used for printing contains following components in ppm levels
- | | | |
|---------------------|-----|-------------------|
| Water | 35% | (CAS # 7732-18-5) |
| Diacetone Alcohol | 3% | (CAS # 123-42-2) |
| Proprietary Mixture | 62% | (CAS # N/A) |
- Contains no substance specified as hazardous, which should be taken into account according to EC directives.

Before Sterilization



After Sterilization



- PRODUCT SPECIFICATIONS** : Each box contains 100 pieces in a sheet and 1000 indicators in a pack. The indicator is temperature, time and steam concentration sensitive.
- STORAGE AND DISPOSAL** : The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization colour for 6 month.
- PACKAGING** : 50 packs in cardboard.
- MANUFACTURER** : **DORT-A TIP MALZEMELERI SAN. ITH. IHR. TIC. LTD. STI.**
Balıkhisar Mahallesi, Köyiçi Serpmeleri, No:795 Akyurt, Ankara, Turkey
Tel: +90 312 363 50 52 / 53 Fax: +90 312 363 52 10
Web: <http://www.4amedical.com> E-mail: export@4amedical.com