

/ You can find more information about this product in our web-shop.




TUBE HOLDERS

Simple handling and safety are of the highest importance during daily routine procedures.

In the development of tube holders, special attention is paid to ensuring that the products meet the requirements of the users. Tube Holders should be used

exclusively with VACUETTE® Multiple Use Drawing Needles or VACUETTE® VISIO PLUS Needles, as both components are designed for optimal compatibility with one another. The use of these products with needles from other manufacturers is the responsibility of the user.

- / Standard Tube Holders
- / Speedy Quick Release Holder
- / HOLDEX® Single-Use Holders
- / SAFELINK
- / Blood Transfer Unit
- / Blood Culture Holder




Standard Tube Holders

Holders are used in the daily blood collection routine. Tube Holders should be used exclusively with VACUETTE® Multiple Use Drawing Needles or VACUETTE® VISIO PLUS Needles, as both components are designed for optimal compatibility with one another.



Item no.	Qty. Inner / outer
450201	10 / 1,200
450209	100 / 2,000



Speedy Quick Release Holder


The Speedy Quick Release Holder enables safe and reliable ejection of the used blood collection needle into a suitable sharps container, tested for use up 100 times. Nevertheless, Greiner Bio-One recommends single use.



Item no.	Qty. Inner / outer
450212	10 / 1,000


CONFIDENTIAL

SECTION A-A



- Notes:**
- ① Cavity identification
 - ② Two gates – opposite one another
 - ③ 2x engraving according master, 0.03mm embossed / matt finish (Rz 10-12)

polished
Ra0.1



1:1

Drawn	G. Kofler	Date	10.03.2004	Scale	2:1	Standards:	DIN ISO 8015 DIN ISO 1101	<p>greiner bio-one</p> <small>Greiner Bio-One GmbH, Postfach 101563, 42699 Solingen, Germany</small>
Checked	TE		May 2013					
Eng. No.	sk450201	Rev. No.	01	Sub	A3	Dimensions in millimeters		
Subst.	sk450201.00			Material	PP	Title		
<small>This drawing is the property of Greiner Bio-One and is not to be distributed outside the company without our written permission.</small>		<small>CAD generated drawing and not a final drawing. Do not scale drawing!</small>		Greiner Bio-One Standard tube holder		60_EUR_SK450201_10_038 60_EUR_SK450201_03_08 Sheet 1 of 1		