

Serum Clot Activator Tube



Serum clot activator tube is to obtain serum specimen used for serum determinations in chemistry, serology, and immunohematology.

Vacuum Blood Collection Tube, Pro-coagulation Tube

Material: PET, Glass

Color: Red

Cat. No.	Specification	Volume	Additive	Qty/case (Glass)	Qty/case (PET)
630301	13×75mm	3ml	Clot Activator	1800	1800
630302	13×75mm	4ml	Clot Activator	1800	1800
630303	13×75mm	5ml	Clot Activator	1800	1800
630304	13×100mm	5ml	Clot Activator	1200	1200
630305	13×100mm	6ml	Clot Activator	1200	1200
630306	13×100mm	7ml	Clot Activator	1200	1200
630307	16×100mm	8ml	Clot Activator	900	900
630308	16×100mm	9ml	Clot Activator	900	900
630309	16×100mm	10ml	Clot Activator	900	900