
ORDIN DE PLATA NR.: 1920
TIP.DOC. 1
DATA EMITERII:7 martie 2023
=====

PLATITI: 50-00
LEI: Cincizeci lei 00 bani

=====

PLATITOR: (R) "BIOSISTEM
MLD" S.R.L.
CONTUL DE PLATI/CODUL IBAN
MD95ML000000002251429243
CODUL FISCAL :1010600028048 /

=====

PRESTATORUL PLATITOR
BC"Moldindconbank"S.A. suc."Invest" Chisinau
CODUL BANCII:
:MOLDMD2X329:

=====

BENEFICIAR (R) Institutul
de Medicina Urgenta IMSP
CONTUL DE PLATI/CODUL IBAN
MD55VI022510300000002MDL
CODUL FISCAL :1003600152606 /

=====

PRESTATORUL BENEFICIAR
B.C."VICTORIABANK"S.A.
CODUL BANCII:
:VICBMD2X

=====

DESTINATIA PLATII:Pentru garantia pentru: TIPUL TRANSFERULUI
oferta la procedura de achizi?ie public: NORMAL/URGENT :N
a nr. ocds-b3wdp1-MD-1677658749400 din 0:
9.03.2023

=====

CODUL TRANZACTIEI:001:
DATA PRIMIRII:07/03/2023
DATA EXECUTARII:
SEMNAURILE
EMITENTULUI

CONDUCATOR:Web Poiata Vitalie
MIIGYwYJKoZIhvcNAQcCoIIGVDCCBlACAQExCzAJBgUrDgMCGGUAMAsGCSqGSib:
DQEHAaCCBGwwggRoMIIDUKADAgECAhNHAACjbi1rgFksQ0G4AAAAAKNuMA0GCSq:
Sib3DQEBcWUAMCixIDAeBgNVBAMTF0NFULQxLUNBLU1vbGRpbmRjb25iYW5rMB4:
DTIxMDEyODExMzgwNVVoXDTI0MDEyODExNDgwNVowgZ8xCzAJBgNVBAYTA1EMRA:
gYDVQQIEwdNb2xkb3ZhMREwDwYDVQQHEwhDaGlzaw5hdTEWMBQGA1UEChMNQm1

(semnatura electronica)
CONTABIL-SEF:Web Nasedchin Alexandr
MIIGZwYJKoZIhvcNAQcCoIIGWDCCBlQCAQExCzAJBgUrDgMCGGUAMAsGCSqGSib3:
DQEHAaCCBHAWggRsMIIDVKADAgECAhNHAACjcahRKqbJeg8QAAAAAKNxMA0GCSqG:
Sib3DQEBcWUAMCixIDAeBgNVBAMTF0NFULQxLUNBLU1vbGRpbmRjb25iYW5rMB4X:
DTIxMDEyODExMzkwOFoXDTI0MDEyODExNDkxOFowgaMxCzAJBgNVBAYTA1EMRAW:
YDVQQIEwdNb2xkb3ZhMREwDwYDVQQHEwhDaGlzaw5hdTEWMBQGA1UEChMNQm1v

L.S.
CONDUCATOR:
(semnatura electronica)
(semnatura manuala)
CONTABIL-SEF:
(semnatura manuala)
SEMNAURA PRESTATORUL L.S.

MOTIVUL REFUZULUI L.S.



BC "MOLDINDCONBANK" S.A.

Filiala "Invest"

Republica Moldova, MD-2068
mun. Chișinău, bd. Moscovei, 14/1
Tel. : (373-22) 43-44-81, 43-46-24
Fax : (373-22) 43-44-22
cod: MOLDMD2X329

Data 14. IAN. 2016
Nr. 03/2 - 19/23

Республика Молдова, MD-2068
мун. Кишинэу, бул. Московской, 14/1
Тел. : (373-22) 43-44-81, 43-46-24
Факс : (373-22) 43-44-22
код: MOLDMD2X329

Filiala „Invest” BC „Moldindconbank” SA confirmă existența contului curent în moneda națională al **“BIOSISTEM MLD” S.R.L. (c/f 1010600028048)**, cu **IBAN MD95ML000000002251429243**.

Codul băncii MOLDMD2X329.

Director

Nina Turcan

Director financiar



Nina Balmuş

Ex. Diana Brinza
Tel. 43-45-96

REPUBLICA



MOLDOVA

CERTIFICAT DE ÎNREGISTRARE

Societatea cu Răspundere Limitată "BIOSISTEM MLD"
— ESTE ÎNREGISTRATĂ LA CAMERA ÎNREGISTRĂRII DE STAT —

Numărul de identificare de stat - codul fiscal
1010600028048

Data înregistrării

12.08.2010

Data eliberării

12.08.2010

Svirepova Ludmila, registrator

*Funcția, numele, prenumele persoanei
care a eliberat certificatul*

L. Svirepova
semnătura

MD 0101250





I.P. "AGENȚIA SERVICII PUBLICE"

Departamentul înregistrare și licențiere a unităților de drept

EXTRAS

din Registrul de stat al persoanelor juridice

nr. 8506 din 28.04.2021

Denumirea completă: **Societatea cu Răspundere Limitată «BIOSISTEM MLD».**

Denumirea prescurtată: **«BIOSISTEM MLD» S.R.L.**

Forma juridică de organizare: **Societate cu Răspundere Limitată.**

Numărul de identificare de stat și codul fiscal: **1010600028048.**

Data înregistrării de stat: **12.08.2010.**

Sediul: **MD-2001, str. Albișoara, 16/1, ap.(of.) 7, mun. Chișinău, Republica Moldova.**

Obiectul principal de activitate:

- 1 Activitatea farmaceutică;**
- 2 Importul, fabricarea, comercializarea, asistența tehnică și (sau) reparația dispozitivelor medicale și (sau) a opticii;**
- 3 Acordarea asistenței medicale de către instituțiile medico-sanitare private;**
- 4 Comerțul cu ridicata al calculatoarelor, echipamentelor periferice și software-ului;**
- 5 Întreținerea și repararea mașinilor de birou și a tehnicii de calcul;**
- 6 Consultații în domeniul sistemelor de calcul.**

Capitalul social: **5400 lei.**

Administrator: POIATA VITALIE,

Asociați:

- 1. POIATA VITALIE 33,40 %**
- 2. NASEDCHIN ALEXANDR 33,30 %**
- 3. KOJEVNIKOV DMITRII 33,30 %.**

Prezentul extras este eliberat în temeiul art. 34 al Legii nr. 220-XVI din 19 octombrie 2007 privind înregistrarea de stat a persoanelor juridice și a întreprinzătorilor individuali și confirmă datele din Registrul de stat la data de: 28.04.2021.

Specialist coordonator
tel. 022-207-840



Lazari Aliona



EB 0358735

Lista fondatorilor Biosistem-mld SRL

Nr.	Nume, Prenume	IDNP
1.	Vitalie Poiata	0983103892591
2.	Alexandr Nasedchin	2002001070747
3.	Dmitrii Kojevnikov	0972305012362

Staplers | Hemostats | Energy Devices | Hernia Repair Meshes | Ligation Clips | Adhesives | Access Devices
IUDs | Umbilical Cotton Tape | Bonewax

Meril Infrastructure

Situated on 180 acres of land and built over a 300,000 sq.ft. area, Meril's ISO 13485 & cGMP certified ultra-modern manufacturing facility is constructed to satisfy the stringent needs of our Cardiovascular, Orthopedics, Diagnostics and Endo-Surgery businesses, as well as our R&D endeavors.

Fully integrated manufacturing systems ensure backward integration, flawless man-material movement and complete control over processes to eliminate production errors. All manufacturing and sterilization processes are conducted in-house, in addition to the analytical and microbiological QA/QC tests required to meet the world-class production standards.

Meril Park, Vapi



The Meril Organization

Meril's 4000+ personnel are a strong, experienced team comprising designers, engineers, chemists, microbiologists, regulatory affair experts, R&D scientists, clinical affairs experts, legal, finance, sales and marketing professionals with innovative capabilities. The team is continuously striving towards improving and saving lives every day by providing revolutionary solutions for diagnosis, prevention and treatment.

Currently Meril has direct sales & representative office in India, Germany, Russia, USA, Brazil, UK, Turkey, South Africa, China and Australia with distribution network spread over an additional 102 countries in 6 continents.



Meril Academy

Meril has a strong commitment towards spreading awareness, imparting education and training of all our stake holders which range from patients to their care providers.

Meril Training Academy has a large campus spread over 100,000 sq. ft. at Vapi which allows sufficient room for doctors and para-medical personnel to be trained by global KOL faculty while using innovative and modern training modalities.

Meril Academy is dedicated to dissemination of scientific knowledge amongst Physicians, Surgeons, Technicians and Paramedic Staff across the world with varying training requirements using state-of-the-art and best in class teaching tools and programs.

The training programs are comprised of various reflective, interactive and experiential activities based workshops comprising of advanced simulators based learning, live cases, and distance learning, modular learning courses and includes Wet labs (human anatomy phantom) etc



MERIL STUDIO

Meril studio spans a broad operational canvas ranging from Cardiovascular Devices to Orthopaedics, In-Vitro Diagnostics, Endo-Surgery and ENT.

EXECUTIVE LOUNGE

The executive lounge is a well-lighted seating space inside the academy for executives.

MERIL SIMULATION ROOM

Meril Simulation room spread up in 3000 sq. ft. area is a dedicated area for Cardiology, Orthopaedic and Endo-Surgery practitioners.

TAKSHASHILA AUDITORIUM

The Auditorium is an acoustically designed venue of international quality. Equipped with state-of-the-art infrastructure enabling direct telecast, the auditorium is ideal for large conferences, seminars and symposiums.

SUSHRUTA HALL

It is a well resourced, well lighted, wi-fi enabled and spacious conference hall equipped with state-of-the-art visual and interactive tools for smoother management of a wide range of events.

OPERATION THEATRE

The operation theatre provides healthcare professionals with opportunities to make mistakes, refine techniques and ensure theatre effectiveness.



Meril Endo-Surgery

In line with the objective to design and develop novel, clinically relevant, state-of-art surgical devices, Meril is expanding its Endo-Surgery franchise. At Meril, we strive to create and develop products through innovative thinking, cutting edge technologies and collaboration with medical professionals and healthcare organizations, with an aim to identify clinical needs and translate them into market-leading solutions.

Meril's vision is to become a 'one-stop' solution for surgeons by providing them with an extensive product portfolio comprising absorbable & non-absorbable surgical sutures, tissue sealants, absorbable hemostats, hernia repair meshes, intra uterine devices, energy devices and mechanical closure devices.



Contents

Synthetic Absorbable Sutures	
MITSU	15
MITSU C+	47
MITSU FST	53
MEGASORB	59
MEGASORB PLUS	69
MEGASORB RAPID	71
FILAPRON	75
FILAXYN	79
Natural Absorbable Sutures	
MERISOFT Plain Catgut	89
MERISOFT Chromic Catgut	93
Natural Non-Absorbable Sutures	
FILASILK	99
Synthetic Non-Absorbable Sutures	
MERISTEEL	111
FILAMIDE	115
FILAPROP / P	127
MERICRON XL / P	149



Wound Closure

Wounds which do not heal naturally are treated with wound closure devices. The affected tissues need to be kept in apposition until the healing process provides the wound with proper strength to stand against the stress without any assistance. A wound can be treated with the help of various devices such as sutures, staplers, topical adhesives or clips. Different tissues have their own unique functions and requirements. When one decides on a particular wound closure device, one must take care of the healing needs of the particular tissue which is being treated.

The most common device used for wound closure is suture. Sutures, in general term 'stitches', are sterilized surgical threads, which are used to ligate blood vessels or to approximate tissues.

Sutures are generally categorized as:

- Absorbable and Non-Absorbable
- Natural and Synthetic
- Monofilament and Braided

Absorbable sutures have a tendency to degrade as the wound heals. Non-Absorbable sutures oppose the body's efforts to dissolve them. On complete healing of the wound, surgeons remove the non-absorbable sutures. Meril Endo-Surgery has a wide range of synthetic absorbable sutures which are made from Polyglactin 910 (MITSU family). Polydioxanone (FILAXYN), Poliglecaprone 25 (FILAPRON) Polyglycolic Acid (MEGASORB family) as well as Barbed Sutures with two variants (Polydioxanone and Poly Glycolide-Co-Caprolacrone). Among non absorbable sutures Meril portfolio covers Polypropylene (FILAPROP), Polyamide (FILAMIDE) and Polyester (MERICRON XL). In natural absorbable sutures, Meril portfolio covers Chromic and Plain Catgut (MERISOFT) and in natural non-absorbable we have silk (FILASILK).


For placement of sutures, surgical needles are used, which are made of stainless steel. Different types and sizes of needles are used depending on the area of application. We are dedicated to assist surgeons in tissue repair by giving them a broad gamut of products with our unmatched quality and service.

Meril Endo-Surgery Needles

In pursuit of offering best quality products, Meril sutures are paired with the superior stainless steel needles for an unparalleled surgical experience.

Meril sutures come with sharper and stronger needles (300 series Stainless Steel) specifically designed to meet the ever growing needs in all the surgical disciplines.

With every needle swaged using Multi Stroke Technology (MST) and tested for a needle pull check, Meril needles are distinguished with a high penetration coating and a square body design for a firm grip.



Round body needle, recommended for blood vessels



Round body/ Taper needle with black finish for enhanced visibility under glaring OT lights



Taper needle, for tough, difficult-to-penetrate tissues



Round body / Spatulated needle for finer surgeries



For easy penetrations pass after pass



Needle tip with cross sectional geometry giving improved penetration and surgical outcome

Needle Shapes



Needle Specifications

Type	Taper	Reverse cutting	Cutting	Flash Point	Trocar Point	Spatula	Blunt	Taper cut
Design								
Point								
Shape	Typical Applications							
Straight	Skin/subcutaneous layer closure; Purse string							
1/4 Circle	Ophthalmic/Microsurgical procedures							
3/8 Circle	Skin/subcutaneous layer/abdominal wall closure; Ophthalmic/Vascular/Plastic procedures							
1/2 Circle	Peritoneum closure; Gastrointestinal/Respiratory/Urological procedures							
5/8 Circle	Pelvic & Oral cavity; Urological procedures							
J-Shape	Laparotomy closure							
Ski	Laparoscopic procedures							

Understanding the Suture Codes

Meril code logic represents the necessary information needed to identify the required suture-needle combination. The suture codes are alphanumeric in nature to indicate the suture material, the suture size, and the reference code. The alphabetical suffixes are a key to understand the additional features about the suture and/or the needle.

Material Prefixes

ABBREVIATION	MEANING
PGN	Polyglactin 910
PGT	Polyglactin 910 (with Triclosan)
PGP	Polyglactin 910 (with Chlorhexidine)
PGF	Polyglactin 910 Fast
MS	Polyglycolic Acid
PCL	Poliglecaprone 25
PDX	Polydioxanone
PNDX	Barbed Poly (P-Dioxanone)
PNCL	Barbed Polyglycolic Acid-poly (Glycolide-Co-Caprolacrone)
ME	Plain Catgut
ME	Chromic Catgut
SLK	Silk
SLS	Silk Non Needled
NYL	Polyamide
NYS	Polyamide Non Needled
ME	Polyester
PME	Polyester with Pledgets
PPL	Polypropylene
PPG	Polypropylene with Pledgets
STC	Steel Wire

*In case of

- Polyester (Mericon XL), codes start with ME6XX, ME6XXX, ME4XXX, MEXXXX; and polyester with pledgets (Mericon XL P), codes start with PMEXXX or MMEXXXX
- Chromic Catgut (Merisoft Chromic), codes start with ME40XX, ME42XX, ME50XX
- Chromic Catgut non needled (Merisoft Chromic), codes start with SS22XX
- Plain Catgut (Merisoft Plain), codes start with ME50XX
- Plain Catgut non needled (Merisoft Plain), codes start with SS20XX


Suture Sizes

- For sutures having a single-digit size, the size is written as: zero followed by the size digit.
E.g.: Size 01 (read as size 1 according to USP), Size 03 (read as size 3 according to USP).
- For finer sutures, the size is written as: the size digit followed by a zero.
E.g.: Size 60 (read as size 6-0 according to USP), Size 70 (read as size 7-0 according to USP), Size 100 (read as size 10-0 according to USP)

Suffixes

ABBREVIATION	MEANING
B/BL	Black Needle
BX	Black Needle; Extra Length Suture
DN	Double Needle
DNL	Double Needle; Long Length Suture
DS	Doctor Special Code
G	Green Suture
ML	Multi Strands
SL	Short Length Suture
U	Undyed Suture
V	Violet Suture
W	White Suture
X	Extra Pair of Suture
XF	Extra Fine Needle
XL	Extra Length Suture

Decoding the Codes



Abbreviation for Specification

ABBREVIATION	MEANING
HC	1/2 CIRCLE
QC	1/4 CIRCLE
CU	3/8 CIRCLE
FE	5/8 CIRCLE
RC	REVERSE CUTTING
CT	CUTTING
RB	ROUND BODIED
TC	TAPER CUT
ST	STRAIGHT
SM	SPATULA
TRO	TROCAR
BP	BLUNT POINT
TP	PLEDGET
VIO	VIOLET
UND	UNDYED
BD	BIDIRECTIONAL

Needle Profiles	Point type	Symbol
	Round bodied/Taper point	
	Conventional cutting	
	Reverse cutting	
	Flash Point	
	Tapercutting	
	Trocac Point	
	UltraPoint spatula curved	
	Straight	

MITSUTM

Polyglactin 910 Suture

Properties

- 75% of the original is retained until initial 14 Days
- 50% of the original is retained until initial 21 Days
- There is a subsequent loss between four to five weeks post implantation
- Complete absorption usually takes place between 56 to 70 Days

Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Ophthalmic procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction



Product Overview

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Composition	: Poly (glycolide-co-l-lactide 90/10)
Coating	: Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate
Sizes	: Dyed - USP 8/0 (metric 0.4) - USP 2 (metric 5)
	: Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Meril Endo-Surgery Product Catalogue | **Wound Closure** 17

Meril Endo-Surgery Product Catalogue | **Wound Closure** 21

Meril Endo-Surgery Product Catalogue | **Wound Closure** 25

Meril Endo-Surgery Product Catalogue | **Wound Closure** 27

Meril Endo-Surgery Product Catalogue | **Wound Closure** 29

Meril Endo-Surgery Product Catalogue | **Wound Closure** 31

Meril Endo-Surgery Product Catalogue | **Wound Closure** 37

Meril Endo-Surgery Product Catalogue | **Wound Closure** 39

MITSU C+™

Polyglactin 910 Suture with Chlorhexidine

Properties

- Retention of 75% knot tensile strength after 14 days
- Retention of 50% knot tensile strength after 21 days
- Complete mass absorption within 56 to 70 days

Indications


- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction



Product Overview

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Composition	: Poly(glycolide-co-l-lactide 90/10)
Coating	: Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate + antibacterial agent Chlorhexidine
Sizes	: USP 5/0 (metric 1) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU C+ (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 WITH CHLORHEXIDINE)								MITSU C+ (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 WITH CHLORHEXIDINE)							
PGLA C PLUS VIO 2-0 X 90-40MM HC RC	PGP20 2382	2-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO	PGLA C PLUS UND 3-0 X 70-24MM CU CT	PGP30 936U	3-0		24 mm 3/8 Circle Conventional Cutting	70	12	UND
PGLA C PLUS VIO 2-0 X 70-22MM HC RB	PGP200328	2-0		22 mm 1/2 Circle Round Body	70	12	VIO	PGLA C PLUS VIO 3-0 X 45-17MM HC RB	PGP300306	3-0		17 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS VIO 2-0 X 70-26MM HC RB	PGP200333	2-0		26 mm 1/2 Circle Round Body	70	12	VIO	PGLA C PLUS VIO 3-0 X 70-36MM HC RB	PGP300338	3-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 2-0 X 70-36MM HC RB	PGP200339	2-0		36 mm 1/2 Circle Round Body	70	12	VIO	PGLA C PLUS UND 3-0 X 70-60MM ST CT	PGP300663U	3-0		60 mm Straight Conventional Cutting	70	12	UND
PGLA C PLUS VIO 2-0 X 90-36MM HC TC	PGP200517	2-0		36 mm 1/2 Circle Taper Cut	90	12	VIO	PGLA C PLUS VIO 4-0 X 70-20MM HC RB	PGP40 2304	4-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 2-0 X 90-26MM CU CT	PGP20 2404U	2-0		26 mm 3/8 Circle Conventional Cutting	90	12	UND	PGLA C PLUS VIO 4-0 X 70-17MM HC RB	PGP400304	4-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 2-0 X 90-36MM HC RB	PGP20 945U	2-0		36 mm 1/2 Circle Round Body	90	12	UND	PGLA C PLUS UND 4-0 X 70-60MM ST CT	PGP400662U	4-0		60 mm Straight Conventional Cutting	70	12	UND
PGLA C PLUS VIO 2-0 X 70-36MM FE RB	PGP20375	2-0		36 mm 5/8 Circle Round Body	70	12	VIO	PGLA C PLUS VIO 4-0 X 70-26MM HC RB	PGP400774	4-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 2-0 X 70-26MM FE RB	PGP206023	2-0		26 mm 5/8 Circle Round Body	70	12	VIO	PGLA C PLUS VIO 5-0 X 70-17MM HC RB	PGP500303	5-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 3-0 X 90-25MM CU CT	PGP30 2328U	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	UND	PGLA C PLUS VIO 5-0 X 45-16MM HC RB	PGP50 2303	5-0		16 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS VIO 3-0 X 70-20MM HC RB	PGP30 2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO	PGLA C PLUS VIO 0 X 350 (70X5)	PGP101217	0		Non needled	70	12	VIO
PGLA C PLUS VIO 3-0 X 90-22MM HC CT	PGP30 2472	3-0		22 mm 1/2 Circle Conventional Cutting	90	12	VIO								
PGLA C PLUS UND 3-0 X 70-17MM HC RB	PGP300215	3-0		17 mm 1/2 Circle Round Body	70	12	UND								
PGLA C PLUS VIO 3-0 X 70-22MM HC RB	PGP300327	3-0		22 mm 1/2 Circle Round Body	70	12	VIO								
PGLA C PLUS UND 3-0 X 70-24MM CU RC	PGP304423	3-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND								
PGLA C PLUS VIO 3-0 X 45-26MM HC RB	PGP300773	3-0		26 mm 1/2 Circle Round Body	45	12	VIO								
PGLA C PLUS VIO 3-0 X 70-26MM HC RB	PGP300332	3-0		26 mm 1/2 Circle Round Body	70	12	VIO								
PGLA C PLUS VIO 3-0 X 75-31MM HC RB	PGP300245	3-0		31 mm 1/2 Circle Round Body	75	12	VIO								
PGLA C PLUS VIO 3-0 X 90-36MM HC TC	PGP30 2516	3-0		36 mm 1/2 Circle Taper Cut	90	12	VIO								

MITSU FST™

Polyglactin 910 Fast Suture

Properties

- 50% of the original is retained until initial 5 Days
- approximately 10-14 Days post implantation, the original tensile strength is completely lost
- Complete absorption usually takes place by 42 Days

Indications

- Gynaecology/Obstetric procedures (e.g.: episiotomy repair)
- Urology procedures
- Ophthalmic procedures (e.g.: conjunctiva suturing)
- Odontology procedures (e.g.: oral mucosa)
- Skin closure (intracutaneous, subcutaneous)
- Ligatures






Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Excellent knot security
- Excellent handling properties
- Quick mass absorption



Product Overview

Structure	: Multifilament, Braided
Colour	: Undyed
Composition	: Poly(glycolide-co-l-lactide 90/10)
Coating	: Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate
Sizes	: USP 8/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)								MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 0 X 110-36MM HC TC	PGF10 2763	0		36 mm 1/2 Circle Taper Cut	110	12	UND	PGLA FAST UND 2-0 X 90-36MM HC TC	PGF20 2762	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 0 X 90-40MM HC RB	PGF109977	0		40 mm 1/2 Circle Round Body	90	12	UND	PGLA FAST UND 2-0 X 140-36MM HC RC RB DN	PGF20 2777DN	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	140	12	UND
PGLA FAST UND 0 X 90-48MM HC RB	PGF1016951	0		48 mm 1/2 Circle Round Body	90	12	UND	PGLA FAST UND 2-0 X 75-31MM HC RB	PGF209975	2-0		31 mm 1/2 Circle Round Body	75	12	UND
PGLA FAST UND 0 X 90-36MM HC RB	PGF109036	0		36 mm 1/2 Circle Round Body	90	12	UND	PGLA FAST UND 2-0 X 90-48MM HC RB	PGF209978	2-0		48 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 0 X 60-36MM CU RC	PGF10 2721	0		36 mm 3/8 Circle Reverse Cutting	60	12	UND	PGLA FAST UND 2-0 X 70-26MM HC RB	PGF204171	2-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA FAST UND 0 X 110-40MM HC RB	PGF10 2763RB	0		40 mm 1/2 Circle Round Body	110	12	UND	PGLA FAST UND 2-0 X 90-36MM HC RB	PGF209451	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 0 X 90-36MM HC TC	PGF10 9963	0		36 mm 1/2 Circle Taper Cut	90	12	UND	PGLA FAST UND 2-0 X 135-36MM HC RB RC DN	PGF209321	2-0		36 mm 1/2 Circle Round Body Reverse Cutting	135	12	UND
PGLA FAST UND 0 X 90-30MM CU RC-3D	PGF102287Z	0		30 mm 3/8 Circle Reverse Cutting	90	36	UND	PGLA FAST UND 2-0 X 90-36MM HC RC	PGF209236	2-0		36 mm 1/2 Circle Reverse Cutting	90	12	UND
PGLA FAST UND 1 X 90-40MM HC TC	PGF016190	1		40 mm 1/2 Circle Taper Cut	90	12	UND	PGLA FAST UND 2-0 X 70-19MM CU RC	PGF200428	2-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA FAST UND 1 X 90-40MM HC RB	PGF016191	1		40 mm 1/2 Circle Round Body Heavy Needle	90	12	UND	PGLA FAST UND 2-0 X 60-36MM CU RC	PGF20 2720	2-0		36 mm 3/8 Circle Reverse Cutting	60	12	UND
PGLA FAST UND 1 X 90-36MM HC RB	PGF016192	1		36 mm 1/2 Circle Round Body	90	12	UND	PGLA FAST UND 2-0 X 110-36MM HC TC	PGF20 2762XL	2-0		36 mm 1/2 Circle Taper Cut	110	12	UND
PGLA FAST UND 1 X 90-48MM HC RB	PGF016193	1		48 mm 1/2 Circle Round Body	90	12	UND	PGLA FAST UND 2-0 X 120-36MM HC TC	PGF209948	2-0		36 mm 1/2 Circle Taper Cut	120	12	UND
PGLA FAST UND 1 X 90-36MM HC TC	PGF019190	1		36 mm 1/2 Circle Taper Cut	90	12	UND	PGLA FAST UND 2-0 X 75-26MM CU RC	PGF20 9933	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 0 X 90-45MM HC TC	PGF109947	0		45 mm 1/2 Circle Taper Cut	90	12	UND	PGLA FAST UND 2-0 X 100-36MM HC TC SP	PGF20 2761SP	2-0		36 mm 1/2 Circle Taper Cut	100	12	UND
PGLA FAST UND 2-0 X 90-40MM HC TC	PGF209946	2-0		40 mm 1/2 Circle Taper Cut	90	12	UND	PGLA FAST UND 2-0 X 90-36MM HC TC SP	PGF20 2762SP	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 2-0 X 90-35MM CU RC	PGF209938	2-0		35 mm 3/8 Circle Reverse Cutting	90	12	UND								
PGLA FAST UND 2-0 X 100-36MM HC TC	PGF20 2761	2-0		36 mm 1/2 Circle Taper Cut	100	12	UND								

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)								MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 2-0 X 90-26MM CU RC	PGF209941	2-0		26 mm 3/8 Circle Reverse Cutting	90	12	UND	PGLA FAST UND 3-0 X 75-36MM CU RC	PGF309961	3-0		36 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 2-0 X 90-26MM CU RC-3D	PGF209941Z	2-0		26 mm 3/8 Circle Reverse Cutting	90	36	UND	PGLA FAST UND 3-0 X 70-26MM HC RB EX	PGF304161	3-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA FAST UND 2-0 X 75-24MM CU RC-3D	PGF209275Z	2-0		24 mm 3/8 Circle Reverse Cutting	75	36	UND	PGLA FAST UND 3-0 X 75-25MM FE CT	PGF309930	3-0		25 mm 5/8 Circle Conventional Cutting	75	12	UND
PGLA FAST UND 2-0 X 75-60MM ST CT	PGF209944	2-0		60 mm Straight Conventional Cutting	75	12	UND	PGLA FAST UND 3-0 X 70-26MM HC RB-3D	PGF304161Z	3-0		26 mm 1/2 Circle Round Body	70	36	UND
PGLA FAST UND 3-0 X 75-16MM CU CT	PGF30 9919	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	UND	PGLA FAST UND 3-0 X 75-19MM CU CT	PGF309319	3-0		19 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA FAST UND 3-0 X 90-60MM ST CT	PGF309943	3-0		60 mm Straight Conventional Cutting	90	12	UND	PGLA FAST UND 4-0 X 75-19MM CU RC	PGF409922	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 3-0 X 90-25MM FE RB	PGF309929	3-0		25 mm 5/8 Circle Round Body	90	12	UND	PGLA FAST UND 4-0 X 75-17MM HC RB	PGF409970	4-0		17 mm 1/2 Circle Round Body	75	12	UND
PGLA FAST UND 3-0 X 90-35MM CU RC	PGF309937	3-0		35 mm 3/8 Circle Reverse Cutting	90	12	UND	PGLA FAST UND 4-0 X 75-22MM HC RB	PGF409971	4-0		22 mm 1/2 Circle Round Body	75	12	UND
PGLA FAST UND 3-0 X 70-26MM CU RC	PGF30 2732	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND	PGLA FAST UND 4-0 X 45-13MM CU CT	PGF408351	4-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 45-26MM CU RC	PGF30 2732SL	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	UND	PGLA FAST UND 4-0 X 75-13MM CU RC	PGF409913	4-0		13 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 3-0 X 75-19MM CU RC	PGF309923	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND	PGLA FAST UND 4-0 X 45-16MM CU CT	PGF40 2718	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 75-22MM HC RB	PGF309974	3-0		22 mm 1/2 Circle Round Body	75	12	UND	PGLA FAST UND 4-0 X 75-16MM CU CT	PGF40 9918	4-0		16 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA FAST UND 3-0 X 90-36MM HC RB	PGF309441	3-0		36 mm 1/2 Circle Round Body	90	12	UND	PGLA FAST UND 4-0 X 45-13MM ST CT	PGF409909SL	4-0		13 mm Straight Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 45-22MM HC CT	PGF30 9935	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	UND	PGLA FAST UND 4-0 X 45-16MM CU CT EX	PGF408451	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 70-18MM HC RB	PGF300318	3-0		18 mm 1/2 Circle Round Body	70	12	UND	PGLA FAST UND 4-0 X 45-22MM CU CT	PGF409926	4-0		22 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 45-16MM CU CT	PGF30 2719	3-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND	PGLA FAST UND 4-0 X 45-45MM ST CT	PGF409939	4-0		45 mm Straight Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 70-22MM HC CT	PGF30 2735	3-0		22 mm 1/2 Circle Conventional Cutting	70	12	UND	PGLA FAST UND 4-0 X 45-13MM CU RC	PGF404940	4-0		13 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA FAST UND 3-0 X 70-26MM HC RB	PGF30 2764	3-0		26 mm 1/2 Circle Round Body	70	12	UND	PGLA FAST UND 4-0 X 45-22MM HC CT	PGF409926HC	4-0		22 mm 1/2 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 45-60MM ST CT	PGF309943SL	3-0		60 mm Straight Conventional Cutting	45	12	UND	PGLA FAST UND 5-0 X 75-17MM HC RB	PGF50 9969	5-0		17 mm 1/2 Circle Round Body	75	12	UND
PGLA FAST UND 3-0 X 75-26MM CU RC	PGF309932	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND	PGLA FAST UND 5-0 X 45-13MM CU CT	PGF504931	5-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
								PGLA FAST UND 5-0 X 45-11MM CU CT	PGF504901	5-0		11 mm 3/8 Circle Conventional Cutting	45	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 5-0 X 45-11MM CU RC	PGF50 9915	5-0		11 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA FAST UND 5-0 X 75-16MM CU RC	PGF502296	5-0		16 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 5-0 X 70-12MM CU RC	PGF502493	5-0		12 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA FAST UND 5-0 X 45-16MM CU CT	PGF508441	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 2-0 X 150 (50X3)	PGF200516	2-0		Non needled	50	12	UND
PGLA FAST UND 3-0 X 75-24MM CU RC	PGF300935	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	UND

MEGASORB™

Polyglycolic Acid Suture

Properties

- 70%-80% of the original is retained until initial 14 Days
- 50% of the original is retained until initial 21 Days
- There is a subsequent loss between four to five weeks post implantation
- Complete absorption usually takes place between 60 and 90 Days

Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Ophthalmic procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- High tensile strength
- Excellent knotting ability
- Smooth passage through tissue
- Easy handling



Product Overview

Structure : Multifilament, Braided
Colour : Dyed, Undyed



Composition : Polyglycolic acid
Coating : Polycaprolactone and Calcium Stearate
Sizes : Violet - USP 8/0 (metric 0.4) - USP 2 (metric 5)
: Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)

Type of absorption : Hydrolysis
Sterilisation : Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)								MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)							
PGA VIO 0 X 90-30MM HC RB	MS2338	0		30 mm 1/2 Circle Round Body	90	12	VIO	PGA VIO 0 X 100-48MM HC RB	MS109250XL	0		48 mm 1/2 Circle Round Body	100	12	VIO
PGA VIO 0 X 90-40MM HC RB	MS2346	0		40 mm 1/2 Circle Round Body	90	12	VIO	PGA VIO 0 X 75-50MM HC RB	MS109251	0		50 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 140-40MM HC RB TC DN	MS2346DN	0		40 mm 1/2 Circle Round Body Taper Cut	140	12	VIO	PGA VIO 0 X 45-40MM HC TC	MS2518SL	0		40 mm 1/2 Circle Taper Cut	45	12	VIO
PGA VIO 0 X 180-40MM HC RB TC DN	MS2346DNL	0		40 mm 1/2 Circle Round Body Taper Cut	180	12	VIO	PGA UND 0 X 90-40MM HC RB	MS2346U	0		40 mm 1/2 Circle Round Body	90	12	UND
PGA VIO 0 X 45-40MM HC RB	MS2346S	0		40 mm 1/2 Circle Round Body	45	12	VIO	PGA VIO 0 X 70-30MM HC RB	MS2342	0		30 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 0 X 90-40MM HC TC	MS2518	0		40 mm 1/2 Circle Taper Cut	90	12	VIO	PGA VIO 0 X 75-45MM HC RB	MS109244	0		45 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 90-36MM HC RC	MS2534	0		36 mm 1/2 Circle Reverse Cutting	90	12	VIO	PGA VIO 0 X 75-48MM HC RB	MS109250	0		48 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 75-40MM HC RB	MS100352	0		40 mm 1/2 Circle Round Body	75	12	VIO	PGA VIO 1 X 90-40MM HC RB TN	MS2347	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 0 X 70-30MM HC RB EX	MS100808	0		30 mm 1/2 Circle Round Body	70	12	VIO	PGA VIO 1 X 140-40MM HC RB TC DN	MS2347DN	1		40 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGA VIO 0 X 70-22MM HC RB	MS100329	0		22 mm 1/2 Circle Round Body	70	12	VIO	PGA VIO 1 X 180-40MM HC RB TC DN	MS2347DNL	1		40 mm 1/2 Circle Round Body Taper Cut	180	12	VIO
PGA VIO 0 X 70-30MM CU RC	MS100587	0		30 mm 3/8 Circle Reverse Cutting	70	12	VIO	PGA VIO 1 X 45-40MM HC RB	MS2347S	1		40 mm 1/2 Circle Round Body	45	12	VIO
PGA VIO 0 X 100-40MM HC RB	MS2346XL	0		40 mm 1/2 Circle Round Body	100	12	VIO	PGA VIO 1 X 70-40MM HC RB	MS2359	1		40 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 0 X 70-36MM HC RB	MS100260	0		36 mm 1/2 Circle Round Body	70	12	VIO	PGA VIO 1 X 90-40MM HC RB TH	MS2360	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 0 X 75-26MM HC RB	MS100334	0		26 mm 1/2 Circle Round Body	75	12	VIO	PGA VIO 1 X 90-40MM HC RC	MS2421	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 0 X 100-48MM HC RB	MS109250XL	0		48 mm 1/2 Circle Round Body	100	12	VIO								

Meril Endo-Surgery Product Catalogue | **Wound Closure** 65

Meril Endo-Surgery Product Catalogue | **Wound Closure** 67

MEGASORB PLUS™

Polyglycolic Acid Suture with Chlorhexidine

Properties

- Retention of 70% knot tensile strength after 14 days
- Complete mass absorption within 60 to 90 days

Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Smooth knot run-down
- Gentle knot repositioning
- High knot holding capacity
- Minimized inflammatory reaction



Also available **MEGASORB T+** Polyglycolic Acid with Triclosan

Product Overview

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
Chemical composition	: Polyglycolic acid
Coating	: Polycaprolactone + calcium stearate + Antibacterial agent Chlorhexidine
Sizes	: USP 8/0 (metric 0.4) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB PLUS (POLYGLYCOLIC ACID SUTURE WITH CHLORHEXIDINE)							
PGA C PLUS VIO 0 X 90-40MM HC RB	MSP2346	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 1 X 90-40MM HC RB TH	MSP2360	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 1 X 90-40MM HC RB TN	MSP2347	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 1 X 90-40MM HC RC	MSP2421	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA C PLUS VIO 2-0 X 90-30MM HC RB	MSP2317	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 3-0 X 70-20MM HC RB	MSP2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO

MEGASORB RAPID™

Polyglycolic Acid Fast Suture

Properties

- 50% of the original is retained until initial 6 to 8 Days
- There is a subsequent loss by 17 days post implantation
- Complete absorption usually takes place by 42 Days

Indications

- Gynaecology/Obstetric procedures (e.g.: episiotomy repair)
- Urology procedures
- Ophthalmic procedures (e.g.: conjunctiva suturing)
- Odontology procedures (e.g.: oral mucosa)
- Paediatric procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Good knot security
- Excellent handling properties
- Quick mass absorption



Product Overview

Structure	: Multifilament, Braided
Colour	: Undyed
	
Chemical composition	: Polyglycolic acid
Sizes	: USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB RAPID (POLYGLYCOLIC ACID FAST SUTURE)								MEGASORB RAPID (POLYGLYCOLIC ACID FAST SUTURE)							
PGA RAPID UND 0 X 90-36MM HC RB	MS109037	0		36 mm 1/2 Circle Round Body	90	12	UND	PGA RAPID UND 4-0 X 75-19MM CU RC	MS409922	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA RAPID UND 1 X 90-40MM HC RB	MS019140	1		40 mm 1/2 Circle Round Body	90	12	UND	PGA RAPID UND 4-0 X 75-24MM CU RC	MS409924	4-0		24 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA RPID UND 2-0 X 140-36MM HC RC RB DN	MS2777	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	140	12	UND	PGA RAP UND 4-0 X 75-26MM CU CT	MS9918	4-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGA RPID UND 2-0 X 90-36MM HC TC	MS2762	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND	PGA RAPID UND 4-0 X 70-19MM CU CT	MS406529	4-0		19 mm 3/8 Circle Conventional Cutting	70	12	UND
PGA RAPID UND 2-0 X 90-26MM CU RC	MS209933	2-0		26 mm 3/8 Circle Reverse Cutting	90	12	UND	PGA RAPID UND 4-0 X 75-17MM HC RC	MS409418	4-0		17 mm 1/2 Circle Reverse Cutting	17	12	UND
PGA RAPID UND 2-0 X 90-36MM HC RB	MS209451	2-0		36 mm 1/2 Circle Round Body	90	12	UND	PGA RAPID UND 5-0 X 45-12MM CU RC	MS9915	5-0		12 mm 3/8 Circle Reverse Cutting	45	12	UND
PGA RAPID UND 2-0 X 75-30MM HC RB	MS209975	2-0		30 mm 1/2 Circle Round Body	75	12	UND	PGA RAPID UND 2-0 X 45 (45X3)	MS200616	2-0		Non needled	45	12	UND
PGA RAPID UND 2-0 X 90-36MM HC RC	MS209236	2-0		36 mm 1/2 Circle Reverse Cutting	36	12	UND	PGA RAPID UND 3-0 X 45 (45X3)	MS300285	3-0		Non needled	45	12	UND
PGA RAPID UND 3-0 X 70-26MM CU RC	MS302732	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND								
PGA RAPID UND 3-0 X 70-22MM HC RB	MS309974	3-0		22 mm 1/2 Circle Round Body	70	12	UND								
PGA RAPID UND 3-0 X 75-19MM CU RC	MS309923	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND								
PGA RAP UND 3-0 X 75-26MM CU RC	MS9932	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND								
PGA RAPID UND 3-0 X 45-22MM HC RC	MS306230	3-0		22 mm 1/2 Circle Reverse Cutting	45	12	UND								
PGA RAPID UND 3-0 X 70-26MM HC RB	MS304161	3-0		26 mm 1/2 Circle Round Body	70	12	UND								
PGA RAPID UND 3-0 X 45-22MM HC CT	MS309927	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	UND								
PGA RAPID UND 4-0 X 45-16MM CU CT	MS408451	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND								
PGA RAPID UND 4-0 X 75-17MM HC RB	MS409970	4-0		17 mm 1/2 Circle Round Body	75	12	UND								

FILAPRON™

Poliglecaprone 25 Suture

Properties

- 80%-60% of the original is retained until initial 7 Days (DYED)
- 60%-50% of the original is retained until initial 7 Days (UNDYED)
- 40%-30% of the original is retained until initial 14 Days (DYED)
- 30%-20% of the original is retained until initial 14 Days (UNDYED)
- There is a subsequent loss 28 days post implantation
- Complete absorption usually takes place between 90 to 120 Days

Indications

When short wound support is desirable

- Skinclosure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Short-term absorption profile
- High knot pull tensile strength
- Superior handling features



Product Overview

Structure	: Monofilament
Colour	: Dyed, Undyed
Composition	: Poly(glycolide-co-caprolactone)PGCL
Sizes	: USP 6/0 (metric 0.7) - USP 1 (metric 4)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Meril Endo-Surgery Product Catalogue | **Wound Closure** 73

Meril Endo-Surgery Product Catalogue | **Wound Closure** 75

FILAXYN™

Polydioxanone Suture

Properties

- 70% of the original is retained until initial 14 Days
- 50%-60% of the original is retained until initial 28 Days
- There is a subsequent minimal absorption until about 90 th post implantation day
- Complete absorption usually takes place between 180 - 210 days

Indications

- When prolonged wound support is desirable
- Abdominal wall closure
- Orthopedic procedures
- Oncological procedures

Advantages

- High initial tensile strength
- Predictable absorption
- Prolonged wound support
- Smooth passage through tissues
- Minimal tissue reaction
- Good handling characteristics
- Standard monofilament knot for security



Product Overview

Structure	: Monofilament
Colour	: Violet
Composition	: Poly (p-dioxanone)
Sizes	: USP 6/0 (metric 1.5) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Meril Endo-Surgery Product Catalogue | **Wound Closure** 85

Meril Endo-Surgery Product Catalogue | **Wound Closure** 87

MERISOFT PLAIN™

Plain Catgut Suture

Properties

The tensile strength is maintained for not more than 7 to 10 days post implantation, and its absorption is complete within 60 to 90 days.

Indications

MERISOFT PLAIN™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability
- Rigidly controlled spinning & drying to ensure higher tensile strength



Product Overview

Primary packaging	: Packed in an aluminum foil pouch
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Yellow-tan
Sterilization	: Gamma radiation

Meril Endo-Surgery Product Catalogue | **Wound Closure** 91

MERISOFT CHROMIC™

Chromic Catgut Suture

Properties

MERISOFT CHROMIC™ gut suture is treated with a chromium salt solution to resist the enzymatic action, prolonging the absorption time to 90 to 120 days.

Indications

MERISOFT CHROMIC™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability



Product Overview

Primary packaging	: Packed in an aluminum foil pouch
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Brown dyed
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)								MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)							
CHROMIC CATGUT 0 X 76-30MM CU RB	ME4217	0		30 mm 3/8 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 1 X 75-30MM HC RB	ME014995	1		30 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM CU RB	ME4221	0		40 mm 3/8 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 1 X 100-40MM HC TRO	ME4239	1		40 mm 1/2 Circle Trocar	100	12	CHROMIC
CHROMIC CATGUT 0 X 100-45MM HC RB	ME4229	0		45 mm 1/2 Circle Round Body	100	12	CHROMIC	CHROMIC CATGUT 1 X 76-63MM CU RB	Me4267	1		63 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-30MM HC RB	ME4242	0		30 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 1 X 76-50MM HC RC	ME4273	1		50 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC RB TN	ME4246	0		40 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 1 X 100-40MM HC TC	ME4297	1		40 mm 1/2 Circle Taper Cut	100	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC RB TH	ME4258	0		40 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 2 X 100-45MM HC RB	ME4228	2		45 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 0 X 76-36MM HC RB	ME104036	0		36 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 2 X 100-50MM HC RB	ME024250	2		50 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 0 X 100-40MM HC TC	ME4279	0		40 mm 1/2 Circle Taper Cut	100	12	CHROMIC	CHROMIC CATGUT 2 X 100-40MM HC TRO	ME4244	2		40 mm 1/2 Circle Trocar	100	12	CHROMIC
CHROMIC CATGUT 0 X 76-45MM ST RB	ME4205	0		45 mm Straight Round Body	76	12	CHROMIC	CHROMIC CATGUT 2 X 100-65MM HC RB	ME4278	2		65 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 0 X 76-45MM CU RB	ME4225	0		45 mm 3/8 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 2 X 75-40MM HC RC	ME024244	2		40 mm 1/2 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 0 X 76-45MM HC RB	ME4250	0		45 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 2-0 X 76-30MM CU RB	ME4216	2-0		30 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC TC	ME4257	0		40 mm 1/2 Circle Taper Cut	76	12	CHROMIC	CHROMIC CATGUT 2-0 X 76-30MM HC RB	ME4241	2-0		30 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC RC	ME104040	0		40 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC	CHROMIC CATGUT 2-0 X 152-36MM HC RCRB DN	ME4277	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	152	12	CHROMIC
CHROMIC CATGUT 1 X 76-45MM CU RB	ME4226	1		45 mm 3/8 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 2-0 X 75-40MM HC RB	ME200991	2-0		40 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 1 X 100-45MM HC RB	ME4227	1		45 mm 1/2 Circle Round Body	100	12	CHROMIC	CHROMIC CATGUT 2-0 X 76-26MM HC RB	ME204225	2-0		26 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 1 X 100-50MM HC RB	ME4248	1		50 mm 1/2 Circle Round Body	100	12	CHROMIC								
CHROMIC CATGUT 1 X 76-40MM HC RB	ME4259	1		40 mm 1/2 Circle Round Body	76	12	CHROMIC								

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)								MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)							
CHROMIC CATGUT 2-0 X 76-40MM HC RB	ME4245	2-0		40 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 4-0 X 75-19MM CU RC	ME400419	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-45MM HC RB	ME4249	2-0		45 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 4-0 X 38-12MM CU RC	ME4266	4-0		12 mm 3/8 Circle Reverse Cutting	38	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM HC CT	ME4255	2-0		30 mm 1/2 Circle Conventional Cutting	76	12	CHROMIC	CHROMIC CATGUT 4-0 X 75-26MM HC RB	ME400121	4-0		26 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-50MM ST RB	ME4208	2-0		50 mm Straight Round Body	76	12	CHROMIC	CHROMIC CATGUT 5-0 X 90-12MM CU RC	ME500268	5-0		12 mm 3/8 Circle Reverse Cutting	90	12	CHROMIC
CHROMIC CATGUT 2-0 X 75-45MM HC CT	ME200863	2-0		45 mm 1/2 Circle Conventional Cutting	75	12	CHROMIC	CHROMIC CATGUT 5-0 X 30-10MM FE TC	ME4263	5-0		10 mm 5/8 Circle Taper Cut	30	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-22MM CU CT	ME4201	3-0		22 mm 3/8 Circle Conventional Cutting	76	12	CHROMIC	CHROMIC CATGUT 5-0 X 76-13MM CU RC	ME4276	5-0		13 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-20MM HC RB	ME4237	3-0		20 mm 1/2 Circle Round Body	76	12	CHROMIC	CHROMIC CATGUT 5-0 X 76-13MM CU RB	ME4289	5-0		13 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM HC RB	ME300882	3-0		26 mm 1/2 Circle Round Body	75	12	CHROMIC	CHROMIC CATGUT 1 X 75-60MM HC RB	ME016075	1		60 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-30MM HC RB	ME300172	3-0		30 mm 1/2 Circle Round Body	75	12	CHROMIC	CHROMIC CATGUT 2 X 75-80MM HC RB	ME028075	2		80 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM CU RC	ME304326	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	CHROMIC	CHROMIC CATGUT 0 X 152	SS2214	0		Non needed	152	12	CHROMIC
CHROMIC CATGUT 3-0 X 38-16MM FE CT	ME4265	3-0		16 mm 5/8 Circle Conventional Cutting	38	12	CHROMIC	CHROMIC CATGUT 1 X 152	SS2215	1		Non needed	152	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-16MM CU CT	ME4274	3-0		16 mm 3/8 Circle Conventional Cutting	76	12	CHROMIC	CHROMIC CATGUT 2 X 152	SS2216	2		Non needed	152	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-12MM HC RC	ME4286	3-0		12 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC	CHROMIC CATGUT 2-0 X 152	SS2213	2-0		Non needed	152	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-26MM CU RC	ME4828	3-0		26 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC								
CHROMIC CATGUT 3-0 X 75-26MM HC RC	ME300886	3-0		26 mm 1/2 Circle Reverse Cutting	75	12	CHROMIC								
CHROMIC CATGUT 4-0 X 76-16MM CU RB	ME4048	4-0		16 mm 3/8 Circle Round Body	76	12	CHROMIC								
CHROMIC CATGUT 4-0 X 76-16MM CU RC	ME4287	4-0		16 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC								
CHROMIC CATGUT 4-0 X 76-17MM CU TC	ME5047	4-0		17 mm 3/8 Circle Taper Cut	76	12	CHROMIC								
CHROMIC CATGUT 4-0 X 75-20MM HC RB	ME400181	4-0		20 mm 1/2 Circle Round Body	75	12	CHROMIC								
CHROMIC CATGUT 4-0 X 75-17MM HC RB	ME400117	4-0		17 mm 1/2 Circle Round Body	75	12	CHROMIC								

FILASILK™

Silk Suture

Indications

- General surgery procedures
- Skin closure
- Oral surgery procedures
- Ophthalmic procedures
- Ligatures

Advantages

- Excellent handling properties
- Good knot security



Product Overview

Structure	: Braided
Colour	: Black
Origin	: Silkworm larval cocoons (Bombyx mori)
Composition	: Natural silk filaments
Sizes	: USP 8/0 (metric 0.4) - USP 2 (metric 5)
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE								FILASILK SILK SUTURE							
SILK BLK 0 X 90-60MM ST CT	SLK100789	0		60 mm Straight Conventional Cutting	90	12	BLACK	SILK BLK 0 X 76-30MM HC RB	SLK10 5334	0		30 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 0 X 90-50MM HC RC	SLK100552	0		50 mm 1/2 Circle Reverse Cutting	90	12	BLACK	SILK BLK 0 X 90-40MM HC RB	SLK105041	0		40 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 0 X 90-35MM CU RC	SLK100371	0		35 mm 3/8 Circle Reverse Cutting	90	12	BLACK	SILK BLK 0 X 90-40MM HC RC	SLK102224	0		40 mm 1/2 Circle Reverse Cutting	90	12	BLACK
SILK BLK 0 X 76-45MM CU RC	SLK10 5037	0		45 mm 3/8 Circle Reverse Cutting	76	12	BLACK	SILK BLK 0 X 90-75MM ST CT	SLK100791	0		75 mm Straight Conventional Cutting	90	12	BLACK
SILK BLK 0 X 90-30MM HC RB	SLK10 5332	0		30 mm 1/2 Circle Round Body	90	12	BLACK	SILK BLK 0 X 70-39MM CU RC	SLK101039	0		39 mm 3/8 Circle Reverse Cutting	70	12	BLACK
SILK BLK 0 X 75-30MM HC RC	SLK102223	0		30 mm 1/2 Circle Reverse Cutting	75	12	BLACK	SILK BLK 1 X 76-60MM CU CT	SLK01 5062	1		60 mm 3/8 Circle Conventional Cutting	76	12	BLACK
SILK BLK 0 X 45-23MM HC RC	SLK100373	0		23 mm 1/2 Circle Reverse Cutting	45	12	BLACK	SILK BLK 1 X 75-30MM CU CT	SLK012243	1		30 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 0 X 75-39MM CU CT	SLK105039	0		39 mm 3/8 Circle Conventional Cutting	75	12	BLACK	SILK BLK 1 X 75-30MM HC RB	SLK010723	1		30 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 0 X 75-30MM HC RB	SLK100334	0		30 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 1 X 90-40MM HC RB	SLK010339	1		40 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 0 X 75-24MM CU RC	SLK105024	0		24 mm 3/8 Circle Reverse Cutting	75	12	BLACK	SILK BLK 1 X 75-27MM HC RB	SLK0115175	1		27 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 0 X 75-37MM HC RB	SLK100337	0		37 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 1 X 100-90MM CU CT	SLK012793	1		90 mm 3/8 Circle Conventional Cutting	100	12	BLACK
SILK BLK 0 X 75-26MM HC RB	SLK100383	0		26 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 1 X 90-24MM CU RB	SLK015124	1		24 mm 3/8 Circle Round Body	90	12	BLACK
SILK BLK 0 X 75-45MM CU RC	SLK105076	0		45 mm 3/8 Circle Reverse Cutting	75	12	BLACK	SILK BLK 1 X 90-26MM CU CT	SLK015126	1		26 mm 3/8 Circle Conventional Cutting	90	12	BLACK
SILK BLK 0 X 90-30MM CU RC	SLK102508	0		30 mm 3/8 Circle Reverse Cutting	90	12	BLACK	SILK BLK 1 X 75-30MM HC RC	SLK015130	1		30 mm 1/2 Circle Reverse Cutting	75	12	BLACK
SILK BLK 0 X 90-30MM CU RB	SLK105030	0		30 mm 3/8 Circle Round Body	90	12	BLACK	SILK BLK 1 X 75-38MM CU RC	SLK018138	1		38 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 0 X 76-30MM CU RB	SLK10 5017	0		30 mm 3/8 Circle Round Body	76	12	BLACK	SILK BLK 1 X 75-37MM HC RB	SLK01 5137N	1		37 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 0 X 75-60MM ST CT	SLK1013113	0		60 mm Straight Conventional Cutting	75	12	BLACK	SILK BLK 1 X 90-40MM HC RC	SLK01 5005	1		40 mm 1/2 Circle Reverse Cutting	90	12	BLACK
SILK BLK 0 X 100-90MM CU CT	SLK1013114	0		90 mm 3/8 Circle Conventional Cutting	100	12	BLACK	SILK BLK 1 X 75-40MM CU CT	SLK0113116	1		40 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 0 X 75-31MM HC RB	SLK103341	0		31 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 1 X 100-60MM ST CT	SLK010391	1		60 mm Straight Conventional Cutting	100	12	BLACK
SILK BLK 0 X 75-40MM CU RC	SLK105040	0		40 mm 3/8 Circle Reverse Cutting	75	12	BLACK								

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE								FILASILK SILK SUTURE							
SILK BLK 1 X 75-30MM CU RC	SLK010689	1		30 mm 3/8 Circle Reverse Cutting	75	12	BLACK	SILK BLK 2-0 X 75-26MM HC RB	SLK200327	2-0		26 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2 X 90-90MM CU RC	SLK020794	2		90 mm 3/8 Circle Reverse Cutting	90	12	BLACK	SILK BLK 2-0 X 75-17MM HC RB	SLK205610	2-0		17 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2 X 90-40MM HC RB	SLK020936	2		40 mm 1/2 Circle Round Body	90	12	BLACK	SILK BLK 2-0 X 75-25MM HC CT	SLK205225	2-0		25 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2 X 75-30MM HC RC	SLK025230	2		30 mm 1/2 Circle Reverse Cutting	75	12	BLACK	SILK BLK 2-0 X 76-50MM ST RB	SLK20 5208	2-0		50 mm Straight Round Body	76	12	BLACK
SILK BLK 2 X 75-40MM HC RC	SLK025340	2		40 mm 1/2 Circle Reverse Cutting	75	12	BLACK	SILK BLK 2-0 X 90-15MM CU RC	SLK200503	2-0		15 mm 3/8 Circle Reverse Cutting	90	12	BLACK
SILK BLK 2 X 75-40MM HC RB	SLK020936SL	2		40 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 2-0 X 75-22MM HC CT	SLK2052221	2-0		22 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2 X 100-90MM CU CT	SLK020792	2		90 mm 3/8 Circle Conventional Cutting	100	12	BLACK	SILK BLK 2-0 X 76-60MM ST CT	SLK20 5290	2-0		60 mm Straight Conventional Cutting	76	12	BLACK
SILK BLK 2 X 70-36MM CU RC	SLK021236	2		36 mm 3/8 Circle Reverse Cutting	70	12	BLACK	SILK BLK 2-0 X 100-90MM CU CT	SLK200791	2-0		90 mm 3/8 Circle Conventional Cutting	100	12	BLACK
SILK BLK 2 X 70-39MM CU RC	SLK021239	2		39 mm 3/8 Circle Reverse Cutting	70	12	BLACK	SILK BLK 2-0 X 45-26MM CU RC	SLK200685	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2 X 75-35MM HC RB	SLK028235	2		35 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 2-0 X 75-30MM HC CT	SLK205231	2-0		30 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 1 X 50-90MM CU CT (2 suture per foil)	SLK012797	1		90 mm 3/8 Circle Conventional Cutting	50	12	BLACK	SILK BLK 2-0 X 90-26MM HC RB	SLK200327XL	2-0		26 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 2-0 X 90-20MM CU RC	SLK200367	2-0		20 mm 3/8 Circle Reverse Cutting	90	12	BLACK	SILK BLK 2-0 X 75-26MM CU RC	SLK200685XL	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 90-60MM ST CT	SLK200790	2-0		60 mm Straight Conventional Cutting	90	12	BLACK	SILK BLK 2-0 X 90-24MM CU RB	SLK205224	2-0		24 mm 3/8 Circle Round Body	90	12	BLACK
SILK BLK 2-0 X 76-45MM CU RC	SLK20 5036	2-0		45 mm 3/8 Circle Reverse Cutting	76	12	BLACK	SILK BLK 2-0 X 90-26MM CU CT	SLK205226	2-0		26 mm 3/8 Circle Conventional Cutting	90	12	BLACK
SILK BLK 2-0 X 76-50MM HC RC	SLK20 5052	2-0		50 mm 1/2 Circle Reverse Cutting	76	12	BLACK	SILK BLK 2-0 X 45-30MM HC RC	SLK200530	2-0		30 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 90-30MM HC RB	SLK20 5331	2-0		30 mm 1/2 Circle Round Body	90	12	BLACK	SILK BLK 2-0 X 75-30MM CU RC	SLK205230	2-0		30 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 76-25MM HC TC	SLK20 5670	2-0		25 mm 1/2 Circle Taper Cut	76	12	BLACK	SILK BLK 2-0 X 75-16MM CU CT	SLK200509	2-0		16 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 45-23MM HC RC	SLK200737	2-0		23 mm 1/2 Circle Reverse Cutting	45	12	BLACK	SILK BLK 2-0 X 76-25MM HC RB	SLK205671	2-0		25 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 2-0 X 75-22MM HC RB	SLK205222	2-0		22 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 2-0 X 100-55MM ST CT	SLK200775	2-0		55 mm Straight Conventional Cutting	100	12	BLACK
								SILK BLK 2-0 X 76-36MM CU RC	SLK205236	2-0		36 mm 3/8 Circle Reverse Cutting	76	12	BLACK
								SILK BLK 2-0 X 76-26MM HC RB VB	SLK20 5065	2-0		26 mm 1/2 Circle Round Body	76	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE								FILASILK SILK SUTURE							
SILK BLK 2-0 X 90-17MM HC TC	SLK20 5611	2-0		17 mm 1/2 Circle Taper Cut	90	12	BLACK	SILK BLK 2-0 X 70-36MM HC RB	SLK201336	2-0		36 mm 1/2 Circle Round Body	70	12	BLACK
SILK BLK 2-0 X 76-30MM HC RB VB	SLK20 533VB	2-0		30 mm 1/2 Circle Round Body	76	12	BLACK	SILK BLK 3-0 X 90-60MM ST CT	SLK300622	3-0		60 mm Straight Conventional Cutting	90	12	BLACK
SILK BLK 2-0 X 100-75MM ST CT	SLK2013102	2-0		75 mm Straight Conventional Cutting	100	12	BLACK	SILK BLK 3-0 X 90-25MM HC RC	SLK300379	3-0		25 mm 1/2 Circle Reverse Cutting	90	12	BLACK
SILK BLK 2-0 X 90-35MM CU RC	SLK200667	2-0		35 mm 3/8 Circle Reverse Cutting	90	12	BLACK	SILK BLK 3-0 X 90-16MM CU CT	SLK30 5003	3-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLACK
SILK BLK 2-0 X 75-36MM CU RC	SLK205275	2-0		36 mm 3/8 Circle Reverse Cutting	75	12	BLACK	SILK BLK 3-0 X 76-26MM CU RC	SLK30 5028	3-0		26 mm 3/8 Circle Reverse Cutting	76	12	BLACK
SILK BLK 2-0 X 75-35MM CU CT	SLK202429	2-0		35 mm 3/8 Circle Conventional Cutting	75	12	BLACK	SILK BLK 3-0 X 90-25MM HC RB	SLK30 5085	3-0		25 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 2-0 X 75-80MM CU CT	SLK200261	2-0		80 mm 3/8 Circle Conventional Cutting	75	12	BLACK	SILK BLK 3-0 X 76-20MM HC RB	SLK30 5087	3-0		20 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 2-0 X 75-30MM HC RB	SLK200333	2-0		30 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 3-0 X 75-24MM CU RC	SLK300684	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 75-48MM HC CT	SLK200552	2-0		48 mm 1/2 Circle Conventional Cutting	75	12	BLACK	SILK BLK 3-0 X 75-17MM HC RB	SLK3054601	3-0		17 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 75-60MM ST CT	SLK200790SL	2-0		60 mm Straight Conventional Cutting	75	12	BLACK	SILK BLK 3-0 X 75-22MM HC RB	SLK300587	3-0		22 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 45-24MM CU RC	SLK200695	2-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLACK	SILK BLK 3-0 X 75-30MM CU RC	SLK300675	3-0		30 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 100-40MM CU CT	SLK203959	2-0		40 mm 3/8 Circle Conventional Cutting	100	12	BLACK	SILK BLK 3-0 X 75-26MM HC RB	SLK300832	3-0		26 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 45-30MM HC RC C1	SLK205227	2-0		30 mm 1/2 Circle Reverse Cutting	45	12	BLACK	SILK BLK 3-0 X 75-30MM HC RB	SLK300571	3-0		30 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 75-75MM ST CT	SLK200766	2-0		75 mm Straight Conventional Cutting	75	12	BLACK	SILK BLK 3-0 X 45-21MM HC RC	SLK300632	3-0		21 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 76-30MM HC RB	SLK20 5333	2-0		30 mm 1/2 Circle Round Body	76	12	BLACK	SILK BLK 3-0 X 75-16MM CU RC	SLK300882	3-0		16 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 90-35MM HC CT	SLK20 5335	2-0		35 mm 1/2 Circle Conventional Cutting	90	12	BLACK	SILK BLK 3-0 X 75-22MM HC CT	SLK302508	3-0		22 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 75-39MM CU RC	SLK200648	2-0		39 mm 3/8 Circle Reverse Cutting	75	12	BLACK	SILK BLK 3-0 X 45-16MM CU RC	SLK300883	3-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 100-60MM ST RC	SLK200723	2-0		60 mm Straight Reverse Cutting	100	12	BLACK	SILK BLK 3-0 X 75-19MM CU RB	SLK305319	3-0		19 mm 3/8 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 75-40MM HC RB	SLK200584	2-0		40 mm 1/2 Circle Round Body	75	12	BLACK	SILK BLK 3-0 X 90-20MM HC RB	SLK300588	3-0		20 mm 1/2 Circle Round Body	90	12	BLACK
								SILK BLK 3-0 X 90-19MM CU RC	SLK302231	3-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLACK
								SILK BLK 3-0 X 75-30MM HC RC	SLK305330	3-0		30 mm 1/2 Circle Reverse Cutting	75	12	BLACK

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE						
SILK BLK 1 X 45 (45X12)	SLS014512	1	Non needed	45	12	BLACK
SILK BLK 2 X 76 (76X2)	SLS02 216	2	Non needed	76	12	BLACK
SILK BLK 2 X 180	SLS02180	2	Non needed	180	12	BLACK
SILK BLK 2 X 270 (45X6)	SLS02364	2	Non needed	45	12	BLACK
SILK BLK 2 X 250	SLS02252	2	Non needed	250	12	BLACK
SILK BLK 2 X 45 (45X12)	SLS024512	2	Non needed	45	12	BLACK
SILK BLK 3 X 76 (76X2)	SLS03 217	3	Non needed	76	12	BLACK
SILK BLK 2-0 X 76 (76X2)	SLS20 213	2-0	Non needed	76	12	BLACK
SILK BLK 2-0 X 75 (75X12)	SLS20305	2-0	Non needed	75	12	BLACK
SILK BLK 2-0 X 180	SLS20180	2-0	Non needed	180	12	BLACK
SILK BLK 2-0 X 75 (75X5)	SLS20375	2-0	Non needed	75	12	BLACK
SILK BLK 2-0 X 35 (35X10)	SLS20 223	2-0	Non needed	35	12	BLACK
SILK BLK 2-0 X 45 (45X17)	SLS20203	2-0	Non needed	45	12	BLACK
SILK BLK 2-0 X 75 (75X10)	SLS201075	2-0	Non needed	75	12	BLACK
SILK BLK 2-0 X 250-2D	SLS20253Y	2-0	Non needed	250	24	BLACK
SILK BLK 2-0 X 150 (150X2)	SLS20605	2-0	Non needed	150	12	BLACK
SILK BLK 2-0 X 250	SLS20253	2-0	Non needed	250	12	BLACK
SILK BLK 2-0 X 45 (45X12)	SLS204512	2-0	Non needed	45	12	BLACK
SILK BLK 3-0 X 76 (76X2)	SLS30 212	3-0	Non needed	76	12	BLACK
SILK BLK 3-0 X 35 (35X10)	SLS30 222	3-0	Non needed	35	12	BLACK
SILK BLK 3-0 X 75 (75X10)	SLS30750	3-0	Non needed	75	12	BLACK
SILK BLK 3-0 X 180	SLS30180	3-0	Non needed	180	12	BLACK
SILK BLK 3-0 X 45 (45X17)	SLS30202	3-0	Non needed	45	12	BLACK
SILK BLK 3-0 X 45 (45X12)	SLS304512	3-0	Non needed	45	12	BLACK
SILK BLK 4-0 X 76 (76X2)	SLS40 211	4-0	Non needed	76	12	BLACK
SILK BLK 0 X 25MTR (6 EA PER BOX)	SLK10 824R	0	Non needed	25 Mtr	6	BLACK
SILK BLK 1 X 25MTR (6 EA PER BOX)	SLK01 825R	1	Non needed	25 Mtr	6	BLACK
SILK BLK 2 X 25MTR (6 EA PER BOX)	SLK02 826R	2	Non needed	25 Mtr	6	BLACK
SILK BLK 2-0 X 25MTR (6 EA PER BOX)	SLK20 823R	2-0	Non needed	25 Mtr	6	BLACK
SILK BLK 3-0 X 25MTR (6 EA PER BOX)	SLK30 822R	3-0	Non needed	25 Mtr	6	BLACK

MERISTEELTM

Steel Suture

Indications

- Orthopedic procedures
- Sternum closure




Advantages

- Exceptional tensile strength
- Excellent tissue compatibility



Product Overview

Structure	: Monofilament
Colour	: Silver/Steel
	
Composition	: Stainless steel
Sizes	: USP 1, 2, 4, 5, 6 and 7
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISTEEL (STEEL SUTURE)								MERISTEEL (STEEL SUTURE)							
STEEL WIRE 0 X 45-17MM HC BP (45X2)	STC600	0		17 mm 1/2 Circle Blunt Point	45	12	-	STEEL WIRE 5 X 45-48MM HC TC (45X4)	STC05650	5		48 mm 1/2 Circle Taper Cut	45	12	-
STEEL WIRE 1 X 45-40MM HC TC (45X4)	STC660	1		40 mm 1/2 Circle Taper Cut	45	12	-	STEEL WIRE 6 X 45-48MM HC CT (45X4)	STC654	6		48 mm 1/2 Circle Conventional Cutting	45	12	-
STEEL WIRE 1 X 45-40MM HC BP (45X4)	STC661	1		40 mm 1/2 Circle Blunt Point	45	12	-	STEEL WIRE 6 X 45-44MM HC RB BP (45X4)	STC9456	6		44 mm 1/2 Circle Blunt Point	45	6	-
STEEL WIRE 1 X 45-35MM HC RB BLUNT	STC9435	1		35 mm 1/2 Circle Blunt Point	45	12	-	STEEL WIRE 6 X 75-44MM HC RB BP (75X2)	STC9494	6		44 mm 1/2 Circle Blunt Point	75	6	-
STEEL WIRE 2 X 45-40MM HC TC (45X4)	STC650	2		40 mm 1/2 Circle Taper Cut	45	12	-	STEEL WIRE 6 X 45-48MM HC CT (45X1)	STC9654	6		48 mm 1/2 Circle Conventional Cutting	45	12	-
STEEL WIRE 2 X 45-40MM HC BP (45X4)	STC662	2		40 mm 1/2 Circle Blunt Point	45	12	-	STEEL WIRE 6 X 45-48MM CU CT (45X6)	STC9654CU	6		48 mm 3/8 Circle Conventional Cutting	45	12	-
STEEL WIRE 2 X 45-40MM HC RB	STC9540	2		40 mm 1/2 Circle Round Body	45	12	-	SS WIRE 6X45-48 HCCT (45X1)44 HCBP(45X3)	STC700	6		48/44 mm 1/2 Circle Conventional Cutting Blunt Point	45	12	-
STEEL WIRE 4 X 45-40MM HC BP (45X4)	STC664	4		40 mm 1/2 Circle Blunt Point	45	12	-	STEEL WIRE 6 X 75-65MM HC RB BLUNT	STC9656	6		65 mm 1/2 Circle Blunt Point	75	12	-
STEEL WIRE 4 X 45-48MM HC CT (45X4)	STC652	4		48 mm 1/2 Circle Conventional Cutting	45	12	-	STEEL WIRE 7 X 75-50MM HC RC (75X2)	STC070624	7		50 mm 1/2 Circle Reverse Cutting	75	12	-
STEEL WIRE 4 X 75-48MM HC TC (75X2)	STC040651	4		48 mm 1/2 Circle Taper Cut	75	12	-	STEEL WIRE 6 X 45-48MM HC TC (45X4)	STC06649	6		48 mm 1/2 Circle Taper Cutting	45	12	-
STEEL WIRE 5 X 75-55MM HC RC(75X2)	STC050945	5		55 mm 1/2 Circle Reverse Cutting	75	12	-								
STEEL WIRE 5 X 45-48MM HC CT (45X4)	STC653	5		48 mm 1/2 Circle Conventional Cutting	45	12	-								
STEEL WIRE 5 X 45-44MM HC RB BP (45X4)	STC9455	5		44 mm 1/2 Circle Blunt Point	45	6	-								
STEEL WIRE 5 X 75-44MM HC RB BP (75X2)	STC9495	5		44 mm 1/2 Circle Blunt Point	75	6	-								
STEEL WIRE 5 X 45-48MM HC RB BP (45X4)	STC9497	5		48 mm 1/2 Circle Blunt Point	45	12	-								
STEEL WIRE 5 X 45-48MM HC CT (45X1)	STC9653	5		48 mm 1/2 Circle Conventional Cutting	45	12	-								

FILAMIDE™

Polyamide Suture

Indications

- Skin closure (intracutaneous, subcutaneous)
- Plastic and reconstructive procedures
- Microsurgery procedures
- Ophthalmic procedures

Advantages

- Flexible, easy to handle and tie
- Smooth passage through tissue
- Excellent histocompatibility




Product Overview

Structure	: Monofilament
Colour	: Black
Composition	: Polyamide 6 - 6.6
Sizes	: USP 10/0 (metric 0.2) - USP 2 (metric 5)
Sterilization	: Gamma radiation

Meril Endo-Surgery Product Catalogue | **Wound Closure** 117

Meril Endo-Surgery Product Catalogue | **Wound Closure** 119

Meril Endo-Surgery Product Catalogue | **Wound Closure** 121

Description								Description							
FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour				
FILAMIDE (POLYAMIDE SUTURE)								FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BU 3-0 X 45-19MM CU RC	NYL3013498	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYAMIDE BLK 4-0 X 60-19MM CU CT	NYL40 3357	4-0		19 mm 3/8 Circle Conventional Cutting	60	12	BLACK
POLYAMIDE BU 3-0 X 75-26MM CU RC	NYL301685	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	BLUE	POLYAMIDE BLK 4-0 X 38-10MM CU RC	NYL40 3326	4-0		10 mm 3/8 Circle Reverse Cutting	38	12	BLACK
POLYAMIDE BU 3-0 X 75-60MM ST CT	NYL306441	3-0		60 mm Straight Conventional Cutting	75	12	BLUE	POLYAMIDE BLK 4-0 X 70-19MM CU CT	NYL40 3567	4-0		19 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BU 3-0 X 45-26MM CU RC	NYL301625BU	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYAMIDE BU 4-0 X 45-11MM CU RC	NYL401612BU	4-0		11 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BU 3-0 X 45-24MM CU RC	NYL301857	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYAMIDE BU 4-0 X 45-19MM CU RC	NYL401620BU	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BLK 3-0 X 100-30MM CU RC-3D	NYL303330Z	3-0		30 mm 3/8 Circle Reverse Cutting	100	36	BLACK	POLYAMIDE BLK 4-0 X 70-16MM CU CT	NYL40 3318	4-0		16 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 3-0 X 90-22MM HC RB	NYL309022	3-0		22 mm 1/2 Circle Round Body	90	12	BLACK	POLYAMIDE BLK 4-0 X 70-10MM CU RC	NYL40 3326XL	4-0		10 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 3-0 X 45-26MM CU RC	NYL301625SL	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK	POLYAMIDE BLK 4-0 X 45-19MM CU RC-3D	NYL401620Z	4-0		19 mm 3/8 Circle Reverse Cutting	45	36	BLACK
POLYAMIDE BLK 3-0 X 75-24MM CU RC	NYL305669	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLACK	POLYAMIDE BLK 4-0 X 90-24MM CU RC	NYL403424XL	4-0		24 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 75-60MM ST CT	NYL300627	3-0		60 mm Straight Conventional Cutting	75	12	Black	POLYAMIDE BLK 4-0 X 75-18MM CU RB	NYL403418	4-0		18 mm 3/8 Circle Round Body	75	12	BLACK
POLYAMIDE BU 3-0 X 100-26MM CU CT	NYL3013500	3-0		26 mm 3/8 Circle Conventional Cutting	100	12	Blue	POLYAMIDE BU 5-0 X 90-16MM CU CT	NYL501854	5-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYAMIDE BLK 3-0 X 100-30MM CU RC	NYL303330	3-0		30 mm 3/8 Circle Reverse Cutting	100	12	Black	POLYAMIDE BLK 5-0 X 90-13MM CU RC	NYL50 3316	5-0		13 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 75-20MM HC RC	NYL303370	3-0		20 mm 1/2 Circle Reverse Cutting	75	12	Black	POLYAMIDE BLK 5-0 X 45-10MM CU RC	NYL501601	5-0		10 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 90-16MM CU CT	NYL40 3319	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLACK	POLYAMIDE BLK 5-0 X 45-19MM CU RC	NYL501618	5-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 45-11MM CU RC	NYL401612	4-0		11 mm 3/8 Circle Reverse Cutting	45	12	BLACK	POLYAMIDE BLK 5-0 X 45-13MM CU RC	NYL501852	5-0		13 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 75-19MM CU RC	NYL401619	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLACK	POLYAMIDE BLK 5-0 X 45-15MM CU CT	NYL5018951	5-0		15 mm 3/8 Circle Conventional Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 45-19MM CU RC	NYL401620	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLACK	POLYAMIDE BLK 5-0 X 75-28MM HC RB	NYL505528	5-0		28 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 4-0 X 45-16MM CU RC	NYL401861	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLACK	POLYAMIDE BLK 5-0 X 60-13MM CU RC	NYL50 3353	5-0		13 mm 3/8 Circle Reverse Cutting	60	12	BLACK
POLYAMIDE BLK 4-0 X 90-19MM CU RC	NYL401618	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLACK	POLYAMIDE BLK 5-0 X 35-16MM CU CT	NYL50 3315	5-0		16 mm 3/8 Circle Conventional Cutting	35	12	BLACK
POLYAMIDE BLK 4-0 X 75-19MM CU CT	NYL403419	4-0		19 mm 3/8 Circle Conventional Cutting	75	12	BLACK	POLYAMIDE BLK 5-0 X 70-16MM CU SB	NYL50 3726	5-0		16 mm 3/8 Circle Slim Blade	70	12	BLACK
POLYAMIDE BLK 4-0 X 75-28MM HC RB	NYL405428	4-0		28 mm 1/2 Circle Round Body	75	12	BLACK	POLYAMIDE BLK 5-0 X 70-13MM CU RC	NYL50 3752	5-0		13 mm 3/8 Circle Reverse Cutting	70	12	BLACK

122

Wound Closure | Meril Endo-Surgery Product Catalogue

Meril Endo-Surgery Product Catalogue |

Wound Closure

123

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)						
POLYAMIDE BLK 0 X 76 (76X2)	NYS10 904	0	Non needled	76	12	BLACK
POLYAMIDE BLK 1 X 76 (76X2)	NYS01 905	1	Non needled	76	12	BLACK
POLYAMIDE BLK 1 X 750 (75X10)	NYS011075	1	Non needled	75	12	BLACK
POLYAMIDE BLK 1 X 35 (35X5)	NYS01 705	1	Non needled	35	12	BLACK
POLYAMIDE BLK 1 X 100 (100X2)	NYS01185	1	Non needled	100	12	BLACK
POLYAMIDE BLK 2 X 76 (76X2)	NYS02 906	2	Non needled	76	12	BLACK
POLYAMIDE BLK 2-0 X 76 (76X2)	NYS20 903	2-0	Non needled	76	12	BLACK
POLYAMIDE BLK 2-0 X 35 (35X5)	NYS20 703	2-0	Non needled	35	12	BLACK

FILAPROP™

Polypropylene Suture



Indications


- Skin closure (intracutaneous, subcutaneous)
- Cardiovascular procedures
- Neurological procedures
- Ophthalmic procedures










Advantages





- Smooth passage through tissue
- Minimal acute inflammatory reaction in tissues
- Unsurpassed strength
- Optimal elasticity and elongation properties
- Available in Less Memory Tray (LMT)

Product Overview

Structure	: Monofilament
Colour	: Blue
Composition	: Polypropylene
Sizes	: USP 10/0 (metric 0.2) - USP 2 (metric 5)
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)								FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 0 X 90-30MM CU RC	PPL102815	0		30 mm 3/8 Circle Reverse Cutting	90	12	BLUE	POLYPROPYLENE BU 0 X 70-30MM HC RB HB	PPL10 842	0		30 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 0 X 70-30MM HC RB	PPL10 830	0		30 mm 1/2 Circle Round Body	70	12	BLUE	POLYPROPYLENE BU 0 X 100-26MM HC RB DN	PPL108026XDN	0		26 mm 1/2 Circle Round Body Round Body	100	12	BLUE
POLYPROPYLENE BU 0 X 90-30MM HC RB	PPL10 846	0		30 mm 1/2 Circle Round Body	90	12	BLUE	POLYPROPYLENE BU 0 X 75-37MM HC TC	PPL108037	0		37 mm 1/2 Circle Taper Cut	75	12	BLUE
POLYPROPYLENE BU 0 X 90-30MM HC TC	PPL10 894	0		30 mm 1/2 Circle Taper Cut	90	12	BLUE	POLYPROPYLENE BU 0 X 75-38MM CU CT	PPL108038	0		38 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 0 X 75-40MM HC RB	PPL108434	0		40 mm 1/2 Circle Round Body	75	12	BLUE	POLYPROPYLENE BU 1 X 90-45MM HC TC	PPL018450	1		45 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 0 X 75-37MM HC RB	PPL108424	0		37 mm 1/2 Circle Round Body	75	12	BLUE	POLYPROPYLENE BU 1 X 90-50MM HC RC	PPL018475	1		50 mm 1/2 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 0 X 75-30MM HC RB	PPL100830	0		30 mm 1/2 Circle Round Body	75	12	BLUE	POLYPROPYLENE BU 1 X 150-48MM HC RB LOOP	PPL01 834	1		48 mm 1/2 Circle Round Body	150	12	BLUE
POLYPROPYLENE BU 0 X 75-24MM CU CT	PPL108024	0		24 mm 3/8 Circle Conventional Cutting	75	12	BLUE	POLYPROPYLENE BU 1 X 70-30MM HC RB	PPL01 840	1		30 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 0 X 75-26MM HC RB	PPL108026	0		26 mm 1/2 Circle Round Body	75	12	BLUE	POLYPROPYLENE BU 1 X 70-40MM HC RB	PPL01 843	1		40 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 0 X 75-45MM HC RC	PPL108045	0		45 mm 1/2 Circle Reverse Cutting	75	12	BLUE	POLYPROPYLENE BU 1 X 100-45MM HC RC	PPL01 883	1		45 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 0 X 75-45MM HC RB	PPL108075	0		45 mm 1/2 Circle Round Body	75	12	BLUE	POLYPROPYLENE BU 1 X 90-90MM CU CT	PPL012989	1		90 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 0 X 90-26MM HC RB DN	PPL10 824MIPH	0		26 mm 1/2 Circle Round Body Round Body	90	1	BLUE	POLYPROPYLENE BU 1 X 100-40MM HC RC	PPL01884	1		40 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 0 X 100-40MM HC RB	PPL108434XL	0		40 mm 1/2 Circle Round Body	100	12	BLUE	POLYPROPYLENE BU 1 X 100-36MM HC RB	PPL01842	1		36 mm 1/2 Circle Round Body	100	12	BLUE
POLYPROPYLENE BU 0 X 100-31MM J TYPE TC	PPL10982J	0		31 mm J Shape Taper Cut	100	12	BLUE	POLYPROPYLENE BU 1 X 90-90MM CU RC	PPL018190	1		90 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 0 X 100-31MM HC RB	PPL108430	0		31 mm 1/2 Circle Round Body	100	12	BLUE	POLYPROPYLENE BU 1 X 75-30MM HCRB TRODN	PPL018130	1		30 mm 1/2 Circle Round Body TRO	75	12	BLUE
POLYPROPYLENE BU 0 X 90-40MM HC RC	PPL108470	0		40 mm 1/2 Circle Reverse Cutting	90	12	BLUE	POLYPROPYLENE BU 1 X 100-40MM HC RB	PPL01 843DS	1		40 mm 1/2 Circle Round Body	100	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)								FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 2-0 X 90-31MM HC RB DN	PPL208526	2-0		31 mm 1/2 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 2-0 X 75-70MM ST RB DN	PPL208400	2-0		70 mm Straight Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-31MM HC TC DN	PPL208850	2-0		31 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE	POLYPROPYLENE BU 2-0 X 90-30MM HC RB DN	PPL20844DN	2-0		30 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 45-26MM CU RC	PPL208026	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYPROPYLENE BU 2-0 X 90-17MM HC TC	PPL208938	2-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC TC DN	PPL208852DN	2-0		40 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE	POLYPROPYLENE BU 2-0 X 100-40MM CU CT	PPL208690	2-0		40 mm 3/8 Circle Conventional Cutting	100	12	BLUE
POLYPROPYLENE BU 2-0 X 100-75MM ST CT	PPL200752	2-0		75 mm Straight Conventional Cutting	100	12	BLUE	POLYPROPYLENE BU 3-0 X 70-22MM CU CT	PPL30 018	3-0		22 mm 3/8 Circle Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 2-0 X 90-55MM ST CT	PPL208571	2-0		55 mm Straight Conventional Cutting	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-25MM CU CT	PPL30 800	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC RB DN	PPL20 2040DN	2-0		40 mm 1/2 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-25MM HC RB	PPL30 825	3-0		25 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 70-30MM HC RB	PPL20 841	2-0		30 mm 1/2 Circle Round Body	70	12	BLUE	POLYPROPYLENE BU 3-0 X 70-25MM HC RB	PPL30 838	3-0		25 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 2-0 X 90-90MM CU CT	PPL200988	2-0		90 mm 3/8 Circle Conventional Cutting	90	12	BLUE	POLYPROPYLENE BU 3-0 X 70-20MM HC RB	PPL30 887	3-0		20 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RB-3D	PPL20844Z	2-0		30 mm 1/2 Circle Round Body	90	36	BLUE	POLYPROPYLENE BU 3-0 X 70-60MM ST CT	PPL30 888	3-0		60 mm Straight Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HCTCDN-3D	PPL208977Z	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	36	BLUE	POLYPROPYLENE BU 3-0 X 90-26MM HC RB DN	PPL30 8522	3-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X90-31MM HCRBVB-3D	PPL20295Z	2-0		31 mm 1/2 Circle Round Body	90	36	BLUE	POLYPROPYLENE BU 3-0 X 90-17.5MM HCTC DN	PPL30 8936	3-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC TC	PPL208926	2-0		26 mm 1/2 Circle Taper Cut	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-17MM HC RB	PPL308558	3-0		17 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RC	PPL208931	2-0		30 mm 1/2 Circle Reverse Cutting	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-16MM HC TC DN	PPL308667	3-0		16 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
								POLYPROPYLENE BU 3-0 X 75-22MM HC RB	PPL308770	3-0		22 mm 1/2 Circle Round Body	75	12	BLUE
								POLYPROPYLENE BU 3-0 X 90-26MM HC TC DN	PPL308976	3-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
								POLYPROPYLENE BU 3-0 X 45-19MM CU RC	PPL308665	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)								FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 3-0 X 45-24MM CU RC	PPL308684	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYPROPYLENE BU 3-0 X 90-31MM HC RB DN	PPL308525DN	3-0		31 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 75-26MM HC RB	PPL308832	3-0		26 mm 1/2 Circle Round Body	75	12	BLUE								
POLYPROPYLENE BU 3-0 X 75-19MM CU RC	PPL308687	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLUE	POLYPROPYLENE BU 3-0 X 45-22MM HC CT	PPL308731	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM HC RB	PPL308559	3-0		16 mm 1/2 Circle Round Body	90	12	BLUE	POLYPROPYLENE BU 3-0 X 45-19MM CU CT	PPL308878	3-0		19 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM HC TC	PPL308977	3-0		25 mm 1/2 Circle Taper Cut	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-26MM HC RB DNV	PPL308355	3-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-30MM HC RB	PPL308525	3-0		30 mm 1/2 Circle Round Body	90	12	BLUE								
POLYPROPYLENE BU 3-0 X 90-16MM CU RC	PPL308549	3-0		16 mm 3/8 Circle Reverse Cutting	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-31MM HC TC DN	PPL308849DN	3-0		31 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM CU RC	PPL308685	3-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLUE								
POLYPROPYLENE BU 3-0 X 90-17MM HC TC	PPL308317	3-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-40MM HC TC DN	PPL308851	3-0		40 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-36MM HC RB	PPL308336	3-0		36 mm 1/2 Circle Round Body	90	12	BLUE								
POLYPROPYLENE BU 3-0 X 75-22MM HC RB DN	PPL308771	3-0		22 mm 1/2 Circle Round Body Round Body	75	12	BLUE	POLYPROPYLENE BU 3-0 X90-26MM HCRB TRODN	PPL308833	3-0		26 mm 1/2 Circle Round Body TRO	90	12	BLUE
POLYPROPYLENE BU 3-0 X 75-60MM ST CT	PPL3013413	3-0		60 mm Straight Conventional Cutting	75	12	BLUE								
POLYPROPYLENE BU 3-0 X 75-16MM CU CT	PPL3013453	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	BLUE	POLYPROPYLENE BU 3-0 X 90-20MM HC RBTC DN	PPL309557	3-0		20 mm 1/2 Circle Round Body Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 120-26MM HC RB	PPL308534	3-0		26 mm 1/2 Circle Round Body	120	12	BLUE								
POLYPROPYLENE BU 3-0 X 75-24MM CU RC	PPL307694	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLUE	POLYPROPYLENE BU 3-0 X 90-17MM HC RB DN	PPL308558DN	3-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-36MM HC RB DN	PPL307499	3-0		36 mm 1/2 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 3-0 X 90-16MM CU RC-3D	PPL308549Z	3-0		16 mm 3/8 Circle Reverse Cutting	90	36	BLUE
POLYPROPYLENE BU 3-0 X 120-26MM HC TC DN	PPL307778	3-0		26 mm 1/2 Circle Taper Cut Taper Cut	120	12	BLUE	POLYPROPYLENE BU 3-0 X 90-60MM ST CT	PPL300621	3-0		60 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-26MM CU RC	PPL308025	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYPROPYLENE BU 3-0 X 45-26MM CU CT EX	PPL308021	3-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
								POLYPROPYLENE BU 3-0 X 75-36MM CU CT	PPL300626	3-0		36 mm 3/8 Circle Conventional Cutting	75	12	BLUE
								POLYPROPYLENE BU 4-0 X 90-19MM CU RC	PPL408015	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLUE
								POLYPROPYLENE BU 4-0 X 90-16MM CU CT	PPL408875	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLUE



134

Wound Closure | Meril Endo-Surgery Product Catalogue

Meril Endo-Surgery Product Catalogue |

Wound Closure

135


Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)								FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 4-0 X 90-25MM CU CT	PPL408879	4-0		25 mm 3/8 Circle Conventional Cutting	90	12	BLUE	POLYPROPYLENE BU 4-0 X 90-24MM CU RC	PPL408684	4-0		24 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 70-16MM HC RB	PPL40 849	4-0		16 mm 1/2 Circle Round Body	70	12	BLUE	POLYPROPYLENE BU 4-0 X 70-19MM CU RC	PPL408682XL	4-0		19 mm 3/8 Circle Reverse Cutting	70	12	BLUE
POLYPROPYLENE BU 4-0 X 70-16MM CU CT	PPL40 870	4-0		16 mm 3/8 Circle Conventional Cutting	70	12	BLUE	POLYPROPYLENE BU 4-0 X90-25MM HC TC	PPL408425	4-0		25 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X90-25MM HC TC DNC1	PPL40 867	4-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE	POLYPROPYLENE BU 4-0 X 90-16MM CU RC	PPL408416	4-0		16 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17.5MM HCTC DN	PPL40 8935	4-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE	POLYPROPYLENE BU 4-0 X 75-22MM HC RB DN	PPL40 8526DN	4-0		22 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB DN	PPL40 8521	4-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 70-17MM HC TC	PPL40 835	4-0		17 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 4-0 X 90-16MM HC RB DN	PPL40 8557	4-0		16 mm 1/2 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 45-16MM CU RC	PPL408884	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB	PPL408355	4-0		26 mm 1/2 Circle Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 45-16MM CU CT	PPL408007	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC RB	PPL408558	4-0		17 mm 1/2 Circle Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 45-26MM CU CT	PPL408020SL	4-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 75-13MM HC CT DN	PPL408666	4-0		13 mm 1/2 Circle Conventional Cutting Conventional Cutting	75	12	BLUE	POLYPROPYLENE BU 4-0 X90-26MM HC TC DNC2	PPL40 8845	4-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC RB	PPL408761	4-0		20 mm 1/2 Circle Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 45-19MM CU CT	PPL408011	4-0		19 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC CT DN	PPL408840	4-0		20 mm 1/2 Circle Conventional Cutting Conventional Cutting	90	12	BLUE	POLYPROPYLENE BU 4-0 X 90-17MM HC TC	PPL40 857	4-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC CT DN	PPL408846	4-0		26 mm 1/2 Circle Conventional Cutting Conventional Cutting	90	12	BLUE	POLYPROPYLENE BU 4-0 X 45-12MM CU RC	PPL40 879	4-0		12 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-19MM CU RC	PPL408682	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYPROPYLENE BU 4-0 X 90-16MM HC TC DN	PPL408417	4-0		16 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 45-20MM CU RC	PPL408683	4-0		20 mm 3/8 Circle Reverse Cutting	45	12	BLUE	POLYPROPYLENE BU 4-0 X 90-17MM HC TC EX	PPL400835	4-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-13MM HC RB DN	PPL408204	4-0		13 mm 1/2 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 70-17MM HC TC C1	PPL40 899	4-0		17 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 4-0 X 90-18MM CU RB DN	PPL408861	4-0		18 mm 3/8 Circle Round Body Round Body	90	12	BLUE	POLYPROPYLENE BU 4-0 X 90-26MM HC RB EX	PPL408831	4-0		26 mm 1/2 Circle Round Body	90	12	BLUE
								POLYPROPYLENE BU 4-0 X 90-20MM HC TC DN	PPL408841DN	4-0		20 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
								POLYPROPYLENE BU 4-0 X 70-16MM CU CT-3D	PPL40870Z	4-0		16 mm 3/8 Circle Conventional Cutting	70	36	BLUE

Meril Endo-Surgery Product Catalogue | **Wound Closure** 139

Meril Endo-Surgery Product Catalogue | **Wound Closure** 141

Meril Endo-Surgery Product Catalogue | **Wound Closure** 143

Meril Endo-Surgery Product Catalogue | **Wound Closure** 145

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)								FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 7-0 X 75-9.3MM CU TC DN	PPL707237	7-0		9.3 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE	POLYPROPYLENE BU 8-0 X 75-6.5MM CU RB DN	PPL808732XL	8-0		6.5 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9.5MM CU RB DN	PPL70 8704RB	7-0		9.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 8-0 X 60-6.5MM CURBDNEX	PPL805662	8-0		6.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU TC DNVB	PPL707237BL	7-0		9 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE	POLYPROPYLENE BU 8-0 X 75-6.5MM CURB DNV	PPL808732BXL	8-0		6.5 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 75-6MM CU RB DN	PPL702076XL	7-0		6 mm 3/8 Circle Round Body Round Body	75	12	BLUE	POLYPROPYLENE BU 9-0 X 13-5.1MM CU RB DN	PPL902783	9-0		5 mm 3/8 Circle Round Body Round Body	13	12	BLUE
POLYPROPYLENE BU 2 X 8-0X60-8MM CU RB DN	PPL80 876ML	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 9-0 X 90-5MM CU RB DN	PPL902780	9-0		5 mm 3/8 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6MM CU RB DNTN	PPL80 812	8-0		6 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 0 X 100-40MM HC RC	PPL108470XL	0		40 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 8-0 X 60-9MM CU RB DN	PPL80 813	8-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 1 X 100-48MM HC RC	PPL018576	1		48 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU RB DNTN	PPL80 876	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 4-0 X90-17MM HC RB DNTP	PPG40 846	4-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6.5MM CU RB DN	PPL80 8732	8-0		6.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 4-0X90-17.5MM HC TC DTP3	PPG40 8935	4-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU RB DNTH	PPL80 8741	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 5-0 X 75-13MM HC RB DN TP	PPG50 810	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU TC VBDN	PPL80 8101BL	8-0		8 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE	POLYPROPYLENE BU 5-0 X60-10MM HC RB DN TP3	PPG50 8803	5-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU RB DN	PPL80 8088	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 5-0 X 75-13MM HC RB DN TP3	PPG50 8710	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6MM CU RB DNTH	PPL80 8086XF	8-0		6 mm 3/8 Circle Round Body Round Body	60	12	BLUE	POLYPROPYLENE BU 6-0 X60-10MM HC RB DN TP3	PPG60 8207	6-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 13-06MM CU RB	PPL802775	8-0		6 mm 3/8 Circle Round Body	13	12	BLUE	POLYPROPYLENE BU 6-0X60-10MM HC RB DN TP3	PPG60 8260SR	6-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 75-8MM CU RB DNTH	PPL808741XL	8-0		8 mm 3/8 Circle Round Body Round Body	75	12	BLUE	POLYPROPYLENE BU 7-0X60-09MM HC RB DN TP3	PPG70 8270SR	7-0		9 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 45-6MM CU RB DN	PPL802777	8-0		6 mm 3/8 Circle Round Body Round Body	45	12	BLUE	POLYPROPYLENE BU 7-0X75-8MM CU RB DN (2S)	PPL708733DN	7-0		8 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 8-0 X 60-9MM CU TC DNV	PPL808703	8-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE								

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 3-0 X 90-22MM HC RB DN	PPL308771XL	3-0		22 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 75-17MM HC RB	PPL408558SL	4-0		17 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 60-9MM HC RB	PPL608029	6-0		9 mm 1/2 Circle Round Body	60	12	BLUE

MERICRON XLTM

Polyester Suture




Indications



- Orthopedics procedures
 - Pediatric sternum closure
- Cardiovascular procedures
 - Ophthalmic procedures
 - Neurological procedures



Advantages

- Smoother suture surface
- Less resistance to knot sliding
- Better strength
- Pliable structure

Product Overview

Structure	: Braided
Colour	: Green, White
	
Composition	: Polyethylene terephthalate (polyester)
Coating	: Silicon
Sizes	: USP 6/0 (metric 0.7) - USP 5 (metric 7)
Sterilization	: Gamma Radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)								MERICRON XL (POLYESTER SUTURE)							
POLYESTER GR 0 X 90-50MM HC TC	Me100979	0		50 mm 1/2 Circle Taper Cut	90	12	GREEN	POLYESTER GREEN 1 X 90-36MM HC TC	ME014407	1		36 mm 1/2 Circle Taper Cut	90	12	GREEN
POLYESTER GR 0 X 90-26MM HC TC DN	ME6978	0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN	POLYESTER GR 2 X 100-45MM HC TC	ME643	2		45 mm 1/2 Circle Taper Cut	100	12	GREEN
POLYESTER GR 0 X 75-22MM HC RC	ME517G	0		22 mm 1/2 Circle Reverse Cutting	75	12	GREEN	POLYESTER GR 2 X 75-45MM HC TC (4 in 1)	ME4843	2		45 mm 1/2 Circle Taper Cut	75	6	GREEN
POLYESTER GR 0 X 90-36MM HC RB	ME844G	0		36 mm 1/2 Circle Round Body	90	12	GREEN	POLYESTER GR 2 X 75-45MM HC TC	ME4842	2		45 mm 1/2 Circle Taper Cut	75	12	GREEN
POLYESTER GR 0 X 75-30MM HC RB	ME106030	0		30 mm 1/2 Circle Round Body	75	12	GREEN	POLYESTER GR 2 X 90-37MM HC RB	ME026237	2		37 mm 1/2 Circle Round Body	90	12	GREEN
POLYESTER GR 0 X 75-26MM HC RB	ME10361	0		26 mm 1/2 Circle Round Body	75	12	GREEN	POLYESTER GR 5 X 75-55MM HC TC	ME646	5		55 mm 1/2 Circle Taper Cut	75	12	GREEN
POLYESTER GR 0 X 75-31MM HC RB	ME100975	0		31 mm 1/2 Circle Round Body	75	12	GREEN	POLYESTER GR 5 X 75-55MM HC TC (4 in 1)	ME4846	5		55 mm 1/2 Circle Taper Cut	75	6	GREEN
POLYESTER GR 0 X 90-36MM HC RC	ME1032049	0		36 mm 1/2 Circle Reverse Cutting	90	12	GREEN	POLYESTER GR 2-0 X 90-17MM HC TC DN C2	ME617	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GR 0 X 90-30MM CU RC	ME106680	0		30 mm 3/8 Circle Reverse Cutting	90	12	GREEN	POLYESTER WH 2-0 X 76-17MM HC TC	ME673W	2-0		17 mm 1/2 Circle Taper Cut	76	12	WHITE
POLYESTER GR 0 X 75-36MM CU RC	ME107562	0		36 mm 3/8 Circle Reverse Cutting	75	12	GREEN	POLYESTER WH 2-0 X 90-25MM HC TC DN	ME687W	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GREEN 1 X 90-30MM HC RC	ME0132047	1		30 mm 1/2 Circle Reverse Cutting	90	12	GREEN	POLYESTER GR 2-0 X 90-17MM HC TC DN C1	ME6937	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GREEN 1 X 75-22MM HC RC	ME010518	1		22 mm 1/2 Circle Reverse Cutting	75	12	GREEN	POLYESTER GR 2-0 X 90-26MM HC TC DN	ME6977	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GREEN X 90-37MM HC RB	ME016137	1		37 mm 1/2 Circle Round Body	90	12	GREEN	POLYESTER WH 2-0 X 90-26MM HC TC DN	ME6987	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GREEN 1 X 90-30MM HC RB	ME016130	1		30 mm 1/2 Circle Round Body	90	12	GREEN	POLYESTER WH 2-0 X 90-17MM HC TC DN	ME617W	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GREEN 1 X 75-31MM J TYPE RB H	ME010969J	1		31 mm J Shape Round Body	75	12	GREEN								
POLYESTER GREEN 1 X 90-30MM HC RB-3D	ME016130Z	1		30 mm 1/2 Circle Round Body	90	36	GREEN								

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour	Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)								MERICRON XL (POLYESTER SUTURE)							
POLYESTER GR 4-0 X 76-8MM QC SM	ME682	4-0		8 mm 1/4 Circle Spatula Micropoint	76	12	GREEN	POLYESTER GREEN 3-0 X 75-17MM HC TC-3D	ME306935Z	3-0		17 mm 3/8 Circle Reverse Cutting	75	36	GREEN
POLYESTER GREEN 4-0 X 75-26MM HC RB	ME406831	4-0		26 mm 1/2 Circle Round Body	75	12	GREEN	POLYESTER GREEN 4-0 X 45-13MM CU RC-3D	ME407573Z	4-0		13 mm 3/8 Circle Reverse Cutting	45	36	GREEN
POLYESTER GR 4-0 X 90-17MM HC TC DN	ME4013395	4-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN	POLYESTER GREEN 4-0 X 45-19MM CU RC-3D	ME407142Z	4-0		19 mm 3/8 Circle Reverse Cutting	45	36	GREEN
POLYESTER GREEN 4-0 X 75-20MM HC RB	ME406761	4-0		20 mm 1/2 Circle Round Body	75	12	GREEN	POLYESTER WHGR 2-0X75-25MM HC TC DN TP6	MME6577X	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GREEN 4-0 X 75-13MM CU RB DN	ME6891G	4-0		13 mm 3/8 Circle Round Body Round Body	75	12	GREEN	POLYESTER WHGR 2-0X75-17MM HC TC DN TP6	MME6555X	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GREEN 4-0 X 90-17.5MM HC RC DN	ME6935G	4-0		17.5 mm 1/2 Circle Reverse Cutting Reverse Cutting	90	12	GREEN	POLYESTER WHGR 2-0X75-25MM HC TC DN TP5	MME6577	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GR 4-0 X 75-20MM HC RB DN	ME406761DN	4-0		20 mm 1/2 Circle Round Body Round Body	75	12	GREEN	POLYESTER WHGR 2-0X75-17MM HC TC DN TP5	MME6555	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GR 5-0 X 45-8MM QC SM DN	ME683	5-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	GREEN	POLYESTER GR 2-0 X 90-17MM HC TC DN TP	PME617	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GREEN 5-0 X 60-17MM HC RB DN	ME506757	5-0		17 mm 1/2 Circle Round Body Round Body	60	12	GREEN	POLYESTER WH 2-0 X 90-17MM HC TC DN TP	PME617W	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GREEN 5-0 X 45-8MM HC SM DN	ME500883DN	5-0		8 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	GREEN	POLYESTER GR 2-0 X 90-25MM HC TC DN TP	PME687	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GREEN 5-0 X 75-13MM HC RB DN	ME6889	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	GREEN	POLYESTER WH 2-0 X 90-25MM HC TC DN TP	PME687W	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GREEN 5-0 X 75-13MM CU RB DN	ME6890G	5-0		13 mm 3/8 Circle Round Body Round Body	75	12	GREEN	POLYSTER XL GR 2-0 X 75-26MM HC TC DN TP	ME202613	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	75	12	GREEN
POLYESTER WH 5-0 X 45-8MM QC SM DN	ME50893DN	5-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	WHITE	POLYSTER XL GR 2-0 X 75-17MM HC TC DN TP	ME200553	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	12	GREEN
POLYESTER WH 5-0 X 45-8MM HC SM DN	ME50883DN	5-0		8 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	WHITE	POLYESTER XL WHGR 2-0X75-17MM HCTCDN TP3	MME6555MP	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GR 6-0 X 75-11MM CU RB DN	ME6597	6-0		11 mm 3/8 Circle Round Body Round Body	75	12	GREEN								
POLYESTER GREEN 1 X 90-36MM HC RC - 3D	ME012520Z	1		36 mm 1/2 Circle Reverse Cutting	90	36	GREEN								

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER XL WHGR 2-0X75-25MM HCTCDN TP3	MME6577MP	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER XL GR 2-0 X 90-17MM HC TCDN P3	PME617MP	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER XL WH 2-0 X 90-17MM HC TCDN P3	PME617WMP	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GR 2-0 X 90-17MM HC TC DN TP6	PME6937	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GR 2-0 X 90-26MM HC TC DN TP6	PME6977G	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 90-26MM HC TC DN TP6	PME6987W	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER XL WH/GR 2-0X90-25MM HC TC DNTP5	SM10B77	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	6	GREEN/ WHITE
POLYESTER GR 5-0 X 75-13MM HC RB DN TP3	PME6889G	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	GREEN
POLYESTER GR 2 X 75-40MM HC TC	ME02869	2		40 mm 1/2 Circle Taper Cut	75	12	GREEN

Note: _____

158 **Wound Closure** | Meril Endo-Surgery Product Catalogue

Note: _____

Meril Endo-Surgery Product Catalogue | **Wound Closure** 159